

Abbott LibreLink App User Guide

[Home](#) » [Abbott](#) » Abbott LibreLink App User Guide 

Contents

- 1 [Abbott LibreLink App](#)
- 2 [Downloading the FreeStyle LibreLink app1](#)
- 3 [How to start the FreeStyle Libre 2 sensor](#)
- 4 [Documents / Resources](#)
 - 4.1 [References](#)



Abbott LibreLink App



Downloading the FreeStyle LibreLink app¹

1. Go to App Store or Play Store
2. Search FreeStyle LibreLink
3. Look for the FreeStyle LibreLink – ZA app
4. Press Get or Install (the app is free to download and use)²
5. Enter your ID information if requested
6. Open the app and follow the on screen instructions to get set up

How to start the FreeStyle Libre 2 sensor


1. Open the FreeStyle LibreLink app¹ on your smartphone.
iPhone: Tap the scan button.
2. iPhone: Hold the top of the iPhone near the sensor.
Android: Hold the back of the Android smartphone near the sensor and position the Near Field Communication (NFC) antenna over the sensor.
3. iPhone: Do not move your iPhone until you hear a tone and/or feel a vibration.
Android: Continue holding your smartphone near the sensor until you hear a second tone and/or feel a vibration.
This completes the scan.



- 1** Open the FreeStyle LibreLink app¹ on your smartphone.
iPhone: Tap the scan button.



- 2** iPhone: Hold the top of the iPhone near the sensor.
Android: Hold the back of the Android smartphone near the sensor and position the Near Field Communication (NFC) antenna over the sensor.



- 3** iPhone: Do not move your iPhone until you hear a tone and/or feel a vibration.
Android: Continue holding your smartphone near the sensor until you hear a second tone and/or feel a vibration.
This completes the scan.

Once you have set up your FreeStyle Libre 2 sensor your glucose readings are automatically³ updated in real-time and sent directly to the FreeStyle LibreLink app¹

QUICK TIPS

- Scan your sensor with your smartphone¹ for glucose readings during a signal loss
- Remove bulky or metallic phone cases as they can interfere with the NFC signal
- A 60-minute warm-up required when applying the sensor

Images are for illustrative purposes only. Not actual patient or data.

1. The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check the website for more information about device compatibility before using the app. Use of FreeStyle LibreLink may require registration with LibreView.
2. The mobile app is free to download and use. Internet connection is required to download the app, charges may apply.
3. Glucose readings are automatically displayed in the FreeStyle LibreLink app only when your smartphone and sensor are connected and in range.

It's easy to set alarms on your FreeStyle LibreLink app¹

1 Tap Alarms in the menu

2 Touch Low Glucose Alarm and turn on alarm (alarms are off by default)

Signal Loss Alarm is automatically turned on when glucose alarm is turned ON

3 Scroll to select a Low Glucose Value²

4 Touch Alarm Tone and make your tone choice

5 Set your High Glucose Alarm using the same steps³

Visit www.FreeStyleLibre.za.com for more information


Images are for illustrative purposes only. Not actual patient data.

The FreeStyle LibreLink app is only compatible with certain mobile devices and operating systems. Please check the website for more information about device compatibility before using the app. Use of FreeStyle LibreLink may require registration with LibreView.

1. To receive alarms on your phone, you should start your FreeStyle Libre 2 sensor with your compatible phone. Notifications will only be received when alarms are turned on and the sensor is within 6 meters of the reading device. You must have override do not disturb settings enabled to receive alarms and alerts on your smartphone.
2. 3.9 mmol/L is the default Low Glucose Alarm level and can be set between 3.3-5.6 mmol/L.
3. 13.3 mmol/L is the default High Glucose Alarm level and can be set between 6.7-22.2 mmol/L.

© 2024 Abbott. The sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott. Abbott Laboratories S.A. (Pty) Ltd., Abbott Place, 219 Golf Club Terrace, Constantia Kloof, 1709. Tel: (011) 858-2000. Fax: (011) 858 2137. ADC-88751 V1 Feb 2024.

Documents / Resources



[Abbott LibreLink App](#) [pdf] User Guide

LibreLink App, App

References

- [S life.to - This website is for sale! - life Resources and Information.](#)
- [User Manual](#)

Manuals+ Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.