



AbleNet 10000045 iTalk2 Speech Device with Levels Dual-Message Communicator User Guide

[Home](#) » [AbleNet](#) » AbleNet 10000045 iTalk2 Speech Device with Levels Dual-Message Communicator User Guide 

AbleNet[®]

10000045 iTalk2 Speech Device with Levels Dual-Message Communicator User Guide

iTalk2™

speech device

WARRANTY

AbleNet manufactured products include a 2-year limited warranty. This warranty is against defects in materials and manufacturing for 2-years from the date of purchase. Full warranty details are available at www.ablenetinc.com.

AbleNet, Inc.

2625 Patton Road Roseville,

MN 55113

United States of America

(651) 294-3101

ablecare@ablenetinc.com

www.ablenetinc.com



Made with recycled content

Contents

- [1 Product Registration](#)
- [2 Getting Started](#)
- [3 Device Overview](#)
- [4 Documents / Resources](#)
 - [4.1 References](#)
- [5 Related Posts](#)

Product Registration

Registering your product provides you access to AppleCare, product updates, and resources for your product. Scan the QR code below to register your product.



<https://www.ablenetinc.com/product-registration/>

Getting Started

Scan the QR code below to watch a short getting started video or follow the instructions listed.



<https://ablenetinc.zendesk.com/hc/en-us/articles/360060500011>

To get started:

1. Add one alkaline 9-volt battery to the device (battery not included).
2. Turn the device on.
3. Select the desired level.
4. Press and hold the record button until the recording indicator light turns on.

5. Press and hold one colored switch top. Begin speaking into the microphone after you hear a beep.
6. When done recording message, release colored switch top.
7. Repeat steps 5 and 6 for another colored switch top.
8. Press the record button to turn the recording indicator light off.
9. Optional: Select a different recording level to record a second and third set of messages.
10. Adjust volume as needed.
11. Your device is ready to use.

This device can do more! Full Instructions for Use are available at www.ablenetinc.com.



ableCARE™

Download the AppleCare app to access live support from the AppleCare Product Success team, an online knowledge base filled with videos and getting started information, and other resources. Scan the QR code below to download the free AppleCare app to your phone or tablet.




Download on the
App Store



<https://apps.apple.com/us/app/ablecare/id1564779986?ign-mpt=uo%3D2>

<https://play.google.co>

Device Overview



Create printable picture symbol overlays for your device with the FREE AbleNet Symbol Overlay Maker app. Scan the QR code to download today!





Documents / Resources

	AbleNet 10000045 iTalk2 Speech Device with Levels Dual-Message Communicator [pdf] User Guide 10000045, iTalk2 Speech Device with Levels Dual-Message Communicator
---	--

References

- [A Home - AbleNet](#)
- [A Home - AbleNet](#)