

Amico APC-IFU-P02-01EN Push Button Lock Instructions


Home » Amico » Amico APC-IFU-P02-01EN Push Button Lock Instructions

Contents

- 1 Amico APC-IFU-P02-01EN Push Button Lock
- **2 Product Information**
- **3 Product Usage Instructions**
- **4 Frequently Asked Questions**
- **5 Required Tools**
- **6 Instructions to Setup/Change Combination**
- 7 Documents / Resources
 - 7.1 References



Amico APC-IFU-P02-01EN Push Button Lock



Product Information

Specifications

• Product: Push Button Lock

• Manufacturer: Amico Patient Care Corporation

• Required Tool: Phillips screwdriver Phillips #2 (PH-2)

Product Usage Instructions

Setup/Change Combination

- Step 1: Fully extend the drawer with the push button lock installed. Keep the drawer open throughout the process.
- 2. Step 2: Locate the lock cover behind the drawer front with two holes.
- 3. **Step 3:** Insert a Phillips screwdriver in each hole and loosen the screws by turning them counter-clockwise for 1-1.5 turns. Remove the cover by lifting it.
- 4. **Step 4:** To change the combination, turn the knob left (counter-clockwise) until it stops and release it. Do not force the knob.
- 5. **Step 5:** Follow the specific combination instructions provided using the reference numbers assigned to each button in the image.
- 6. Step 6: Push the combination change slide on the back side of the lock inward using a slotted screwdriver.
- 7. Step 7: Turn the knob left (counter-clockwise) to clear the existing combination.
- 8. **Step 8:** Select a new combination, depress each button fully, and release it. Ensure a distinctive click is felt for each button press.
- 9. Step 9: Turn the knob right (clockwise) to activate the new combination.

Frequently Asked Questions


Q: What should I do if I enter the wrong combination?

A: Turn the knob left (counter-clockwise) to the stop position and release to clear any previously depressed buttons. Enter the correct combination afterward.


Q: Does clearing the existing combination reset it to the factory default combination?

A: No, clearing the existing combination does not reset it to the factory default combination.

Required Tools



Phillips screwdriver Phillips #2 (PH-2)




Slotted screwdriver Tip Size – 13/64"

Instructions to Setup/Change Combination

| Procedure | Images | Remarks |
|---|--------|---|
| Step 1 Fully extend the drawer with the pu sh button lock installed. | | The drawer should be kept o pen throughout the process. |
| Step 2 The lock cover can be found behind the drawer front. The lock cover has two holes as shown. | Holes | |
| Step 3 Insert a Phillips screwdriver in each hole until the screwdriver engages with the screw head. Loosen each screw by turning the screwdriver counter-clockwise 1-1 .5 turns. | | Do not remove the screws c ompletely. |
| Step 4 Remove the cover by lifting it. | | Place the cover on a clean s urface to avoid any scratche s to its surface. |


Step 5 To change the default combination or existing combination of the lock, turn the knob left (counter-clockwise) until it stops, and then release the knob. Do not force the knob at any time.



It is recommended to chang e the default combination of the lock to make sure the co mbination is private.


The number assigned to eac h button in the image will be used as a reference for any combination mentioned in th e following instructions.

Step 6 Enter the existing combinat ion. The factory default combination for a new lock is "2"&"4" – "3": Press buttons "2" and "4" simultaneously, release, then press button "3" and release. Make sure each button is fully depressed before releasing.




A distinctive click should be f elt when the buttons are pre ssed correctly. If a wrong co mbination is entered, turn th e knob left (counter-clockwis e) to the stop position and re lease. This will clear any pre viously depressed buttons. Now enter in the correct combination. The existing co mbination must be used to s et a new one.

Step 7 Push the combination chan ge slide on the back side of the loc k inward using a slotted screwdrive r. (Slide is located at the end oppo site to the knob) .



A slotted screwdriver tip size 13/64" or lower is recommen ded.

Step 8 Turn the knob left (counterclockwise) to the stop position and release. This clears the existing combination.



NOTE: Clearing the existing combination does NOT reset it to the factory default combination.


| Procedure | Images | Remarks |
|--|--|--|
| Step 9 Select a <u>new combination</u> and <u>write it down</u> . | Some or all of the buttons may be used for your new combination. A group of buttons in a combination can be pressed individually or simultaneously to generate a combination. A combination CANNOT be having the same button pressed twice (e.g 2-2-3-4 or 2-3-2-4 will not work). | IMPORTANT: Make sure the correct buttons for the new combination are depressed in the following step. The existing combination must be used to set a new one. |
| Step 10 Enter the new combination. Depress each button fully and release it. (E .g .: "1" – "2" – "3" & "4") | (Example only for combination "1"-"2"-"3"& "4") | A distinctive click must be fel t each time a button is press ed to know that the button is fully depressed. |
| Step 11 Turn the knob right (clock wise) to the stop position and relea se to activate the new combination . | 00000 | Once the knob is turned cloc kwise fully, the new combina tion for the lock is set . |
| Procedure | Images | Remarks |

Step 12

Test the new combination:

- a) Turn the knob left (counter-cloc kwise) until it stops and then relea se. Do not force the knob.
- b) Enter the new combination. (E .g .: "1" "2" "3" & "4")
- c) Turn the knob right (clockwise). The latch bolt should retract, com pressing the spring.
- d) If the knob does not turn, it ca n be due to an incorrect combinati on being entered in Step 10.


Try a few possible combinations to find the right combination entered in Step 10 then go back to Step 6 to set the desired combination again.



Step 13

To make sure the new combination works correctly:

- a) Turn the knob left (counter-cloc kwise) until it stops and then relea se. Do not force the knob.
- b) Enter any wrong combination. (E .g . "1" - "4" - "3" - "2")
- c) Try to turn the knob right (cloc kwise). The knob should not turn.
- d) If the knob turns, it means no combination is set. In this case, tur n the knob left (counter-clockwise), release, and go back to Step 6 to s et the combination.



Procedure Images Remarks Step 14 After setting up the new lock combination and testing the functio nality, align the two slots on the co ver with the two screws and reassemble the cover. Step 15 After making sure the back cover i s properly placed, tighten the two s crews using a Phillips screwdriver.

- www.amico.com
- Amico Patient Care Corporation
- 122A East Beaver Creek Road, Richmond Hill, ON, L4B 1G6 Canada

Toll Free Tel: 1.877.462.6426Toll Free Fax: 1.866.440.4986

• Tel: 905.764.0800 | Fax: 905.764.0862

• Email: apc-csr@amico.com

• www.amico.com

• APC-IFU-P02-01EN 02.22.2024

Documents / Resources



Amico APC-IFU-P02-01EN Push Button Lock [pdf] Instructions

APC-IFU-P02-01EN Push Button Lock, APC-IFU-P02-01EN, Push Button Lock, Button Lock, Lock

References

- A Amico Group of Companies | Home Page
- User Manual

Manuals+, Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.