




## BARCO MXRT-7500 Display Controller User Guide

[Home](#) » [BARCO](#) » BARCO MXRT-7500 Display Controller User Guide 



Display Controller



User Guide

Windows 7, Windows 8.1 and Windows 10

Barco NV

Beneluxpark 21, 8500 Kortrijk, Belgium

[www.barco.com/en/support](http://www.barco.com/en/support)

[www.barco.com](http://www.barco.com)

## Contents

- 1 Welcome!
- 2 Display Controller installation
- 3 Driver and software installation
- 4 Configuring Barco displays in Windows
- 5 Driver and software features
- 6 Important information
- 7 Documents / Resources
  - 7.1 References

## Welcome!

### 1.1 About the product

#### Display Controller

Thank you for choosing this Barco Display Controller!

Barco's state-of-the-art Display Controllers deliver the performance, quality and stability required for today's advanced medical imaging applications. The powerful boards ensure ultra-fast and smooth image loading, and graphics processing of images in every resolution.

Use the instructions in this guide to install your Barco Display Controller.

### 1.2 What's in the box

#### Contents

Your Barco Display Controller comes with:

- This Barco Display Controller User Guide
- 1 extender bracket is included with the MXRT-7600.
- 2 Single-Link Dongles are included with the MXRT-4500, MXRT-5500 & MXRT-7500.
- 1 low-profile bracket is included with the MXRT-2400, MXRT-2500 & MXRT-2600.
- 1 DMS-59-to-DVI adapter cable and 1 low-profile bracket are included with the MXRT-1450 & MXRT1451.



Keep your original packaging. It is designed for this Display Controller and is the ideal protection during transport and storage.

## Display Controller installation

### 2.1 Which Display Controller ?

#### Display Controller range

Your Barco medical display is compatible with a large range of Barco Display Controller boards. Depending on the order details, the display can be delivered with or without a Display Controller.


If you are installing a Barco Display Controller, please follow the installation instructions in this section.

Barco displays are compatible with a limited set of non-Barco display controllers. If you are installing a non-Barco display controller, please consult its corresponding documentation.

### 2.2 Installing a Barco Display Controller

#### Guidance

This chapter will guide you through the physical installation of a Barco Display Controller for your display system.

 **WARNING:** Wear a grounded, protective ESD strap when handling or during installation of the Display Controller. Electrostatic charges can damage the Display Controller.

## Overview

Prior to installing the Barco Display Controller (s) for your Barco Display System in your workstation, please take a few minutes to familiarize yourself with the Display Controller(s) and the PCIe slots.

### Types of Display Controllers for Barco Display Systems

The following models of Barco Display Controllers are available for your display system. Please check which of the following models is delivered with your system and follow the corresponding installation instructions:

| Barco Model      | Compatible PCIe Slot |
|------------------|----------------------|
| Barco MXRT-1450  | x11, x8, x16         |
| Barco MXRT-1451  | x11, x8, x16         |
| Barco MXRT-2400  | x16                  |
| Barco MXRT-2500  | x16                  |
| Barco MXRT-2600  | x16                  |
| Barco MXRT-45002 | x16                  |
| Barco MXRT-5400  | x16                  |
| Barco MXRT-5450  | x16                  |
| Barco MXRT-5500  | x16                  |
| Barco MXRT-5550  | x16                  |
| Barco MXRT-5600  | x16                  |
| Barco MXRT-7400  | x16                  |
| Barco MXRT-7500  | x16                  |
| Barco MXRT-7600  | x16                  |

### Which PCIe slot to use

The table above lists the various Display Controller model(s) available for your Barco Display System and the recommended PCIe slot to use for optimum performance. The figure below shows the different types of PCIe slots that can be used.

1. Recommended PCIe slot. You can also use x16 & x8 slots for x1 boards.
2. Available in limited markets




Image 2-1: Examples of PCIe slots

1. x16 slot
2. x8 slot
3. x1 slot

## 2.3 Installation procedure

**⚡ WARNING:** Wear a grounded, protective ESD strap when handling or during installation of the Display Controller. Electrostatic charges can damage the Display Controller.

**👉** If you are using a motherboard containing an on-board graphics solution and do not intend to use it as part of a multiple-display setup, disable it either in the computer's System Set-up utility (BIOS) or the Windows device manager.

### How to install

The following instructions will take you step by step through the installation of the Barco Display Controller(s) for your Barco Display System.

1. If you are not going to use your old Display Controller, uninstall its drivers and software.
2. Turn off the computer, display(s), and other peripheral devices.
3. Unplug the computer's power cord and disconnect all cables from the back of your computer.
 

**⚡ Warning:** Wait approximately 20 seconds after unplugging the power cord before disconnecting a peripheral or removing a component from the motherboard to avoid possible damage to the motherboard.
4. Remove the computer cover. If necessary, consult your computer's manual for instructions.
5. If necessary, unscrew or unfasten and remove any existing display controller(s) from your computer.
6. Locate the appropriate slot and, if necessary, remove the metal back-plate cover(s).
7. Align the Barco Display Controller(s) for your Barco Display System with the slot(s) and press it(them) in firmly until the card(s) is(are) fully seated.
8. Connect the power cable to the 6-pin power connection on the Display Controller. Make sure the cables are not interfering with anything inside the computer (for example, a cooling fan).



Tip: This step only applies to the MXRT-7400, MXRT-7500 & MXRT-7600.





Image 2-2: Power connection for the MXRT-7400, MXRT-7500 & MXRT-7600 controllers

9. Screw in or firmly fasten the Display Controller. Replace and secure the computer cover.

## 2.4 Connecting your Barco Displays

 For a detailed description of the display installation and signal connection, please refer to the Display User Guide.

### IO-Panel for the Barco MXRT-1450 & MXRT-1451




Image 2-3: MXRT-1450 & MXRT-1451

1. DMS-59 connector provides DVI-I / Head 1 & Head 2 output connections through included the Y-adapter cable.

### IO-Panel for the Barco MXRT-2400 & MXRT-2500




Image 2-4: MXRT-2400 & MXRT-2500

1. DisplayPort Connection
2. DVI-I Connection

### IO-Panel for the Barco MXRT-2600




Image 2-5: MXRT-2600

1. DisplayPort #1
2. DisplayPort #2

**IO-Panel for the Barco MXRT-5450 & MXRT-5550**




Image 2-6: MXRT-5450 & MXRT-5550

1. Head 1- DVI-I Connection
2. Head 2- DVI-I Connection

**IO-Panel for the Barco MXRT-4500, MXRT-5400, MXRT-5500 & MXRT-7400**




Image 2-7: MXRT-4500, MXRT-5400, MXRT-5500 & MXRT-7400

1. DisplayPort #1
2. DisplayPort #2
3. DVI-I

**IO-Panel for the Barco MXRT-5600, MXRT-7500 & MXRT-7600**




Image 2-8: MXRT-5600, MXRT-7500 & MXRT-7600

1. DisplayPort #1
2. DisplayPort #2
3. DisplayPort #3
4. DisplayPort #4

For displays with native DisplayPort input, use a native DisplayPort cable to connect the DisplayPort output of the Display Controller to the DisplayPort input of the display. You cannot connect the DVI output of a Display Controller to the DisplayPort input of a display.

## 2.5 Dongles

### About

Barco dongles are designed to allow Barco Display Controllers with a DisplayPort connector to a display with only a DVI input. Most current Barco displays support direct DisplayPort connection.

If video cable conversion is not required, you may bypass this section.

Barco dongles are available for purchase independently.

### 2.5.1 Passive Single-Link Dongle

#### About

The passive Single-Link Dongle converts DisplayPort input signals to single-link DVI output signals. It is compatible with all Barco grayscale displays and up to 2MP color models. For color displays of 3MP and greater resolutions, the Dual-Link Dongle is necessary.



Image 2-9: Single-Link dongle

1. To display Single-Link DVI cable
2. To DisplayPort connector on Barco Display Controller

### 2.5.2 Active Single-Link Dongle

#### About

Barco Display Controllers are not compatible with third-party active single-link dongles. Please use Barco passive Single-Link dongles.

### 2.5.3 Active Dual-Link Dongle

#### About

The active Dual-Link Dongle converts DisplayPort input signals to dual-link DVI output signals. Unlike the passive Single-Link dongle, the Dual-Link dongle allows higher resolutions (greater than 1920×1200) on color displays.

 USB calibration must be used when connecting a display through a Dual-Link dongle.




Image 2-10: Dual-Link dongle

1. To computer's USB port
2. To display via DVI Dual-Link cable
3. To DisplayPort connector on Barco Display Controller

## 2.6 Shipping the Barco Display Controller

### About

After the installation and validation of software components, Barco recommends removing the display controllers from the workstation and returning them to their original packaging prior to shipment.

 Barco does not recommend shipping display controllers installed in the workstation. If it is necessary to ship a controller installed in a workstation, the MXRT-7600 requires an extender bracket to protect against shock and vibration. Assemble the extender bracket as shown below. Refer to the workstation user documentation on proper installation to its card guide.




Image 2-11: MXRT-7600 Extender Bracket

1. Bracket
2. Extender
3. MXRT-7600

## Driver and software installation

### 3.1 Introduction



## About

This chapter will guide you through the installation of the drivers, software and documentation associated with your Barco Display System or Barco Display Controller(s).

### Prerequisites

Before starting the installation of the Barco drivers, software and documentation following prerequisites must be adhered:

- Your operating system must be installed and running. Following versions of Windows are supported:
  - Windows 7 (32-bit or 64-bit)
  - Windows 8.1 (64-bit)
  - Windows 10 (64-bit)
- You must be logged on as a user with administrator privileges.
- All Barco displays must be connected to the appropriate Display Controller(s) in your system.



For optimal system performance, Barco recommends installing no more than two drivers on a system at one time. If the configuration will require three drivers, the Barco Driver Installer will alert the user to replace one board to eliminate one of the drivers.

- When there is a non-Barco board in the system, you must first install the driver for the non-Barco display controller before installing the Barco driver.



After each driver installation, you should reboot the system before proceeding with the installation of another driver.

You will need to install the Barco Display Controller system drivers and software in the following cases:

- After you have installed the Barco Display Controller(s) for your Barco Display System in your system for the first time.
- After you have reinstalled or upgraded your operating system.
- When upgrading to a newer version of the MXRT driver and software, manual uninstallation of the prior version is not necessary. The Barco Product Installation Wizard will detect any prior installations and start the uninstallation process automatically.



The installation dialog will display in English if your operating system's language is not supported.

## 3.2 Installation procedure

### Installation procedure

1. Start your system.

If you have a fresh OS installation, or you have uninstalled an existing driver, the OS may automatically install an inbox driver from the Windows driver store, either an AMD driver or a standard VGA driver, for the Barco Display Controller(s). If this occurs, the OS prompts you to restart your computer, click Yes to allow the automatic driver installation to complete and reboot the system.

2. Launch the Barco Product Installation Wizard. The installation Wizard should start automatically when you insert the Barco Display System Installation DVD into your computer's DVD drive. If the installation does not start, in the "auto play" window, click Run setup.exe.
3. The first page of the installation wizard is the license agreement. You must accept the license agreement in order to proceed.
4. The second page of the installation wizard will show you the display controller driver and software components

that will be installed. To accept the installation of default software components, click Next. To custom select software components, click on specific software components to unselect them.




Image 3-1: Barco Product Installation Wizard

- MXRT Driver and Intuitive Workflow Tools: Barco MXRT Display Controller drivers and accompanying software for supporting Intuitive Workflow tools
  - BarcoMed Self Exam: Barco diagnostic tool
  - MediCal QAWeb Agent: Barco calibration software
5. If a previous installation of an MXRT driver exists, the installation wizard will detect it and guide you through the uninstallation process if necessary.
  6. During installation, the desktop may flash, and the Installation Wizard window may appear on different displays. This is expected behavior.
  7. When the installation of all components has completed, the system must reboot to complete the changes, and an automatic reboot window will be displayed.

### 3.3 Silent installation

#### Installation procedure

Navigate to the Barco installation folder and execute the command `setup.exe -silent`.

This can be done from the command shell, from the Run command, or from a command shortcut. The setup program will automatically install the drivers for any MXRT boards that are present, the BarcoMed Self Exam program, and QAWeb (if part of the installation package).

#### Configure silent installation options

You can modify the `setup.ini` file at the Barco root folder to customize certain silent install behavior. The configurable options are listed in the [Custom] section of the `setup.ini` file.

#### Reboot

- Locate the [Custom] section of `setup.ini`.
- If set to Yes (default), the installer will prompt the user or launch a timer to reboot following software installation.

If No, the prompt/timer is not shown.

### **3.4 Installation options**

#### **Configure installation options**

You can modify the setup. ini file at the Barco root folder to change the default setting of MXRT driver after either installation or silent installation.

#### **Install driver in 24-bit**

- Locate the [MXRT\_WDDM] section of setup.ini
- Remove -30bit command line parameter from both Install and SilentInstall lines
- Modify the default.ini file, found in:  
Barco\_MXRT\_Driver\_SoftwarePackage\_xxxx\Setup\_Barco\_Productivity\_Tools.x.x.x\Change  
Color%20Depth\30BitDesktop to false.

#### **Install driver with Coronis Fusion displays in DualView mode**

- Locate the [MXRT\_WDDM] section of setup.ini
- Remove -singleview command line parameter from both Install and SilentInstall lines

#### **Install driver with VirtualView enabled**

- Locate the [MXRT\_WDDM] section of setup.ini
- Add -virtualview parameter to both Install and SilentInstall lines
- Modify the default.ini file, found in:  
Barco\_MXRT\_Driver\_SoftwarePackage\_xxxx\Setup\_Barco\_Productivity\_Tools.x.x.x\ Change  
VirtualView\FeatureEnabled to true.

### **3.5 After installation**

#### **Installation verification**

To verify that the driver was installed, go to the Windows Control Panel, select System, select Device Manager, then select Display Adapters. Verify that Barco Display Controllers are properly identified, as shown below:




Image 3-2: Verifying driver installation

### Automatic display configuration

Once the drivers, software and documentation have been installed and the system has been rebooted, the computer should automatically detect your Barco displays and attach them to the desktop with the correct resolution. If the computer fails to detect your Barco displays or fails to attach them to the desktop correctly, please use the Windows Screen Resolution to set the correct resolution.

### Upgrade drivers

When performing a driver upgrade, the Barco System Settings Control Panel default profile will be applied. Any user profile that was saved on the system previously is still available and can be selected through the Barco System Settings Control Panel.

## 3.6 Uninstallation

### Uninstalling the drivers and software

To uninstall the Barco drivers, software or documentation for your Barco Display System, please use the Windows Add/Remove Programs. This can be found in the Windows Control Panel under Programs & Features.





Image 3-3: Windows 7 Add/Remove Programs

## Barco System Cleaner

The Barco System Cleaner is a tool that will remove all Barco software components from your workstation. This includes the display driver, the accompanying software for supporting Intuitive Workflow features, BMSE, and QAWeb calibration software. The application can be found at C:\Program Files\Barco.

 The Barco System Cleaner will remove all Barco components from your system. It is recommend to only use the System Cleaner under direction of Barco Customer Support.

## Configuring Barco displays in Windows

### 4.1 Display resolution

#### Changing resolutions

1. Click on the Configure Displays button on the Barco System Settings Control Panel or right click on the desktop and select Screen resolution in Windows 7 and 8.1. This will open the Windows Control Panel shown in Image 4-1. In Windows 10, right click on the desktop and select Display settings to launch a control panel with similar functions but a different look and feel.
2. Click on the Resolution drop down box to show the list of resolutions.
3. Select the desired resolution, and click Apply.





Image 4-1: Windows Display Control Panel

 In Windows 7, the maximum horizontal resolution of the desktop is 8192 pixels for a single display controller. Your calculation should include the VirtualView display, if used. Please refer to Microsoft Knowledge Base article 2724530.

### 4.2 Software rotation

#### Configuring screen orientation

Software rotation is only necessary for displays that do not support hardware rotation, such as MDRC displays and some third party displays.

1. Right click on the Desktop, and select Screen Resolution in the context menu.
2. Select a Display.
3. In the Orientation drop down list, these options are available:

- Landscape
- Portrait
- Landscape (flipped)
- Portrait (flipped)

4. Select the desired setting, and click Apply.

## Driver and software features

### 5.1 Barco System Settings Control Panel

#### 5.1.1 Description

##### Overview

The Barco System Settings Control Panel provides a centralized configuration interface for users to personalize their Barco Display System environment.




Image 5-1: Barco System Settings Control Panel



#### 5.1.2 Operation

Accessing the Barco System Settings Control Panel

- Right click on the desktop and select Barco System Settings, or
- Click on the Barco icon in the System Tray, or
- Press the hotkey Control+ALT+O

### Layout of the Barco System Settings Control Panel

The Barco System Settings Control Panel consists of three sections: the Control Bar at the top, the Navigation pane on the left, and the configuration tiles on the right. The Control Bar shows the current active profile and provides four additional buttons.

- Open User Manual: This opens a PDF version of this user guide if a PDF viewer is installed in the system.
- Identify Displays: Some settings in the Barco System Settings Control Panel require selection of individual displays. Press this button to see which display corresponds to which configuration numbers.  
 The identification numbers used by the Barco System Settings Control Panel do not reflect the display ID number assigned by Windows as shown in the Windows Screen Resolution configuration page.
- Configure Displays: This opens the Windows Screen Resolution Control Panel to allow changes to the location and resolution of displays in the Windows desktop.
-  : Clicking the gear symbol opens the Control Panel Settings Box and allows the user to change settings of the Barco System Settings Control Panel.
  - Hotkey to Show This Control Panel: Click Clear to remove the hotkey. Click on the hotkey box to enter a new one.
  - Disable Navigation Bar Popups: This button can enable popups on the Navigation Pane, making the icons larger and more legible.

Through the Navigation Pane, the user can get quick access to the configuration tile of a feature by clicking on its icon. The Navigation Pane also allows customization of the Barco System Settings Control Panel. Dragging icons in the Navigation Pane will change the order of the tiles, allowing commonly used tiles to appear at the top. In the middle of the Navigation Pane is the visibility line. The tile for any feature can be hidden by dragging its corresponding Navigation Pane icon below the visibility line, reducing control panel clutter. If the icons in the Navigation Pane appear too small, popups with larger icons can be enabled in the Control Panel Settings Box.

### Profile management

The Profile Management tile is an exception to the rules above; it cannot be reordered or hidden. Selected options in the Barco System Settings Control Panel can be saved in profiles, and those profiles are managed with this tile. The profiles are specific to the current user, and they will be automatically applied when the user logs into the system.

The profile can be exported to and imported from a remote disk, so it can be deployed to multiple systems.

- Profile: The current profile is selected with this drop-down menu. If the profile is edited, it will be marked here as “modified”.
- Save: This saves changes to the current profile. Changes to the Default profile cannot be saved.
- Save as...: This saves the current configuration setting as a profile with a new name. Provide a profile name in the popup dialog box.
- Undo: This reverts to the saved version of the current profile.
- Import...: This imports a profile from file and makes it available in the Profile list to be selected.
- Export...: This saves the current profile into a file. To apply a profile system wide:

1. Export the profile, and name it Default.ini



2. Edit the new file:
    - Change line 3 from names=<profileName> to names=Default
    - Change line 5 from [<profileName>] to [Default]
  3. Replace the default profile in c:\ProgramData\Barco\ProductivityTools\Default.ini with the new one
- Delete: This removes the current profile from the profile list and restores the system to the Default profile.

## 5.2 Application Appearance Manager

### 5.2.1 Description

#### Overview

The high luminance of Barco diagnostic displays may not be necessary when using the displays to view text documents, emails, or other non-diagnostic applications. The Application Appearance Manager (AAM) feature allows the user to set all windows of specific applications to a lower desired luminance, while retaining the full diagnostic luminance for all other applications.

For the SteadyColor™ displays, like Coronis Uniti and the MDNC-6121, AAM can also change the output color profile of the windows of specific applications to match their expectations. This may be useful, for instance, on a web browser used for non-diagnostic purposes like viewing artwork and photos.



Image 5-2

#### Supported display controllers

When driven by an MXRT-x500 or later display controller, AAM can be enabled on selected Barco Coronis, Nio, and Mammography displays. These include the MDCC-6430, MDCC-6330, MDCC-4330, MDCC-4230, MDNC-6121, and MDMC-12133.

### 5.2.2 Operation

#### Using Application Appearance Manager

Select each desired application and add it to the list of managed applications. Independently edit the luminance and color profile for each managed application.


Certain applications can be blocked from being managed by AAM. In the directory C:\Program Files\Barco\ProductivityTools, open the AAMBlackList.txt file. Enter the name of the application executable, and after rebooting the system, the application will not appear in AAM. If the name of the executable isn't known, launch the application, enable AAM, and note the name of the executable from the unmanaged list. Administrative rights are required to edit the file.

### 5.2.3 Configuration

#### Configuring color management

- Enabled: This slider will enable or disable the Application Appearance Manager.
- Manage Selected App Name: Check the box of an application to control its appearance.
- Color Profile: Select the desired color profile for each controlled application.



 Color profile selection is only available on Coronis Uniti displays that are calibrated with QAWeb 1.13.01 or later and MDNC-6121 displays that are calibrated with QAWeb 1.13.10 or later.

- **Luminance:** Change the luminance of the selected managed application with this slider. Although the slider will allow entry up to 1000 candela, the maximum luminance is the calibration luminance of the display showing the application window. The minimum luminance is 250 candela.

## **5.3 Conference CloneView™**

### **5.3.1 Description**

#### **Overview**

This feature allows the user to clone the images sent to one or more displays to other displays or projectors attached to the same Barco Display Controller. Conference CloneView supports zooming and panning on the cloned image(s) for ease of viewing.



**Image 5-3**

#### **Supported display controllers**


Conference CloneView is supported by all systems with an MXRT-x400 or later display controller.

### **5.3.2 Operation**


#### **Using Conference CloneView**

Create a new clone session and select up to 3 source displays and up to 3 target displays for the clone session. The source and target displays must be attached to the same display controller. The cloned image can be scaled to fit the resolution of the target display. With a Barco Display Controller with 4 outputs, it is possible to have two independent cloning sessions.

When the cursor is over to the cloned image, it will change to the Barco cursor. The user can zoom in on the cloned image by rolling the mouse wheel and pan the image with the left mouse button.

 When zoomed in, the cloned image may be bigger than the clone display and part of the image may be off screen. Click and drag the left mouse button and move the cursor to pan to the portion of the image that is off screen.

Right clicking on the cloned image will bring up the Conference CloneView context menu. The menu options are described below.



**Image 5-4: Conference CloneView Context Menu**

- **Minimize:** This minimizes the cloned image to show the desktop.
- **Zoom to 1:1:** Changes the scaling to 1 target pixel per source pixel. If the source resolution is larger than the target resolution, panning will be necessary to see the entire image. If the source resolution is smaller than the target resolution, black borders will show around the image. “Stretch to fill” has priority over this option.
- **Reset Zoom Level:** When “Stretch to fill” is selected in the Barco System Settings Control Panel, this will reset to the minimal zoom possible to make the stretched aspect ratio possible. When “Stretch to fill” is not selected, this option will reset the zoom to 1:1. When the image is already at the minimum zoom level, this option is grayed out.
- **Pause Cloning:** This suspends updates to the clone image; the source can continue to change while the target image remains static. Select it a second time to resume active cloning.
- **Stop Cloning:** Select this option to stop cloning on the current Display Controller. It has the same effect as clicking on the Stop Cloning button in the Barco System Settings Control Panel.

### **5.3.3 Configuration**

#### **Configuring Conference CloneView**

- **Add clone configuration:** Click this button to define a new clone configuration.
- **Clone From and Clone To:** Click on one or more source displays in the Clone From list and one or more target displays from the Clone To list for the session. An active source cannot be used as a target, and an active target cannot be used as a source.
- **Stretch to fill:** If this is unchecked, the aspect ratio of the source displays will be preserved. If this is checked, the cloned image will be stretched to fill the target display(s) with the cloned image.
- **Start/Stop:** Clicking the Start button activates the clone session, and the button will change to Stop, which would end the session.
- **Remove this configuration:** Deletes this configuration.

## **5.4 DimView™**

### **5.4.1 Description**

#### **Overview**

The DimView feature reduces ambient light during diagnostic readings by dimming navigational displays when the cursor is moved off those displays. While the feature is intended for use with navigational heads, it can be enabled on any display.



Image 5-5

#### **Supported display controllers**

DimView is supported on the MXRT-1450 , MXRT-1451, and all MXRT-x500 or later display controllers.

#### **5.4.2 Operation**

##### **Using DimView**

DimView can be individually enabled on each display. All DimView-enabled displays can operate independently, or they can be configured to dim and brighten together.

#### **5.4.3 Configuration**

##### **Configuring DimView**

- Enable DimView for these displays: All displays that support DimView are listed in the configuration section. Click on the checkbox to enable the feature on that display.
- Dim/brighten selected displays together: When this box is checked, all DimView-enabled displays will brighten when the cursor is moved onto any one of those displays, and dim only when the cursor is off all of them.
- Luminance of dimmed displays: This slider bar sets the luminance while dimmed.

### **5.5 Film Clip**

#### **5.5.1 Description**

##### **Overview**

Film Clip allows the user to view a physical radiological film by using the I-Luminate™ feature of the display as a virtual light box.



Image 5-6

#### **Supported display controllers**

When driven by an MXRT-x400 or later display controller, Film Clip can be enabled on the MDMG-5221 and MDMC-12133 displays.

#### **5.5.2 Operation**

##### **Using Film Clip**

The size and location of the film clip light box are programmable. It has an automatic time-out with a programmable duration. An optional hotkey can quickly turn the light box on or off. Supported film sizes for both MDMG-5221 and MDMC-12133 include the 18 cm x 24 cm and 24 cm x 30 cm; the MDMC-12133 also supports the larger 34 cm x 43 cm size.

### 5.5.3 Configuration

#### Configuring Film Clip

- Enable Film Clip hotkey for these displays: All displays that support Film Clip are listed in the configuration section. Click on the checkbox to launch the feature on that display with the hotkey.
- Timeout: The slider bar sets the time-out period for Film Clip mode.
- Hotkey to enabled Film Clip: Click Clear to remove a hotkey. Click on the hotkey box to enter a new one.
- Position: This sets the location of the light box image on enabled displays.
- Size: This sets the size of the light box image to match the physical film dimensions.

## 5.6 FindCursor™

### 5.6.1 Description

#### Overview

The FindCursor feature provides a method to quickly locate the cursor on a system with multiple displays.



Image 5-7

#### Supported display controllers

FindCursor is supported on all MXRT-x400 or later display controllers.

#### Default behavior

FindCursor is enabled by default.

### 5.6.2 Operation

#### Using FindCursor

To quickly locate the cursor, hold down the hotkey (default: Control +Shift +F). The cursor location will be highlighted by a circle, which appears yellow on color displays and gray on grayscale displays.

### 5.6.3 Configuration

#### Configuring FindCursor

- Enable/Disable: Use the checkbox to enable or disable FindCursor.
- Hotkey: The currently selected hotkey is displayed in the edit box. To program a new hotkey, highlight the edit box, and input the new keystrokes. The change will take affect right away.

## 5.7 I-Luminate™

### 5.7.1 Description

#### Overview

This feature boosts the luminance of supported displays.



Image 5-8

### **Supported display controllers**

When driven by an MXRT-x400 or later display controller, I-Luminate can be enabled on supported displays, including MDMC-12133, MDMG-5221, and MDCG-5221.

#### **5.7.2 Operation**

##### **Using I-Luminate**

A hotkey will boost the luminance on all selected displays.

#### **5.7.3 Configuration**

##### **Configuring I-Luminate**

- Enable I-Luminate for these displays: All displays that support I-Luminate are listed in the configuration section. Click on the checkbox to enable the feature on that display.
- Hotkey to enabled I-Luminate: Click Clear to remove a hotkey. Click on the hotkey box to enter a new one.
- Timeout: The slider bar sets the time-out period for the I-Luminate mode.

### **5.8 Reading Environment**

#### **5.8.1 Description**

##### **Overview**


The reading environment settings for diagnostic displays allow the user to specify the color temperature and luminance of the MDMC-12133 and MDNC-6121. The reading environment for non-diagnostic displays allows the user to specify the maximum luminance on Barco clinical displays and third party displays, however color temperature selection is not available.



Image 5-9

### **Supported display controllers**


When driven by an MXRT-x500 or later display controller, the Reading Environment can be configured on MDMC-12133, MDNC-6121, Barco MDRC displays, and third party displays.

 Diagnostic reading environment configuration is available on Coronis Uniti displays that are calibrated with QAWeb 1.13.01 or later and on MDNC-6121 displays that are configured with QAWeb 1.13.10 or later. Additionally, the system must not be connected to the QAWeb Server. If connected to the server, set the reading environment using the QAWeb Agent. There is no QAWeb requirement for non-diagnostic reading environment configuration.

## 5.8.2 Operation

### Using Reading Environment for SteadyColor Diagnostic Displays

Configure the white point chromaticity, ambient light, color calibration model, and calibrated luminance to user preference. After the settings have been changed, the QAWeb Agent will check if it has calibration files for the new settings. If it does, it will upload the calibration data to the display, which may take 20-60 seconds. If it does not have the calibration files, it will calibrate the display to create them, which may take up to 10 minutes.

 When changes to the Reading Environment for Diagnostic Displays are made, the Barco System Settings Control Panel will indicate that changes are still in progress. Please let the system apply its changes before proceeding with readings.


### Using Reading Environment for Non-Diagnostic Displays

Reading Environment can be individually enabled on each non-diagnostic display to reduce its luminance. Unlike DimView, the luminance of these displays will not change with cursor movement. The Reading Environment feature can be used in conjunction with DimView.

## 5.8.3 Configuration

### Configuring Reading Environment for Diagnostic Displays

- Use the settings QAWeb already has: When selected, reading environment settings for diagnostic displays are disabled in the Barco System Settings Control Panel.
- Suggest the following settings to QAWeb: When selected, the reading environments are modifiable.
- White point chromaticity: This selects between clearbase, bluebase, and native white points.
- Ambient light condition: This selects the expected ambient light condition based on reading room class.
- SteadyColor calibration: This selects the color calibration model of the SteadyColor display.
- Luminance of white: This slider sets the calibration luminance of the display.

 In some installations with QAWeb, Ambient light condition, SteadyColor calibration, and Luminance of White do not appear in the Barco System Settings Control Panel. These can be configured in QAWeb with an appropriate image quality policy.

### Configuring Reading Environment for Non-Diagnostic Displays

- Select the displays: All supported non-diagnostic displays are listed in the configuration section. Click on the checkbox to enable the feature on that display.
- Display luminance: This slider bar sets the luminance reduction for the selected displays.

## 5.9 Screen Capture

### 5.9.1 Description

#### Overview

The Screen Capture feature captures the desktop into an image, including the Intuitive Workflow features, such as

SpotView.



Image 5-10

### **Supported display controllers**

This feature is supported by MXRT-x400 and above display controllers.

### **5.9.2 Operation**

#### **Using Screen Capture**

The screen capture is triggered through a hotkey. The user can choose to capture the image of the display with the cursor to the clipboard. The user can also choose to capture an image of each display to a file in either PNG or PPM format.

### **5.9.3 Configuration**

#### **Configuring Reading Environment**

- Hotkey: Click Clear to remove a hotkey. Click on the hotkey box to enter a new one.
- Copy captured image to Windows clipboard: Check this box to copy the image of the display with the cursor to the clipboard.
- Output directory: To save the capture images as files, enter a target location here.
- File format: Use these radio buttons to choose the file format as either PNG or PPM.

## **5.10 SingleView™**

### **5.10.1 Description**

#### **Overview**

SingleView enables the use of a Coronis Fusion display as a single display in the Windows desktop and eliminates any tearing down the center. SingleView is enabled by default on Coronis Fusion displays.



Image 5-11

### **Supported displays**

When driven by an MXRT-x400 or later display controller, SingleView is supported on all Coronis Fusion displays.

### 5.10.2 Operation

#### Using SingleView

A button will switch all Coronis Fusion displays into SingleView, and another will switch all into DualView.

Hotkeys are available for both actions.

When in SingleView, it is possible to have the left and right halves out of position, and a button will swap the two.

### 5.10.3 Configuration

#### Configuring SingleView

- Enable SingleView: This check box will toggle between SingleView and DualView.
- Swap Left/Right: This can correct misaligned SingleView displays one at a time.
- Hotkeys: SingleView and DualView activations have individual hotkeys. Click Clear to remove a hotkey.  
Click on the hotkey box to enter a new one.

## 5.11 SmartCursor™

### 5.11.1 Description

#### Overview

The Barco SmartCursor feature prevents the cursor from becoming stuck on edges of adjacent displays of different sizes.



Image 5-12

#### Supported display controllers

SmartCursor is supported by all systems with an MXRT-x400 or later display controller.

### 5.11.2 Operation

#### Using SmartCursor

The SmartCursor operation is illustrated in Image 5-13.

Consider two points, A and B, on two displays of different sizes. Without SmartCursor, the cursor cannot move left from point A because it will be stuck on that edge. With SmartCursor, moving the cursor left from point A will move the cursor to point B. For symmetry, with SmartCursor moving to the right from point B, the cursor will appear back at point A.





Image 5-13: *SmartCursor* moving left from point A and right from point B.

### 5.11.3 Configuration

#### Configuring SmartCursor

Enable: Click on this checkbox to turn SmartCursor on or off.

### 5.12 SoftGlow™

#### 5.12.1 Description

##### Overview

The Barco Coronis Uniti display supports SoftGlow. It consists of a task light, which sheds a light on the desktop, and a wall light, which provides ambient lighting for the reading room to reduce eye fatigue. The brightness of each is configurable.



Image 5-14

#### Supported display controllers

When driven by an MXRT-x400 or later display controller, SoftGlow can be configured on MDMC-12133 displays.

#### 5.12.2 Operation

##### Using SoftGlow

The Uniti lights will be set to the SoftGlow settings when the user has logged on to the Barco display system.

### 5.12.3 Configuration

#### Configuring SoftGlow

Task Light and Wall Light: Use these sliders to set the brightness of the lights. Click and enter 0 to shut the light off.

### 5.13 SpotView™

#### 5.13.1 Description

##### Overview

The SpotView feature allows focused observation during readings by dimming images outside a region of interest and optionally enhancing the contrast in the region of interest.

The SpotView Mag feature offers 2x zoom within the SpotView region of interest. SpotView Invert inverts the pixels in the region of interest. SpotView Align creates a bar-shaped region of interest. SpotView Align has two modes of operation, a straight bar, or a V-shaped bar.



Image 5-15


### Supported displays and display controllers

When driven by an MXRT-x500 or later display controller, SpotView can be enabled on selected Barco Coronis, Nio, and Mammography displays, including MDCC-6430, MDCC-6330, MDCC-4330, MDCC-4230, MDCG-5221, MDCG-3221, MDNC-6121, MDNC-3421, MDNC-3321, MDNC-2221, MDNG-5221, MDMG5121, MDMG-5221, MDMC-12133.

### 5.13.2 Operation

#### Using SpotView

SpotView highlights a region of interest. The region of interest is selected by the use of the Barco Touchpad or by the mouse and a hotkey (default: Control+Shift+X). To control SpotView with the touchpad, hold and move one finger. The highlighted region of interest is always bound to displays that support SpotView.

 To show SpotView when the Barco Touchpad is in Mouse Emulation Mode, hold one finger then tap a second. To enhance viewing on Coronis and Mammography displays, the SpotView feature boosts the luminance of the display, if supported. The boost feature will turn off after one minute of continuous use. To further enhance viewing, the SpotView feature will optionally enhance the contrast of the region of interest. Contrast enhancement is not available with SpotView Align.

SpotView Mag highlights a region of interest, boosts luminance, and applies 2x zoom to the area. It is controlled through the Barco Touchpad or by the mouse and a hotkey (default: Control +Shift +Z). For the touchpad, while holding one finger to show SpotView, tap a second finger to toggle on (or off) SpotView Mag.

SpotView Invert inverts pixels in the region of interest. It is controlled through the Barco Touchpad or by a mouse and a hotkey (default: Control +Shift +S). For the touchpad, while holding one finger to show SpotView, double-tap a second finger to toggle on (or off) SpotView Invert. SpotView Invert can be used simultaneously with SpotView Mag with the Barco Touchpad or by a mouse and hotkey (default: Control +Shift+A).

When using the hotkey to turn on SpotView, the Spot will appear with the current cursor position being the center. When the Spot is moved onto a display that does not support the feature, it will not show the spot.

SpotView Align implements SpotView technology in different shapes to enable alternative uses. The two alternate shapes are a bar and a vee. They are only available through use of the Barco Touchpad, and both allow custom angles specified by the user.



Image 5-16

To enable the SpotView Align bar shape, first hold one finger on the touchpad to show SpotView then hold two fingers on it to show SpotView Align. Rotate two fingers for the desired angle, and keep just one finger to lock the angle. The bar can be moved by dragging the one finger. The angle of SpotView Align can be adjusted again by rotating two fingers on the touchpad. While holding one finger to show SpotView Align, tap or doubletap a second finger to toggle on (or off) SpotView Mag or SpotView Invert.



Image 5-17

To enable the SpotView Align vee shape, first enable the bar shape, then hold three fingers to show the mirrorimage vee shape. Rotating two fingers will rotate the right bar, and the left image will follow. While holding one finger to show SpotView Align, tap or double-tap a second finger to toggle on (or off) SpotView Mag or SpotView Invert.

### 5.13.3 Configuration

#### Configuring SpotView

- Enable SpotView for these displays: All displays that support SpotView are listed in the configuration section. Click on the checkbox to enable the feature on that display.
- Bar: Check this box to enable the SpotView Align bar shape by holding down two fingers on the touchpad.
- Vee: Check this box to enable the SpotView Align vee shape by holding down three fingers on the touchpad.
- Luminance outside of spot: This slider bar sets the luminance on SpotView enabled displays outside of the spot.
- SpotView size, SpotView Mag size, and SpotView Align Width: These sliders control the diameter of their spots and the width of the SpotView Align bar.  
**Note:** These sizes are described in centimeters. The apparent size will vary due to parallax, rounding of the display size, and the shaded penumbra of the spot. It is not intended to be used alone to make exact measurements of body parts.
- Enable Dynamic Contrast Enhancement: Click on this checkbox to turn on the contrast enhancement feature. This is only available for the standard SpotView circle shape.
- Hotkeys: SpotView and SpotView Mag, both with and without SpotView Invert, have individual hotkeys. Click Clear to remove a hotkey. Click on the hotkey box to enter a new one.

## 5.14 Touchpad gestures

### 5.14.1 Description

#### Overview

In addition to controlling SpotView, the Barco Touchpad can control the cursor and allows the user to control PACS and other applications with multi-touch gestures and touchpad buttons. The user can program the gestures and buttons to send shortcut keystrokes recognized by the PACS system.



Image 5-18

### Supported display controllers

The touchpad gestures are supported by all systems with an MXRT-x400 or later display controller and the Barco Touchpad CTH-480 or later.

### 5.14.2 Operation







#### Using the Barco Touchpad Mouse Emulation

When Mouse Emulation Mode is enabled, the user can control the cursor using the Barco Touchpad instead of a mouse. The gestures recognized in Mouse Emulation Mode are listed in the table below. To show SpotView with Mouse Emulation Mode enabled, hold one finger and tap a second.

| Touchpad gesture  | Emulated Mouse action                |
|---|--------------------------------------|
| Move one finger   | Move cursor                          |
| Tap one finger  | Left click                           |
| Double-tap one finger   | Left double-click                    |
| Tap two fingers   | Right click                          |
| Double-tap two fingers  | Right double-click                   |
| Double-tap one finger, holding it the second time and moving  | Left click and drag                  |
| Double-tap two finger, holding one the second time and moving | Right click and drag                 |
| Pinch two fingers together                                    | Control+scroll wheel down (Zoom out) |
| Pinch two fingers out   | Control+scroll wheel up (Zoom in)    |
| Drag two fingers up   | Scroll wheel up                      |
| Drag two fingers down   | Scroll wheel down                    |
| Drag two fingers left   | Click scroll Wheel and drag left     |
| Drag two fingers right  | Click scroll Wheel and drag right    |

#### Using the Barco Touchpad Gesture Recognition

The Barco Touchpad recognizes ten two- and three-finger gestures: Two- and Three-finger Swipe Left, Two and Three-finger Swipe Right, Two- and Three-finger Swipe Up, Two- and Three-finger Swipe Down, Two-finger Pinch In, and Two-finger Pinch Out. These are shown in the table below.

|            |   |             |   |
|------------|---|-------------|---|
| Swipe Left |   | Swipe Right |   |
| Swipe Up   |  | Swipe Down  |  |
| Pinch In   |  | Pinch Out   |  |

Each gesture can be programmed to send keystroke shortcuts, as if those keys were pressed on the keyboard. The gesture can send a single keystroke or send the keystroke continuously until the gesture ceases. Some gestures are used by default in Mouse Emulation Mode; those defaults can be overwritten in the Barco System Settings Control Panel.

### 5.14.3 Configuration

#### Configuring the Barco Touchpad

From the Barco System Settings Control Panel, navigate to the Touchpad Gestures section to associate the Barco Touchpad gestures with desired shortcut keys.

- **Enable Mouse Emulation:** Click on this checkbox to control the cursor with the touchpad.
- **Speed:** This slider controls the responsiveness of the cursor in Mouse Emulation Mode.
- **Emitted Shortcut:** This field shows the currently defined keyboard shortcut for each gesture. Click on the field to define a new shortcut.
- **Clear:** Click the Clear button to remove the hotkey for a given gesture or button.
- **Continuous:** This checkbox enables continuous shortcut emission.
- **Rate:** This slider controls the rate of continuous shortcut emission.
- **Shortcut Description:** Click on this free text field to describe the intention of the shortcut.

## 5.15 VirtualView™

### 5.15.1 Description

#### Overview

VirtualView gives the user additional real-estate on the screen by creating a virtual display in Windows without the need for an additional physical display on the desk. A virtual display is created for the user to use as a navigational head, or for other software, such as dictation. The user can set the location of the virtual display, and when the cursor moves into that virtual area, or hotkey is triggered, the virtual display appears on the Windows desktop.



Image 5-19

### **Supported display controllers**

VirtualView is supported by all systems with an MXRT-2500, MXRT-4500, MXRT-5500, MXRT-7500, MXRT5600, or MXRT-7600 display controller.

### **5.15.2 Operation**

#### **Using VirtualView**

VirtualView creates a virtual display within the Windows desktop. In the Windows control panel, it appears as a normal display with a display number, and its resolution and location can be changed similarly to a physical display.

By moving the cursor to the virtual display location on the desktop, the virtual display appears on a physical display. VirtualView can also be activated and hidden by a hotkey (default: Control +Shift +V). It can be stretched to make the contents larger or shrunk to make the window take less desktop space. Windows and applications can be dragged and dropped onto the virtual display, and they will show only when VirtualView is showing.

To promote usability with PACS applications, VirtualView has a visibility setting. In the Diagnostic visibility mode (default), it will appear on the top of all other windows, and in the Administrative visibility mode, it may be hidden by other windows, including the PACS software.

The virtual display window has Minimize, Maximize, Restore, and Close buttons.

- Clicking the Minimize button hides VirtualView and creates an icon on the Task Bar, and clicking on the icon will restore it to its previous size and location.
- Clicking the Maximize button resizes VirtualView to the largest size possible on that physical display while maintaining aspect ratio, and that button will change to Restore.
- Clicking on the Restore button returns the maximized window to its previous size and location.
- Clicking the Close button does not terminate the VirtualView feature. It hides the window, and when shown again, it will appear in the default location and size.

### **5.15.3 Configuration**

#### **Configuring VirtualView through Barco System Settings Control Panel**

- Connect: VirtualView is disabled by default. Click on this checkbox to connect a VirtualView display to the desktop.

Once VirtualView has been connected to the desktop through the Barco System Settings Control Panel, the VirtualView display will show up in the Windows Display Control Panel as Barco Virtual display, and can be controlled as a regular display on the desktop.

- Reset window size: This resets the virtual display window to 1-to-1 scaling.

- Configure Visibility: Selects between Diagnostic and Administrative visibility modes.
- Configure Hotkey: Specify a hotkey to show/hide virtual display window.

## Configuring VirtualView through Windows Display Control Panel

The Windows Display Control Panel can be opened through the Barco System Settings Control Panel. Or, from the Control Panel, select Display, and select Screen Resolution. Or, click on the desktop and select Screen Resolution.

The Barco virtual display resolution and location on the Windows desktop can also be changed in the control panel in the same method as normal displays. VirtualView supports portrait and landscape resolutions up to 1200 x 1920 portrait (1920 x 1200 for landscape).

## Important information

### 6.1 Safety information

#### General recommendations

Read the safety and operating instructions before operating the device.

Retain safety and operating instructions for future reference.

Adhere to all warnings on the device and in the operating instructions manual.

Follow all instructions for operation and use.

#### Electrical Shock or Fire Hazard

To prevent electric shock or fire hazard, do not remove cover.

No serviceable parts inside. Refer servicing to qualified personnel.

Do not expose this apparatus to rain or moisture.

#### Modifications to the unit

Do not modify this equipment without authorization of the manufacturer.

#### Type of protection (Electrical)

Device with external power supply: Class I equipment.

#### Degree of safety (flammable anesthetic mixture)

Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

### 6.2 Environmental information

#### Disposal Information

#### Waste Electrical and Electronic Equipment



This symbol on the product indicates that, under the European Directive (EU) 2015/863 governing waste from electrical and electronic equipment, this product must not be disposed of with other municipal waste.

Please dispose of your waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate these items from other types of waste and recycle them responsibly to promote the sustainable reuse of material resources.

For more information about recycling of this product, please contact your local city office or your municipal waste disposal service.

For details, please visit the Barco website at: <http://www.barco.com/en/AboutBarco/weee>

#### Turkey RoHS compliance



Türkiye Cumhuriyeti: AEEE Yönetmeliğine Uygundur.

[Republic of Turkey: In conformity with the WEEE Regulation] Chinese Mainland RoHS

According to the "Management Methods for the Restriction of the Use of Hazardous Substances in Electrical and Electronic Products" (Also called RoHS of Chinese Mainland), the table below lists the names and contents of toxic and/or hazardous substances that Barco's product may contain. The RoHS of Chinese Mainland is included in the MCV standard of the Ministry of Information Industry of China, in the section "Limit Requirements of toxic



substances in Electronic Information Products”.

| Component name             | Hazardous substances and elements |    |    |      |     |      |
|----------------------------|-----------------------------------|----|----|------|-----|------|
|                            | Pb                                | Hg | Cd | Cr6+ | PBB | PBDE |
| Printed Circuit Assemblies | x                                 | o  | o  | o    | o   | o    |

This table is prepared in accordance with the provisions of SJ/T 11364.

o: Indicates that this toxic or hazardous substance contained in all of the homogeneous materials for this part is below the limit requirement in GB/T 26572.

x: Indicates that this toxic or hazardous substance contained in at least one of the homogeneous materials used for this part is above the limit requirement in GB/T 26572.

All Electronic Information Products (EIP) that are sold within Chinese Mainland must comply with the “Marking for the restriction of the use of hazardous substances in electrical and electronic product” of Chinese Mainland, marked with the Environmental Friendly Use Period (EFUP) logo. The number inside the EFUP logo that Barco uses (please refer to the photo) is based on the “General guidelines of environment-friendly use period of electronic information products” of Chinese Mainland.



Taiwan RoHS

#### Declaration of the Presence Condition of the Restricted Substances Marking

Equipment name 102-c58708-01; MXRT-5600; MXRT-7600 Type designation (Type)

| Unit          | Restricted substances and its chemical symbols |      |      |        |       |        |
|---------------|--|------|------|--------|-------|--------|
|               | (Pb)   | (Hg) | (Cd) | (Cr+6) | (PBB) | (PBDE) |
| Fansink       |  | o    | o    | o      | o     | o      |
| Circuit Board |  | o    | o    | o      | o     | o      |
| Bracket       | o  | o    | o    | o      | o     | o      |
| Screw         | o  | o    | o    | o      | o     | o      |

**Note 1** “Exceeding 0.1 wt %” and “exceeding 0.01 wt %” indicate that the percentage content of the restricted substance exceeds the reference percentage value of presence condition

**Note 2** “o” indicates that the percentage content of the restricted substance does not exceed the percentage of reference value of presence.

**Note 3** The “-” indicates that the restricted substance corresponds to the exemption

#### 6.3 Regulatory compliance information

##### FCC class B

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This device generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this device does cause harmful



interference to radio or television reception, which can be determined by turning the device off and on, the user is encouraged to try to correct the interference by one or more of the following measures:





- Reorient or relocate the receiving antenna.
- Increase the separation between the device and receiver.
- Connect the device into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.





## 6.4 Explanation of symbols













### Symbols on the device









On the device or power supply, you may find the following symbols (nonrestrictive list):

|   |  |
|---|--|
|    | Indicates compliance with Part 15 of the FCC rules (Class A or Class B)            |
|    | Indicates the device is approved according to the UL regulations                   |
|   | Indicates the device is approved according to the UL regulations for Canada and US |
|  | Indicates the device is approved according to the UL regulations for Canada and US |

|   |  |
|---|--|
|  | Indicates the device is approved according to the UL Demko regulations |
|  | Indicates the device is approved according to the CCC regulations      |
|  | Indicates the device is approved according to the VCCI regulations     |
|  | Indicates the device is approved according to the KC regulations       |

|   |   |
|---|---|
|     | Indicates the device is approved according to the BSMI regulations  |
|    | Indicates the device is approved according to the PSE regulations   |
|    | Indicates the device is approved according to the EAC regulations   |
|    | Caution: Federal law (United States of America) restricts this device to sale by or on the order of a licensed healthcare practitioner. |
|    | Indicates the USB connectors on the device  |
|   | Indicates the DisplayPort connectors on the device  |
|  | Indicates the legal manufacturer  |
|  | Indicates the manufacturing date  |
|  | Indicates the temperature limitations <sup>3</sup> for the device to safely operate within specs  |
|  | Indicates the device serial number  |
|  | Indicates the device part number or catalogue number  |
|  | <b>Warning:</b> dangerous voltage   |

|   |  |
|---|--|
|   | Caution  |
|  | Consult the operating instructions   |
|  | Indicates this device must not be thrown in the trash but must be recycled, according to the European WEEE (Waste Electrical and Electronic Equipment) directive |
|  | Indicates Direct Current (DC)  |
|  | Indicates Alternating Current (AC)   |
|  | Stand-by   |
|  | Equipotentiality   |
|  | Protective earth (ground)  |

## 6.5 Legal disclaimer

### Disclaimer notice

Although every attempt has been made to achieve technical accuracy in this document, we assume no responsibility for errors that may be found. Our goal is to provide you with the most accurate and usable documentation possible; if you discover errors, please let us know.

Barco software products are the property of Barco. They are distributed under copyright by Barco NV or Barco Inc., for use only under the specific terms of a software license agreement between Barco NV or Barco Inc. and the licensee. No other use, duplication, or disclosure of a Barco software product, in any form, is authorized.

The specifications of Barco products are subject to change without notice.

### Trademarks

All trademarks and registered trademarks are property of their respective owners.

### Copyright notice

This document is copyrighted. All rights are reserved. Neither this document, nor any part of it, may be reproduced or copied in any form or by any means – graphical, electronic, or mechanical including photocopying, taping or information storage and retrieval systems – without written permission of Barco.

© 2018 Barco NV all rights reserved.

### Privacy policy

Barco is concerned about respecting the privacy of our users. Please see our privacy policy on our website at <http://www.barco.com/en/about-barco/legal/privacy-policy> to learn more.

## 6.6 Technical specifications

### MXRT-1450

|                          |   |
|--------------------------|---|
| Product acronym          | MXRT-1450   |
| Bus compatibility        | Could be installed within PCI Express x1, x8, x16 mechanical slots Works in PCI Express x1, x4, x8, x16 electrical slots / operates at x1 speed   |
| Power consumption        | 17 Watt   |
| Form factor              | Low profile, half-length, 2.3"x 6.6"  |
| Operating system         | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit  |
| Platforms                | Intel® and AMD architectures  |
| Graphics accelerator     | ATI FirePro   |
| Graphics memory          | 512 MB DDR3   |
| Memory interface         | 64-bit  |
| Memory bandwidth         | 9.6 GB/s  |
| Pixel depth              | 32-bit color  |
| Electrical standard      | Single-link DVI complying to v1.0 specification   |
| Connectors               | DMS-59  |
| Connectivity             | One DMS-59-to-DVI adaptor included  |
| Supported resolutions    | 2560×1600 DisplayPort, 1920×1200 DVI/Analog   |
| DirectX support          | Microsoft® DirectX v11.0  |
| OpenGL support           | OpenGL 4.0  |
| OpenCL support           | OpenCL 1.0  |
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI C NS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2002/95/EC), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 60°C (32° to 140° F)  |

## MXRT-1451

|                       |   |
|-----------------------|---|
| Product acronym       | MXRT-1451   |
| Bus compatibility     | Could be installed within PCI Express x1, x8, x16 mechanical slots Works in PCI Express x1, x4, x8, x16 electrical slots / operates at x1 speed |
| Power consumption     | 17 Watt   |
| Form factor           | Low profile, half-length, 2.3"x 6.6"  |
| Operating system      | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit  |
| Platforms             | Intel® and AMD architectures  |
| Graphics accelerator  | ATI FirePro   |
| Graphics memory       | 1 GB DDR3   |
| Memory interface      | 64-bit  |
| Memory bandwidth      | 9.6 GB/s  |
| Pixel depth           | 32-bit color  |
| Electrical standard   | Single-link DVI complying to v1.0 specification   |
| Connectors            | DMS-59  |
| Connectivity          | One DMS-59-to-DVI adaptor included  |
| Supported resolutions | 2560×1600 DisplayPort, 1920×1200 DVI/Analog   |
| DirectX support       | Microsoft® DirectX v11.0  |
| OpenGL support        | OpenGL 4.1  |
| OpenCL support        | OpenCL 1.0  |

|                          |   |
|--------------------------|---|
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC60950-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EC), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 60°C (32° to 140° F)  |

## MXRT-2400

|                          |   |
|--------------------------|---|
| Product acronym          | MXRT-2400   |
| Bus compatibility        | PCIe Gen2 x16   |
| Power consumption        | 43 Watt   |
| Form factor              | 169.67 mm (L) x 64.46 mm (H) single PCIe slot wide  |
| Operating system         | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit  |
| Platforms                | Intel® and AMD architectures  |
| Graphics accelerator     | ATI FirePro   |
| Graphics memory          | 512 MB DDR3   |
| Memory interface         | 64-bit  |
| Memory bandwidth         | 14.4 GB/s   |
| Pixel depth              | 32-bit pixels (supports 8-bit and 10-bit per color channel)   |
| Electrical standard      | Dual-link DVI complying to v1.0, DisplayPort complying to v1.1a   |
| Connectors               | 1- DVI-I, 1- DisplayPort  |
| Connectivity             | One DisplayPort to Single-Link DVI-I dongle included  |
| Supported resolutions    | Up to 6MP color   |
| DirectX support          | Microsoft® DirectX v11.0, Vertex Shader 5.0, Pixel Shader 5.0   |
| OpenGL support           | OpenGL 4.0  |
| OpenCL support           | OpenCL 1.0  |
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI C NS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2002/95/EC), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 60°C (32° to 140° F)  |

## MXRT-2500

|                       |  |
|-----------------------|--|
| Product acronym       | MXRT-2500  |
| Bus compatibility     | PCIe Gen2 x16  |
| Power consumption     | 50 Watt  |
| Form factor           | 168 mm (L) x 68 mm (H) single PCIe slot wide                     |
| Operating system      | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit |
| Platforms             | Intel® and AMD architectures                                     |
| Graphics accelerator  | ATI FirePro  |
| Graphics memory       | 1 GB DDR3  |
| Memory interface      | 128-bit  |
| Memory bandwidth      | 28.8 GB/s  |
| Pixel depth           | 32-bit pixels (supports 8-bit and 10-bit per color channel)      |
| Electrical standard   | Dual-link DVI complying to v1.0, DisplayPort complying to v1.2   |
| Connectors            | 1- DVI-I, 1- DisplayPort   |
| Connectivity          | One DisplayPort to Single-Link DVI-I dongle included             |
| Supported resolutions | Up to 6MP color  |
| DirectX support       | Microsoft® DirectX v11.0, Vertex Shader 5.0, Pixel Shader 5.0    |
| OpenGL support        | OpenGL 4.2   |

|                          |  |
|--------------------------|--|
| OpenCL support           | OpenCL 1.1   |
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EC), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 60°C (32° to 140° F)   |

## MXRT-2600

|                          |   |
|--------------------------|---|
| Product acronym          | MXRT-2600   |
| Bus compatibility        | PCIe Gen3 x16 (wired x8)  |
| Power consumption        | 26 Watt   |
| Form factor              | 168 mm (L) x 68 mm (H) single PCIe slot wide  |
| Operating system         | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit  |
| Platforms                | Intel® and AMD architectures  |
| Graphics accelerator     | ATI FirePro   |
| Graphics memory          | 2 GB DDR3   |
| Memory interface         | 128-bit   |
| Memory bandwidth         | 28.8 GB/s   |
| Pixel depth              | 32-bit pixels (supports 8-bit and 10-bit per color channel)   |
| Electrical standard      | DisplayPort complying to v1.2a  |
| Connectors               | 2- DisplayPort  |
| Supported resolutions    | Up to 6MP color   |
| DirectX support          | Microsoft® DirectX v11.2, Vertex Shader 5.0, Pixel Shader 5.0   |
| OpenGL support           | OpenGL 4.4  |
| OpenCL support           | OpenCL 1.2  |
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI C NS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EU), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 55°C (32° to 131° F)  |

## MXRT-4500



|                       |   |
|-----------------------|---|
| Product acronym       | MXRT-4500   |
| Bus compatibility     | PCIe Gen2.1 x16   |
| Power consumption     | 75 Watt   |
| Form factor           | 163 mm (L) x 97 mm (H) single PCIe slot wide  |
| Operating system      | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit  |
| Platforms             | Intel® and AMD architectures  |
| Graphics accelerator  | ATI FirePro   |
| Graphics memory       | 1 GB GDDR5  |
| Memory interface      | 128-bit   |
| Memory bandwidth      | 64 GB/s   |
| Pixel depth           | 32-bit pixels (supports 8-bit and 10-bit per color channel)   |
| Electrical standard   | Dual-link DVI complying to v1.0, DisplayPort complying to v1.2  |
| Connectors            | 1- DVI-I, 2- DisplayPort  |
| Connectivity          | Two DisplayPort to Single-Link DVI-I dongles included DisplayPort to Dual- Link DVI-I dongle available from Barco; Part Number K9305104 |
| Supported resolutions | Up to 6MP color and 10MP grayscale  |
| DirectX support       | Microsoft® DirectX v11.0, Vertex Shader 5.0, Pixel Shader 5.0   |

|                          |  |
|--------------------------|--|
| OpenGL support           | OpenGL 4.1   |
| OpenCL support           | OpenCL 1.1   |
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EC), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 60°C (32° to 140° F)   |

## MXRT-5400

|                          |  |
|--------------------------|--|
| Product acronym          | MXRT-5400  |
| Bus compatibility        | PCIe Gen2 x16  |
| Power consumption        | 72 Watt  |
| Form factor              | 230.53 mm (L) x 98.34 mm (H) single PCIe slot wide   |
| Operating system         | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit   |
| Platforms                | Intel® and AMD architectures   |
| Graphics accelerator     | ATI FirePro  |
| Graphics memory          | 1 GB GDDR5   |
| Memory interface         | 128-bit  |
| Memory bandwidth         | 64 GB/s  |
| Pixel depth              | 32-bit pixels (supports 8-bit and 10-bit per color channel)  |
| Electrical standard      | Dual-link DVI complying to v1.0, DisplayPort complying to v1.1a  |
| Connectors               | 1- DVI-I, 2- DisplayPort   |
| Connectivity             | Two DisplayPort to Single-Link DVI-I dongles included DisplayPort to Dual- Link DVI-I dongle available from Barco; Part Number K9305104  |
| Supported resolutions    | Up to 6MP color and 10MP grayscale   |
| DirectX support          | Microsoft® DirectX v11.0, Vertex Shader 5.0, Pixel Shader 5.0  |
| OpenGL support           | OpenGL 4.0   |
| OpenCL support           | OpenCL 1.0   |
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2002/95/EC), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 60°C (32° to 140° F)   |

## MXRT-5450

|                      |  |
|----------------------|--|
| Product acronym      | MXRT-5450  |
| Bus compatibility    | PCIe Gen2 x16  |
| Power consumption    | 72 Watt  |
| Form factor          | 230.53 mm (L) x 98.34 mm (H) single PCIe slot wide               |
| Operating system     | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit |
| Platforms            | Intel® and AMD architectures                                     |
| Graphics accelerator | ATI FirePro  |
| Graphics memory      | 1 GB GDDR5   |
| Memory interface     | 128-bit  |
| Memory bandwidth     | 64 GB/s  |
| Pixel depth          | 32-bit pixels (supports 8-bit and 10-bit per color channel)      |
| Electrical standard  | Dual-link DVI complying to v1.0                                  |
| Connectors           | 2- DVI-I   |

|                          |   |
|--------------------------|---|
| Connectivity             | Dual-Link DVI Cables  |
| Supported resolutions    | Up to 6MP color and 10MP grayscale  |
| DirectX support          | Microsoft® DirectX v11.0, Vertex Shader 5.0, Pixel Shader 5.0   |
| OpenGL support           | OpenGL 4.0  |
| OpenCL support           | OpenCL 1.0  |
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI C NS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2002/95/EC), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 60°C (32° to 140° F)  |

## MXRT-5500

|                          |  |
|--------------------------|--|
| Product acronym          | MXRT-5500  |
| Bus compatibility        | PCIe Gen3 x16  |
| Power consumption        | 75 Watt  |
| Form factor              | 184.15 mm (L) x 111 mm (H) single PCIe slot wide   |
| Operating system         | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit   |
| Platforms                | Intel® and AMD architectures   |
| Graphics accelerator     | ATI FirePro  |
| Graphics memory          | 2 GB GDDR5   |
| Memory interface         | 256-bit  |
| Memory bandwidth         | 102.3 GB/s   |
| Pixel depth              | 32-bit pixels (supports 8-bit and 10-bit per color channel)  |
| Electrical standard      | Dual-link DVI complying to v1.0, DisplayPort complying to v1.2   |
| Connectors               | 1- DVI-I, 2- DisplayPort   |
| Connectivity             | Two DisplayPort to Single-Link DVI-I dongles included DisplayPort to Dual-Link DVI-I dongle available from Barco; Part Number K9305104   |
| Supported resolutions    | Up to 6MP color and 10MP grayscale   |
| DirectX support          | Microsoft® DirectX v11.1, Vertex Shader 5.0, Pixel Shader 5.0  |
| OpenGL support           | OpenGL 4.2   |
| OpenCL support           | OpenCL 1.2   |
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EC), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 60°C (32° to 140° F)   |

## MXRT-5550

|                      |  |
|----------------------|--|
| Product acronym      | MXRT-5550  |
| Bus compatibility    | PCIe Gen3 x16  |
| Power consumption    | 75 Watt  |
| Form factor          | 184.15 mm (L) x 111 mm (H) single PCIe slot wide                 |
| Operating system     | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit |
| Platforms            | Intel® and AMD architectures                                     |
| Graphics accelerator | ATI FirePro  |
| Graphics memory      | 2 GB GDDR5   |
| Memory interface     | 256-bit  |
| Memory bandwidth     | 102.3 GB/s   |

|                          |  |
|--------------------------|--|
| Pixel depth              | 32-bit pixels (supports 8-bit and 10-bit per color channel)  |
| Electrical standard      | Dual-link DVI complying to v1.0  |
| Connectors               | 2- DVI-I   |
| Connectivity             | Dual-Link DVI Cables   |
| Supported resolutions    | Up to 6MP color and 10MP grayscale   |
| DirectX support          | Microsoft® DirectX v11.1, Vertex Shader 5.0, Pixel Shader 5.0  |
| OpenGL support           | OpenGL 4.2   |
| OpenCL support           | OpenCL 1.2   |
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EC), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 60°C (32° to 140° F)   |

## MXRT-5600

|                          |   |
|--------------------------|---|
| Product acronym          | MXRT-5600   |
| Bus compatibility        | PCIe Gen3 x16   |
| Power consumption        | 75 Watt   |
| Form factor              | 172 mm (L) x 110 mm (H) single PCIe slot wide   |
| Operating system         | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit  |
| Platforms                | Intel® and AMD architectures  |
| Graphics accelerator     | ATI FirePro   |
| Graphics memory          | 4 GB GDDR5  |
| Memory interface         | 128-bit   |
| Memory bandwidth         | 96 GB/s   |
| Pixel depth              | 32-bit pixels (supports 8-bit and 10-bit per color channel)   |
| Electrical standard      | DisplayPort complying to v1.2a  |
| Connectors               | 4- DisplayPort  |
| Supported resolutions    | Up to 12MP color and 10MP grayscale   |
| DirectX support          | Microsoft® DirectX v11.2, Vertex Shader 5.0, Pixel Shader 5.0   |
| OpenGL support           | OpenGL 4.4  |
| OpenCL support           | OpenCL 1.2  |
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI C NS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EU), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 45°C (32° to 113° F)  |

## MXRT-7400

|                      |  |
|----------------------|--|
| Product acronym      | MXRT-7400  |
| Bus compatibility    | PCIe Gen2 x16  |
| Power consumption    | 138 Watt   |
| Power connector      | One 2×3 power connector  |
| Form factor          | 281.29 mm (L) x 98.53 mm (H) single PCIe slot wide               |
| Operating system     | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit |
| Platforms            | Intel® and AMD architectures                                     |
| Graphics accelerator | ATI FirePro  |

|                          |  |
|--------------------------|--|
| Graphics memory          | 2 GB GDDR5   |
| Memory interface         | 256-bit  |
| Memory bandwidth         | 128 GB/s   |
| Pixel depth              | 32-bit pixels (supports 8-bit and 10-bit per color channel)  |
| Electrical standard      | Dual-link DVI complying to v1.0, DisplayPort complying to v1.1a  |
| Connectors               | 1- DVI-I, 2- DisplayPort   |
| Connectivity             | Two DisplayPort to Single-Link DVI-I dongles included DisplayPort to Dual- Link DVI-I dongle available from Barco; Part Number K9305104  |
| Supported resolutions    | Up to 6MP color and 10MP grayscale   |
| DirectX support          | Microsoft® DirectX v11.0, Vertex Shader 5.0, Pixel Shader 5.0  |
| OpenGL support           | OpenGL 4.0   |
| OpenCL support           | OpenCL 1.0   |
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2002/95/EC), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 60°C (32° to 140° F)   |

## MXRT-7500

|                          |  |
|--------------------------|--|
| Product acronym          | MXRT-7500  |
| Bus compatibility        | PCIe Gen3 x16  |
| Power consumption        | 140 Watt   |
| Power connector          | One 2×3 power connector  |
| Form factor              | 242 mm (L) x 98.53 mm (H) single PCIe slot wide  |
| Operating system         | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit   |
| Platforms                | Intel® and AMD architectures   |
| Graphics accelerator     | ATI FirePro  |
| Graphics memory          | 4 GB GDDR5   |
| Memory interface         | 256-bit  |
| Memory bandwidth         | 154 GB/s   |
| Pixel depth              | 32-bit pixels (supports 8-bit and 10-bit per color channel)  |
| Electrical standard      | DisplayPort complying to v1.2  |
| Connectors               | 4- DisplayPort   |
| Connectivity             | Two DisplayPort to Single-Link DVI-I dongles included DisplayPort to Dual- Link DVI-I dongle available from Barco; Part Number K9305104  |
| Supported resolutions    | Up to 6MP color and 10MP grayscale   |
| DirectX support          | Microsoft® DirectX v11.1, Vertex Shader 5.0, Pixel Shader 5.0  |
| OpenGL support           | OpenGL 4.2   |
| OpenCL support           | OpenCL 1.2   |
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI CNS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EC), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 60°C (32° to 140° F)   |

## MXRT-7600

|                   |               |
|-------------------|---------------|
| Product acronym   | MXRT-7600     |
| Bus compatibility | PCIe Gen3 x16 |



|                          |   |
|--------------------------|---|
| Power consumption        | 150 Watt  |
| Power connector          | One 2x3 power connector   |
| Form factor              | 248 mm (L) x 110 mm (H) single PCIe slot wide   |
| Operating system         | Windows 7 – 32/64-bit, Windows 8.1 – 64-bit, Windows 10 – 64-bit  |
| Platforms                | Intel® and AMD architectures  |
| Graphics accelerator     | ATI FirePro   |
| Graphics memory          | 8 GB GDDR5  |
| Memory interface         | 256-bit   |
| Memory bandwidth         | 160 GB/s  |
| Pixel depth              | 32-bit pixels (supports 8-bit and 10-bit per color channel)   |
| Electrical standard      | DisplayPort complying to v1.2   |
| Connectors               | 4- DisplayPort  |
| Supported resolutions    | Up to 12MP color and 10MP grayscale   |
| DirectX support          | Microsoft® DirectX v11.2, Vertex Shader 5.0, Pixel Shader 5.0   |
| OpenGL support           | OpenGL 4.4  |
| OpenCL support           | OpenCL 1.2  |
| Approvals and compliance | FCC Part 15 Class B, CE EN 55022 Limit B, EN 55024, UL-60950-1, BMSI C NS, CISPR- 22/24, IEC609050-1, VCCI, CSA C22.2, EU RoHS directive (2011/65/EU), Certificate of Information & Communication Equipment (Republic of Korea) |
| Operating temperature    | 0° to 45°C (32° to 113° F)  |



**K5905271-13 | 19/06/2018**


Registered address: Barco NV | President Kennedypark 35, 8500 Kortrijk, Belgium

Barco NV | Beneluxpark 21, 8500 Kortrijk, Belgium




[www.barco.com](http://www.barco.com)

---

## Documents / Resources

|   |  |
|---|--|
|  | <p><b><a href="#">BARCO MXRT-7500 Display Controller</a></b> [pdf] User Guide<br/> MXRT-7500 Display Controller, MXRT-7500, Display Controller, Controller</p> |
|---|--|

## References

-  [Manual-Hub.com - Free PDF manuals!](#)
-  [Inspired sight and sharing solutions - Barco](#)
-  [Privacy policy - Barco](#)
-  [Support - Barco](#)
-  [Manual-Hub.com - Free PDF manuals!](#)
- [User Manual](#)

### Manuals+. [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.