



## Bellman Classic Alarm Clock BE1350 User Manual

[Home](#) » [Bellman](#) » Bellman Classic Alarm Clock BE1350 User Manual 

### Contents

- 1 User Manual
  - 1.1 Bellman Classic Alarm Clock BE1350
- 2 Buttons and Controls
- 3 Technical Specifications
- 4 Getting Started
- 5 Setting the Time
- 6 Setting the Alarm
- 7 Activating the Alarm
- 8 Using the Snooze
- 9 Settings
- 10 Using Battery Backup
- 11 Troubleshooting
- 12 Related Posts


**User Manual**



**Bellman Classic Alarm Clock BE1350**

## **Buttons and Controls**





## Technical Specifications

### In the box

- BE1350 Classic alarm clock
- BE1272 Bed shaker with sound
- Power supply
- User manual

### Power and battery

- Mains power  
7.5 V DC / 1000 mA  
External power supply unit
- Backup batteries  
4 x 1.2 V AAA NiMH rechargeable  
batteries (not included)
- Battery backup operating time  
~ 24 h when fully charged
- Battery backup charging time  
~ 10 h from fully depleted

### Dimensions and weight

- Height: 108 mm, 4.3"
- Width: 121 mm, 4.7"
- Depth: 92 mm, 3.6"
- Weight: 300 g, 10.6 oz. excl. batteries

### Alarm and snooze

- Alarm timeout: 20 minutes
- Decreasing snooze time interval  
from 9 to 2 minutes
- Snooze timeout: 20 minutes

### Display

- Type: Twisted Nematic LCD
- Adjustable backlight in 3 steps
- Display timeout: 8 seconds

### Output signals

- Sound  
100 dBA @ 10 cm, 94 dBA @ 30 cm  
Multiple frequencies: 950 Hz - 3 kHz
- Bed shaker  
Vibrator power: 2.0 - 4.0 V DC


### Environmental requirements

- For indoor use only
- Temperature: 32° to 95° F, 0°-35° C
- Relative humidity: 5% -95%  
Non-condensing

### Accessories

- BE1270 Bed shaker


## Getting Started



1. Connect the bed shaker on the back of the alarm clock.
2. Connect the power supply to the alarm clock and the mains outlet to turn it on.
3. Tuck the bed shaker under the pillow or mattress and place the alarm clock on the bedside table. Setting


## Setting the Time

- Press the left dial and turn it to set hours.
- To set minutes, press and turn the dial again.
- Press once again to save your settings.





## Setting the Alarm

- Press the right dial and turn it to set hours.
- To set minutes, press and turn the dial again.
- Press once again to save your settings.




### Activating the Alarm

- Press the alarm on/off button to activate the alarm.
- The alarm symbol  appears on the clock face.
- To turn off the alarm, press the button again.



### Using the Snooze

- Press the snooze button briefly to snooze the alarm.
- The snooze time interval will automatically decrease from 9 to 2 minutes in 2 minute steps.



## Settings

### Sound on/off

Press the sound on/off button marked with  on the back of the alarm clock to turn on/off the sound. A  symbol will appear on the clock face when the sound is muted.

### Display backlight


Press the button marked with  repeatedly to adjust the intensity in 4 steps. Step 5 results in permanent backlight.

### Time format

Press the time format button marked with **12/24h** on the back of the alarm clock to toggle between a 24h and a 12h setting.

## Using Battery Backup

The alarm clock features a battery backup function that steps in when the power goes out. It keeps the alarm clock and bed shaker going for about 24 hours (when the batteries are fully charged).



Here is how you install backup batteries:






1. Remove the battery cover on the bottom of the alarm clock using a Philips screwdriver.
2. Insert 4 x 1.2 V AAA NiMH rechargeable batteries and put the cover back on.



**Warning!** ONLY USE RECHARGEABLE NiMH BATTERIES. Non-rechargeable batteries will start to leak and the battery acid will damage the electronics. The resulting damage is not covered by warranty.

## Troubleshooting

Most problems with the alarm clock can be solved quickly by following the advice below.

Problem	Try this
The alarm clock seems to be turned off	<ul style="list-style-type: none"> <li>▪ Check that the power supply is connected correctly, see <b>Getting started</b>.</li> <li>▪ Charge the backup batteries (not included) for a couple of hours.</li> </ul>
The  symbol on the clock face starts to blink	<ul style="list-style-type: none"> <li>▪ The power supply is disconnected and the backup batteries are nearly depleted. Connect the power supply and charge the batteries for a couple of hours.</li> </ul>
A  symbol appears on the clock face	<ul style="list-style-type: none"> <li>▪ The alarm clock detects no backup batteries. Open the battery cover and fit four 1.2 V AAA NiMH rechargeable batteries.</li> </ul>
A  symbol appears on the clock face	<ul style="list-style-type: none"> <li>▪ The sound is muted. Press the sound on/off button marked with  until the symbol disappears if you wish to turn it back on.</li> </ul>
The display is too bright or too dark or I want it to be permanently lit	<ul style="list-style-type: none"> <li>▪ Press the button marked with  repeatedly to adjust the intensity in 4 steps.</li> <li>▪ Choose step 3 that results in permanent backlight.</li> </ul>
The bed shaker does not vibrate	<ul style="list-style-type: none"> <li>▪ Make sure it's connected properly to the alarm clock, see <b>Getting started</b>.</li> </ul>
The alarm volume is too low	<ul style="list-style-type: none"> <li>▪ The volume increases gradually and reaches over 100 dB.</li> </ul>

## Read More About This Manual & Download PDF:

Bellman Classic Alarm Clock BE1350 User Manual – [Optimized PDF](#)  
 Bellman Classic Alarm Clock BE1350 User Manual – [Original PDF](#)

**Questions about your Manual? Post in the comments!**