



Corsano App BP Cuff Device User Guide

[Home](#) » [Corsano](#) » Corsano App BP Cuff Device User Guide 

Contents

- 1 Corsano App BP Cuff Device
- 2 NIBP – Pair a BP cuff device
- 3 NIBP – BP cuff measurement
- 4 NIBP – Messages
- 5 NIBP – Messages
- 6 Specifications
- 7 Frequently Asked Questions
- 8 Documents / Resources
 - 8.1 References


Corsano 

Corsano App BP Cuff Device



NIBP – Pair a BP cuff device


Add a Corsano BP cuff device to the Corsano App



1. Press the right top corner to go to Settings
2. Select BP Cuff
3. Press [Start Pairing] Don't press the Start button of the CuffB monitor
4. Follow instructions. Press "Cancel" to abort the pairing process and try again.
5. BP cuff pairing success. Make sure this is your device by visually checking the identifier number on the Cuff monitor.


NIBP – BP cuff measurement

- Follow these steps to perform a BP cuff measurement – Instructions



1. Press [Start BP Measurement] from the Dashboard.
2. Tick the box if you take medication to lower Blood Pressure.
3. Follow the instructions to properly adjust the tightness of the bracelet.
4. Follow the instructions to properly adjust the position of the bracelet.
5. Prepare for the Blood Pressure Cuff measurement. Sit and follow the instructions.
6. Rest your arm on a table at the same level as your heart.
7. Check the tightness of the Cuff. Sit still, don't move and relax.


Follow these steps to perform a BP cuff measurement – Measurement




1. Press the [START] button on the Cuff monitor. The Cuff will inflate and measure your Blood Pressure. Wait for the measurement to finish. This takes about 30s.
2. After completion of the third measurement, the Non-Invasive Blood Pressure is calibrated. The result of the calibration is available in the APP.

NIBP – Messages

Initial setup, day 1 – No BP cuff init





**Good Morning,
Matthias**




- Alert message on top of the dashboard
- Detail message with init status in the BP cell

NIBP – Messages


Initial setup – One or more BP cuff init done



**Good Morning,
Matthias**




- Alert message on top of the dashboard, with details



- Detail message with init status in the BP cell. First init measurement done. Continuous Blood Pressure will start coming.

Initial setup,- One or more BP cuff init done, with error.



- Alert message on top of the dashboard, With details



- Detail message with init status in the BP cell. First init measurement done. Second init measurement with error. Continuous Blood Pressure start coming.

Specifications

- **Device:** Corsano BP Cuff
- **App:** Corsano App
- **Measurement Type:** Blood Pressure
- **Measurement Method:** Non-Invasive

Frequently Asked Questions

Initial Setup, Day 1 – No BP Cuff Init

- Alert message on top of the dashboard.
- Detailed message with initialization status in the BP cell.

Initial Setup – One or More BP Cuff Init Done

- Alert message on top of the dashboard with details.
- Detailed message with initialization status in the BP cell. First init measurement done.
- Continuous Blood Pressure readings will start appearing.

Initial Setup – One or More BP Cuff Init Done, with Error

- Alert message on top of the dashboard with details.
- Detailed message with initialization status in the BP cell. First init measurement done.
- Second init measurement with error. Continuous Blood Pressure readings will start appearing.

Documents / Resources



[Corsano Corsano App BP Cuff Device](#) [pdf] User Guide
Corsano App BP Cuff Device, App BP Cuff Device, Cuff Device, Device

References

- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.