

# Dino-Lite AF7115MZTW Edge Digital Microscope User Guide

Home » Dino-Lite » Dino-Lite AF7115MZTW Edge Digital Microscope User Guide The Dino-Lite » Dino-Lite »



### **Contents**

- 1 Dino-Lite AF7115MZTW Edge Digital Microscope
- **2 Product Usage Instructions**
- 3 AF7115MZTW Edge
- 4 Documents / Resources
- 4.1 References
- **5 Related Posts**

# Dino-Lite

**Dino-Lite AF7115MZTW Edge Digital Microscope** 



# **Specifications**

• Brand: Dino-Lite

• Product Code: AF7115MZTW – Edge

Lighting: White LEDNumber of LEDs: 8

• LED on/off switchable: Yes

• Diffuser available: Yes

Emission filter: NoPolarizer: Yes, linear

Magnification (T Zone): 15X~50X
 Magnification (W Zone): 5X~30X

• Working distance: 15 mm (T Zone), 30 mm (W Zone)

Lens type: Macro zoomSensor type: CMOS

• Sensor resolution: 5 Megapixel (2592×1944)

• Maximum frame rate: 30 fps

• Compatibility: Windows 7, 8, 10 & 11, MacOS 10.14 and up

• Interface: USB 2.0

• Included software: DinoCapture 2.0 (Windows), DinoXcope (Mac OS)

• Supported image formats (Windows): BMP, GIF, PNG, MNG, TIF, TGA, PCX, WBMP, JP2, JPC, JPG, PGX, RAS, PNM

• Supported video formats (Windows): WMV, FLV, SWF

• Supported image formats (MacOS): JPEG, PNG

• Supported video formats (MacOS): MOV (max 1.3MP)

• Imaging standards: DirectShow, UVC

• Wifi: Wireless-ready, requires the WF-20 adapter (optional)

· Housing material: Metal housing

• Magnification lock: Yes

• **Dimensions**: 15 mm (H) x 34 mm (D)

• Weight: 97g

• Cable length: 1.8m

• Features: Flexible LED Control (FLC), Microtouch sensor, ESD safe

• Special feature: Measurement Calibration

Microscope package contents: Microscope, carry pouch, software CD, calibration target, user manual, N3C-O-Open cap, N3C-C Closed cap, N3C-D- diffuser cap, N3C-E extension cap, N3C-L- Long cap, N3C-S- Side light cap

• Warranty information: 2 years European warranty

• Regulatory approval: CE, FCC, ROHS

• Price range: €900.00 – €1100.00

# **Product Usage Instructions**

### **Setup and Installation**

- 1. Unbox the Dino-Lite AF7115MZTW Edge microscope and remove all packaging materials.
- 2. Connect the microscope to your computer using the provided USB cable.
- 3. Install the included software, DinoCapture 2.0 for Windows or DinoXcope for Mac OS, by following the onscreen instructions.
- 4. Once the software is installed, launch it and connect the microscope by selecting the appropriate device from the software's device list.
- 5. Ensure that the microscope is securely connected and recognized by the software.

### **Basic Operation**

To use the Dino-Lite AF7115MZTW – Edge microscope, follow these steps

- 1. Place the object you want to examine under the microscope lens.
- 2. Adjust the focus by rotating the focus wheel located on the microscope body.
- 3. Use the zoom control to adjust the magnification level according to your needs.
- 4. Toggle the LED lights on or off using the LED on/off switch.
- 5. If necessary, attach the diffuser cap to diffuse the light for better illumination.
- 6. If reflections are causing issues, adjust the polarizer to reduce them.

# **Advanced Features**

The Dino-Lite AF7115MZTW - Edge microscope offers several advanced features

- Flexible LED Control (FLC): Allows you to adjust the brightness of the LED lights for optimal illumination.
- **Measurement Calibration:** The microscope supports measurement calibration, enabling you to measure dimensions of objects directly on the captured images.
- **Microtouch Sensor:** The microscope is equipped with a Micro touch sensor that allows you to capture images with a single touch.

• **ESD Safe:** The microscope is designed to be electrostatic discharge (ESD) safe, ensuring protection against static electricity.

# Frequently Asked Questions (FAQ)

• Q: Can I use the Dino-Lite AF7115MZTW – Edge microscope with my Mac computer?

A: Yes, the microscope is compatible with Mac computers running MacOS 10.14 and up. You can use the DinoXcope software provided for Mac OS.

Q: What are the supported image and video formats for Windows?

A: The microscope supports various image formats, including BMP, GIF, PNG, MNG, TIF, TGA, PCX, WBMP, JP2, JPC, JPG, PGX, RAS, PNM. For video formats, it supports WMV, FLV, and SWF.

Q: Does the microscope have wireless capabilities?

A: Yes, the Dino-Lite AF7115MZTW – Edge is wireless-ready. However, to enable wireless functionality, you need to purchase and use the WF-20 adapter (optional).

Q: What is the warranty period for the microscope?

A: The microscope comes with a 2-year European warranty.

# AF7115MZTW - Edge

• Brand: Dino-Lite


• Product Code: AF7115MZTW - Edge

# **Short Description**

- · High Optical Resolution
- 5.0 Megapixels
- Flexible LED Control (FLC)
- Adjustable Polarizer
- · Macro zoom optics

# **Description**

The 5MP Edge AM7115MZTW – Edge provides macro-level imaging with zoomin capabilities. This microscope includes a polarizer to reduce reflections and features FLC for adjustable illumination. AF7115MZTW is compatible with the optional WF-20 Wi-Fi streamer.



# The main features of the AM7115MZTW - Edge are

- 1. 5 Megapixel sensor
- 2. Macro zoom optics
- 3. Adjustable polarizer
- 4. Robust Metal housing
- 5. Flexible LED Control (FLC)
- 6. High Optical Resolution
- · Adaptable Interface
- · And more...

# Working distance/field of view/depth of field

MAGNIFICATI WORKING FIELD OF FIELD OF DEPTH OF
ON RATE DISTANCE VIEW(X) VIEW(Y) FIELD
Listed values may differ slightly Without front Unit = mm
cap

# T Zone

50	23	8	6	4
40	43	10	7	6
30	43	13	10	11
20	128	20	15	27
15	183	26	20	48
MAGNIFICATI	WORKING	FIELD OF	IELD OF	DEPTH OF
ON RATE	DISTANCE*	VIEW(X)	VIEW(Y)	FIELD

MAGNIFICATI WORKING FIELD OF FIELD OF DEPTH OF

сар

# W Zone

30	19		13	10	7	
20	33		20	15	7	
15	47		26	20	25	
10	76		40	30	57	
5						
Specificatio n	163		79	59	297	
Lighting	Lighting					
Light/ LED type		White				
Number of LEDs		8				
LED on/off switchable		Yes				
Diffuser available		Yes (N3C-D included)				
Emission filter		No				
Polarizer		Yes, linear				
Optics						
Magnification		T 15X~50X, W 5X~30X				
Macro zoom		No				
Working distance		Standard				

Lens type	Glass with anti-reflection coating			
Sensor				
Sensor type	CMOS			
Resolution	5 Megapixel (2592×1944)			
Maximum frame rate	30 fps			
Compatibility				
Interface	USB 2.0			
Operating system	Windows 7, 8, 10 & 11, MacOS 10.14 and up			
Included software	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)			

Supported image formats (Windows)	BMP, GIF, PNG, MNG, TIF, TGA, PCX, WBMP, JP2, JPC, JPG, PGX, R AS, PNM			
Supported video formats (Windows)	WMV, FLV, SWF			
Supported image formats (MacOS)	JPEG, PNG			
Supported video formats (MacOS)	MOV (max 1.3MP)			
Imaging standards	DirectShow, UVC			
Wifi	Wireless-ready, requires the WF-20 adapter (optional)			
Housing				
Housing material	Metal housing			
Magnification lock	Yes			
Dimensions	15 mm (H) x 34 mm (D)			
Weight	97g			
Cable length	1.8m			
Features				
Special feature	Flexible LED Control (FLC)			
Measurement	Yes			
Calibration	Yes			
Micro touch sensor	Yes			
ESD safe	Yes			
Information				
Package contents	Microscope, carry pouch, software CD, calibration target, user manual, N3C-O- Open cap, N3C-C- Closed cap, N3C-D- diffuser cap, N3C-E- extension cap, N3C-L- Long cap, N3C-S- Side light cap			
Warranty information	2 years European warranty			
Regulatory approval	CE, FCC, ROHS			
Price range	€900,00 - €1100,00			



# <u>Dino-Lite AF7115MZTW Edge Digital Microscope</u> [pdf] User Guide

AF7115MZTW Edge Digital Microscope, AF7115MZTW, Edge Digital Microscope, Digital Microscope, Microscope

# References

- TCPDF
- User Manual

Manuals+, Privacy Policy