

## DRYEASY Bedwetting Alarm with Volume Control



# DRYEASY Bedwetting Alarm with Volume Control Instruction Manual

[Home](#) » [DRYEASY](#) » DRYEASY Bedwetting Alarm with Volume Control Instruction Manual 

### Contents

- 1 DRYEASY Bedwetting Alarm with Volume Control
- 2 Specifications
- 3 Product Usage Instructions
- 4 Steps for Usage
- 5 Bedwetting Alarm Instructions
- 6 How to use
- 7 Cleaning
- 8 One-year warranty
- 9 Instructions for using the armband
- 10 FAQ
- 11 Documents / Resources
  - 11.1 References

# DRYEASY™

DRYEASY Bedwetting Alarm with Volume Control



## Specifications

- **Age Recommendation:** Not intended for children under 5 years of age
- **Usage:** Should be used under adult supervision
- **Cord Safety:** Ensure the cord does not become tangled around the user

## Product Usage Instructions


### Safety Precautions

1. This product is not suitable for children under 5 years of age. Always use it under adult supervision.
2. Make sure the cord is not tangled around the user to prevent accidents.
3. Read the enclosed manual thoroughly before using the product.
4. If you have any doubts about using the product, consult your doctor.

### Steps for Usage

1. Unpack the product and read the manual carefully.
2. Ensure the cord is untangled and in good condition.
3. Follow the instructions provided in the manual for operating the product.
4. If you experience any discomfort or unusual symptoms, stop using the product immediately and seek medical advice.

## Bedwetting Alarm Instructions



## How to use

### A. To set up the alarm:

1. Remove the insulation tag from the battery compartment (for first-time use only).
2. Insert the plug of the cord into the alarm socket.
3. Test the alarm by touching the sensor tip (metal prong) with a wet finger. Press the "Reset" button to stop the alarm from sounding.

### B. To adjust the volume level, select operation mode and sound:

1. Open the battery compartment door and press the "+" or "-" key for volume adjustment.
2. Slide the switch inside the battery compartment for a change of operation mode from "vibration only" to "vibration and sound" and "sound only"
3. To select a desired sound, choose a sound group (A or B) by sliding the sound group selection switch and then choose a sound (1, 2 or 3) from the sound selection switch (6 sounds available for selection. A1, A2, A3, B1, B2 & B3).
4. To play the sounds at random, slide the sound group selection switch to R.

### C. To attach the alarm and the sensor:

1. Attach the alarm to the pyjamas by following the steps set out below (see Figure A):
  1. Lift up the lever
  2. Place the pyjamas into the mouth of the alarm clip
  3. Lock the alarm clip by pressing down the alarm lever

**(Note:** For best results, fix the alarm close to the collar bone)
2. Run the cord under the pyjamas as it can prevent the sensor from being dislodged (see Figure B).




figure A



figure B

3. Attach the sensor to the outside of close-fitting underpants in the following manner:


(**Tip:** Pull down the underpants slightly for easy fixing of sensor )



A) Grasp the area which is most likely to get wet first and hook it onto the cover.



B) Push the sensor towards the cover, so that the underpants are caught in between the cover and the sensor.



C) Press down the tail of the cover to lock.

**Important:** The sensor should be attached with the cover facing up (as shown below)



(**Note:** To detach the sensor from the underpants, simply lift up the tail of the cover)

#### **D. To stop the alarm from sounding when bedwetting occurs:**

1. Remove the wet material from the sensor.
2. Press the “Reset” button on the alarm.

#### **Cleaning**

<b>Alarm</b>	<b>Sensor</b>
1. Use a dry cloth to wipe the surface. 2. DO NOT immerse in water.	1. Detach the cord from the alarm before cleaning the sensor. 2. Wash the sensor with soapy water, rinse it and dry with cloth. ( Note: The sensor should be cleaned regularly so as to maintain its sensitivity. )

#### **Useful tips**

1. Parents are recommended to have a trial run of all steps for the use of the alarm during the time with the child before the first time use.
2. The child must wear a pair of close-fitting underwear as the urine can be detected faster. DO NOT wear loose underwear or boxer shorts.
3. Rapid response to the alarm alert can enhance the success rate in overcoming bedwetting. Parents play an important role, in particular at the initial stage, in awakening the child as soon as the alarm goes off and are therefore recommended to sleep in the same room with the child for the first 2 weeks of using the alarm.
4. Use the stickers and progress chart to track the progress.
5. Typical treatment takes 6-12 weeks but it varies from person to person. Motivation and persistence are very important for success.
6. Keep using the alarm for 2 more weeks after dryness to lower the chance of relapse.

#### **One-year warranty**

We offer you a one-year warranty for manufacturer’s defects. For details, please visit our website or contact us by email.

#### **Instructions for using the armband**

For small to medium arm size (from 15 cm to 27 cm)

1. Clip the alarm near the buckle
2. Thread the strap through the buckle to form a loop
3. Slip onto the arm and tighten the strap until it is comfortable (it should be snug enough to keep the alarm in place, but not so tight that it restricts arm circulation).

**For instructions in other languages, please visit the link below or scan the QR code as shown:**

[www.bedwettingsolving.com/en/page.php?kev=info](http://www.bedwettingsolving.com/en/page.php?kev=info)



If you have any questions about our products or warranty, please email to: [info@bedwettingsolving.com](mailto:info@bedwettingsolving.com)  
We will give you a reply within 24 hours.

### **Satisfied with our product?**

Please take a couple of minutes to give us a product review on Amazon.

### **To submit a review, please follow the steps below:**

- Log in to your Amazon account
- Click “Your Orders”
- Click “Write a product review”

Thank you for your support!

### **Attention**

1. The alarm can be attached to the pyjama top or to the arm with the use of an armband.



DO NOT attach the alarm to the underwear



2. The proper way to stop the alarm from going off:

- Remove the sensor from the wet underwear
- Press the “Reset” button on the alarm

To watch the video, visit the link below or scan the QR code as shown

[www.bedwettingsolving.com/en/video2.mp4](http://www.bedwettingsolving.com/en/video2.mp4)

3. Do not unplug the cord from the alarm except for sensor- cleaning. To unplug the cord, press the tab and pull it out. Never pull the cord out by force as it may get damaged.
4. Do not wet the plug when cleaning the sensor. If it gets wet, dry it before use.
5. The alarm is NOT waterproof! Never immerse it in water or urine.

### **Caution**

1. This product is NOT intended to be used for children under 5 years of age and should be used under adult supervision.

## FAQ

Q: Can children under 5 years old use this product?

A: No, this product is not intended for children under 5 years of age. Always supervise children when using this product.

Q: What should I do if the cord becomes tangled?


A: Stop using the product immediately and carefully untangle the cord to avoid any accidents or injuries.

Q: How often should I consult a doctor regarding the usage of this product?

A: If you have any doubts or concerns about using this product, it is recommended to consult your doctor for guidance and advice.

---

## Documents / Resources

	<p><a href="#">DRYEASY Bedwetting Alarm with Volume Control</a> [pdf] Instruction Manual Bedwetting Alarm with Volume Control, Alarm with Volume Control, Volume Control, Control</p>
---	---

## References

- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.