



# dynamic BIOSENSORS AS-1-Rb v5.1 Adapter Strand User Manual

[Home](#) » [dynamic BIOSENSORS](#) » dynamic BIOSENSORS AS-1-Rb v5.1 Adapter Strand User Manual 

## Contents

- [1 dynamic BIOSENSORS AS-1-Rb v5.1 Adapter Strand](#)
- [2 Product Usage Instructions](#)
- [3 Key Features](#)
- [4 Product Description](#)
- [5 Preparation | MIX&RUN](#)
- [6 Documents / Resources](#)
  - [6.1 References](#)
- [7 Related Posts](#)

**dynamic**  
BIOSENSORS

**dynamic BIOSENSORS AS-1-Rb v5.1 Adapter Strand**



## Specifications

- **Product Name:** heliX+ ADAPTER STRAND 1 with red dye Rb
- **Manufacturer:** Dynamic Biosensors GmbH & Inc.
- **Order Number:** AS-1-Rb
- **Key Features:** 2 spots with 2 different anchor sequences for DNA-encoded addressing
- **Concentration:** 400 nM
- **Storage:** Black cap, limited shelf life, avoid freeze-thaw cycles, aliquot the nanolever

## Product Usage Instructions

### Preparation | MIX&RUN

1. Mix Adapter strand 1 – Rb (400 nM) and conjugated Ligand strand (500 nM) at 1:1 ratio (v/v).
2. Incubate the solution from step 1 at room temperature at 600 rpm for 30 minutes for complete hybridization.
3. Mix the solution from step 2 with Adapter strand 1 – Rb – lfs (200 nM) at 1:1 ratio (v/v).
4. The solution is now ready for biochip functionalization.

**Note:** Stability of the solution is dependent on the stability of the ligand molecule.

## FAQs

- **Q: What is the concentration of Adapter strand 1 – Rb – lfs?**  
A: The concentration is 200/250 nM as per the product description.
- **Q: Can Adapter strand 2 be used interchangeably with Adapter strand 1?**  
A: No, Adapter strand 2 with red dye Rb is prehybridized with ligand-free strand and has a different order number (HK-NHS-1 AS-2-Rb-lfs).

## ADAPTER STRAND 1

with red dye Rb

Dynamic Biosensors GmbH & Inc.


AS-1-Rb v5.1

## Key Features

- Adapter strand 1 – Rb for functionalization of heliX® Adapter Chip Spot 1.
- Compatible with heliX® Adapter Chip.
- Includes Adapter strands for 50 regenerations.
- Ideal for MIX&RUN sample preparation.
- Adapter strand 1 carries a moderately hydrophilic red dye (Rb) with a neutral net charge.

## heliX® Adapter Chip Overview

2 spots with 2 different anchor sequences for DNA-encoded addressing.



## Product Description

Order Number: AS-1-Rb

**Table 1. Contents and Storage Information**

Material	Cap	Concentration	Amount	Buffer	Storage
<b>Adapter strand 1 – Rb</b>	Black	400 nM	5 x 100 µL	TE40 <a href="#">[1]</a>	-20°C


**For research use only.**

- This product has a limited shelf life, please see expiry date on label.
- To avoid many freeze thaw cycles please aliquot the nanolever

## Preparation | MIX&RUN

**In-solution hybridization of adapter and ligand strands:**

- 1. Mix Adapter strand 1 – Rb (400 nM) and conjugated Ligand strand (500 nM) at 1:1 ratio (v/v).
- 2. Incubate the solution of step 1 at RT at 600 rpm for 30 min to ensure complete hybridization.
- 3. Mix solution of step 2 and Adapter strand 1 – Rb – lfs (200 nM) at 1:1 ratio (v/v). Solution is ready to use for biochip functionalization. Stability of the solution is related to the stability of the ligand molecule



**Table 2. Additional material for functionalization of spot 1 and reference spot 2.**

Material	Concentration	Buffer	Related Product Name	Order No
<b>Ligand strand</b> carrying the conjugated ligand	500 nM	P2] E40 [	<b>heliX®</b> Amine Coupling Kit 1	HK-NHS-1
<b>Adapter strand 1 – Rb – lfs</b>	200/250 nM	TE40 [1 ]	<b>Adapter strand 2</b> with red dye <b>Rb</b> pr ehybridized with <i>ligand-free strand</i>	AS-2-Rb-lf s

**Example**

Required volume for 3 functionalizations: 100 µL with a final concentration of 100 nM.

Vial 1		Vial 2
<b>Adapter strand 1 – Rb</b> (400 nM)	Conjugated <b>Ligand strand</b> (500 nM)	<b>Adapter strand 1 – Rb – lfs</b> (200/250 nM)
25 µL	25 µL	50 µL

After incubation time, mix vial 1 and vial 2 to obtain 100 µL of ready-to-use DNA solution.

## Contact

### Dynamic Biosensors GmbH

- Perchtinger Str. 8/10
- 81379 Munich
- Germany

### Dynamic Biosensors, Inc.

- 300 Trade Center, Suite 1400
- Woburn, MA 01801
- USA

Order Information [order@dynamic-biosensors.com](mailto:order@dynamic-biosensors.com)  
Technical Support [support@dynamic-biosensors.com](mailto:support@dynamic-biosensors.com)


### [www.dynamic-biosensors.com](http://www.dynamic-biosensors.com)

Instruments and chips are engineered and manufactured in Germany. ©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.

1. TE40: 10 mM Tris, 40 mM NaCl, 0.05 % Tween20, 50 µM EDTA, 50 µM EGTA
2. If the protein is not stable in PE40 (TE40, HE40), please check buffer compatibility with the switchSENSE® compatibility sheet.

### [www.dynamic-biosensors.com](http://www.dynamic-biosensors.com)

## Documents / Resources

	<p><a href="#">dynamic BIOSENSORS AS-1-Rb v5.1 Adapter Strand</a> [pdf] User Manual AS-1-Rb, AS-2-Rb-lfs, AS-1-Rb v5.1 Adapter Strand, AS-1-Rb v5.1, Adapter Strand, Strand</p>
---	---

## References

- [📄 HomePage | Biosensors International Ltd](#)
- [🏠 Home - Dynamic Biosensors](#)
- [🏠 Home - Dynamic Biosensors](#)
- [📖 User Manual](#)

### [Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.