



dynamic BIOSENSORS AS-2-Gb-Ifs Adapter Strand User Manual

[Home](#) » [dynamic BIOSENSORS](#) » dynamic BIOSENSORS AS-2-Gb-Ifs Adapter Strand User Manual 

dynamic BIOSENSORS AS-2-Gb-Ifs Adapter Strand



Contents


- [1 Key Features](#)
- [2 Adapter Chip Overview](#)
- [3 Product Description](#)
- [4 Preparation | MIX&RUN](#)
- [5 Example](#)
- [6 Customer Support](#)
- [7 Documents / Resources](#)
 - [7.1 References](#)
- [8 Related Posts](#)

Key Features

- Adapter strand 1 – Gc for functionalization of helix ® Adapter Chip Spot 1.
- Compatible with helix ® Adapter Chip.
- Includes Adapter strands for 50 regenerations.
- Ideal for MIX&RUN sample preparation.
- Adapter strand 1 carries a hydrophobic green dye (Gc) with a neutral net charge.

Adapter Chip Overview

2 spots with 2 different anchor sequences for DNA-encoded addressing.



Product Description

Order Number: AS-1-Gc

Table 1. Contents and Storage Information

Material	Cap	Concentration	Amount	Buffer	Storage
Adapter strand 1 – Gc	Black	400 NM	5 x 100 µL	TE40 [1]	-20°C

For research use only.

This product has a limited shelf life, please see expiry date on label.

To avoid many freeze thaw cycles please aliquot the nanolayer.

Preparation | MIX&RUN

n-solution hybridization of adapter and ligand strands:

1. Mix Adapter strand 1 – Gc (400 NM) and conjugated Ligand strand (500 NM) at 1:1 ratio (v/v).
2. Incubate the solution of step 1 at RT at 600 rpm for 30 min to ensure complete hybridization.
3. Mix solution of step 2 and Adapter strand 1 – Gc – lfs (200 nM) at 1:1 ratio (v/v).

Solution is ready to use for biochip functionalization.

Stability of the solution is related to the stability of the ligand molecule.




Table 2. Additional material for functionalization of spot 1 and reference spot 2.

Material	Concentration	Buffer	Related Product Name	Order No
Ligand strand carrying the conjugated ligand	500 NM	PE40 [helix ® Amine Coupling Kit	HK-NHS-1
Adapter strand 1 – Gc – lfs.	200/250 NM	TE40 [1]	Adapter strand 2 with green dye Gc prehybridized with ligand-free strand	AS-2-Gb-lfs

Example

Required volume for 3 functionalization's: 100 µL with a final concentration of 100 nM.

Vial 1		Vial 2
Adapter strand 1 – Gc (400 NM)	Conjugated Ligand strand (500 NM)	Adapter strand 1 – Gc – lfs (200/250 NM)
25 µ	25 µL	50 µL

After incubation time, mix vial 1 and vial 2 to obtain 100 µL of ready-to-use DNA solution.

Customer Support

Contact

Dynamic Biosensors GmbH
Perchtinger Str. 8/10
81379 Munich

Germany

Dynamic Biosensors, Inc.
300 Trade Center, Suite 1400
Woburn, MA 01801
USA


Order Information order@dynamic-biosensors.com

Technical Support support@dynamic-biosensors.com
www.dynamic-biosensors.com

Instruments and chips are engineered and manufactured in Germany. ©2024 Dynamic Biosensors GmbH |
Dynamic Biosensors, Inc. All rights reserved.



Documents / Resources

	<p>dynamic BIOSENSORS AS-2-Gb-lfs Adapter Strand [pdf] User Manual AS-2-Gb-lfs, AS-2-Gb-lfs Adapter Strand, Adapter Strand, Strand</p>
--	--

References

- [📄 HomePage | Biosensors International Ltd](#)
- [🔗 Home - Dynamic Biosensors](#)
- [🔗 Home - Dynamic Biosensors](#)
- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.