



dynamic BIOSENSORS AS-2-Gc Adapter Strand 2 User Manual

[Home](#) » [dynamic BIOSENSORS](#) » dynamic BIOSENSORS AS-2-Gc Adapter Strand 2 User Manual 



Contents

- [1 dynamic BIOSENSORS AS-2-Gc Adapter Strand 2 User Manual](#)
- [2 Key Features](#)
- [3 Product Description](#)
- [4 Preparation | MIX&RUN](#)
- [5 Contact](#)
- [6 Documents / Resources](#)
 - [6.1 References](#)

dynamic BIOSENSORS AS-2-Gc Adapter Strand 2 User Manual



ADAPTER STRAND 2

with green dye Gc


Dynamic Biosensors GmbH & Inc.
AS-2-Gc v5.1

Key Features

- Adapter strand 2 for functionalization of heliX® Adapter Chip Spot 2.
- Compatible with heliX® Adapter Chip.
- Includes Adapter strands for 50 regenerations.
- Ideal for MIX&RUN sample preparation.
- Adapter strand 2 carries a hydrophobic green dye (Gc) with a neutral net charge.

heliX® Adapter Chip Overview

2 spots with 2 different anchor sequences for DNA-encoded addressing.



Product Description

Order Number: AS-2-Gc

Table 1. Contents and Storage Information

Material	Cap	Concentration	Amount	Buffer	Storage
Adapter strand 2 - Gc	White	400 nM	5 x 100 µL	TE40 ^[1]	-20°C

For research use only.

This product has a limited shelf life, please see expiry date on label.

To avoid many freeze thaw cycles please aliquot the nanolever.

Preparation | MIX&RUN

In-solution hybridization of adapter and ligand strands:


1. Mix **Adapter strand 1 - Gc** (400 nM) and conjugated **Ligand strand** with ligand 1 (500 nM) at 1:1 ratio (v/v).
2. Mix **Adapter strand 2 - Gc** (400 nM) and conjugated **Ligand strand** with ligand 2 (500 nM) at 1:1 ratio (v/v).
3. Incubate separately the two solutions of step 1 and 2 at **RT** at **600 rpm** for **30 min** to ensure complete hybridization.
4. Mix solution of step 1 and 2 at 1:1 ratio (v/v).

Solution is ready to use for biochip functionalization.

Stability of the solution is related to the stability of the ligand molecules.

Table 2. Additional material for functionalization of spot 1 and reference spot 2.

Material	Concentration	Buffer	Related Product Name	Order No
Adapter strand 1 - Gc	400 nM	TE40 ^[1]	Adapter strand 2 with green dye Gc	AS-1-Rb
Ligand strand carrying the conjugated ligand 1	500 nM	PE40 ^[2]	helix ® Amine Coupling Kit 1	HK-NHS-1
Ligand strand carrying the conjugated ligand 2	500 nM	PE40 ^[2]	helix ® Amine Coupling Kit 1	HK-NHS-1



Example

Required volume for 3 functionalizations: 100 µL with a final concentration of 100 nM.

Vial 1		Vial 2	
<i>Adapter strand 1 - Gc</i> (400 nM)	Conjugated <i>Ligand strand</i> with ligand 1 (500 nM)	<i>Adapter strand 2 - Gc</i> (400 nM)	Conjugated <i>Ligand strand</i> + with ligand 2 (500 nM)
25 µL	25 µL	25 µL	25 µL

After incubation time, mix vial 1 and vial 2 to obtain 100 µL of ready-to-use DNA solution.

Contact

Dynamic Biosensors GmbH

Perchtinger Str. 8/10
81379 Munich
Germany

Dynamic Biosensors, Inc.

300 Trade Center, Suite 1400
Woburn, MA 01801
USA

Order Information order@dynamic-biosensors.com
Technical Support support@dynamic-biosensors.com

www.dynamic-biosensors.com

Instruments and chips are engineered and manufactured in Germany.
©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.

[1] TE40: 10 mM Tris, 40 mM NaCl, 0.05 % Tween20, 50 µM EDTA, 50 µM EGTA

[2] If the protein is not stable in PE40 (TE40, HE40), please check buffer compatibility with the switchSENSE® compatibility sheet.

www.dynamic-biosensors.com

Read More About This Manual & Download PDF:

Documents / Resources



[dynamic BIOSENSORS AS-2-Gc Adapter Strand 2](#) [pdf] User Manual
AS-2-Gc v5.1, AS-2-Gc Adapter Strand 2, AS-2-Gc, Adapter Strand 2, Strand 2

References

- [📄 HomePage | Biosensors International Ltd](#)
- [📄 Home - Dynamic Biosensors](#)
- [📄 Home - Dynamic Biosensors](#)
- [User Manual](#)

Manuals+, Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.