

dynamic
BIOSENSORS

**AS-2-Rc Fully
Automated
Laboratory
Analysis System**



dynamic BIOSENSORS AS-2-Rc Fully Automated Laboratory Analysis System User Manual

[Home](#) » [dynamic BIOSENSORS](#) » dynamic BIOSENSORS AS-2-Rc Fully Automated Laboratory Analysis System User Manual 

Contents

- 1 [dynamic BIOSENSORS AS-2-Rc Fully Automated Laboratory Analysis System](#)
- 2 [Product Usage Instructions](#)
- 3 [Key Features](#)
- 4 [heliX® Adapter Chip Overview](#)
- 5 [Product Description](#)
- 6 [Preparation | MIX&RUN](#)
- 7 [FAQ](#)
- 8 [Documents / Resources](#)
 - 8.1 [References](#)
- 9 [Related Posts](#)

dynamic
BIOSENSORS

dynamic BIOSENSORS AS-2-Rc Fully Automated Laboratory Analysis System



Specifications

- **Product Name:** heliX+ ADAPTER STRAND 2 with red dye Rc
- **Manufacturer:** Dynamic Biosensors GmbH & Inc.
- **Model:** AS-2-Rc v5.1
- **Key Features:**
 - 2 spots with 2 different anchor sequences for DNA-encoded addressing.
- **Order Number:** AS-2-Rc
- **Concentration:** 400 nM
- **Storage:** White cap, please see expiry date on label. Avoid many freeze-thaw cycles by aliquoting the nanolever.

Product Usage Instructions

Preparation | MIX&RUN

1. Mix Adapter strand 1 – Rc (400 nM) and conjugated Ligand strand with ligand 1 (500 nM) at 1:1 ratio (v/v).
2. Mix Adapter strand 2 – Rc (400 nM) and conjugated Ligand strand with ligand 2 (500 nM) at 1:1 ratio (v/v).
3. Incubate separately the two solutions of step 1 and 2 at room temperature at 600 rpm for 30 min to ensure complete hybridization.
4. Mix solutions of step 1 and 2 at 1:1 ratio (v/v). Solution is ready for biochip functionalization.


Key Features

- Adapter strand 2 for functionalization of heliX® Adapter Chip Spot 2.
- Compatible with heliX® Adapter Chip.
- Includes Adapter strands for 50 regenerations.
- Ideal for MIX&RUN sample preparation.

- Adapter strand 2 carries a hydrophobic red dye (Rc) with a neutral net charge.

heliX® Adapter Chip Overview

2 spots with 2 different anchor sequences for DNA-encoded addressing.



Product Description

Order Number: AS-2-Rc


Table 1. Contents and Storage Information

Material	Cap	Concentration	Amount	Buffer	Storage
<i>Adapter strand 2 - Rc</i>	White	400 nM	5 x 100 µL	TE40 ^[1]	-20°C

- For research use only.
- This product has a limited shelf life, please see expiry date on label.
- To avoid many freeze thaw cycles please aliquot the nanolever

Preparation | MIX&RUN

In-solution hybridization of adapter and ligand strands:



1. Mix Adapter strand 1 – Rc (400 nM) and conjugated Ligand strand with ligand 1 (500 nM) at 1:1 ratio (v/v).
2. Mix Adapter strand 2 – Rc (400 nM) and conjugated Ligand strand with ligand 2 (500 nM) at 1:1 ratio (v/v).
3. Incubate separately the two solutions of step 1 and 2 at RT at 600 rpm for 30 min to ensure complete hybridization.
4. Mix solution of step 1 and 2 at 1:1 ratio (v/v). Solution is ready to use for biochip functionalization. Stability of the solution is related to the stability of the ligand molecules.

Table 2. Additional material for functionalization of spot 1 and reference spot 2.

Material	Concentration	Buffer	Related Product Name	Order No
Adapter strand 1 - Rc	400 nM	TE40 ^[1]	Adapter strand 2 with red dye Rc	AS-1-Rb
Ligand strand carrying the conjugated ligand 1	500 nM	PE40 ^[2]	heliX [®] Amine Coupling Kit 1	HK-NHS-1
Ligand strand carrying the conjugated ligand 2	500 nM	PE40 ^[2]	heliX [®] Amine Coupling Kit 1	HK-NHS-1

Example

Required volume for 3 functionalizations: 100 µL with a final concentration of 100 nM.

Vial 1		Vial 2	
Adapter strand 1 - Rc (400 nM)	Conjugated Ligand strand with ligand 1 (500 nM)	Adapter strand 2 - Rc (400 nM)	Conjugated Ligand strand + with ligand 2 (500 nM)
25 µL	25 µL	25 µL	25 µL

After incubation time, mix vial 1 and vial 2 to obtain 100 µL of ready-to-use DNA solution.

Contact

Dynamic Biosensors GmbH Perchtinger Str. 8/10 81379 Munich Germany

Dynamic Biosensors, Inc.

300 Trade Center, Suite 1400 Woburn, MA 01801 USA

- Order Information order@dynamic-biosensors.com
- Technical Support support@dynamic-biosensors.com
- www.dynamic-biosensors.com

Instruments and chips are engineered and manufactured in Germany.

©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.

1. **TE40:** 10 mM Tris, 40 mM NaCl, 0.05 % Tween20, 50 µM EDTA, 50 µM EGTA
2. If the protein is not stable in PE40 (TE40, HE40), please check buffer compatibility with the switchSENSE® compatibility sheet.

www.dynamic-biosensors.com

FAQ

Q: What is the shelf life of the product?

A: The product has a limited shelf life, please refer to the expiry date on the label.


Q: How should I store the product to maintain stability?

A: Store the product with a white cap to avoid exposure to light. Avoid freeze-thaw cycles by aliquoting the nanolever.

Q: How should I prepare the solution for biochip functionalization?

A: Follow the MIX&RUN preparation steps provided in the user manual to prepare the solution for biochip functionalization.

Documents / Resources

	<p>dynamic BIOSENSORS AS-2-Rc Fully Automated Laboratory Analysis System [pdf] User Manual</p> <p>AS-2-Rc, AS-2-Rc Fully Automated Laboratory Analysis System, Fully Automated Laboratory Analysis System, Automated Laboratory Analysis System, Laboratory Analysis System, Analysis System, System</p>
---	--

References

- [User Manual](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.