

**dynamic**  
BIOSENSORS  
**BIOSENSORS HK-  
EA-2 RNA Enzyme  
Activity Kit**



# dynamic BIOSENSORS HK-EA-2 RNA Enzyme Activity Kit User Manual

[Home](#) » [dynamic BIOSENSORS](#) » dynamic BIOSENSORS HK-EA-2 RNA Enzyme Activity Kit User Manual 

## Contents

- [1 dynamic BIOSENSORS HK-EA-2 RNA Enzyme Activity Kit](#)
- [2 Key Features](#)
- [3 Helix Adapter Chip Overview](#)
- [4 Product Description](#)
- [5 Preparation](#)
- [6 Documents / Resources](#)
  - [6.1 References](#)
- [7 Related Posts](#)

**dynamic**  
BIOSENSORS

**dynamic BIOSENSORS HK-EA-2 RNA Enzyme Activity Kit**




## Key Features

- ssRNA template (48 + 32 bp) for functionalization of heliX
- Adapter Chip on Spot 1.
- dsDNA-RNA template (48 + 32 bp) for functionalization of heliX
- Adapter Chip on Spot 2.
- Compatible with heliX Adapter Chip.
- Includes Adapter strands for 100 regenerations.
- This RNA template carries a moderately hydrophilic red dye (Ra) with a single positive net charge.

## Helix Adapter Chip Overview

- 2 spots with 2 different anchor sequences for RNA-encoded addressing.



## Product Description

**Order Number: HK-EA-2**

**Table 1.** Contents and Storage Information

Material	Cap	Concentration	Amount	Buffer	Storage
<i>Adapter strand 1 - nf / Template-32-P-Ra RNA</i>	White	200/250 nM	10 x 100 µL	TE40 <sup>[1]</sup>	-20°C
<i>cTemplate Adapter strand 2 - nf / Template-32-P-Ra RNA</i>	Black	200/250 nM	10 x 100 µL	TE40 <sup>[1]</sup>	-20°C

For research use only.

This product has a limited shelf life, please see the expiry date on the label. To avoid many freeze-thaw cycles please aliquot the nanolever.

## Preparation

### IMPORTANT

Both Adapter strands are already pre-hybridized. Adapter strand 1 to the Template-32-P-Ra RNA strand, leaving the upper part as ssDNA, and cTemplate-Adapter strand 2 to the Template-32-P-Ra RNA strand, leaving Spot 2 completely as dsDNA-RNA. Next, simply mix in the same vial the Adapter strand 1 – nf / Template-32-P-Ra RNA (200/250 nM) and the template Adapter strand 2 – nf / Template-32-P-Ra RNA (200/250 nM) at a 1:1 ratio (v/v). The solution is ready to use for biochip functionalization.

## Useful Order Numbers

**Table 2.** Order Numbers

Product Name	Comment	Order No
heliX <sup>®</sup> Adapter Chip	Chip with 2 detection spots	ADP-48-2-0
10x Passivation solution	For passivation of chip surface	SOL-PAS-1-5
Regeneration solution	For regeneration of chip surface	SOL-REG-1-5

## Assay Setup in heliOS

For studying enzymatic activity of a nucleic acid modifying enzyme (e.g., polymerase, reverse transcriptase, helicase, etc.). Go to heliOS > create a New Assay Workflow > add Custom Assay > load Enzyme Binding and Activity > modify the parameters based on your needs and run the assay. Suggested assay parameters (e.g., flow rates, functionalization time, LED power, etc.) are within the heliOS assay.

### TIP

Is strongly recommended to perform binding kinetics of the enzyme beforehand. The obtained K<sub>d</sub> during enzyme kinetics can be the initial test concentration for the association of the enzyme during enzyme kinetics. This concentration is a good compromise to not overcrowd the surface and avoid multiple enzymes binding to the same template.

### TIP

For an initial scouting of the substrate, choose a broad concentration splitting spanning the low nanomolar to high micromolar, and a blank. A minimum of 6 concentrations of the substrate are recommended to obtain a reliable sigmoidal fit during the extraction of the K<sub>M</sub>. For inhibition assay, analysis, or any further questions, please contact the support team at [support@dynamicbiosensors.com](mailto:support@dynamicbiosensors.com).

## Contact

Dynamic Biosensors GmbH  
Perchtinger Str. 8/10

81379 Munich

## Germany


Dynamic Biosensors, Inc.  
300 Trade Center, Suite 1400  
Woburn, MA 01801

## USA

Order Information [order@dynamic-biosensors.com](mailto:order@dynamic-biosensors.com)  
Technical Support [support@dynamic-biosensors.com](mailto:support@dynamic-biosensors.com)  
[www.dynamic-biosensors.com](http://www.dynamic-biosensors.com)

Instruments and chips are engineered and manufactured in Germany.  
2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.

## Documents / Resources

	<p><a href="#">dynamic BIOSENSORS HK-EA-2 RNA Enzyme Activity Kit</a> [pdf] User Manual v1.1, HK-EA-2 RNA Enzyme Activity Kit, HK-EA-2, RNA Enzyme Activity Kit, Enzyme Activity Kit , Activity Kit, Kit</p>
---	--

## References

- [HomePage | Biosensors International Ltd](#)
- [Home - Dynamic Biosensors](#)
- [Home - Dynamic Biosensors](#)
- [User Manual](#)

### Manuals+, Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.