

[Home](#) » [ELEMENTAL MACHINES](#) » **ELEMENTAL MACHINES Element-D WiFi Data and Metadata Collecting Installation Guide** 

ELEMENTAL MACHINES Element-D WiFi Data and Metadata Collecting



Contents [[hide](#)]

- [1 Applicable Asset Models](#)
- [2 Package Contents](#)
- [3 Safety Information](#)
- [4 Prerequisites](#)
- [5 Installation Instructions](#)
- [6 Customer Support](#)
- [7 Documents / Resources](#)
 - [7.1 References](#)

Applicable Asset Models

- Five Easy Plus FP20

Package Contents

- Element-D
- Power Supply
- Mounting Plate
- RS-232 to USB connector cable

Note: You may receive multiple types of cables if you are connecting Element-D's to different types of assets.

Safety Information



Review and abide by the safety information found in the Element-D and Mettler Toledo Five Easy Plus FP20 pH Meter user manuals. Tidy connection cable to prevent tripping hazards.

Prerequisites

The Element-D only captures data for operational use. The capture of calibration data is not supported. Unplug the Element-D from the pH meter during the calibration process.

Installation Instructions

Receiving the Shipment

1. If this is your first time setting up the Elemental Machines system, you will receive an email for your dashboard account verification. Save this email for when the devices arrive.
2. When devices are shipped, they will be added to your Elemental Insights™ Dashboard with default names.
When you first log in you will see all your devices with a 'disconnected' status.
3. When devices arrive, carefully remove the contents from the packaging. If any of the contents look damaged, please send a picture to help@elementalmachines.com.
4. Confirm the Element-D's in your package are labeled with the assets to which they are to connect. If the labels do not match the intended assets for connection, please

contact help@elementalmachines.com.


Element-D Initial Power-up

5. Plug the Element-D in using the provided power supply.
6. If the device has not been preloaded with your company's Wi-Fi, complete Wi-Fi Connection following instructions in the provided 771-00010 Wi-Fi Connection Guide – Gateway GW2/3, Element-C/D.
 - a. 771-00010 is included with your Element-D as a printed copy. Digital versions of this are available in the Support section of the Elemental Insights™ Dashboard.

Attaching Connection Cable to Asset

7. Locate the RS-232 data port on the asset (shown in Figure 1) to which the Element-D will connect.

Figure 1 Data port on asset



8. Plug the RS-232 end of the Connection Cable (Figure 2) into the data port on the pH meter.

Figure 2 Connector cable USB end (left) and RS-232 end (right)



Configuring Asset Menus

9. To receive automatic readings from the pH meter, if the print icon is showing (circled below) simultaneously hold down the STO/RCL and Mode/Setup buttons until the data transfer icon disappears.


Figure 3 pH meter print icon (circled in blue) and menu buttons



Attaching Connection Cable to Element

10. Plug the other end of the Connection Cable into the data port on the Element-D. This will be in the USB port which are any of the four on the left as pictured in Figure 4.

Figure 4 Data ports on Element-D



Positioning the Element-D

11. Determine where the Element-D should be placed to allow for power access and a tidy connection cable.
12. Securely position the Element-D in the chosen location. The Element-D has magnets which make it easily mountable to magnetic surfaces.
 - a. **Note:** If applicable, an adhesive metal mounting plate available from EM may be used.

Finalization

13. To send a pH measurement to the Elemental Insights™ Dashboard, press Read.
14. Make sure the end point format is Manual by waiting until the \sqrt{M} icon (see Figure 5) appears on the screen the press Read again.

Figure 5 Icon circled in blue



15. Once configured, press the Read button to begin receiving measurements at one second intervals. To stop measurements, press the Read button again.
16. Confirm that the value appears on the dashboard.
17. If the dashboard does not display any value, check cable connection, menu settings

and power cycle the Element-D. Contact help@elementalmachines.com if problem persists.

18. Contact any third-party vendor for next steps regarding communication to the external platform (if applicable)

Customer Support

772-00079 rev 01 Element-D, ED1 Mettler Toledo Five Easy Plus

FP20 pH Meter Serial Configuration Installation Guide

Copyright © 2024 Elemental Machines

help@elementalmachines.com




Documents / Resources

	<p>ELEMENTAL MACHINES Element-D WiFi Data and Metadata Collecting [pdf] Installation Guide</p> <p>ED1, FiveEasy Plus FP20, Element-D WiFi Data and Metadata Collecting , Element-D, WiFi Data and Metadata Collecting, and Metadata Collecting , Metadata Collecting</p>
--	--

References

- [User Manual](#)

 ELEMENTAL
MACHINES

 and Metadata Collecting, ED1, Element-D, Element-D WiFi Data and Metadata Collecting, ELEMENTAL MACHINES, FiveEasy Plus FP20, Metadata Collecting, WiFi Data and Metadata Collecting

—Previous Post

[ELEMENTAL MACHINES ED1 Element IoT Device Installation Guide](#)

Leave a comment

Your email address will not be published. Required fields are marked *

Comment *

Name

Email

Website

☐ Save my name, email, and website in this browser for the next time I comment.

Post Comment

Search:

Search

[Manuals+](#) | [Upload](#) | [Deep Search](#) | [Privacy Policy](#) | [@manuals.plus](#) | [YouTube](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.