



euromex Delphi-X Inverso Microscope User Manual

[Home](#) » [euromex](#) » euromex Delphi-X Inverso Microscope User Manual 



quick start manual for use of DIC



**User Manual
additional
v.425271**

Contents

1 Quickstart

2 Documents / Resources

2.1 References

3 Related Posts

Quickstart

Preface

The Delphi-X Inverso has been designed for easy use of DIC imaging. The following steps will quickly guide you through the set-up process of the microscope

To work properly the setup of the DIC system needs a combination objective, condenser prisms, rotatable analyzer and the objective DIC prisms

The DIC prisms for objective and objectives work in specific pairs. Therefore, the DIC prism sliders are marked with the magnification of the matching objective

DIC prisms for condenser

The DIC prisms are already pre-mounted in the factory, there is no need to adjust anything Setting up

1. Insert the objective prisms into the slots of the condenser positioned under the objective. The DIC prism sliders are marked with the magnification of the matching objective (photo 1)
2. Insert the rotatable analyzer into its slot (photo 2)
3. Set the condenser position to BF (photo 3)
4. Focus on your specimen (start with lowest magnification)
5. Rotate the analyzer till the image turns dark
6. Set the condenser position to DIC or DICII (photo 4)
7. Fine-tune the DIC image by rotate the analyzer





Phase contrast, fluorescence and brightfield microscopy

Remove the DIC objective prism, rotatable analyzer, change the condenser position from DIC or DICII to brightfield or phase contrast. before starting microscopy. Doing so will increase the brightness of the image and overall image quality



euromex.academy



All information mbyte changed with . out prior notice v.425271

Documents / Resources



[euromex Delphi-X Inverso Microscope](#) [pdf] User Manual

DIC_manual_EN_2, Delphi-X Inverso Microscope, Delphi-X, Inverso Microscope, Microscope

References

- [User Manual](#)

Manuals+, Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.