



**Flow
Headset
Medically
Approved
Brain**



FLOW NEUROSCIENCE Flow Headset Medically Approved Brain Stimulation Instructions

[Home](#) » [FLOW NEUROSCIENCE](#) » FLOW NEUROSCIENCE Flow Headset Medically Approved Brain Stimulation Instructions 

Contents

- [1 FLOW NEUROSCIENCE Flow Headset Medically Approved Brain](#)
- [2 Product Information](#)
- [3 Product Usage Instructions](#)
- [4 Frequently Asked Questions](#)
- [5 Mechanism of action](#)
- [6 Treatment schedule](#)
- [7 Warnings and Precautions](#)
- [8 Support and contact information](#)
- [9 Documents / Resources](#)
 - [9.1 References](#)



FLOW NEUROSCIENCE Flow Headset Medically Approved Brain



Product Information

Specifications:

- Product Name: Flow Clinician
- Classification: Class IIa medical device
- Intended Use: Treatment of unipolar major depressive disorder (MDD) in adults
- Mechanism: Transcranial Direct Current Stimulation (tDCS)
- Current Strength: 2mA
- Treatment Sessions: 30 minutes each

Product Usage Instructions

Mechanism of Action:

The Flow headset delivers a 2mA electric current to the left dorsolateral prefrontal cortex, enhancing neuronal activity to improve mood and emotional regulation.

Efficacy:

Flow has shown significant results in treating depression, achieving a 57.5% remission rate in 10 weeks. It can be used as a standalone therapy or with antidepressants.

Treatment Schedule:

The standard treatment includes an activation phase for weeks 1-3 and a strengthening phase from week 4 onwards, with each session lasting 30 minutes.

Behavioural Therapy:

The Flow app offers optional behavioral therapy courses such as activation, exercise, meditation, diet, and sleep hygiene, complementing the stimulation protocol.

Contraindications:

No absolute contraindications reported for Flow Clinician.

Warnings and Precautions:

Flow Clinician can be used by individuals with various chronic conditions. It is not licensed for use in pregnancy but appears safe for breastfeeding and postpartum depression.

Adverse Reactions:

Common side effects may include skin irritation, headaches, and rare cases of worsening depression or anxiety symptoms.

Frequently Asked Questions

- **Q: Can Flow Clinician be used by children and adolescents?**
 - A: Flow Clinician is not licensed for use in children and adolescents, but studies are ongoing for ages 14+.
- **Q: Are there any known interactions with medications?**
 - A: Flow Clinician has no known interactions with medications, making it suitable for individuals with various chronic conditions.

Indication

Class IIa medical device approved for the treatment of unipolar major depressive disorder (MDD) in adults, either as monotherapy or as an adjunct with antidepressants and psychological therapies.


Mechanism of action

The Flow headset delivers a gentle 2mA electric current (transcranial Direct Current Stimulation, tDCS) to the left dorsolateral prefrontal cortex. This region is crucial for emotional regulation and is often hypoactive in individuals with depression. tDCS enhances neuronal activity in this area, counteracting hypoactivity, and improving mood and emotional regulation.

Efficacy

Results at 10 weeks from a multisite, double-blind placebo-controlled randomised superiority trial were statistically significant for both primary and secondary endpoints. Flow achieved a 57.5% remission rate using the MADRS scale and was effective as both standalone and adjunctive therapy.¹


Remission rate when used as standalone and adjunct:



Open-label patient cohort study in an NHS primary care general practice found at 6 weeks: 2


- 58.1% of patients showed reliable improvement, and 32.3% achieved remission using the PHQ-9
- Significant improvements were observed in functioning (WSAS) and health-related quality of life (EQ-5D-5L)

Improvements in Depression Severity and Functional Impairment



Treatment schedule

The standard treatment is split into 2 phases – activation and strengthening. Each session (termed stimulation) is 30 minutes long. Patients complete a MADRS-s survey at the beginning of each week to monitor progress



Activation Phase

Weeks 1-3

- Includes 5 stimulations a week
- Introduces optional behavioral therapy courses via the app



Strengthening Phase:


Week 4 onwards

- Includes up to 2 stimulations a week

- Continues to strengthen and preserve results
- Continue behavioral therapy courses

Behavioural therapy

The Flow app includes 7 behavioral therapy courses including behavioral activation, exercise, meditation, diet, and sleep hygiene, totaling over 50 short therapy sessions. These courses are entirely optional and do not interfere with the Flow stimulation protocol.



Contraindications

No absolute contraindications.

Warnings and Precautions

- Broken/inflamed/infected skin (including, for example, psoriasis) at the electrode site
- Cranial or intracranial implant (e.g. brain clips, deep brain stimulators)
- Craniofacial abnormalities (e.g., congenital deformities, severe trauma, or reconstructive surgeries) which may affect electrode placement
- Epilepsy or history of seizures
- Active suicidal ideation (requires closer monitoring)
- History of hypomania or mania (may require closer monitoring)

Special populations

- Other chronic conditions: No known interactions with medications; use by individuals with diabetes, heart disease, hypertension, asthma, co-occurring mental health conditions, neurodivergence, and brain injuries/ disorders without reported safety concerns.
- Pregnancy: Not licensed for use; no safety concerns identified in existing studies, research ongoing.
- Breastfeeding: Safe to use; no effect on breast milk production.

- Postpartum: Appears safe for postpartum depression and successfully used in NHS pilots.
- Children and Adolescents: Not licensed for use; studies ongoing for ages 14+

Adverse reactions

Based on real world evidence from >20,000 users, the incidence rate of adverse reactions is 4.5%, which is favorable against antidepressants

Side effect	Incidence	Probable duration	Recommendations to reduce intensity and/or duration
Skin irritation (including tingling, stinging redness)	$\geq 10\%$ to $< 20\%$	Less than a day. Rarer cases up to 3 days	Avoid stimulating on irritated skin <ul style="list-style-type: none"> • new pads with each stimulation • well-moisturized skin • well-hydrated • ice pack following stimulation • correct electrode placement
Headaches	$\geq 10\%$ to $< 20\%$	Less than a day	Adequate hydration • relaxed environment • correct electrode and headset placement • post-stimulation rest
Worsening depression symptoms and/or anxiety	$\geq 1\%$ to $< 10\%$	Up to several days. Usually only occurs in the first few weeks of treatment	Relaxed environment • correct electrode placement
Fatigue, malaise and sleep disturbances	$\geq 1\%$ to $< 10\%$	Up to several days	Changing/fixing time of day of stimulation • relaxed environment
Tinnitus	$\geq 0.1\%$ to $< 1\%$	Less than a day. In rare cases persisting for a few weeks	Correct electrode placement • monitoring during session and stop use if persisting
1st degree skin burns	$\geq 0.1\%$ to $< 1\%$	Up to 3 days	New pads with each session • avoid stimulating on irritated skin • well-moisturised skin • Well-hydrated
Photopsia (Flashes of Light)	$\geq 0.1\%$ to $< 1\%$	Up to a minute	Ensuring correct electrode contact and placement • not removing headset without powering down
Hypomania/ Mania	$< 0.01\%$	Up to several days	Closer monitoring if history of bipolar disorder


- Adverse events should be reported. Report forms and information can be found at www.mhra.gov.uk/yellowcard
- Adverse events should also be reported to Flow Neuroscience by emailing support@flowneuroscience.com.

Clinician review and customization

Review your patients' progress and adherence through the Flow Clinician Platform (CPP). Accessible anytime and anywhere through a secure web browser, this platform streamlines patient visits by providing informed insights. It also allows for the customisation of treatment protocols (maximum 1 per day, 7 per week) based on clinical judgment to address individual patient needs. This setup ensures personalized treatment adjustments to help optimize patient outcomes.



The suggested initial treatment duration is 10 weeks. Following 10 weeks, patients' response should be assessed and length of continued treatment determined, taking into consideration the effectiveness to date and the risk of relapse.



Support and contact information

For technical support and assistance, patients can contact support@flowneuroscience.com. Clinicians can

enquire about setup, pricing, and training support by contacting clinical@flowneuroscience.com.

- www.flowneuroscience.com

References:

1. Woodham, R.D., Selvaraj, S., Lajmi, N. et al. Nat Med (2024).
2. Griffiths, C. et al. (2024). Open Journal of Depression. 13(02), 25–39.

Class IIa medical device


CE 673331

UKCA 776047

Holder: Flow Neuroscience AB, Södra Tullgatan 3, 21140 Malmö

www.flowneuroscience.com

Documents / Resources



[FLOW NEUROSCIENCE Flow Headset Medically Approved Brain Stimulation](#) [pdf] Instructions

Flow Headset Medically Approved Brain Stimulation, Medically Approved Brain Stimulation, Approved Brain Stimulation, Brain Stimulation, Stimulation

References

- [🕒 tDCS depression treatment - Flow Neuroscience](#)
- [🕒 tDCS depression treatment - Flow Neuroscience](#)
- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.