

GALENVS UPX0100 Universal Pathogen Complex Module



GALENVS UPX0100 Universal Pathogen Complex Module User Guide


[Home](#) » [GALENVS](#) » **GALENVS UPX0100 Universal Pathogen Complex Module User Guide** 

Contents

- [1 GALENVS UPX0100 Universal Pathogen Complex Module](#)
- [2 Product Usage Instructions](#)
- [3 Universal Pathogen Complex Matrix Module](#)
- [4 Milk Fecal](#)
- [5 Other matrices](#)
- [6 Documents / Resources](#)
 - [6.1 References](#)
- [7 Related Posts](#)



GALENVS UPX0100 Universal Pathogen Complex Module



Specifications:

- **Product Name:** Universal Pathogen Complex Matrix Module
- **Sample Types:** Milk, Fecal, oral fluid, viscous, swab, stool, exudate, urine, etc.
- **Components:** Pre-Treatment Beads (small and large), Proteinase K, Universal Pathogen DNA/RNA Extraction Kit
- **Speed:** 13,000 RPM

Product Usage Instructions

Step A: For Milk and Fecal Samples

1. Ensure Pre-Treatment Beads (small) are used for milk samples and Pre-Treatment Beads (large) for fecal samples.
2. Add the sample to the Universal Pathogen Complex Matrix Module.
3. Include Proteinase K in the sample.
4. Run the module at 13,000 RPM.

Step B: For Other Matrices

1. Proceed with the Universal Pathogen DNA/RNA Extraction Kit.
2. Add the sample to the module.
3. Include Proteinase K in the sample.
4. Operate the module at a speed of 13,000 RPM.

FAQ:

- **Q: Can I use the same kit for all sample types?**

A: Yes, the Universal Pathogen Complex Matrix Module is designed to handle various sample types efficiently.

Just follow the specific instructions for each type.

- **Q: What is the purpose of Proteinase K in the process?**

A: Proteinase K is used to digest proteins in the sample, facilitating the extraction of DNA/RNA for analysis.

- **Q: Is there a difference between manual and automated operation?**

A: The manual operation involves user intervention at different stages, while automated operation streamlines the process with predefined steps.


Universal Pathogen Complex Matrix Module

Quick Start Guide

- For milk and fecal samples, proceed to step A
- For all other matrices, eg. oral fluid, viscous, swab, stool, exudate, urine, etc., proceed to step B


Milk Fecal

- A To the 2mL Bead Beating tube provided, add 200μL of sample, 300μL Pre-Treatment Buffer (PTB) and 40μL Proteinase K. Shake well before use.
- Vortex for 30 seconds and heat at 70°C for 5 mins. To pellet debris and solids, briefly centrifuge at 13000rpm for 1 min.
- **Note:** use Pre-Treatment Beads (small) for milk and Pre-Treatment Beads (large) for fecal samples.



Other matrices

- B To a 1.5mL centrifuge tube, add 200μL of sample, 300μL Pre-Treatment Buffer (PTB) and 40μL Proteinase K. Shake well before use.
- Vortex for 30 seconds and heat at 70°C for 5 mins. To pellet debris and solids, briefly centrifuge at 13000rpm for 1 min.




Collect 300µL and proceed to
Universal Pathogen DNA/RNA Extraction Kit

- **Manual:** UP0050 UP0020
- **Automated:** UP0016 UP0096

UPX Universal Pathogen Complex Matrix Module

Documents / Resources

	GALENVS UPX0100 Universal Pathogen Complex Module [pdf] User Guide UPX0100 Universal Pathogen Complex Module, UPX0100, Universal Pathogen Complex Modul e, Pathogen Complex Module, Complex Module
---	--

References

- [User Manual](#)

[Manuals+.](#) [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.