




# heliXcyto CCK-1-1 Dynamic Biosensors User Manual

[Home](#) » [heliXcyto](#) » heliXcyto CCK-1-1 Dynamic Biosensors User Manual 

## Contents

- [1 heliXcyto CCK-1-1 Dynamic Biosensors](#)
- [2 Specifications](#)
- [3 Product Usage Instructions](#)
- [4 FAQ](#)
- [5 Product Description](#)
- [6 Reconstitution](#)
- [7 CONTACT](#)
- [8 Documents / Resources](#)
  - [8.1 References](#)
- [9 Related Posts](#)

heliXcyto

**heliXcyto CCK-1-1 Dynamic Biosensors**



## Specifications

- **Product Name:** heliXcyto Chip Cleaning Kit
- **Order Number:** CCK-1-1
- **Contents:**
  - Cleaning Solution 1 (detergent-based) – 1 mL
  - Lyophilized Enzyme Powder – 5 mg
  - Tris-HCl Buffer (reconstitution buffer) – 1 mL
- **Storage Information:**
  - Cleaning Solution 1: Transparent Cap
  - Lyophilized Enzyme Powder: Red Cap

## Product Usage Instructions

### Reconstitution

- Follow the reconstitution instructions provided in the user manual for each component of the kit.

### Application Note

Manually cleaning helicity chips extends their lifetime. To clean the chips, inject the cleaning solutions into the microfluidic channel of the chip using a pipette. To ensure effective cleaning:

1. Wear gloves when handling chips.
2. Use freshly filtered buffers only.
3. Inject solutions into the entry port on the left side of the flipped chip.

## Chip Cleaning Protocol

- Refer to the detailed chip cleaning protocol provided in the user manual for step-by-step instructions on how to clean helicity chips using Cleaning Solutions 1 and 2.

## Contact Information

- If you have any questions or need assistance, you can contact Dynamic Biosensors GmbH in Germany or Dynamic Biosensors, Inc. in the USA using the provided contact information.

## FAQ

- **Q:** Can I use other cleaning solutions with the heliXcyto Chip Cleaning Kit?
- **A:** It is recommended to use only the provided Cleaning Solution 1 and Cleaning Solution 2 for optimal cleaning results and to maintain chip performance.
- **Q:** How often should I clean my helicity chips?
- **A:** The frequency of chip cleaning depends on usage. It is recommended to clean the chips regularly to ensure accurate and reliable results.

## CHIP CLEANING KIT

- Detergent- and enzyme-based solutions for manual cleaning of chips
- Dynamic Biosensors GmbH & Inc.
- CCK-1-1 v1.1

## Product Description

- Order Number: CCK-1-1

**Table 1. Contents and Storage Information**

Material	Cap	Amount	Storage
<b>Cleaning Solution 1</b> (detergent-based)	Transparent	1 mL	RT
<b>Lyophilized Enzyme Powder</b>	Red	5 mg	2-8°C
<b>Tris-HCl Buffer</b> (reconstitution buffer)	Transparent	1 mL	2-8°C

- The heliXcyto Chip Cleaning Kit contains Cleaning Solution 1 and Cleaning Solution 2 for manual cleaning of the heliXcyto chip outside of the heliXcyto device.
- For research use only.
- This product has a limited shelf life, please see the expiry date on the label.


## Reconstitution

- The Cleaning Solution 2 has to be reconstituted before first use by dissolving the 5 mg lyophilized enzyme powder in 1 mL of Tris-HCl reconstitution buffer.
- Reconstituted Cleaning solution 2 should be aliquoted (10 – 100 µL volumes) and stored at -20 °C.
- Cleaning solution 1 is ready-to-use and does not require any preparation.

### Application Note

Manually cleaning heliXcyto chips extends their lifetime. This can be performed by injecting cleaning solutions into the microfluidic channel of the chip with a pipette. To reduce the introduction of dust particles into heliXcyto chips during manual handling and optimally wash out cell remnants it is recommended to:

- wear gloves when handling chips
- only use freshly filtered buffers
- always inject solutions into the entry port on the left side of the flipped chip



### heliXcyto chip cleaning protocol

- Fill the microfluidic channel with 10 µl fresh Cleaning solution 2
- Seal chip in- and outlet with clear sticky tape (e.g. Sellotape™, do not use tape that leaves glue remnants)
- Incubate for at least 20 min at RT (or longer, e.g. overnight)
- Wash microfluidic channel 3 x with 10 µl Cleaning solution 1
- Wash the microfluidic channel 3x with 10 µl water (Milli-Q or similar)
- Empty the chip by pipetting 10 µl air.

## CONTACT

### Dynamic Biosensors GmbH


- Perchtinger Str. 8/10
- 81379 Munich
- Germany

### Dynamic Biosensors, Inc.

- 300 Trade Center, Suite 1400
- Woburn, MA 01801
- USA
- Order Information [order@dynamic-biosensors.com](mailto:order@dynamic-biosensors.com)
- Technical Support [support@dynamic-biosensors.com](mailto:support@dynamic-biosensors.com)
- [www.dynamic-biosensors.com](http://www.dynamic-biosensors.com)

Instruments and chips are engineered and manufactured in Germany. ©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.

## Documents / Resources

	<p><a href="#">heliXcyto CCK-1-1 Dynamic Biosensors</a> [pdf] User Manual CCK-1-1 Dynamic Biosensors, CCK-1-1, Dynamic Biosensors, Biosensors</p>
---	---

## References

- [HomePage | Biosensors International Ltd](#)
- [Home - Dynamic Biosensors](#)
- [Home - Dynamic Biosensors](#)
- [User Manual](#)

### Manuals+ Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.