

# ibidi 12170 Objective Heater Universal Instruction Manual


[Home](#) » [ibidi](#) » ibidi 12170 Objective Heater Universal Instruction Manual 

## Contents

- [1 ibidi 12170 Objective Heater Universal](#)
- [2 Product Information](#)
- [3 Product Usage Instructions](#)
- [4 Objective Heater Universal – Silver Line Version](#)
  - [1.0](#)
- [5 Preamble](#)
- [6 Intended Use](#)
- [7 Principle](#)
- [8 Equipment](#)
- [9 Operation](#)
- [10 Documents / Resources](#)
  - [10.1 References](#)



**ibidi 12170 Objective Heater Universal**



## Product Information


The product is a flexible heating unit with a cable. It operates on DC 24V, 2A, and has a power rating of 48W. The product is designed to be used in an operating site with specific temperature and humidity conditions. It is compliant with EU legislation as indicated by the CE marking on the product identification label.

## Product Usage Instructions

1. Before using the device, read the instruction manual carefully.
2. Ensure that the device is connected to a power source with a voltage of DC 24V.
3. Place the device in an appropriate operating site with the recommended temperature and humidity conditions.
4. Make sure to dispose of the product separately from other household waste. Follow the instructions on page 9 for proper disposal.
5. Use the device according to the intended purpose, which is to provide flexible heating.
6. In case of an emergency, refer to the emergency release mechanism described in section 5.2 of the instruction manual.

## Objective Heater Universal – Silver Line Version 1.0

for use with ibidi Heating System Slide/Dish – Silver Line and ibidi Heating System 4 Slides – Silver Line



## Preamble

## Introduction

This manual is your guide to using the Objective Heater Universal – Silver Line for cell culture experiments on an optical microscope. It instructs first-time users on how to use the instrument, and serves as a reference for experienced users. Before using the Objective Heater Universal – Silver Line, please read this instruction manual carefully and make sure that the contents are fully understood. This manual should be easily accessible to the operator at all times during instrument operation. If this manual gets lost, see the online version on [ibidi.com](http://ibidi.com) or request a replacement via [techsupport@ibidi.com](mailto:techsupport@ibidi.com). To ensure safe operation, the Objective Heater Universal – Silver Line must only be operated with the supplied components and according to the instruction manual. For Research Use Only! Not for use in diagnostic procedures. The Objective Heater Universal – Silver Line (#12170) can be used with the Heating System Slide/Dish – Silver Line (product numbers 12110, 12720, 12722, 12111) and the Heating System 4 Slides – Silver Line (product numbers 12130, 12131).

## Safety Symbols

Note that the signal words WARNING, CAUTION and NOTE have specific meanings in this manual. Do not proceed beyond a signal word until you have performed the indicated actions. Warning messages in the text are displayed in a gray shaded box. Please see section ?? for general safety considerations.

- **WARNING** – A potentially hazardous situation which, if not avoided, could result in serious injury or even death.
- **CAUTION** – A potentially hazardous situation which, if not avoided, could result in minor or moderate injury. It is also used to alert against damaging the equipment or the instrument.
- **NOTE** – Additional information to help achieve optimal instrument and assay performance.


Symbols on the product identification label and back panel of the device:

- CE Marking: This symbol indicates the product's compliance with EU legislation.
- This label is positioned on the back of the device and prompts you to read the manual before using the device.

- Product disposal: The symbol indicates that this product must be recycled/disposed of separately from other household waste. See page 9 for details.

### Nomenclature

For better readability, the Objective Heater Universal – Silver Line is hereafter abbreviated and re-ferred to as “Objective Heater” in body text and figure captions. The Objective Heater consists of the following three components:



### Specifications

**Table 1:** Specifications of the Objective Heater Universal – Silver Line

|  |   |
|--|---|
| Electrical Specifications              |   |
| Overvoltage category                   | III   |
| Input line voltage                     | DC 24 V, 2 A, 48 W  |
| Operating Conditions                   |   |
| Operating site                         | Indoor use only   |
| Operating temperature                  | 18–30°C/64–86°F (min 5°C/9°F less than set temperature)                                   |
| Humidity                               | 80% relative humidity (RH) up to 31°C/87°F<br>30% relative humidity (RH) up to 40°C/104°F |
| Operating altitude                     | max. 2000 m (atmospheric pressure 800–1060 hPa/11.6–15.4 psi)                             |
| Storage conditions                     | –5–50°C/23–122°F, humidity <60% relative humidity (RH)                                    |
| Outer Dimensions and Characteristics   |   |
| Length of flexible heating unit        | 12 cm   |
| Height of flexible heating unit        | 2 cm  |
| Length of cable                        | 232.5 cm  |
| Weight of heating unit including cable | 110 g   |
| Fit With Microscope Objectives         |   |
| Diameter of Objectives                 | ~ 2.5–4.5 cm  |

**Table 1:** (continued)

---

## Compatibility With Temperature Controllers and Temperature Control Range

---

|                                  |  |
|----------------------------------|--|
| Compatibility                    | ibidi Temperature Controller – Silver Line         |
| Temperature control range        | Ambient temperature (min. 18°C/64°F) to 45°C/113°F |
| Recommended temperature settings | See Section 5.3                                    |

---

**NOTE** – The Objective Heater is only compatible with the ibidi Temperature Controller –Silver Line, hereafter abbreviated and referred to as "Temperature Controller".

### Disclaimer

- ibidi shall not be held liable, either directly or indirectly, for any damage incurred as a result of product use.
- The contents of this manual are subject to change without notice for product improvement.
- This manual is considered complete and accurate at publication.
- This manual does not guarantee the validity of any patent rights or other rights.
- If an ibidi software program doesn't function properly, this may be caused by a conflict from another program operating on the computer. In this case, take corrective action by uninstalling the conflicting product(s).
- ibidi is a registered trademark of ibidi GmbH in Germany and other countries.

### Safety Considerations

#### WARNING

- Only operate the Objective Heater with the supplied components and the Temperature Controller.
- Only use the cables and plugs delivered with the Objective Heater. The power plug of the control unit must be inserted in an outlet with a ground (earth) contact.
- Do not replace detachable power cables with inadequate specifications power cables. By violating these instructions you risk electric shock and fire.
- Do not operate the Objective Heater under conditions that pose a risk of explosion, implosion, or the release of gases.
- Do not place flammable solids, liquids, gases, or gas outlets near the Objective Heater (e.g., matches, ethanol, propane gas, solvents). Do not bring these products in contact with any other component of the system either.
- Do not operate a damaged Objective Heater. If the Objective Heater seems to be damaged, contact [techsupport@ibidi.com](mailto:techsupport@ibidi.com).
- Only operate the Objective Heater properly installed to a microscope objective.
- Some accessible parts of the Objective Heater can reach temperatures up to 65°C. Avoid touching the temperature-controlled parts of the system when you have set the Temperature Controller to high temperatures.

#### CAUTION

- Not all objectives/lenses tolerate heat! Please check with the manufacturer of your objective if it can be heated.
- Ensure that the external power supply is easily accessible. The Objective Heater must be installed in a manner such that none of its components hinders access to the external power supply.
- Immediately replace damaged cords, plugs, or cables to avoid the risk of personal injury or damage to the

instrument.

- Only ibidi technical staff and technical staff instructed by ibidi are permitted to open and service the Objective Heater.
- The Objective Heater should not come in contact with moisture. If the housing is damaged, the external power supply should not be used.
- Avoid strong magnetic fields and sources of high frequency. The Objective Heater might not function properly when located near a strong magnetic field or high frequency source.
- Avoid vibrations from vacuum pumps, centrifuges, electric motors, processing equipment, and machine tools.
- Avoid dust and corrosive gas. Do not install the Objective Heater where it could be exposed to high levels of dust or outside air or ventilation outlets.
- Install the Objective Heater in a location that enables easy access for maintenance.
- Do not place heavy objects on the Objective Heater and do not fold the Objective Heater.

### **Limited Warranty**

Products manufactured by ibidi, unless otherwise specified, are warrantied for a period of one year from the date of shipment to be free of defects in materials and workmanship. If any defects in the product are found during this warranty period, ibidi will repair or replace the defective part(s) or product free of charge.

This warranty does not apply to defects resulting from the following:

1. Improper or inadequate installation.
2. Improper or inadequate operation, maintenance, adjustment, or calibration.
3. Unauthorized modification or misuse.
4. Use of consumables, disposables, and parts not supplied by an authorized ibidi distributor.
5. Corrosion due to the use of improper solvents, samples, or due to surrounding gases.
6. Accidents beyond ibidi's control, including natural disasters.

This warranty does not cover consumables, such as cell culture chambers and dishes, tubes, fluidic connectors, reagents, etc.

The warranty for all parts supplied and repairs provided under this warranty expires on the warranty expiration date of the original product.

### **Repairing the Objective Heater Universal – Silver Line**

For inquiries concerning repair service contact ibidi GmbH technical support (techsupport@ibidi.com) and provide the model name and serial number of your system.

**CAUTION** Do not try to repair the Objective Heater by yourself. Disassembly of the Objective Heater is not allowed. Disassembly poses a risk of personal injury or damage to the devices. Contact ibidi technical support if there is a need to disassemble a device.

### **Waste Disposal – WEEE/RoHS Compliance Statement**

The European Union (EU) has enacted two directives, the first on product recycling (Waste Electrical and Electronic Equipment, WEEE) and the second on limiting the use of certain substances (Restriction on the use of Hazardous Substances, RoHS).

### **EU Directive WEEE**

The Objective Heater Universal – Silver Line must be disposed of in compliance with the WEEE Directive 2012/19/EC. This symbol on the product is in accordance with the European Union's Waste Electrical and Electronic Equipment (WEEE) Directive. The symbol indicates that this product must be recycled/disposed of

separately from other household waste. It is the end user's responsibility to dispose of this product by taking it to a designated WEEE collection facility for the proper collection and recycling of the waste equipment. The separate collection and recycling of waste equipment will help to conserve natural resources and protect human health and the environment. For more information about recycling, please contact your local environmental office, an electrical/electronic waste disposal company or the distributor where you purchased the product.

### **EU Directive RoHS**

RoHS conformity is declared in the EU conformity in Section 1.11.

### **Regulatory Statement**

### **EC Declaration of Conformity**

- We declare the compliance of the product with the requirements of Directive 2014/35/EU – Low Voltage Directive and with Directive 2014/30/EU about Electromagnetic Compatibility. Any modification to the product, not authorized by us, will invalidate this declaration.
- The object of the declaration described above is in conformity with Directive 2011/65/EU of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.
- The product meets the requirements of the following standards:

### **DIN EN 61010-1:2020**

Teil 1: Allgemeine Anforderungen Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1: General requirements

### **DIN EN 61326-1:2013**

Electrical equipment for measurement, control and laboratory use. MC requirements. General requirements This declaration certifies the conformity to the specified directives but not includes any warranted quality of the instrument. The safety documentation of the product shall be considered in detail

### **Intended Use**

The Objective Heater Universal – Silver Line is designed for use in combination with a stage top incubation system to keep cells on a microscope stage at 37°C. The Objective Heater, operated with the ibidi Temperature Controller – Silver Line, is a heating device for microscope objectives used for live cell imaging. Heating of the microscope objective with the Objective Heater prevents strong temperature gradients between the sample and the objective. This is especially important for high magnification objectives that are in close contact with the sample (e.g., water or oil immersion objectives). Using the Objective Heater improves keeping a defined and stable temperature of the sample. Please note that the Objective Heater is designed for research purposes only.

### **Principle**

Microscope objectives that require close approximation to the sample (e.g., water and oil immersion objectives) might cause cooling of the sample by thermal transfer from the sample to the objective, when the objective is not heated. Adequate heating of the objective by the Objective Heater prevents cooling of the sample. The Objective Heater includes a flexible heating element that can be easily wrapped around microscope objectives and fixed by a Velcro® fastener to enable efficient heating by a large contact area between the microscope objective and the Objective Heater.









Figure 1: Schematic illustration of a live cell imaging setup with an ibidi Stage Top Incubator and the Objective Heater for optimal temperature stability of the sample.

## Equipment

This section gives an overview of the components and their functionality and characteristics. The components of the Objective Heater Universal – Silver Line (#12170) are listed below.

| Description   | Drawing  |
|---|--|
| Flexible heating unit (red) with piece of cable and part (female) of the break-away connector   |  |
| Cable for connecting the Objective Heater to the Temperature Controller, with part (male) of the break-away connector, converter box and plug |  |
| Velcro® fastener for holding the Objective Heater on the objective  |  |

## Operation



Before starting an experiment, check that the Objective Heater fits your microscope objectives (see Section 1.5).

- **CAUTION** Not all objectives/lenses tolerate heat! Please check with the manufacturer of your objective if it can be heated.
- **NOTE** – The Temperature Controller controls the temperature of the Objective Heater according to the sensor inside the heating unit of the Objective Heater. As the Objective Heater dissipates heat to the objective and the sample, the temperature of the Objective Heater is typically set higher than the target temperature of the sample. In order to adjust the temperature within the sample to a specific target temperature follow the instructions in Section 5.3.

### **Installation and Connection of the Objective Heater Universal – Silver Line**

Objectives come in different sizes and lengths, depending on manufacturer and magnification. Hence, there are two different ways to mount the Objective Heater to a microscope objective (Figure 2).

1. Test if the heating unit creates a significant overlap when wrapped around the objective. If there is no overlap, proceed with steps 2A and 3A. If there is an overlap, which is not optimal for homogenous heating, proceed with steps 2B and 3B.
  - 2A. Thread the blue end of the Velcro® fastener through the metal buckle and create a loose loop that can be easily placed over the objective.
  - 3A. Bend the heating unit of the Objective Heater around the objective, the cable directing towards the revolving nosepiece (turret) of the microscope. Next, place the Velcro® fastener around the heating unit while holding it in place, and secure the Objective Heater in the desired position by closing the Velcro®.
  - 2B. Thread the blue end of the Velcro® fastener as well as one end of the heating unit carefully through the metal buckle and create a loose loop that can be easily placed over the objective.
  - 3B. Place the Objective Heater around the objective, the cable directing towards the revolving nose-piece (turret) of the microscope. Fasten the Velcro® fastener and heating unit around the objective so that there is enough contact between heating unit and objective. Close the Velcro®.
2. Connect the electrical cable of the Objective Heater to port 3 or 6 at the back of the Temperature Controller.




Figure 2: Mounting of the Objective Heater to an objective.

**Figure 2:** Mounting of the Objective Heater to an objective. To control the temperature of the Objective Heater, please follow the instructions of the Temperature Controller (see Instructions). Optional: To setup communication with the IncuControl Software, use a USB cable to connect the Temperature Controller to the computer.

### Emergency Release Mechanism

A break-away connector in the cable between the Objective Heater and the Temperature Controller allows for emergency release of the Objective Heater from the Temperature Controller (see Figure 3). If a force above about 14 N acts on the cable of the Objective Heater, the break-away connector releases open to prevent any force from being transferred to the microscope's objective.




Figure 3: Principle of the emergency release mechanism by a break-away connector in the cable between the Objective Heater and the Temperature Controller.

### Temperature Adjustment of the Objective Heater Universal – Silver Line

ibidi recommends a temperature adjustment before starting the first experiment to ensure that the sample reaches the correct temperature. The temperature in the sample is affected primarily by the temperature of the Stage Top Incubator. In addition, the temperature of the microscope objective can contribute to the sample temperature by heat transfer from the sample to the objective and vice versa. The temperature of the Objective Heater is regulated by the Temperature Controller. Please follow the instructions for temperature adjustment of the sample as described in the instructions of the Temperature Controller (see Instructions).

1. Adjust temperatures of Heated Plate and Heated Lid in a way that the target temperature for the sample is reached (e.g., 37°C), while the objective is not yet in contact with the sample.
2. As soon as the target temperature of the sample is reached, bring the cold objective in contact with the sample. The temperature of the sample will most likely drop. Wait until the temperature of the sample reaches a stable plateau.
3. In order to bring back the sample temperature to the target temperature, heat up the Objective Heater to an initial temperature between 38°C and 45°C. The optimal temperature depends on the diameter and material of your objective. Wait until the temperature is stable again.
4. Adjust the temperature of the Objective Heater, if necessary. We recommend beginning the temperature adjustment with a low temperature and increasing the temperature stepwise by 0.5–1 °C.




Figure 4: Temperature adjustment procedure to adjust the sample temperature to the target temperature. The temperature adjustment must be repeated from time to time (at least once a year), especially if one of the following conditions has been changed:

- Room temperature
- Air conditioning
- Chamber type or objective lens
- Humidification and gas flow
- Use of an XL-incubator

#### **ibidi GmbH**

Lochhamer Schlag 11 82166 Gräfelfing Germany


**Toll free within Germany:** Phone: 0800 / 00 11 11 28 Fax: 0800 / 00 11 11 29

#### **International calls:**

- **Phone:** +49 89 / 520 46 17 – 0 Fax: +49 89 / 520 46 17 – 59
- **E-Mail:** [info@ibidi.com](mailto:info@ibidi.com)
- [ibidi.com](http://ibidi.com)

---

#### **Documents / Resources**

|   |  |
|---|--|
|  | <p><b><a href="#">ibidi 12170 Objective Heater Universal</a></b> [pdf] Instruction Manual</p> <p>12110, 12720, 12722, 12111, 12170 Objective Heater Universal, Heater Universal, Heater, 12170, 12170 Heater Universal, Objective Heater Universal</p> |
|---|--|

References

- [🌐 12131.com](#)
- [🌐 comibidi.com.com](#)
- [🏠 ibidi – cells in focus](#)
- [🏠 ibidi – cells in focus](#)
- [🏠 ibidi – cells in focus](#)