



IrisVision Wearable Low Vision Glasses for Visually Impaired User Guide

[Home](#) » [IrisVision](#) » IrisVision Wearable Low Vision Glasses for Visually Impaired User Guide 

Contents

- 1 [IrisVision Wearable Low Vision Glasses for Visually Impaired](#)
- 2 [SCHEDULE COACHING](#)
- 3 [CONNECT TO Wi-Fi](#)
- 4 [CONNECT HEADSET TO PHONE](#)
- 5 [BUTTON MAPPING](#)
- 6 [Documents / Resources](#)
 - 6.1 [References](#)
- 7 [Related Posts](#)



IrisVision Wearable Low Vision Glasses for Visually Impaired



Please read this guide before operating your headset and keep it for future reference.

SCHEDULE COACHING


You should receive an email to schedule your coaching sessions. If you have not received an email, please call Customer Support at (855) 207-6665 (Monday- Friday 6 am- 5 pm PST/Saturdays 8 am-1 pm PST).

CONNECT TO Wi-Fi

Connect your Samsung Galaxy S21 phone to a WiFi network using the phone's settings. Please contact Customer Support BEFORE your first coaching session to ensure your headset is connected to WiFi.

CONNECT HEADSET TO PHONE


Plug the cable attached to the headset into the phone.



LAUNCH APP ON PHONE


Power the phone on by pressing the lower button on the right side of the device for 3 seconds Use the folding stand-up charging pad to charge your Inspire device.(Refer to User Guide)

BUTTON MAPPING



1. Press **Button 1** to: navigate forward Close the menu Exit current status
2. Press **Button 2** to: Open a menu item Select a feature
3. Press **Button 3** to: Scroll up in the vertical menu Zoom out
4. Press **Button 4** to: Scroll down in vertical menu Zoom in

Documents / Resources

	IrisVision Wearable Low Vision Glasses for Visually Impaired [pdf] User Guide Wearable Low Vision Glasses for Visually Impaired
---	--

References

- [IrisVision | Wearable Low Vision Glasses for Visually Impaired](#)