

KERN ORA 3AA-AB Beer Refractometer Analogues User Manual

Contents

1 KERN ORA 3AA-AB Beer Refractometer Analogues

2 Product Information

3 Product Usage Instructions

4 Documents / Resources

4.1 References

KERN ORA 3AA-AB Beer Refractometer Analogues

Product Information

Model	Measuring Range	Scales	Dimensions	Net Weight
KERN ORA 3AA/AB	0-32% Brix, 1.000-1.130 SG wort	0.2% Brix, 0.001 SG wort	178x40x40mm	0.145kg
KERN ORA 4AA/AB	Not specified	Example scale: ORA 3AA/AB	189x40x40mm	0.155kg

Product Usage Instructions

General Information

The refractometer is a measuring instrument used to determine the refractive index of transparent substances in the liquid or solid state. It observes the behavior of light as it passes from a prism with known properties to the substance being tested. It is important to use the refractometer only for its intended purpose and not for any other purposes, as it may be hazardous. The manufacturer will not be liable for any damages caused by improper use.

Warranty

The warranty will be void if:

- The instructions in the operating manual are not followed
- The refractometer is used for purposes other than those described
- The device housing is modified or opened
- Mechanical damage and/or damage resulting from media, liquids, or natural wear and tear occur

Basic Safety Information


It is important to follow the instructions in the operating manual. Even if you have prior experience with KERN

refractometer, it is recommended to carefully read through the operating manual. Please note that every language version includes anon-authoritative translation, and the original German document is the definitive version.

Warning

Acids should not come into contact with skin or eyes. If acid comes into contact with skin, flush with copious amounts of water and shower if larger areas of skin are affected. If acid comes into contact with eyes, keep the eyelid open and flush the eye with running lukewarm water from the outer corner to the inner corner for at least 15 minutes. Then consult a doctor or ophthalmologist immediately. Thoroughly clean the refractometer after each use.

Documents / Resources

	<p>KERN ORA 3AA-AB Beer Refractometer Analogues [pdf] User Manual ORA 3AA-AB Beer Refractometer Analogues, ORA 3AA-AB Beer, Refractometer Analogues, R efractometer</p>
---	---

References

-  [Your Name as Your Email | Hover Realnames](#)
-  [KERN & SOHN: Waagen, Mikroskope & Messinstrumente](#)