REF MDS7001B

Blood Pressure Monitor with Bluetooth®

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Records Medical Data via Medline® Health App

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TABLE OF CONTENTS

INTRODUCTION

General Description	1
Intended Use	2
Contraindications (When NOT to Use)	2
Measurement Principle	2
Safety Information	3 – 11

GETTING STARTED

Monitor Components	12
Power Supply and Charger	13 – 15
Setting the Time, Date & Measurement Unit	16 – 17
Using Bluetooth® & the Medline® Health App	18 – 22
Pairing the Blood Pressure Monitor with Your Device	22 – 24
Attaching the Cuff	24 – 25

MEASUREMENT

Starting the Measurement	25 – 26
--------------------------	---------

DATA MANAGEMENT

Recalling Saved Measurements	27
Deleting Saved Measurements	27

INFORMATION FOR THE USER

Tips for Measurement	28
Maintenance	29

ABOUT BLOOD PRESSURE

What are systolic pressure & diastolic pressure?	30
What is the standard blood pressure classification?	30
Irregular Heartbeat Detector	31
Why does my blood pressure fluctuate throughout the day?	32
Why do I get a different BP at home compared to the hospital?	32
Is the result the same if measuring on the dominant arm?	33

TROUBLESHOOTING & FREQUENTLY ASKED QUESTIONS

	34 – 35
--	---------

SPECIFICATIONS

	35
--	----

CONTACT INFORMATION

	36
--	----

FCC STATEMENT

	37 – 38
--	---------

COMPLIED STANDARDS LIST

	39
--	----

EMC GUIDANCE

	40 – 44
--	---------

WARRANTY

	45
--	----

INTRODUCTION

General Description

The Medline Blood Pressure Monitor connects to the Medline Health App via Bluetooth® and features blood pressure measurement, pulse rate measurement and storage of results.

Readings taken are equivalent to those obtained by a trained observer using the cuff and stethoscope auscultation method.

Read this manual thoroughly before using the product.

This manual contains important safety and care information, and step-by-step instructions for using the product.

Features

- 86.1 x 24mm Digital LCD display
- Measure-during-inflating Technology
- Stores 120 measurements total, 60 for each of two users.

Intended Use

This Blood Pressure Monitor is a digital monitor intended for use in measuring blood pressure and heartbeat rate with arm circumference ranging from 22 cm to 42 cm (about 8¾" - 16½"). It is intended for indoor, adult use only.

Contraindications (When NOT to Use)

1. The device should not be used by any person who may be suspected of, or is pregnant.
2. The device is not suitable for use on patients with implanted, electrical devices, such as cardiac pacemakers and defibrillators.

Measurement Principle

This product uses the Oscillometric Measuring method to detect blood pressure. Before every measurement, the unit establishes a "zero pressure" equivalent to the atmospheric pressure. Then it starts inflating the arm cuff while simultaneously detecting pressure oscillations generated by beat-to-beat pulsatile readings. These are used to determine the systolic and diastolic pressure, and also pulse rate.

Safety Information

The images defined in the following chart may be used in this user's manual, labeling or other components to aid in communication of safety and usage information.

	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol for "TYPE BF APPLIED PARTS"
	Symbol for "Class II Equipment"		Symbol for "ENVIRONMENT PROTECTION - Electrical waste products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice."
	Symbol for "MANUFACTURER"		
	Symbol for "SERIAL NUMBER"		
	Symbol for "DIRECT CURRENT"	F1	IH=05A,IT=1A, FUZETEC
	Symbol for "MANUFACTURE DATE"		CAUTION: These notes must be observed to prevent any damage to the device.
	Symbol for "RECYCLE"		For indoor use only

CAUTION:

The device is not intended for patient transport outside a healthcare facility.

The device is not suitable for use on neonatal patients, pregnant women, patients with implanted electronic devices, patients with pre-eclampsia, premature ventricular beats, atrial fibrillation, peripheral arterial disease, arteriovenous shunt, and patients undergoing intravascular therapy or have received a mastectomy. Please consult your doctor prior to using the unit if you suffer from these or other illnesses.

The device is not suitable for measuring the blood pressure of children. Ask your doctor before using it on older children.

This device is intended for non-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement. This unit performs to the accuracy standards of a clinical product.

Do not confuse self-monitoring with self-diagnosis.

This unit allows you to monitor your blood pressure.

Do not begin or end medical treatment without asking a physician for treatment advice.

If you are taking medication, consult your physician to determine the most appropriate time to measure your

blood pressure. Never change a prescribed medication without consulting your physician.

Do not take any therapeutic measures on the basis of a self-measurement. Never alter the dose of a medicine prescribed by a doctor. Consult your doctor if you have any questions about your blood pressure.

When the device is used to measure patients who have common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, the best result may occur with deviation. Please consult your physician about the result.

When using this device, please pay attention to the following situations which may interrupt blood flow, influence blood circulation of the patient, and cause harmful injury:

- The use of the cuff on any arm where intravascular access or therapy, or an arteriovenous (A-V) shunt, is present
- Inflating the cuff on the side of a mastectomy

WARNING: Do not apply the cuff over a wound, or it can cause further injury.

Do not inflate the cuff on the same limb which other monitoring ME equipment (medical electrical equipment) is applied at the same time. This could cause temporary loss of function of simultaneously-used monitoring ME equipment.

On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, open the cuff immediately. Prolonged high pressure (cuff pressure 300mmHg or constant pressure 15mmHg for more than 3 minutes) applied to the arm may lead to an ecchymosis. Always be sure that operation of the device does not result in prolonged impairment of patient blood circulation.

The device cannot be used with HF surgical equipment at the same time.

CAUTION:

This **SPHYGMOMANOMETER** was clinically investigated according to the requirements of ISO 81060-2:2013.

Besides providing inaccurate readings, the effects of this device on the fetus are unknown. Too frequent and consecutive measurements could cause disturbances in blood circulation and injuries.

This unit is not suitable for continuous monitoring during medical emergencies or operations. Such use will cause the patient's arm and fingers to become anaesthetic, swollen and even purple due to a lack of blood.

When not in use, store the device with the adapter in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the storage case.

This device may be used only for the purpose described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application.

This device contains sensitive components and must be treated with caution. Observe the storage and operating conditions described in this booklet.

The user should not have a skin temperature that exceeds 42.8°C (109°F). The temperature of the room in which this device is used should not exceed 40°C (104°F).

This device is not AP/APG equipment and not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

WARNING: Do not service or perform maintenance on this ME equipment while in use.

The patient is an intended operator.

The patient can measure, transmit data and recharge device under normal circumstances and maintain the device and its accessories according to the user's manual.

To avoid measurement errors, please avoid conditions with strong electromagnetic field radiated interference signals or electrical fast transient/burst signals.

The blood pressure monitor, its adapter, and the cuff are suitable for use within the patient environment. If you are allergic to polyester, nylon or plastic, please don't use this device. During use, the patient will be in contact with the cuff.

The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010.

It will not cause any potential sensation or irritation reaction.

Adapter is specified as a part of **ME EQUIPMENT**.

If user experiences discomfort during a measurement, such as pain in the arm or other complaints, press any button to release the air immediately from the cuff. Loosen the cuff and remove it from your arm.

If the cuff pressure reaches 40 kPa (300 mmHg), the unit will automatically deflate. Should the cuff not deflate when pressure reaches 40 kPa (300 mmHg), detach the cuff from the arm and press any button to stop inflation.

Before use, make sure the device functions safely and is in proper working condition. Do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.

Cleaning

Do not wash the cuff in a washing machine or dishwasher!

Dust in the environment may affect the performance of the unit. Please use a soft cloth to clean the whole unit before and after use. Don't use any abrasive or volatile cleaners.

CAUTION:

Please dispose of accessories, detachable parts, and the ME Equipment according to the local guidelines.

The operator shall not touch the adapter and the patient simultaneously.

If you have any problems with this device, such as setting up, maintaining or using, please contact Medline customer service. Don't attempt to or repair the device by yourself in the event of malfunctions.

Please report to Medline if any unexpected operation or events occur.

Keep the unit out of reach of infants, young children or pets to avoid inhalation or swallowing of small parts. It can dangerous or even fatal.

Be careful of strangulation hazards due to cables and hoses, particularly due to excessive length.

ME equipment should warm for at least 30 min. from the minimum storage temperature between uses until it is ready for intended use. ME equipment should cool for at least 30 min. from the maximum storage temperature between uses until it is ready for intended use.

Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, and walkie-talkies can affect this equipment and should be kept at least a distance away from the equipment. The distance is

calculated by the manufacturer from the 80MHz to 5.8 GHz column of Table 4 and Table 9 of IEC 60601-1-2:2014, as appropriate.

Please use accessories and detachable parts specified/ authorized by manufacturer. Otherwise, it may cause damage to the unit or danger to the users/patients.

Please use the device under the environment defined in the user's manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.

LCD Display

SYMBOL	DESCRIPTION	EXPLANATION
SYS	Systolic Blood Pressure	High Blood Pressure
DIA	Diastolic Blood Pressure	Low Blood Pressure
Pul min	Pulse Display	Pulse in Beats Per Minute
+Lo	Low Battery	Low battery and please charge the power.
kPa mmHg	Unit	Measurement Unit of the blood pressure (1mmHg=0.133kPa)
	Irregular Heartbeat	Blood pressure monitor is detecting an irregular heartbeat during measurement.
	Bluetooth Icon	Bluetooth icon blinks when the Bluetooth is working.
LAST 3 AVG.	Average Value	The average value of the latest three groups blood pressure value.
	Memory Query	Indicate it is in the memory mode and which group of memory it is.
	User 1/User 2	Start measurement for User 1/User2
	Current Time	Year/Month/Day (Hour: Minute)
	Motion Indicator	Motion may result in an inaccurate measurement.
	Heartbeat	Blood pressure monitor is detecting a heartbeat during measurement.
	Blood Pressure Level Indicator	Indicates the blood pressure level

GETTING STARTED

Monitor Components

1. Blood Pressure Monitor (MDS7001B)

2. AC Adapter

4. CUFF (22 to 42cm) Type BF Applied Part

(Please use Medline authorized cuff. For sizing, refer to the label on the cuff.)

Power Supply and Charger

1. This device comes with a built-in rechargeable li-polymer battery with a battery current of 1000 mAh.
2. Only use the AC adapter included to charge the battery, as shown below.

When to charge:

 (+Lo) Low power displays on the LCD.

The LCD display is dim.

When powering on the monitor, the LCD doesn't light up.

CAUTION:

The built-in rechargeable lithium-ion battery should not be disassembled by unauthorized maintenance personnel.

If the charger is not charging the battery or the blood pressure monitor is not working normally, please contact Medline customer service.

Store and use the blood pressure monitor in a cool, dry and ventilated environment. Batteries (battery pack or batteries installed) should not be exposed to excessive heat such as fire, hot sunlight, etc.

Only use the Medline authorized AC Adapter to recharge. **DO NOT** use the blood pressure monitor during recharging.

During the recharging, the blood pressure monitor will display .

Unplug the adapter when charging is finished.

Never touch the charging connector and the patient simultaneously.

Do not attempt to replace your blood pressure monitor's battery. It is built-in and not changeable. Only charge the battery in accordance with the user instructions supplied with the blood pressure monitor.

CAUTION:

Avoid charging your blood pressure monitor in extremely high or low temperatures.

Do not use your blood pressure monitor while you are charging it.

Do not attempt to disassemble the blood pressure monitor or force open the built-in battery.

Do not clean the blood pressure monitor when it is being charged. Always unplug the charger first before cleaning the blood pressure monitor.


Do not dispose of your blood pressure monitor in a fire. The battery could explode causing injury or death.

Setting the Time, Date & Measurement Unit

To ensure the stored measurement result has correct time record, please set time and measurement before device is used.

Before use, turn the monitor “ON” by using the switch on top of the unit.

1. When the monitor is “ON”, press and hold  “User 1” button for 5 seconds to start setting the date. The year will appear first.

2. As shown below, the blinking numeral representing [YEAR] will appear. Press  “Query” button to change the numeral. Each press will increase the numeral by one in a cycling manner.

There is no button to go backwards/decrease a value.



3. Press “User 1” button again to confirm [YEAR]. Then the numeral representing [MONTH] blinks.


4. Repeat step 2 and 3 to confirm [MONTH] and [DAY].

NOTE: This is arranged [DAY]/[MONTH]


5. Repeat step 2 and 3 to confirm [HOUR] and [MINUTE]. Time set to 24 hour day. No AM/PM setting.

6. Repeat step 2 and 3 to confirm the unit of measurement: mmHg (0, in US), kPa (0.0, in most other countries)

7. After confirming the unit, the LCD will display “DONE”, then the monitor display will go dark.

NOTE: Changing the ON/RESET switch to “RESET” will clear any previous Date and Time settings. We recommend leaving the switch to “ON”.

Using Bluetooth® & the Medline® Health App

This device is designed with Bluetooth compatibility and use with the Medline Health App on iOS and Android platforms. It is available for consumers to download via the  and  free of charge. Users will be able to record the following values manually or using automation via Bluetooth connectivity with this and other Medline approved medical devices, including, but not limited to:

1. Blood Pressure
2. Heart Rate
3. Weight
4. Temperature
5. Oxygen Saturation
6. Blood Glucose

The data will be stored locally on the user's device and user credentials will be required to log into the app each time the user logs out of the application.

HERE'S HOW TO GET STARTED

User Registration

When a user first downloads the Medline Health App they will be required to register as a user by entering their First Name, Last Name, Email Address and creating

a password. Non-required fields that are available during the registration phase also include DOB, height and profile image. If the user has received a Provider Code from a Healthcare Provider, they will need to enter it when signing up. This will allow the user's health data to be shared with their Healthcare Provider.

Menu Options

Home

When a user logs into the application, they are taken to the main dashboard. A user can also get to this area by selecting the menu icon, , and choosing the **"Home"** option. From the dashboard, the user has the option to manually enter a value by selecting the **"Manual Entry"** link or use a device that has been paired to the app to wirelessly record a value. If a value is manually entered, it will appear on the dashboard once the save option is selected. If a paired device is used, the results will automatically transfer and display on the dashboard once the device completes the reading. **NOTE:** The Medline Health App must be open on the Smart Device for the Blood Pressure measurement to transfer over Bluetooth.

Devices

The **"Devices"** option on the menu will display a list of connected devices. This is where the user will go to add additional devices. Please refer to page 22 of the user's manual for instructions on how to pair the Bluetooth Blood Pressure Monitor to a Smart Device.

Profile

The **“Profile”** option on the menu allows the user to edit their profile. In this area updates can be made to a user’s profile image, name, email address, date of birth, height, and a Provider Code can also be added. Once a new value is entered, the save button must be selected to update the profile.

Settings

The **“Settings”** option on the menu is where the units of measure can be changed based on the user’s preference. For example, the application defaults to blood pressure readings as “mmHg” but there is an option to change readings to “kPa”.

History

The **“History”** option allows a user to review previously recorded results by measurement type. The user can access this information from the Menu or by clicking one of the dashboard icons on the Home screen. A screen is displayed that will show history, charts and fact check details.

When the history icon at the bottom of the screen is selected, a list of previously recorded values will display in order of date. The user can use the drop down menu at the top to view readings from a different measurement type. These results will show the value recorded, date/time of the recording, and whether or not the value was manually

entered. If the value was manually entered, a symbol, , will display at the bottom of the results box. If a user manually inputs a value incorrectly, they have the ability to edit that entry on this page. By selecting the **“Edit”** option on the reading, the user can change the values.

NOTE: Edits can only occur for results that were manually entered. A user also has the ability to delete an erroneous value whether it was manually entered or recorded via Bluetooth by selecting the **“Delete”** option on the reading that needs to be removed.

When the chart icon at the bottom of the screen is selected, a chart will display showing a trend of recorded results. There is a drop down menu to adjust the graph based on the number of results shown.

The last option is the **“Fact Check”** icon. When this icon is selected, a list of one or more charts will display allowing the user to compare their results with suggested averages in the medical industry.

Links

The **“Links”** option gives the user access to additional useful Medline material such as Terms and Conditions of the app, Bluetooth devices that are compatible with the app and how to purchase them.


Support

The “**Support**” option provides helpful content to app users related to usage and measurement questions as well as healthcare resources that provide guidance around results. The content provided is linked to medical organizations that provide additional information.

Pairing the Blood Pressure Monitor with Your Device

1. Turn on the monitor by sliding the power switch on the top of the monitor to the “**ON**” position. You should also have your smart device on and the **Medline Health App** open.
2. In the app, press the menu button, then “**Devices**”, then “**Add Device**”.  **Devices**
3. Make sure the monitor is “**ON**”, then press and hold the **User 2**  button to start pair-up. The symbols shown in the images to the right will blink on the LCD alternatively to indicate pair-up is proceeding.

Pairing Progress

4. When Bluetooth® pairing is completed, the device ID will appear in the app. Select it and the monitor will synch with the app. The app will display information about the device. The monitor will display the “**Successful Pairing**” icon and return to the sleep mode. In sleep mode, the device is “**ON**”, but the screen is dark.

Successful Pairing

Unsuccessful Pairing

NOTE: The Device ID will begin with “0” or “1” + “MDS7001B” + 6 alternating digits of the monitor’s serial number. The serial number can be found on the left side of monitor.

Bluetooth Module No.: LS51802 RF

Frequency Range: 2402 MHz to 2480 MHz

Output Power Range: ≤0 dBm

Supply Voltage: 1.8-3.6 V

Transmitting Distance: 10 meters

List of Compatible Devices

For iOS devices: The operating system must be iOS 8 or more, such as iPhone 4S, iPhone 5/5C/5S, iPhone 6/6 Plus and so on.

For Android devices: The operating system must be 4.3 or later.

⚠ CAUTION:

MDS7001B may interfere with nearby electrical equipment. Interference may occur in the vicinity of equipment marked with the following symbol .

Keep the monitor at least 20 centimeters away from the human body (**especially the head**) when data transmission is proceeding after measurement.

To enable the data transmission function, this product should be paired to Bluetooth end at 2.4 GHz.

To Mitigate Possible Interference

1. The range between the device and BT monitor should be reasonably close, from 1 meter to 10 meters. Please ensure no obstacles between the device and the BT monitor to help obtain quality connection and to lower the RF output range.
2. To avoid interference, other electronic devices (particularly those with wireless transmission/Transmitter) should be kept at least 1 meter away from the monitor.

Attaching the Cuff

1. Remove all accessories (watch, bracelet, etc) from your non-dominant arm. If your physician has diagnosed you with poor circulation in that arm, use the other one.
2. Cuff must be used only over bare skin. Make sure blood flow in arm is not restricted by clothing.
3. Apply cuff to the arm, user's palm up, and the red mark on the cuff is closest to the elbow, lined up with the brachial artery.
4. Position the edge of the cuff about 1-2 cm from elbow.
5. Fasten the cuff around your arm, leaving no extra room between the cuff and your skin. If the cuff is too loose, the measurement will not be accurate.
6. Sit comfortably with your tested arm resting on a flat surface. Place your elbow on a table so that the cuff is at the same level as your heart.

Helpful tips

Rest for 5 minutes before first measurement.


Wait at least 3 minutes between measurements. This allows your blood circulation to recover. The user must relax as much as possible; do not move or talk during the measurement.

Take the measurement in a silent room.

The cuff should be placed and remain at the same level of the heart.

Sit comfortably. Do not cross your legs and keep your feet flat on the ground. Keep your back against the backrest of the chair.

For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.




MEASUREMENT


Starting the Measurement

1. Turn the monitor ON using the switch on top of the unit. **If using the Medline Health app, make sure the application is running on your smart device.**

- Correctly place cuff on arm, then press “User 1” or “User 2” to begin a measurement.


Display & save the results

The data transmission will proceed.

Data Transmission Succeeds

If the data transmission succeeds, the Bluetooth symbol will appear then the screen will go dark.

Data Transmission Fails


If the data transmission fails, the Bluetooth symbol will not appear then the screen will go dark.

DATA MANAGEMENT

Recalling Saved Measurements

- When the monitor is on, even if the display is dark, please press the “Query” button to show the average value of the latest three records. If the memory holds less than three measurements, it will display the latest record first.
- Press the “Query” button again to view each individual record. The unit can store up to 60 measurements, showing the most recent ones first. Once there are more than 60 readings, the oldest ones will be deleted first.


NOTE: Deleted measurements will still be available on the Medline Health App.

- Press “User 1” or “User 2” to view saved measurements for either user.

Deleting Saved Measurements

- Press and hold the “Query” button for 3 seconds when the monitor is in the memory recall mode. The “dEL all + User ID” will flash on the display.
- Press “Query” to confirm deleting and “dEL donE+User ID” will be displayed on the screen. The monitor will then go into sleep mode, it will not turn off.


NOTE: To exit out of the delete mode, please press the “User 1” or “User 2” button.

INFORMATION FOR THE USER

Tips for Measurement

Measurements may be inaccurate if taken in the following circumstances. **Do Not Take Measurements:**

 Within 1 hour after eating or drinking	 After smoking, or drinking coffee or tea
 Within 20 minutes after bathing	 When talking or moving your fingers
 In a very cold environment	 When you have to use the restroom

Maintenance


In order to get the best performance, please follow the instructions below.

 Store unit in a dry place and avoid direct sunlight	 Avoid dusty environments and unstable temperatures
 Use slightly damp cloth to remove dirt	 Avoid getting the cuff or device wet
 Do not submerge to clean	 Avoid intense shaking or dropping the unit

ABOUT BLOOD PRESSURE

What are systolic pressure & diastolic pressure?

When ventricles contract and pump blood out of the heart, blood pressure reaches its maximum value, the highest pressure in the cycle is known as systolic pressure. When the heart relaxes between heartbeats, the lowest blood pressure is diastolic pressure.

What is the standard blood pressure classification?

The blood pressure classification published by World Health Organization (WHO) and International Society of Hypertension (ISH) in 1999 is as follows:

Blood Pressure Classification	Systolic (mmHg)	Diastolic (mmHg)
Optimal	Less than 120	Less than 80
Normal	120 -129	80 -84
High – Normal	130 -139	85 -89
Hypertension Stage 1	140 -159	90 -99
Hypertension Stage 2	160 -179	100 -109
Hypertension Stage 3	180 or Higher	110 or Higher

⚠ CAUTION:

Consult a physician for your normal blood pressure range. Please contact a physician if your measuring result falls out of the range. Please note that only a physician can tell whether your blood pressure value has reached a dangerous point.

Irregular Heartbeat Detector

An irregular heartbeat (IHB) is detected when a heartbeat rhythm varies while the unit is measuring the systolic and diastolic blood pressure. During each measurement, the monitor records the heartbeat intervals and calculates the average. If any average is larger than or equal to 25%, the irregular heartbeat symbol appears on the display when the measurement results appear.

⚠ CAUTION:

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

Why does my blood pressure fluctuate throughout the day?

An individual's blood pressure varies throughout the day. There are many things that can affect blood pressure including the way the cuff is applied, the measurement position, any medicine the individual is taking, etc. Please make sure to take measurements under the same conditions each day for more reliable results.

If the person takes medicine, the pressure will vary more.

Wait at least 3 minutes between measurements.

Why do I get a different blood pressure at home compared to the hospital?

Having your blood pressure taken in a clinical setting may cause you to become nervous, thus artificially raising your reading.

Having the ability to take your blood pressure at home helps you be more relaxed when taking measurements. Be sure to share your information with your clinician.


Is the result the same if measuring on the dominant arm?

The device is approved to be used on both arms, but the results will vary. For the most consistent readings, measure the same arm every time.

What you need to pay attention to when you measure your blood pressure at home:

- Make sure the cuff is applied properly.
- Make sure the cuff is not too tight or too loose.
- Make sure you are relaxed prior to taking your measurement. Try taking 2-3 deep breaths before the measurement.

TROUBLESHOOTING & FREQUENTLY ASKED QUESTIONS

PROBLEM	SYMPTOM	POSSIBLE ISSUE	REMEDY
No power	Display will not light up.	Battery is exhausted.	Charge the monitor
		AC Adaptor is inserted incorrectly.	Make sure AC adapter connection and wall plug are properly plugged in.
Low batteries	Display is dim or shows  +Lo	Power is low.	Charge the monitor
Error message	E 1 shows	Bluetooth pairing timeout.	Make sure that the smart device's Bluetooth is on and measure again; make sure the Medline Health app is open.
	E 3 shows	The cuff is not secure.	Make sure the cuff is applied properly and remeasure. Please refer to page 24.
	E 10 or E 11 shows	The monitor detected movement while measuring.	Relax for a moment and then measure again. Make sure that the cuff is applied properly. Please refer to page 24.
	E 20 shows	The measurement process does not detect the pulse signal.	Loosen the clothing on the arm and then measure again. Make sure that the cuff is applied properly. Please refer to page 24.
	E 21 shows	The measurement was incorrectly performed.	Relax for a moment and then measure again. Make sure that the cuff is applied properly. Please refer to page 24.
	EEx shows	Battery power is too low for proper inflation.	Charge the unit prior to using. Relax for a moment and then measure again.
Warning message	"OUT" shows	Out of measurement range.	The measurement result is out of the measurement range (SYS: 60mmHg to 230mmHg; or DIA: 40mmHg to 130mmHg; or Pulse: 40-199 pulse/minute). Make sure the cuff is applied properly and try taking the measurement again. Please refer to page 24.
Bluetooth Pairing	Can't find device in Bluetooth settings	Bluetooth device not in pairing mode.	Be sure that your Bluetooth device is in "pairing mode." See page 22. Place Bluetooth device closer to the smartphone being used, usually 33 ft (10m) or less.
	Device appears connected but not to correct smartphone	Bluetooth device out of range. Device could be paired with another smartphone.	1) Disconnect Bluetooth on other smartphone 2) Repeat the pairing process to desired smartphone. See page 22.

PROBLEM	SYMPTOM	POSSIBLE ISSUE	REMEDY
Bluetooth Pairing	Blood Pressure recording is not transferred to Smart Device	Bluetooth is turned off on Smart Device; Medline Health app is not open; Pairing was unsuccessful.	Delete the device from the 'Known Devices' list and reconnect the monitor. Please contact our Product Support team at (800) 289-9793

SPECIFICATIONS

Power supply	3.7V 1000mAh Built-in rechargeable li-polymer battery, 6V 1A AC Adaptor
Display mode	Digital LCD display V.A.86.1mm (L) × 24mm (W)
Measurement mode	Oscillographic testing mode
Measurement range	Rated cuff pressure: 0mmHg ~ 299mmHg (8.0kPa ~ 30.7kPa) Measurement pressure: SYS: 60mmHg ~ 230mmHg (8.0Pa ~ 30.7kPa) DIA: 40mmHg ~ 130mmHg (5.3kPa ~ 17.3kPa) Pulse value: (40 ~ 99) beat/minute
Accuracy	Pressure: 5°C to 40°C within ±3kPa (0.4mmHg) pulse value: ±5%
Normal working condition	A temperature range of: +5°C to +40°C. A relative humidity range of: 15% to 90%, non-condensing, but not requiring a water vapor partial pressure greater than 50 hPa An atmospheric pressure range of: 700 hPa to 1060 hPa
Storage & transportation condition	Temperature: -20°C to +60°C. A relative humidity range of ≤ 93%, non-condensing, at a water vapor pressure up to 50hPa
Measurement perimeter of the arm	About 22cm ~ 42cm
Weight	Approx. 290g
External dimensions	Approx. 130.9mm x 73mm x 29.4mm
Attachment	AC Adaptor, user's manual
Mode of operation	Continuous operation
Degree of protection	Type BF applied part
Protection against ingress of water	IP22: The first number 2: Protected against solid foreign objects of 12.5mm Φ and greater. The second number: Protected against vertically falling water drops when enclosure titled up to 15°. Vertically falling drops shall have no harmful effects when the enclosure is titled at any angle up to 15° on either side of the vertical.
Software version	A01
Device classification	Battery Powered Mode: Internally Powered ME Equipment AC Adaptor charged Mode Class II ME Equipment

CONTACT INFORMATION

For more information about our products, please visit
www.medline.com
or call customer service at **1-800-MEDLINE**.

Manufactured for: Medline Industries, Inc.
Three Lakes Drive,
Northfield, IL 60093 USA

FCC STATEMENT

FCC ID: OU9LS808-B-S

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

⚠ CAUTION: The user is cautioned that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

COMPLIED STANDARDS LIST

Risk Management	EN ISO 14971:2012 / ISO 14971:2007 Medical devices – Application of risk management to medical devices
Labeling	ISO 15223-1:2016 Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. Part 1: General requirements
User's Manual	EN 1041:2008 Information supplied by the manufacturer of medical devices
General Requirements for Safety	EN 60601-1:2006+A1:2013/ IEC 60601-1:2005+A1:2012 Medical electrical equipment – Part 1: General requirements for basic safety and essential performance EN 60601-1-11:2015/IEC 60601-1-11:2015 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance – Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
Electromagnetic Compatibility	EN 60601-1-2:2015/ IEC 60601-1-2:2014 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility – Requirements and tests
Performance Requirements	EN ISO 81060-1:2012 Non-invasive sphygmomanometers – Part 1: Requirements and test methods for non-automated measurement type EN 1060-3:1997+A2:2009 Non-invasive sphygmomanometers – Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems IEC 80601-2-30:2013 Medical electrical equipment- Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers
Clinical Investigation	EN 1060-4:2004 Non-invasive sphygmomanometers – Part 4: Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers ISO 81060-2:2013 Non-invasive sphygmomanometers - Part 2: Clinical validation of automated measurement type
Usability	EN 60601-1-6:2010/IEC 60601-1-6:2010+A1:2013 Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability EN 62366:2008/ IEC 62366-1:2015 Medical devices - Application of usability engineering to medical devices
Software Life-Cycle Processes	EN 62304:2006/AC: 2008 / IEC 62304:2006 Medical device software - Software life-cycle processes
Bio-Compatibility	ISO 10993-1:2009 Biological evaluation of medical devices– Part 1: Evaluation and testing within a risk management process ISO 10993-5:2009 Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity ISO 10993-10:2010 Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization

EMC GUIDANCE

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments.

WARNING: Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment MDS7001B, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

- 1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the excepted service life.
- 2. Guidance and manufacturer's declaration electromagnetic emissions and Immunity.

Table 1

Guidance and manufacturer's declaration – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Table 2

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV input/output lines: ±1 kV	Power supply lines: ±2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5	line(s) to line(s): ±1 kV line(s) to earth: ±2 kV 100 kHz repetition frequency	line(s) to line(s): ±1 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0%UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0%UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0° 0% UT; 300 cycle	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0° 0% UT; 300 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE 1: UT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity				
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.				
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance	
Conducted RF IEC 61000-4-6	150 kHz to 80 MHz: 3 Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150 kHz to 80 MHz: 3 Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distances: $d=0.35 \sqrt{P}$ $d=1.2 \sqrt{P}$	
Radiated RF IEC 61000-4-3	10V/m, 80% Am at 1kHz	10V/m, 80% Am at 1kHz	80 MHz to 800 MHz: $d=1.2 \sqrt{P}$ 800 MHz to 2.7 GHz: $d=2.3 \sqrt{P}$	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.				

Table 3

Guidance and manufacturer's declaration – electromagnetic immunity							
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.							
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Fre- quency (MHz)	Band (MHz)	Service	Modulation	Modu- lation (W)	Distance (m)	IMMU- NITY TEST LEVEL (V/m)
	385	380- 390	TETRA 400	Pulse modulation b) 18Hz	1.8	0.3	27
	450	430- 470	GMRS 460, FRS 460	FM c) ±5kHz deviation 1kHz sine	2	0.3	28
	710	704- 787	LTE Band 13, 17	Pulse modulation b) 217Hz	0.2	0.3	9
	745						
	780						
	810	800- 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18Hz	2	0.3	28
	870						
	930						
	1720	1700- 1990	GSM 1800; CDMA 1900; GSM 1900; DECT, LTE Band 1, 3, 4,25; UMTS	Pulse modulation b) 217Hz	2	0.3	28
	1845						
	1970						
	2450	2400- 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5100- 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
	5500						
	5785						

WARRANTY

Medline Industries, Inc. warrants this digital blood pressure monitor against any defects in material and workmanship for a period of three (3) years after the date of purchase. These uses include, but are not limited to, inflation of the bladder and attachment of the cuff hook and loop closure. Excluding cases of freight damage, tampering, clear abuse, misuse, or accidents, Medline Industries, Inc. will, at its discretion, repair or replace this sphygmomanometer and/or its parts during the warranty period without charge. No representative or person is authorized to assume for us any liability in connection with the sale of the products, of Medline Industries, Inc. This warranty gives you specific legal rights, and you may also have other rights that may vary from state to state. Prior to shipping any product, please contact either Medline Industries, Inc. at **1-800-MEDLINE** or an authorized Medline sales representative with any warranty concerns.

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© 2021 Medline Industries, Inc. All rights reserved/Todos los derechos reservados.
Medline is a registered trademark of/es una marca comercial registrada
de Medline Industries, Inc. V1 RE21GTE