

Leading the Way in Patient Care



WE BELIEVE DESIGN IS ABOUT NOT JUST THE SPACE, BUT THE PURPOSE OF THE SPACE.

From the moment Midmark entered the dental industry, helping advance patient care has been our passion. Today, our number one commitment is helping you and your team focus on what you do best—caring for your patients, from the operatory to the mechanical room.

Midmark instrument processing solutions provide the solid foundation you need for worry-free infection prevention measures based on best practices you can trust.

Infection Prevention Matters

Your reputation has been years in the making. Don't let one dirty instrument destroy it. As antibiotic resistant "superbugs" loom ever larger in headlines, the need for infection control in dentistry has reached a new urgency. Sterilization is at the heart of any effective infection control protocol, so an understanding of sterilization and the process for attaining it are truly essential for every dental practice.



INSTRUMENT PROCESSING ESSENTIALS

Follow established instrument processing best practices, outlined here, to create a smooth workflow from dirty to clean.

This workflow will help contain contamination and maximize the efficiency of your instrument cleaning and sterilization process.



YOUR STERILIZATION CENTER

STERILIZATION GOALS MAY BE UNIVERSAL, BUT YOUR PRACTICE IS UNIQUE.

It's one thing to know the best practices in sterilization and instrument processing, but the particular needs and considerations of your own practice can make implementing them a challenge.

Are you focused on getting more done in less time? Managing large teams? Working within space constraints? Or maximizing patient comfort? We have solutions to help you achieve your goals.

Let us help you turn your instrument processing center into your best marketing tool. Our stylish cabinetry options and layout concepts will make exactly the impression on your patients you desire—professionalism, expertise and style. It's time to turn your instrument processing area into a showcase of efficiency, and let it speak for you every time a patient passes by.

FROM OUR CUSTOMER

"My staff says sterilization is
easier than ever and that the
implementation of the new Midmark
steri-center has increased efficiency
and turnover time. The distinct flow
and separation from dirty to clean
have made an enormous impact
on staff confidence in sterilization
protocol. There is a sense of ease that
comes in knowing the instruments are
kept in a sterile environment
until the time of use."

Robert J Stevenson, DDS Ann Arbor, MI











Every Minute Counts

Consider the many ways space and workflow design will drive efficiencies in your instrument processing center. Midmark designers are ready to help you reduce bottlenecks and optimize equipment placement to ensure you are ready to meet the unique demands of your practice.



Galley Layout (overhead view)



EFFICIENCIES WHERE YOU NEED THEM

- **01** With additional upper-cabinet storage, you can keep your instrument processing area clean and organized. Integrated two- or three-slot glove dispensers are designed to maximize the efficiency of storage space.
- **02** Built-in pantry storage makes locating instruments simple with glass-paneled doors. The durable steel construction and asepsissupportive design of our Synthesis Cabinetry let your staff focus on the demands of the day.
- 03 A hands-free gooseneck faucet minimizes surface contamination when cleaning or rinsing instruments. The high arc-shaped spout provides height to clean large items. A pull-out wand provides the maneuverability for cleaning or rinsing.
- **04** VistaCool[™] provides direct-todrain elimination of sterilizer wastewater, so your staff no longer has to deal with condensate bottles and tanks.










GALLEY LAYOUT

Galley layout keeps traffic flowing with entry and exit from both ends for maximum access to workspaces with fewer steps.

Big Team Practice

Is your practice a large, high-volume model of efficiency? Your instrument processing area must be optimized to match, with every aspect finely tuned to help your staff achieve the excellence you expect.

We can help you make smart, sustainable choices in your sterilization center to ensure optimal patient care while maximizing efficiencies. Rely on Midmark ultrasonic cleaners, sterilizers and cabinetry for clockwork functionality, coupled with streamlined workflow, so you can keep your focus on your patients.



U-Shaped Layout (overhead view)




SOLUTIONS THAT HELP YOU GET THE MOST OUT OF YOUR INVESTMENT


- O1 Ample countertop space allows
 easy separation of dirty and clean
 instruments during processing.
 Your staff will feel confident and
 comfortable as they have room to
 do their jobs properly.
- 02 Choose your sterilizer capacity based on the size of your dental facility. You can choose from three models with different capacities, including a small-capacity, fastercycle sterilizer.
- 03 Staff can work at peak efficiency with the QuickClean® Ultrasonic Cleaner for thorough cleaning every time, which is a critical step before sterilization.
- **04** Additional work space can incorporate pass-through cabinets that improve efficiencies by allowing your staff to access instruments from both sides.












U-SHAPED LAYOUT

U-shaped layout accommodates multiple sterilizers with ease and gives your staff the space they need for the preparation and packaging of instruments.

Small Spaces, Big Ideas

For those limited or hard-to-fit places, Midmark can help design a workspace layout that encompasses all of the essential steps in instrument processing. Don't let space constraints compromise your standards.



L-Shaped Layout (overhead view)



SOLUTIONS THAT MAKE THE MOST OF YOUR LIMITED SPACE

- of sterilizers to free up counter space.
- **02** Decentralized storage takes up less room in the sterilization center so your one-direction workflow doesn't have to be interrupted even if the space is limited.
- **03** A pull-out wrap drawer can be located under the sterilizer as a space-saving solution for instrument wrapping, allowing wrapping and loading in the same space.
- **04** Maximize space with sturdy, full-extension slides designed to help you use all available drawer space, including hardto-reach spaces in the back.
- **05** The Midmark M3[™] sterilizer delivers speed and efficiency in a compact space. One-hand operation allows loading and unloading without the need for extra space.











L-SHAPED LAYOUT

L-shaped layout lets you maximize the space you have without compromising efficiency or safety.


15



Patient-Focused Design

The comfort and safety of every patient matters. Nothing builds trust better than a warm and inviting setting—not just in your operatory, but in your sterilization area as well.

Midmark can help you create an open, attractive space that lets your patients see your professionalism extending beyond the operatory. Turn your sterilization center into your best marketing tool by making it the highlight of new patient tours and the feature photo in your promotional materials. Let your patients see your commitment to protecting their health and safety.



Straight-Line Layout (overhead view)

SOLUTIONS THAT PUT YOUR PATIENTS AT EASE

- **01** The heavy-duty pull-out shelves in the sterilization tower cabinet provide easy access to your equipment.
- 02 Interior signal lighting in glassdoor cabinets is an elegant and effective way to let the team know clean instruments are inside, ready to be used, while the clean storage pantry provides ideal storage for cassettes or pouches.
- **03** Multiple trash-drop options minimize cross-contamination and help your staff maintain a professional appearance in your sterilization center, keeping the dirty side of your workflow clean and neat at all times.
- **04** Red-tinted windows serve as a patient-friendly indicator that dirty instruments are waiting to be cleaned.









STRAIGHT-LINE LAYOUT

Open design layout creates an inviting space, showing your patients you take pride in every

Create Your Space

Planning a new sterilization center or equipping an existing room can be more effective when you have a plan in place that considers your space, patient volume and growth needs.

FIRST 1	ГНП	NGS	FIE	ST
FIRSI	ппп	463	ГІГ	73 1

1. What is the project?

Grab a pen and answer these initial questions:

☐ New construction of an office

3. Total number of operatories?

instrument processing area?

□ Galley

□ L-shaped□ U-shaped

☐ Straight-Line

4. Based on the examples on pages 10–17,

what design would work best in your

	 □ Remodel of an existing office □ Remodel of a sterilization area with no structural changes □ Remodel of a sterilization area with structural changes 								
2.	What is you	ur daily patie Hygienist		Total					
	Average								
	Maximum								
	Goal								

5.	Would you like to showcase your instrument processing area to your patients? ☐ Yes ☐ No
6.	Do or will you have a dedicated instrument or sterilization staff member? ☐ Yes ☐ No ☐ Sometimes
7.	What is the maximum number of staff members who might use the instrument processing area at the same time?
8.	How would you describe your most critical space limitation?

DESIGN RECOMMENDATIONS:

- 1. Divide your instrument processing area into sections for: 1) receiving, cleaning, and decontamination; 2) preparation and packaging; 3) sterilization; 4) monitoring; and 5) storage.
- 2. Create one-direction flow in the space that moves from dirty to clean to sterile to storage.
- 3. Anticipate the volume of patients you want to serve and design the space to meet or exceed it.

A LITTLE MORE DETAIL

 1. Do you use cassettes or trays in your sterilizer? □ Cassettes only □ Trays only □ Combination of trays and cassettes 			3. What method of instrument cleaning do you use? ☐ Manual instrument cleaning ☐ Ultrasonic cleaner ☐ Washer-disinfector						
•	rays only, do you anticipate moving to in the future? 🗆 Yes 🗆 No	4.	Will you be using a water purification system to support your instrument processing equipment? ☐ Yes ☐ No						
NEXT STE Now go bac	PS k through the designs we show on page	es 14	1–23. Mark the features and styles that						

Now go back through the designs we show on pages 14–23. Mark the features and styles that appeal to you. What do you like about them?

Every Minute Counts/Galley Layout
Big Team Practice/U-Shaped Layout
Small Spaces, Big Ideas/L-Shaped Layout
Patient-Focused Design/Straight-Line Layout

EQUIPMENT CONSIDERATIONS

	Ultrasonic	Sterilizers
1 Doctor/1 Hygienist	3 gallon	1 chamber autoclave—11" or greater
1 Doctor/2–3 Hygienists	3–6 gallon	2 chamber autoclaves—1, 11" or greater + 1 cassette-style sterilizer or 2, 11" or greater
2 Doctors/2–4 Hygienists	3–6 gallon	2 chamber autoclaves—1, 11" or greater + 1 cassette-style sterilizer
3 Doctors/3–5 Hygienists	6 gallon	3 chamber autoclaves—2, 11" or greater + 1 cassette-style sterilizer

Configure Your Space

From receiving to storage, your instrument processing space needs to work for you. These layout options are designed to follow the 5-step best practices in instrument processing while giving you the options your practice demands.

STEP 1

RECEIVING, CLEANING AND DECONTAMINATION


STEP 2
PREPARATION AND


PACKAGING

STEP 3
STERILIZATION

STEP 4
MONITORING/STERILITY
ASSURANCE


STEP 5






GALLEY

The galley layout consists of work spaces on two opposing walls with a single traffic lane between. This arrangement allows for easy access and efficient workflow, helping your staff keep it all moving with a linear flow while keeping everything within reach.




L-SHAPED

An L-shaped counter arrangement maximizes use of available space where elbow room is limited. The space you have can be all you need to do instrument processing right.



U-SHAPED

Multiple cleaners and sterilizers demand space, and a U-shaped workspace design provides that and more. Ample surface areas allow a larger staff in the room to multi-task and maintain a bustling workflow.



STRAIGHT-LINE

Perfectly suited to the 5-step flow that is instrument processing, a straight-line workspace design is the picture of efficiency. Use it to maximize efficiency, and let it create an open and impressive visual for patients passing by.

Products and Solutions

94% of dental service technicians would recommend a Midmark sterilizer.¹

¹2016 Vennli Survey, Service Technicians. Percentage based on 215 survey responses

"I've been in the industry for 16 years and I only recommend Midmark M9 and M11 sterilizers to my clients.

They rate a '10' for initial quality, ease of installation, product perception/ image of quality, and overall satisfaction. Midmark is always there for us."

Charlie Lenart


Field Service Technician


The Critical First Step

QUICKCLEAN® ULTRASONIC CLEANERS

Providing powerful, effective cleaning and consistent results, the QuickClean Ultrasonic Cleaner eliminates hidden residues manual cleaning may miss. If an instrument is not properly cleaned before being sterilized, the sterilization process can be compromised and ineffective. The advanced technology of the QuickClean can help you create a safer, more efficient work environment by decreasing worker exposure to contaminants and potential sharps injuries while reducing the time and effort needed for cleaning.







- O1 Choose the option that best fits your space configuration and workflow needs. QuickClean comes in three tabletop sizes—1.2, 3.3 and 6.6 gallons.
- **02** QuickClean is available in two recessed options—3.3 or 6.6 gallons.

03 The stainless steel finish and attractive design coordinate well with any décor.

Standard Features

SIMPLE OPERATION AND POWERFUL TECHNOLOGY

QuickClean® is easy to use right out of the box, so your staff can be up and running with minimal training. Clean more thoroughly and consistently throughout the entire bath with advanced Frequency-Leap technology. QuickClean helps ensure all of your instruments are fully cleaned the first time, every time, no matter where they are placed in the basket.





Digital LCD control panel displays prompts



Options for preset or customized settings along with built-in data storage retains the previous cycle settings



Fluid sensor detects insufficient water level in the bath and alerts the

ACCESSORIES AND SERVICE PARTS



Beaker Accessories QC1, 2-Beaker Accessory 9A612001 QC3/QC3R, 4-Beaker Accessory 9A613001



Extra Safety Baskets QC1 Safety Basket 002-10007-00 QC3/QC3R Safety Basket 002-10008-00 QC6/QC6R, 6-Beaker Accessory 9A614001 QC6/QC6R Safety Basket 002-10009-00



Cleaning Solutions Midmark General Purpose Cleaner (32 oz) 9A296001 Midmark Tartar and Stain Remover (32 oz) 9A297001 Midmark Enzymatic Cleaner (32 oz) 9A298001

QUICKCLEAN MODELS

	QC1-01	QC3-01	QC3R-01	QC6-01	QC6R-01
Type of Unit	Table Top Unit	Table Top Unit	Recessed Unit	Table Top Unit	Recessed Unit
Capacity	1.2 gal/4.5 L	3.3 gal/12.5 L	3.3 gal/12.5 L	6.6 gal/25 L	6.6 gal/25 L
Size of Unit	13.5" L x 10.2" W x 10.4" H (34.5 cm x 20 cm x 26.5 cm)	14.7" L x 16.5" W x 10.4" H (37.5 cm x 42 cm x 26.5 cm)	15.55" L x 14.37" W (39.5 cm x 36.5 cm)	21.6" L x 16.5" W x 14.3" H (55 cm x 42 cm x 36.5 cm)	22.48" L x 14.37" W (57.1 x 36.5 cm)
Size of Tank	11.8" L x 5.9" W x 5.9" H (30 cm x 15 cm x 15 cm)	13" L x 11.8" W x 5.9" H (33 cm x 30 cm x 15 cm)	13" L x 11.8" W x 5.9" H (33 cm x 30 cm x 15 cm)	19.9" L x 11.8" W x 7.9" H (50.5 cm x 30 cm x 20 cm)	19.9" L x 11.8" W x 7.9" H (50.5 x 30 x 20 cm)
Limited Warranty	3 Years	3 Years	3 Years	3 Years	3 Years
Average Sound Level (with lid on)	68 dB	61 dB	61 dB	69 dB	66 dB

Absolutely Sterile

MIDMARK M3[™], M9[™], M11[™] STEAM STERILIZERS

No dental practitioner can afford to take chances with sterilizer performance. Buy with confidence from Midmark, the market leader in sterilizers year after year. Our sterilizers are designed to be safe, dependable and easy to use, with every one of our sterilizers inspected and ASME-certified by a third-party-licensed inspector. All of our sterilizers are designed, manufactured and supported in the US with over 100 years of manufacturing experience built into each unit.

SAFETY FEATURES

Pre-Programmed Control: Once a pre-programmed cycle is selected, the unit is designed to sterilize automatically without any operator assistance.

Smart Technology: The sterilizer will automatically stop if the door isn't closed completely or the water level is low, and alert the operator to take proper action.

Open Door Drying: Once sterilization is complete, the door opens automatically and quietly to dissipate steam and dry your instruments.

Pre-Programmed Cycle	Sterilization Temperature	(Fill, He	Cycle 1 eat, and minute	Vent;		zation minute			ry Tim minute		Total H w/o Dr					le Time nutes)
	All Units	M11	M9	M3	M11	M9	M3	M11	M9	M3	M11	M9	МЗ	M11	M9	M3
Unwrapped	270°F (132°C)	15	11	25	3	3	3.5	30	30	25	18	14	6	48	44	31
Pouches	270°F (132°C)	17	12	5	4	4	5.5	30	30	30	22	17	10.5	54	46	40.5
Packs/Low Temp	250°F (121°C)	14	10	4.5	30	30	20	30	30	50	44	40	24.5	74	70	74.5
Handpieces	270°F (132°C)	16	11	N/A	4	4	N/A	30	30	N/A	22	17	N/A	60	49	N/A







WHY EASE OF USE MATTERS

Safety Comes First: Simple operation can make it easy for staff to do the right thing every time, minimizing breaks in safety protocols.

Efficiency Wins the Day: When you need it now, turn to the M3 for the rapid cycle time your busy schedule often demands. Utilize the M11 for maximum capacity to handle your patient flow.

Durability Keeps You Running:

Less intuitive sterilizers can wear out from incorrect usage, but Midmark equipment makes it easier to do the right thing every time and is built for years of reliable service.

MAXIMIZING SAFETY AND EFFICIENCY

UNBEATABLE QUALITY AND STREAMLINED MAINTENANCE



01 Ease of use is a hallmark of Midmark technology. Our instrument processing equipment is designed to be intuitive and simple to use right out of the box.






02 Easy-access drain and fill points make maintenance a breeze.

VISTACOOL™ DIRECT-TO-DRAIN WATER ELIMINATION SYSTEM



Kiss your condensate tanks good-bye. The VistaCool direct-to-drain sterilizer wastewater elimination system makes labor-intensive condensate tanks unnecessary, freeing your staff to spend more time chairside and less time tending to equipment.



Single VistaCool unit 9A586001
Designed to be compatible with
the Midmark M3 sterilizer

Double VistaCool unit 9A586002Designed to be compatible with both the Midmark M9 and M11 sterilizers

STERILIZER ACCESSORIES

Choose accessories to tailor your equipment to the particular demands of your practice

ACCESSORIES FOR THE M9 AND M11 STERILIZERS



M9/M11 Printer 9A599001
Records and prints critical sterilization cycle data including time, temperature and pressure



Cool Hand Tool 9A307001

Designed to help reduce the risk of staff injury by making it safer and simpler to load and unload your sterilizer



Pouch Rack, 6 Slot 002-2108-00
Pouch Rack, 5 Slot 002-2108-01
Separates sterilization pouches for enhanced steam circulation and drying

ACCESSORIES FOR THE M11 STERILIZER



Vertical Cassette Rack 9A215002

Holds up to four 8" x 11" cassettes plus four 3" x 8" cassettes together or twelve 4.5" x 8" cassettes

Horizontal Cassette Rack 9A215001

Holds up to four 8" x 11" cassettes plus four 3" x 8" cassettes together or twelve 4.5" x 8" cassettes



ACCESSORIES FOR THE M3 STERILIZER



Top Cover Protector 9A404001Shields the M3 from scratches or damages



M3 Printer 9A401001
Records and prints critical sterilization cycle data including time, temperature and pressure



Door Tray 9A402001

Allows you to prepare an additional load for sterilization before the previous load is complete, saving valuable time

Midmark Cabinets vs. Local Millwork

Why Midmark? We can think of many reasons. Midmark Artizan® Expressions and Synthesis® Cabinetry lines are carefully designed and constructed with dentists in mind. We know the functional features you appreciate, such as full-extension drawers and asepsis-supportive surfaces that are easy to clean and maintain. We understand wanting to support the shop down the road, but local millwork alternatives are unlikely to match Midmark cabinetry solutions for quality, durability and longevity, creating clinical spaces you will be proud to have your patients see for years to come.

CABINET FRAME









- 3/4" industrial-grade particle board
- Can be configured to specific design requirements

02 Synthesis Cabinetry

- 18-gauge cold-rolled steel
- Modular design

03 Local Millwork

- Often made from 1/2" to 5/8" low-density particle board or various types of plywood
- Basic or limited designs

BASE MATERIAL







01 Artizan Expressions

- Seamless stainless steel
- Integrated adjustable levelers

02 Synthesis Cabinetry

- Thermofoil over electrostatic, powder-coated painted steel (woodgrain and metallic colors)
- Electrostatic, powder-coated painted steel (solid colors)
- Integrated adjustable levelers

03 Local Millwork

- Basic plywood or particle board covered with laminate
- Wood shims used for leveling

PANEL SUBSTRATES







01 Artizan Expressions

- Industrial-grade particle board
- 3/4" panels

02 Synthesis Cabinetry

- Medium-density fiberboard
- 3/4" panels

03 Local Millwork

- Low-density particle board
- Various types of plywood

FINISHES







01 Artizan Expressions

- Factory pre-applied, thermally fused melamine
- 3 mm PVC edgebanding

02 Synthesis Cabinetry

- PVC thermofoil (woodgrain and metallic colors)
- PVC thermofoil and electrostatic, powder-coated paint (solid colors)

03 Local Millwork

- High-pressure laminate
- Do not always follow proper application methods
- May use low-grade substrate material
- Less than 2 mm edgebanding
- Varnished, painted or left unfinished

TYPES OF ASSEMBLY







01 Artizan Expressions

100% glue dowel assembly

02 Synthesis Cabinetry

• Mechanical fasteners (screws, pop rivets) and Tog-L-Loc® sheet metal joining system

03 Local Millwork

• Frequently use staple fasteners

DRAWERS







01 Artizan Expressions

- Wood bottom and backs with factory, pre-applied thermally fused melamine
- Metal sides with front panel adjustment and removal feature

02 Synthesis Cabinetry

• Seamless, polystyrene drawers

03 Local Millwork

• Inside drawers are often varnished, painted or left unfinished

DRAWER SLIDES









HANDLES







01 Artizan Expressions

- Full extension ball-bearing drawer system
- Soft-close feature

02 Synthesis Cabinetry

- Full extension ball-bearing drawer system
- Secured to steel foundation
- Soft-close feature

03 Local Millwork

• Typically use a single, undermounted runner (monorail glide) or slides common in kitchens (less durable)







01 Artizan Expressions

• Brushed nickel handles

02 Synthesis Cabinetry

- Integrated and insert handle option
- Brushed nickel pull options
- Antimicrobial option

03 Local Millwork

• Most use "C" style pulls made of plastic or aluminum

Fabrication of local cabinets may vary. However, materials depicted are typical in local millwork cabinetry.

Synthesis® Sterilization Center Highlights

Aligned with the LEED Green Building Rating System, the steel frame of Synthesis Cabinetry provides the durable structure and asepsis-friendly design required for long product life. Its internal and external modularity allows for easy modification and component replacement as office needs evolve.



Gooseneck faucet



Plexiglass-window doors with red and



Seamless polystyrene drawers and seamless panels



Concealed steel hinges



Full-extension ball-bearing glides



Solid-surface countertops



Drawer dividers



Steel-frame construction

Synthesis Color and Style Options

Synthesis Cabinetry offers different color and styling options to allow you to put your unique signature on a space, whether it is in one office or a multi-facility organization.

HANDLE STYLES



PANEL STYLES



^{*}EPA Registration Numbers 84542-7 and 087753-CT-001

^{**}Available in Pebble Grey or Pebble Grey PVC Free only


CABINET FINISH

Membrane-pressed thermofoils seal the Synthesis MDF door and drawer panels. A special process is used to adhere the thermofoil to the steel structure in woodgrain and metallic colors. Solid colors have a matching painted structure. Thermofoil is a dimensionally flexible material that provides a continuous, seamless surface, eliminating the need for panel edgebanding while offering a high level of scratch and stain resistance.

	.
	Frost
	Pebble Grey and Pebble Grey PVC Free
	Whisper
	Fawn
	Venus Silver (metallic)
	Flat Iron (metallic)
	Pearl Essence (metallic)
	Earthen Bronze (metallic)
	Dune
	Nest
	Flax
	Path
	Sediment
	Hewn
	Radiance
THE PARTY OF THE P	Acorn
	Earth
	Timber
	Storm
AND THE RESERVE OF THE PERSON	Henna

COUNTERTOP COLORS

Solid Surface: Solid surface is a unique material made from an advanced blend of natural minerals and acrylic polymer. It has proven itself to be a remarkably durable, versatile material that's an ideal surface for dental environments. Quartz: Quartz is virtually maintenance free and, because it is one of nature's strongest materials, is also incredibly durable, exceptionally tough and scratch resistant. Laminate: Offering a variety selections, laminate countertops provide function while complementing any style of decor.



L - Green Tigris

ARTIZAN® EXPRESSIONS STERILIZATION CENTER HIGHLIGHTS

Artizan Expressions brings a sleek, modern look to the dental space. We believe "best is standard" and use only quality components and materials on every cabinet we build.



Red and blue plexiglass windows



Glue-dowel construction



Concealed Blum® soft-close steel hinges



Brushed-nickel handles and 3-mm edgebanding



Hands-free systems



Full-extension soft-close drawer systems






- O1 Pull-out chrome faucet and HVE collector
- 02 Clean-swiveling pantry tower


CABINET FINISH

Available in a variety of Arauco® and Formica® colors, Artizan Expressions uses thermally fused melamine panels. These panels are created at high temperature and pressure to form a durable, peel-resistant finish. Matching edgebanding is applied to completely seal the corners from moisture.



COUNTERTOP COLORS

Solid Surface: Solid surface is a unique material made from an advanced blend of natural minerals and acrylic polymer. It has proven itself to be a remarkably durable, versatile material that's an ideal surface for dental environments. Quartz: Quartz is virtually maintenance free and, because it is one of nature's strongest materials, is also incredibly durable, exceptionally tough and scratch resistant. Laminate: Offering a variety selections, laminate countertops provide function while complementing any style of decor.



L - Green Tigris

SS = Solid Surface | **Q** = Quartz | **L** = Laminate





What You Get With Midmark



DURABILITY. CRAFTSMANSHIP. PRECISION.

When you choose Midmark, you get so much more than a product. You get a partner with decades of hands-on manufacturing expertise, inspired design innovation, and landmark thought leadership in workflow efficiency. A pioneer of cabinetry designed specifically for the healthcare market since 1967, Midmark equipment and cabinetry are crafted to fit your space and the way you want to work. And with an extensive warranty and dedicated technical support, we're right by your side the moment you need us.



SERVICE.

We stand by the quality and reliability of our Artizan Expressions and Synthesis Cabinetry lines—which is why we offer a standard five year limited warranty for both lines purchased through an authorized Midmark dealer. In addition to warranty.

Midmark has a large team of expert customer experience and technical support representatives who will assist you with any product questions regarding installation and service.



MANUFACTURING EXPERTISE.

Using the knowledge accrued during our more than 100 years of manufacturing history, Midmark continues to invest in its production systems. We've been innovating and building steam sterilizers in our Ohio facility since 1990. To support our growing cabinetry business, we have expanded the Synthesis manufacturing space by over 50,000 square feet.

We have made significant investments in state-of-the-art machinery at our Ohio facility, as well as at our Kentucky facility where we manufacture Artizan Expressions. In all we do, our goal is to bring you the highest quality and performance in dental equipment and cabinetry.



RESOURCES



HOW IT'S CREATED

Check out our videos on how Artizan Expressions and Synthesis Cabinetry are made. See the quality and craftmanship that go into each and every cabinet we produce.

midmark.com/artizanexpressions_creation midmark.com/synthesiscasework_creation



DENTAL DESIGN TOOL

Our easy-to-use, life-like digital tool helps you find your individual style and visualize what your operatory and sterilization space will look like.


apps.midmark.com/dental-design-tool



WEBSITE

For more information, visit our website. Select "Contact a Sales Rep" to ask about our complete line of specialized cabinetry for your sterilization center.

www.midmark.com/cabinetry





Designing better care."

Arauco is a registered trademark owned by Celulosa Arauco y Constitucion S.A. Formica is a registered trademark owned by The Diller Company. VistaCool is a trademark owned by Vista Research Group, LLC. Tog-L-Loc is a registered trademark owned by BTM Company, LLC. Blum is a registered trademark owned by BTM Company, LLC. Blum is a registered trademark owned by Julius Blum GmbH.

CARB 93120.2 Phase 2 Compliant

The color examples shown are the best representation of the original material. Actual color may vary slightly. We strongly recommend that you contact Midmark Customer Experience at 1.800.MIDMARK to request a sample before placing your order.

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: midmark.com/ISO

For more information or a demonstration, contact your Midmark dealer or call: 1.800.MIDMARK Fax: 1.877.725.6495 Outside the USA call: 1.937.526.3662 Fax: 1.937.526.8214 or visit our website at midmark.com

© 2019 Midmark Corporation, Maimisburg, Ohio USA Products subject to improvement changes without notice. Litho in USA 007-10011-00 Rev. D1 (12/19)

