BIPCARE

DIGITAL TENSIO METER / TENSÍÓMETRO DIGITAL

BIPCARE

DIGITAL TENSIO METER / TENSÍÓMETRO DIGITAL

EN

ENGLISH

4 - 19

ES

ESPAÑOL

20 - 35

PT

PORTUGUÊS

36 - 51

FR

FRANÇAIS

52 - 67

IT

ITALIANO

68 - 83

DE

DEUTSCH

84 - 99

ENGLISH

WELCOME

Thank you for choosing our blood pressure monitor. Before using this device, and to ensure the best use, carefully read these instructions.

The safety precautions enclosed herein reduce the risk of death, injury and electrical shock when correctly adhered to. Keep the manual in a safe place for future reference, along with the completed warranty card, purchase receipt and package. If applicable, pass these instructions on to the next owner. Always follow basic safety precautions and accident-prevention measures when using an electrical appliance. We assume no liability for customer failing to comply with these requirements.

This precision instrument has been designed to measure both blood pressure and pulse rate at home and deliver consistent, reliable results. Your new blood pressure monitor is a fully automatic digital blood pressure measuring device for use on the upper arm. It uses the oscillometric method to provide very fast and reliable measurement of the systolic and diastolic blood pressure as well as the pulse. This device offers clinically proven accuracy whilst having been designed to be user friendly.

IMPORTANT SAFETY PRECAUTIONS

- When using any electrical appliance, basic safety precautions should always be observed.
- This device is not safe for neonates, infants and pregnant women.
- This device cannot be used together with high frequency surgical equipment.
- Contact your physician for specific information about your blood pressure. Self-diagnosis and treatment using measured results may be dangerous. Follow the instructions of your physician or licensed healthcare provider
- DO NOT adjust medication based on measurement results from this blood pressure monitor. Take medication as prescribed by your physician. Only a physician is qualified to diagnose and treat High Blood Pressure.
- The monitor is not intended to be a diagnostic device or replace your regular medical check-ups.
- This device is intended for use in measuring blood pressure and pulse rate in the adult population. Do not use this device on infants or persons who cannot express their intentions.
- In cases of irregular heartbeat (Arrhythmia), measurements obtained with this instrument should only be evaluated after consulting with your doctor.
- Familiarize yourself with the section titled "Important Information on Blood Pressure and its Measurement". It contains important information on the dynamics of blood pressure readings and will help you obtain the best results.
- This product, including its accessories, shall be processed and disposed of, in accordance with local regulations after reaching the end of its life cycle.
- This device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens, etc.) during its use. These can lead to erratic results and temporary impairment of the measuring accuracy.
- Do not use cuffs, AC adapters or batteries other than those included with this product or replacement parts supplied by the manufacturer.
- Do not use the batteries and the AC adapter to provide power at the same time.
- Measurement accuracy may be compromised if the device is operated or stored in temperatures or humidity levels outside the limits stated in the specifications of this manual.
- If this appliance is not going to be used for an extended period of time, please remove the batteries.
- The user must check that the equipment functions safely and see that it is in proper working condition before being used.
- Do not attempt to modify, service or repair this device yourself, there are no user serviceable parts inside. Should a

malfunction occur, please contact your local distributor or the manufacturer. Never service the device while still on the patient.

- The device is not suitable for use in the presence of flammable anaesthetic mixtures with air, oxygen or nitrous oxide.
- The user can safely carry out battery replacements, cleaning the device as well as monitoring blood pressure and pulse rate.
- Do not position the power cord in such a way that interferes and complicates its operation or where it may cause people to trip and fall.
- Keep the appliance and the power adapter away from liquids, heat sources, open flames, humidity, direct sunlight and corrosive gas. Do not use in any of the previous environments.
- Keep this unit away from the reach of children and pets to avoid any possibility of strangulation or other injuries. DO NOT wrap the tube around your neck.
- The standard material used for the bladder and tubing is latex-free.
- This device is intended to keep a record and a control, not to be used as a diagnosis tool or for treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.
- The pulse display is not suitable for checking the frequency of heart pacemakers.
- Using this device too frequently on the same patient can lead to injuries due blood flow interference.
- Do not place the cuff over a wound.
- The pressurization of the cuff can cause temporary loss of function if used simultaneously with other ME equipment on the same limb.

CONTRAINDICATIONS

- Use of this instrument on patients under dialysis therapy or on anticoagulants, antiplatelets, or steroids could cause internal bleeding.

IMPORTANT INFORMATION ON BLOOD PRESSURE AND ITS MEASUREMENT

How does high or low blood pressure arise?


Your level of blood pressure is determined in the circulatory centre of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and speed of the heart (Pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is controlled by fine muscles in the blood vessel walls.

Your level of arterial blood pressure changes periodically during heart activity: During the "blood ejection" (Systole) the value is highest (systolic blood pressure value). At the end of the heart's "rest period" (Diastole) pressure is lowest (diastolic blood pressure value).


Blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

Which values are normal?

Please refer to the diagram below.

There are six grids in the display of the device. Please refer to the picture-01-01. Different grids represent different interval scales of the WHO.

Blood pressure value	WHO grids in device	WHO Classification
DIA<80 & SYS < 120	1	Optimal blood pressure
DIA<85 & SYS < 130	2	Normal blood pressure
DIA<90 & SYS < 140	3	High normal value
DIA<100 & SYS < 160	4	mild hypertension
DIA<110 & SYS < 180	5	moderate hypertension
DIA>=110 or SYS>= 180	6	Severe hypertension

Blood pressure is very high if your diastolic pressure is above 90 mmHg and/or your systolic blood pressure is over 160 mmHg, while at rest. In this case, please consult your physician immediately. Long-term values at this level endanger your health due to continual damage to the blood vessels in your body.

If your systolic blood pressure values are between 140 mmHg and 159 mmHg and/or the diastolic blood pressure values between 90 mmHg and 99 mmHg, consult your physician. Regular self-checks are necessary.

If you have blood pressure values that are too low, (i.e., systolic values under 105 mmHg and/or diastolic values under 60 mmHg), consult your physician.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your physician. Never use the results of your measurements to independently alter the drug doses prescribed by your physician.

Further information:


- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called "labile hypertension." Consult your doctor.
- Correctly measured diastolic blood pressure values above 120mmHg require immediate medical treatment.

What can be done if regular high or low values are obtained?

1. Consult your doctor.
2. Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result from arteriosclerosis. Furthermore, the heart will become structurally damaged with increased blood pressure values.
3. There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organ malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
4. There are other measures (e.g. diet and lifestyle changes) which you can take to reduce and even prevent high blood pressure.


COMPONENTS OF YOUR BLOOD PRESSURE MONITOR

Measuring unit


Pic. 1

Pic. 2

Pic. 3

Pic. 4

The symbols on the LCD display

- | | |
|--|--|
| 1. Systolic blood pressure | 9. Mute symbol |
| 2. Diastolic blood pressure | 10. Average value symbol |
| 3. Irregular heartbeat (arrhythmia) symbol | 11. Low battery symbol |
| 4. Memory symbol | 12. Cuff self-checking function |
| 5. Date/Time display | 13. Pulse display |
| 6. USER A | 14. Heartbeat symbol (Flashes during measurement) |
| 7. USER B | 15. WHO Function symbol (Please, refer to the interval scales table - picture-01-01) |
| 8. Movement error symbol | |

Features

1. Talking function / 2. Two users: 2 x 120 sets memory / 3. Date/time display / 4. Irregular heartbeat checking / 5. Average value function / 6. Low battery display / 7. WHO scale feature / 8. Auto power-off / 9. External power adapter support / 10. Volume adjustment / 11. Cuff self-checking function

Note: The arm circumference should be measured with a measuring tape around the middle of the relaxed upper arm. Do not force the cuff connection into the opening. Make sure the cuff connection is not pushed into the AC adapter port.

USING YOUR MONITOR FOR THE FIRST TIME

Installing the batteries

This device uses four 1.5V "AAA" alkaline batteries. Do not attempt to insert any other type of batteries.

- Press the hook on the bottom of the battery cover and lift the cover off in the direction of the arrow.
- Install the four "AAA" batteries so that the + (positive) and - (negative) polarities match the polarities of the battery compartment, and put the battery cover back on. Make sure that the battery cover is securely in position.

Replacing the batteries

- When the "Low battery" indicator appears on the display, turn the monitor off and remove all of the batteries. Replace them with 4 new batteries at the same time. We recommend using long-life alkaline batteries.
- Batteries may leak. To prevent damaging the monitor, please take the batteries out if the monitor isn't going to be used for a long time (generally for more than three months). Should battery fluid get in your eyes, immediately rinse with plenty of clean water and contact your family doctor.

Dispose of this device and all of its accessories according to local regulations to avoid unnecessary pollution. Do not throw away the batteries with your regular trash.

System settings

Once you have installed the batteries or have connected the device to the mains, long press the SET button for more than 3 seconds to power it on. Use the MEM button to increase the values and the SET button to switch between the different parameters.

Setting the user:

Press the MEM button to select User A or User B. When the screen displays A(/B), press the MEM button to switch between the users. To confirm, press the SET button.

Setting the year:

When the display is flashing.

Setting the Month/Day:

The initial Month/Day displayed is 1/01. When the month is flashing press the MEM button to increase the number by 1. Press the SET button to confirm and repeat the same process for the days.

Setting the time:

When the hours are flashing, press the MEM button to increase the hours by one, press SET to confirm. Repeat the same process for the minutes.

Setting the volume (where applicable):

Some models have an incorporated speaker. When the screen displays SP flashing, press the MEM button to select volume 1, 2, 3 or OFF, then press SET to confirm.

Deleting the record:

To delete existing user data, long press the MEM button.

Note: You can't delete all of the measurement records from the monitor storage at once, if you decide to delete all the recorded data, please make sure to keep all of the recorded information elsewhere for future medical reference. Taking the batteries out won't delete any data.

Connecting the cuff tube

Insert the cuff tube into the port located on the left side of the monitor.

HOW TO TAKE YOUR BLOOD PRESSURE

You should always take your blood pressure seated and in a relaxed state.

Connecting the cuff tube

- Avoid eating, smoking and any form of exertion immediately before taking your blood pressure. These factors can influence the results. Find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before taking a measurement.
- Remove any clothes covering your upper arm.
- Always measure on the same arm (usually your left).
- Always compare measurements taken at the same time of day, since blood pressure changes during the course of the day, as much as 20-40 mmHg.

Common sources of error

Note: Always comparable blood pressure measurements taken under the same conditions.

- Blood pressure should be taken in a quiet environment.
- The arm shall be lying on a surface. Any efforts by the user to support the arm can increase blood pressure. Make sure the patient is in a comfortable, relaxed position and not flexing any of the muscles whilst the measurement is taking course. Use a cushion for support if necessary.
- If the arm lies considerably lower or higher than heart level, you will obtain an erroneous high or low blood pressure. Each 25-30cm difference in height between your heart and the cuff will result in a measurement error of 10 mmHg.

- Cuffs that are too narrow or too short will result in false measurement values. It is extremely important to select the correct cuff size. The cuff size will dependent on the circumference of the arm (measured in the center). The permissible range is printed on the cuff.

The cuff works under a pressure range of 0-299MMHG


- The measurement range of the cuff is: 8.7" – 15.7" (22 - 40 cm)

Note: Only use approved cuffs.

- A loose cuff or a sideways protruding air pocket causes false measurement values.
- With repeated measurements, blood accumulates in the arm leading to false results. Consecutive blood pressure measurements should be repeated after a 1 minute pause or after your arm has been held up in order to allow the accumulated blood to flow back to the rest of your body. If you decide to take your Averaging Mode measurement again, make sure you wait at least one minute beforehand.

Fitting the cuff

- The cuff is pre-formed for easier use. Remove tight or bulky clothing from your upper arm.
- Wrap the cuff around your upper left arm. The rubber tube should be on the inside of your arm extending downward to your hand. Make sure the cuff lies approximately 1/2" to 3/4" (1 to 2 cm) above the elbow. Important! The on the edge of the cuff (Artery Mark) must lie over the artery which runs down the inner side of the arm.
- To secure the cuff, wrap it around your arm and press the hook and loop closure together.
- There should be a little bit of space free between your arm and the cuff. You should be able to fit 2 fingers between your arm and the cuff. Cuffs that don't fit properly result in false measurement values. Measure your arm circumference if you are not sure of proper fit.
- Lay your arm on a table (palm upwards) so that the cuff is at the same height as your heart. Make sure the tube is not kinked.
- Remain seated quietly for at least two minutes before you start the measurement.

Measuring Procedure

The monitor is designed to take and store measurements in its memory for two users.

1. Sit comfortably in a chair with your feet flat on the floor.
2. Select your User ID (A or B).
3. Stretch your arm forward on the desk and keep relaxing, make sure the palm of your hand is turned up. Make sure your arm is in the correct position, to avoid body movement. Sit still and do not talk or move during the measurement.

After the cuff has been appropriately positioned on the arm and connected to the blood pressure monitor, the measurement can begin:

- Press the START/STOP button. The pump will begin to inflate the cuff. In the display, the increasing cuff pressure is continually displayed.
- After automatically reaching an individual pressure, the pump stops and the pressure slowly falls. The cuff pressure is displayed during the measurement.
- When the device has detected your pulse, the heart symbol in the display begins to flash.
- When the measurement has concluded, the measured systolic and diastolic blood pressure values, as well as the pulse will be displayed.
- The results are displayed until you switch the device off. If no button is pressed for 60 seconds, the device will switch off automatically.
- Cuff self-checking symbol ()

The cuff's correct symbol () will be displayed if the cuff position is correct, otherwise the wrong symbol () will be displayed. Please check again the cuff to see if the wrong symbol () is still displayed.

- Movement error symbol ().

The Movement Error Symbol () is displayed if you move your body during the measurement. Please remove the cuff, and wait 2-3 minutes. Reapply the cuff and take another measurement.


NOTE:

The patient should be in the following position:

1. Comfortably seated
2. Legs uncrossed
3. Feet flat on the floor
4. Back and arm supported
5. Middle of the cuff at the level of the right atrium of the heart

HOW TO TAKE YOUR BLOOD PRESSURE

1. We recommend that the patient be as relaxed as possible and in silence whilst taking the blood pressure measurement.
2. We advise that the patient be seated for 5 minutes before the first reading is taken.
3. Any reading can be affected by the environment, the patient's position, exercise, or the patient's physiologic condition.
4. Performance of the automated sphygmomanometer can be affected by extreme temperatures, humidity and altitude.
5. To stop the inflation or measurement, push the START/STOP button. The monitor will stop inflating, start deflating, and will turn off.
6. After the monitor has detected your blood pressure and pulse rate, the cuff automatically deflates and your blood pressure and pulse rate will be displayed.
7. The monitor will automatically turn off after one minute.

IRREGULAR HEARTBEAT DETECTOR

This symbol "-" indicates that certain pulse irregularities were detected during the measurement.

In this case, the result may deviate from your normal basal blood pressure – repeat the measurement.

Please inform your doctor about frequent appearance of the Irregular Heartbeat Symbol.

This instrument is an oscillometric blood pressure monitor device that also analyses pulse frequency during measurement. This instrument has been clinically tested.

If pulse irregularities occur during the measurement, the irregular heartbeat symbol will be displayed afterwards. If the symbol appears more frequently (i.e. several times per week on measurements performed daily) or if it suddenly appears more often than usual, we recommend the patient to seek medical advice. The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

ERROR MESSAGES

SYMBOL	CAUSE	CORRECTION
No display appears	Weak battery or improper placement.	Replace both batteries with new ones. Check the battery installation for proper placement of the battery polarities.
Er1	Sensor abnormality.	Check if the pump is working or not. If it is working, then the problem is an abnormal sensor. Please send it to the local distributor.
Er2	Monitor could not detect pulse wave or cannot calculate the blood pressure data.	Check if the air is releasing too slow or not. If it is too slow, please check if there is any dust in the tube plug of the cuff and the cuff port in the device. If there is, please clean it and start the measurement again. If not, please contact your local distributor.
Er3	Measurement result is abnormal (SYS $\leq$ 45mmHg, DIA $\leq$ 24mmHg).	Occasionally – Measure one more time/ Always – Contact your distributor.
Er4	Too loose cuff or air leakage (Cannot inflate to 30mmHg within 15s).	Fasten the cuff correctly and make sure the air plug is properly inserted in the unit.
Er 5	The air tube is crimped.	Correct it and make the measurement again.
Er 6	The sensor is sensing great fluctuation in the pressure.	Please keep quiet and don't move.
Er 7	The pressure that the sensor sensing is over the limit.	Please send back to the local distributor.
Er 8	The demarcation is incorrect or the device has not been demarcated.	Please send back to the local distributor.

TROUBLESHOOTING

PROBLEM	CAUSE	SOLUTION
It won't power on.	The batteries are flat.	Replace the batteries.
	Battery polarity is reversed.	Insert the batteries correctly according to the polarity.
It doesn't inflate.	Tube is disconnected.	Connect it to the air outlet.
	Tube is broken or has a leak.	Change the tube/cuff.
Err message appears and it stops working.	The patient moved his arm as it was inflating.	Keep your body relaxed.
	The patient was talking as it was inflating.	Please, do not speak while taking a measurement.
The cuff has a leak.	The cuff is fastened too loosely.	Fasten the cuff a little bit tighter.
	The cuff is broken.	Change the cuff.

If the problem persists, please contact your distributor. Do not attempt to repair the product yourself.

MEMORY

At the end of each measurement, this monitor automatically stores each result with date and time. Each unit stores up to 120 sets of measurements for each user and a total of 240 sets.

SYMBOL DESCRIPTIONS

The following symbols may appear in this manual, on the Digital Blood Pressure Monitor B22, or on it's accessories. Some of the symbols represent standards and compliances associated with the Digital Blood Pressure Monitor B22 and it's use.

	Authorized Representative in the European Community
	CE Mark: conforms to essential requirements of the Medical Device Directive 93/42/EEC.
	Date of manufacture.
	Manufacturer
	Specifies serial number
	Type BF applied part
	Direct current
	DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.
	Follow instructions for use
	Put up
	Fragile
	Keep dry
	Avoid the sun
	handle gently
	temperature range
No Sterilize requirement	
Not category AP / APG equipment	
Mode of operation: continuous	

Viewing the stored values

With the unit off, press the Memory button. The display first shows "A", then shows an average of all measurements stored in the unit.

Note: Measurements for each user are averaged and stored separately. Be certain that you are viewing the measurements for the correct user. Pressing the Memory button again displays the previous value. To view a particular stored memory, press and hold the Memory button to scroll to that stored reading.

Discontinuing a Measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the Start/Stop button can be pressed at any time. The device will immediately lower the cuff pressure automatically.


'BATTERY CHANGE' INDICATOR

Batteries discharged- replacements required

When the batteries are discharged, the battery symbol will flash as soon as the instrument is switched on. You cannot take any further measurements and must replace the batteries.

The battery compartment is located on the back of the unit.

1. Remove the cover from the battery compartment, as illustrated below.
2. Insert the batteries (4 x size AAA). Always use AAA long life batteries or alkaline 1.5v batteries.
3. Even though time and date must be reset, the memory retains all values - the year number will automatically flash once the batteries are replaced.
4. To set date and time, follow the procedure described in the previous section.

Which batteries and how to change them?

Use four new, long-life 1.5V AAA batteries. Do not use batteries beyond their expiration date. If the monitor is not going to be used for a prolonged period of time the batteries should be removed.

Using rechargeable batteries

You can also operate this instrument using rechargeable batteries.

- Only use "NiMH" reusable batteries.
- If the battery symbol appears the batteries must be removed and recharged. They must not remain inside the instrument, as they may become damaged and leak. If you do not intend to use the instrument for a week or more, always remove the rechargeable batteries.
- Recharge these batteries using an external charger and carefully follow the manufacturer's instructions.

USING THE AC ADAPTER

You may also operate this monitor using the AC adapter (output 5V DC/1A with Micro USB plug). Only use the AC adapter provided by the manufacturer to avoid damaging the unit.

- Make sure that the AC adapter and cable are not damaged.
- Plug the adapter cable into the AC adapter port on the right side of the blood pressure monitor.
- Plug the adapter into your electrical outlet. When the AC adapter is connected, it doesn't consume any power from the batteries.

Note: If electrical power is interrupted (e.g., by accidentally unplugging the AC adapter from the outlet or the monitor's port) the monitor must be reset by removing the plug from the socket and reinserting the AC adapter connection.

CARE AND MAINTENANCE

Wash hands after each time you use the device.

If the device is going to be used by different patients, please, wash your hands before and after each use.

- Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- The cuff contains a sensitive air-tight bubble. Handle this cuff carefully and avoid all types of stress through twisting or buckling.
- Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and a mild detergent. The cuff with the bladder must not be washed in a dishwasher, a washing machine, or submerged in water.
- Handle the tube carefully. Do not pull from it. Do not allow the tubing to kink and keep it away from sharp edges.
- Do not drop the monitor or treat it roughly in any way. Avoid strong vibrations.
- Never open the monitor! This invalidates the manufacturer's warranty.
- Batteries and electronic instruments must be disposed of in accordance with local regulations, never with domestic waste.

Accuracy test

Sensitive measuring devices must be checked for accuracy from time to time. We recommend a periodical inspection of your unit by an authorized dealer at least once a year. Please contact your local distributor or the manufacturer.

CERTIFICATIONS

Device standard:

This device is manufactured to meet the European blood pressure monitors regulations:

EN1060-1 / 1995 • EN1060-3 / 1997 IEC 80601-2-30/2013• ISO81060-1/2012• IEC60601-1-11/2015• IEC60601-1/2012

Electromagnetic compatibility:

Device fulfils the stipulations of the International standard IEC60601-1-2.

EMC DECLARATION

- This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.
- Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- * Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
- * Caution: this machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used

TECHNICAL SPECIFICATIONS

Weight: 265g (batteries and AC adapter is not included)

Display: 3.74"x2.52" LCD Digital Display

Size: 118 (W) x 110 (L) x 57 (H) mm

Accessories: 1x Main Device, 1x Cuff, 1x User manual, 1x Warranty card

Operating Conditions: Temperature: 5 °C to 40°C; Humidity: 15% to 93% RH;

Storage And Shipping Conditions: Temperature: -25 °C to 70 °C; Humidity: ≤ 93% RH;

Measuring method: Oscillometric

Pressure sensor: Resistive

Measuring range: 0-280mmHg

Pulse: 40 to 199 per minute

Cuff pressure display range: <300mmHg

Memory: Automatically stores the last 120 measurements for 2 users (total 240)

Measuring resolution: 1 mmHg

Accuracy: Pressure within ± 3 mmHg / pulse ± 5 % of the reading

Power source:

- 4*AAA batteries, 1.5 V
- AC adapter INPUT: 100-240VAC 50/60HZ OUTPUT: 5V DC 1A

Accessories: Cuff 8.7" – 15.7" (22 - 40 cm)

Automatically power off: 60 seconds

Users: Adults

Expected service life of the device and accessories: 5 years

Technical alterations reserved.

Guidance and manufacture's declaration – electromagnetic emission

The device is intended for use in the electromagnetic environment specified below. The customer of the user of the device should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The device use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacture's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0 % UT ; 1 cycle 70 % UT; 25/30 cycle 0% UT; 250/300 cycle	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30 A/m 50/60Hz	30 A/m 50/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

**Recommended separation distances between
portable and mobile RF communications equipment and the device**

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter(W)	Separation distance according to frequency of transmitter(m)		
	150 KHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2.7 GHz $d = 2,3\sqrt{P}$
00.1	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (w) according to the transmitter manufacturer.

Note 1 at 80 mhz and 800 mhz, the separation distance for the higher frequency range applies.

Note 2 these guidelines may not apply in all situations. electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacture's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 3 V RMS outside the ISM band, 6 V RMS in the ISM and amateur bands 80% AM at 1kHz	Not applicable	Recommended separation distance $d=0.35\sqrt{p}$ $d=1.2\sqrt{p}$
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1kHz	10 V/m 80 MHz to 2.7 GHz 80% AM at 1kHz	80MHz to 800MHz: $d=1.2\sqrt{p}$ 800MHz to 2.7GHz: $d=2.3\sqrt{p}$ Where, P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

Guidance and manufacturer's declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device, should assure that it is used in such an environment.

Test frequency (MHz)	Band a) (MHz)	Service a)	Modulation a)	Maximum power (w)	Distance (m)	Immunity test level (v/m)
385	380-390	TETRA 400	Pulse Modulation b) 18 Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM c) ±5 kHz deviation 1 kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse Modulation b) 217 Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse Modulation b) 18 Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3 4,25; UMTS	Pulse Modulation b) 217 Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse Modulation b) 217 Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse Modulation b) 217 Hz	0.2	0.3	9
5500						
5785						

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

- For some services, only the uplink frequencies are included.
- The carrier shall be modulated using a 50% duty cycle square wave signal.
- As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E = \frac{6}{d} \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

CONTACT INFORMATION □

Digital Blood Pressure Monitor B22 is legal manufactured by:

Company name: Shenzhen Jamr Technology Co., Ltd.

Address: 2nd Floor, A-building, No.2 Guiyuan Road, Guihua community,
Guanlan town, Longhua new district, Shenzhen, P.R.China

Telephone: 0086 755 85292057

Authorized European Representative:

Company Name: Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg, Germany.

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In compliance with Directives: 2002/95/CE, 2002/96/CE and 2003/108/CE on the restriction of the use of dangerous substances in electric and electronic equipment as well as their waste disposal. The symbol with the crossed dustbin shown on the package indicates that the product at the end of its service life shall be collected as separate waste. Therefore, any products that have reached the end of their useful life must be given to waste disposal centres specialising in separate collection of waste electrical and electronic equipment, or given back to the retailer at the time of purchasing new similar equipment, on a one for one basis. The adequate separate collection for the subsequent start-up of the equipment sent to be recycled, treated and disposed of in an environmentally compatible way contributes to preventing possible negative effects on the environment and health and optimises the recycling and reuse of components making up the apparatus. Abusive disposal of the product by the user involves application of the administrative sanctions according to the laws.

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Name der Firma: Shenzhen Jamr Technology Co.,Ltd
Name der Firma: Shanghai International Holding Corp. GmbH (Europe)