2448-002

ENGLISH

DEUTSCH

ITALIANO

SUOMI

ESPAÑOL

FRANÇAIS

PORTUGUÊS

SVENSKA

TÜRKÇE

NORSK

DANSK

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NEDERLANDS

POLSKI

MiniFESS®

Surgical Instruments

Chirurgische Instrumente

Strumenti chirurgici

Kirurgiset instrumentit

Instrumentos quirúrgicos

Instruments chirurgicaux

Instrumentos cirúrgicos

Kirurgiska instrument

Cerrahi Aletler

Kirurgiske instrumenter

Kirurgiske instrumenter

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chirurgische instrumenten

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INSTRUCTIONS FOR USE.....	3
GEBRAUCHSANWEISUNG.....	4
ISTRUZIONI PER L'USO.....	6
KÄYTTÖOHJEET.....	8
INSTRUCCIONES DE USO.....	10
MODE D'EMPLOI.....	12
INSTRUÇÕES DE UTILIZAÇÃO.....	14
BRUKSANVISNING.....	16
KULLANIM TALİMATLARI.....	18
BRUKSANVISNING.....	20
BRUGSANVISNING.....	22
사용 지침.....	24
GEbruikSAANWIJZING.....	26
INSTRUKCJA UŻYCIA.....	28

stryker

Instructions for Cleaning, Sterilization, and Care

Important: Read all instructions in this manual before using and reprocessing these devices.

Cleaning

Do not allow contaminants to dry on instruments as this makes cleaning more difficult. Immediately after use, place instrument under cold, running water to remove contaminants then pre-soak in enzymatic cleaning solution. Follow enzymatic cleaner instructions for concentration, water temperature, and timeframe. Thoroughly rinse instrument to remove harmful residue from soaking solutions.

Caution: Instruments must not be soaked in caustic or physiological saline solutions as this could cause pitting or rust.

Instruments with joints or lock boxes (scissors, clamps, retractors, etc.) should always be cleaned in the open position. They should be activated several times during the cleaning process to facilitate removal of contaminants in joints.

Use a soft brush to help remove remaining contaminants.

Caution: Do not use metal brushes as this could scratch or deposit metal particles on the instrument which could lead to rusting.

A pressurized water spray can also be used to remove contaminants. Rinse thoroughly under cold, running water followed by a distilled water rinse to remove tap water residue.

Visually inspect instrument to insure it is clean.

Place instrument in an ultrasonic cleaner and follow the ultrasonic cleaner instructions for detergent, temperature, and time. Rinse thoroughly under cold, running water to remove any residue of the detergent once again followed by a distilled water rinse.

Sterilization: Autoclave

Pre-vacuum

Unwrapped

- 132°C (270°F) 3 minutes
- 135°C (275°F) 3 minutes

Wrapped

- 132°C (270°F) 4 minutes
- 135°C (275°F) 3 minutes
- 137°C (279°F) 3 minutes

Gravity

Unwrapped

- 132°C (270°F) 3 minutes
- 135°C (275°F) 3 minutes

Wrapped

- 132°C (270°F) 15 minutes
- 135°C (275°F) 10 minutes

Caring for Stainless Steel Surgical Instruments

For stainless steel to be used for surgical instruments, it must contain a specific amount of carbon. Carbon is required in order to allow the steel to become “hardened.” Hardening strengthens the steel and allows it to retain sharpness. Carbon also allows steel to be more susceptible to tarnishing and rusting. The effect of the carbon is minimized by processing the surface finish of the instrument. As the name implies, the material stainless steel—STAINS LESS, but is not stain-free.

Improper cleaning, disinfection, and sterilizing will contribute more to rusting or staining than the carbon content of the stainless steel.

Common Problems

- 1. Improper drying of the instrument after cleaning, disinfection, or sterilizing especially in joints and lock boxes.
- 2. Using corrosive or caustic cleaning agents. (Do not use any chemical over 10.5 pH)
- 3. Improper rinsing to remove the cleaning or disinfecting solutions;

Using tap water without following with a distilled water rinse (tap water may contain chemicals and minerals such as iron which can leave deposits on the surface).

- 4. Faulty autoclave which may leave deposits on the instrument and attack the instrument’s surface finish.
- 5. Use of stiff metal brushes that roughen the surface of the instrument and leaves it susceptible to rust and stains.

Graphic Symbols Contained on Device Labeling

LOT Lot Number	REF Catalog Number	 Manufacturer
 CE Mark	Rx Only Prescription Use Only	 Device Is Non Sterile
EC REP Authorized Representative in the European Community		 Consult Instructions for use

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MiniFESS jest znakiem towarowym firmy Stryker Corporation z siedzibą w Kalamazoo, MI.
Entellus Medical Inc.
Entellus Medical, Inc., a Stryker Company
ICON (LR) Limited