

**REF** **A0022584, A0025025****Starrsed PC****Supplier:**  
Mechatronics Instruments B.V.  
De Corantijn 13  
1689 AN Zwaag  
The Netherlands  
T +31 229 291 129**Intended purpose**

Reagent for the automatic cleaning of the pipettes in Starrsed ESR analyzers.

**Instructions For Use**

The product is for use on the Starrsed analyzers only. See the instrument manual for detailed instructions.

1. Remove cap from the new bottle.
2. Remove cap with pickup tube from empty bottle.
3. Remove empty bottle.
4. Place the new bottle.
5. Place cap with pickup tube on the new bottle. Avoid touching the pickup tube which goes into the container.

**Chemical/Physical data**

Contains main surfactants (0.5 % by weight in total) as active ingredient, some excipients with high pH to boost the cleaning performance (<0.25 % by weight in total), a coloring agent for distinction and a preservative (< 0.5 ‰ by weight).

**Packaging**

A0022584 bottle 500 ml

A0025025 box with 6x500 ml bottles A0022584

**Safety**

Any serious incident that has occurred in relation to this product must be reported to the distributor and/or manufacturer and the competent local authorities.

The latest version of the Safety Data Sheet is available on the website [www.rrmechatronics.com](http://www.rrmechatronics.com)

**Storage and stability**

Transport stability at 5 - 60 °C (maximum of 1 month)

Shelf life after production: 2 years, but no longer than the expiry date.

Store at 5 - 35 °C in dry conditions in well sealed containers.

In-use stability: Maximum of 2 months at advised temperature of 5 - 30°C.

**Ecological information**

Generally not hazardous for water.

**Waste treatment and cleaning**

Avoid discharge to lakes, streams, sewers, etc. Disposal must be made according to official regulations.

Wherever possible cleaning should be performed with normal cleaning agents. Avoid use of solvents.



RR Mechatronics Manufacturing B.V.  
De Corantijn 13, 1689 AN Zwaag,  
The Netherlands  
[www.rrmechatronics.com](http://www.rrmechatronics.com)



Qserve Group UK Ltd  
282 Farnborough Road,  
Farnborough, GU14 7NA  
United Kingdom



QUNIQUE GmbH  
Bahnhofweg 17  
5610 Wohlen  
Switzerland

rev. A01  
Initial release