

REF A0022584, A0025025**Starrsed PC****Supplier:**
Mechatronics Instruments B.V.
De Corantijn 13
1689 AN Zwaag
The Netherlands
T +31 229 291 129**Intended purpose**

Reagent for the automatic cleaning of the pipettes in Starrsed ESR analyzers.

Instructions For Use

The product is for use on the Starrsed analyzers only. See the instrument manual for detailed instructions.

1. Remove cap from the new bottle.
2. Remove cap with pickup tube from empty bottle.
3. Remove empty bottle.
4. Place the new bottle.
5. Place cap with pickup tube on the new bottle. Avoid touching the pickup tube which goes into the container.

Chemical/Physical data

Contains main surfactants (0.5 % by weight in total) as active ingredient, some excipients with high pH to boost the cleaning performance (<0.25 % by weight in total), a coloring agent for distinction and a preservative (< 0.5 ‰ by weight).

Packaging

A0022584 bottle 500 ml

A0025025 box with 6x500 ml bottles A0022584

Safety

Any serious incident that has occurred in relation to this product must be reported to the distributor and/or manufacturer and the competent local authorities.

The latest version of the Safety Data Sheet is available on the website www.rrmechatronics.com

Storage and stability

Transport stability at 5 - 60 °C (maximum of 1 month)

Shelf life after production: 2 years, but no longer than the expiry date.

Store at 5 - 35 °C in dry conditions in well sealed containers.

In-use stability: Maximum of 2 months at advised temperature of 5 - 30°C.

Ecological information

Generally not hazardous for water.

Waste treatment and cleaning

Avoid discharge to lakes, streams, sewers, etc. Disposal must be made according to official regulations.


Wherever possible cleaning should be performed with normal cleaning agents. Avoid use of solvents.



RR Mechatronics Manufacturing B.V.
De Corantijn 13, 1689 AN Zwaag,
The Netherlands
www.rrmechatronics.com



Qserve Group UK Ltd
282 Farnborough Road,
Farnborough, GU14 7NA
United Kingdom



QUNIQUE GmbH
Bahnhofweg 17
5610 Wohlen
Switzerland

rev. A01
Initial release