



# OWNER'S MANUAL

Arm-type Fully Automatic  
Digital Blood Pressure Monitor



ARM TYPE

**Model DBP-1303 | art. ARM250P**

Document No.: JDBP-1904-002

Version: Z

Date of Issue: 2016.01



The product is in compliance with the requirements of MDD 93/42/EEC, is the identification number of notify body



JOYTECH Healthcare Co., Ltd.

No.365, Wuzhou Road, Yuhang Economic Development Zone, Hangzhou City, 311100 Zhejiang, China



Shanghai International Holding Corp. GmbH (Europe)  
Eiffestrasse 80, 20537 Hamburg, Germany



# CONTENTS

Safety Notice . . . . .	02
Unit Illustration . . . . .	06
Important Testing Guidelines . . . . .	09
Quick Start . . . . .	10
Unit Operation . . . . .	11
Battery Installation . . . . .	11
System Settings . . . . .	12
Applying the Arm Cuff . . . . .	13
Testing . . . . .	14
Power Off . . . . .	16
Memory Check . . . . .	17
Memory Deletion . . . . .	17
Low Battery Indicator . . . . .	18
Troubleshooting . . . . .	19
Blood Pressure Information . . . . .	20
Blood Pressure Q&A . . . . .	22
Maintenance . . . . .	23
Specifications . . . . .	28
Disposal . . . . .	27
Warranty . . . . .	28

# SAFETY NOTICE

Thank you for purchasing the DBP-1303 Blood Pressure Monitor. The unit has been constructed using reliable circuitry and durable materials. Used properly, this unit will provide years of satisfactory use.

This device is intended for non-invasive measuring an adult individual's systolic, diastolic blood pressure and heart rate using the oscillometric method. The device is not intended for use on infants and children. The device is designed for home or clinical use. All functions can be used safely and values can be read out in one LCD DISPLAY. Measurement position is on adult upper arm only.

Please read this manual thoroughly before using the unit. Please retain this manual for future reference.

For specific information about your blood pressure, please  
CONSULT YOUR DOCTOR. The PATIENT is an intended OPERATOR.

To avoid risk and damage follow all warning precautions. Operate unit only as intended. Read all instructions prior to use.

## WARNING SIGNS AND SYMBOLS USED



Caution



Mandatory



Prohibited



Type BF Equipment



Instructions For Use MUST be Consulted



Serial Number



Discard the used product to the recycling collection point according to local regulations



The product conforms to the requirements of the EC Directive MDD (93/42/EEC) on medical devices



Manufacturer



Authorised Representative in the European Community



Keep Dry



Keep off Sunlight



Manufacturing Date

### CAUTION

Individuals with serious circulation problems may experience discomfort. Consult your physician prior to use.

Contact your physician if test results regularly indicate abnormal readings. Do not attempt to self-treat these symptoms without consulting your physician first.



Product is designed for its intended use only. Do not misuse in any way.

Product is not intended for infants or individuals who cannot express their intentions.

Do not disassemble or attempt to repair.

Do not use cell phones and other devices, which generate strong electrical or electromagnetic fields, near the device, as they may cause incorrect readings and interference or become interference source to the device.



Only use a recommended AC adaptor double-insulated complying with EN 60601-1 and EN 60601-1-2 (see page 6). An unauthorized adaptor may cause fire and electric shock.

### BATTERY PRECAUTIONS

Do not mix new and old batteries simultaneously.

Replace batteries when Low Battery Indicator  appears on screen.

Be sure battery polarity is correct.



Do not mix battery types. Long-life alkaline batteries are recommended.

Remove batteries from device when not in operation for more than 3 months.

Dispose batteries properly; observe local laws and regulations.

## IMPORTANT INSTRUCTIONS BEFORE USE


- 01.** Do not confuse self-monitoring with self-diagnosis. Blood pressure measurements should only be interpreted by a health professional who is familiar with your medical history.
- 02.** Contact your physician if test results regularly indicate abnormal readings.
- 03.** If you are taking medication, consult with your physician to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting with your physician.
- 04.** Individuals with serious circulation problems may experience discomfort. Consult your physician prior to use.
- 05.** For persons with irregular or unstable circulation resulting from diabetes, liver disease, arteriosclerosis or other medical conditions, there may be variations in blood pressure values measured at the wrist versus at the upper arm. Monitoring the trends in your blood pressure taken at either the arm or the wrist is nevertheless useful and important.
- 06.** People suffering from vascular constriction, liver disorders or diabetes, people with cardiac pacemakers or a weak pulse, and women who are pregnant should consult their physician before measuring their blood pressure themselves. Different values may be obtained due to their condition.
- 07.** People suffering from arrhythmias such as atrial or ventricular premature beats or atrial fibrillation only use this blood pressure monitor in consultation with your doctor. In certain cases oscillometric measurement method can produce incorrect readings.
- 08.** Too frequent measurements can cause injury to the patient due to blood flow interference.
- 09.** The cuff should not be applied over a wound as this can cause further injury.
- 10.** DO NOT attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.
- 11.** The cuff should not be placed on the arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm.
- 12.** Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring equipment on the same limb.
- 13.** A compressed or kinked connection hose may cause continuous cuff pressure resulting in blood flow interference and potentially harmful injury to the patient.
- 14.** Check that operation of the unit does not result in prolonged impairment of the circulation of the patient.

15. Product is designed for its intended use only. Do not misuse in any way.
16. Product is not intended for infants or individuals who cannot express their intentions.
17. Prolonged over-inflation of the bladder may cause ecchymoma of your arm.
18. Do not disassemble the unit or arm cuff. Do not attempt to repair.
19. Use only the approved arm cuff for this unit. Use of other arm cuffs may result in incorrect measurement results.
20. The system might produce incorrect readings if stored or used outside the manufacturer's specified temperature and humidity ranges. Make sure to store the blood pressure monitor, children, pets and pests are outside of accessible range.
21. Do not use the device near strong electrical or electromagnetic fields generated by cell phones or other devices, they may cause incorrect readings and interference or become interference source to the device.
22. Do not mix new and old batteries simultaneously
23. Replace batteries when Low Battery Indicator  appears on screen. Replace both batteries at the same time.
24. Do not mix battery types. Long-life alkaline batteries are recommended.
25. Remove batteries from device when not in operation for more than 3 months.
26. Do not insert the batteries with their polarities incorrectly aligned.
27. Dispose batteries properly; observe local laws and regulations.
28. Only use a recommended AC adaptor double-insulated complying with EN 60601-1 and EN 60601-1-2. An unauthorized adapter may cause fire and electric shock.
29.  Advising operator that Instruction manual/ Booklet must be consulted .
30. The PC with connection to the device with USB shall meet the requirements of standard IEC 60601-1 or IEC 60950-1.
31. Do not use the device during transport vehicles for influencing measurement accuracy, such as patient transport in an ambulance or helicopter.
32. Contains small parts that may cause a choking hazard if swallowed by infants.
33. Please align the polarities of each battery with the +ve and -ve signs imprinted on the battery housing when you replace the batteries .
34. The time required for the device to warm from the minimum storage temperature (-25°C) between use until the device is ready for use at Ambient Temperature 20°C: about 2 hours.
35. .The time required for the device to cool from the maximum storage temperature (70°C) between use until the device is ready for use at ambient temperature (20 °C) : about 2 hours

# UNIT ILLUSTRATION




## CONTENTS




MONITOR UNIT



OWNER'S MANUAL



STORAGE CASE



ARM CUFF




2MOP Medical AC Adapter  
(DC6.0 V, 600mA)  
(recommended, not provided)






MONITOR  
UNIT






ARM  
CUFF



DISPLAY



# IMPORTANT TESTING GUIDELINES

- 01.** Avoid eating, exercising, and bathing for 30 minutes prior to testing.
- 02.** Sit in a calm environment for at least 5 minutes prior to testing.
- 03.** Do not use cell phones and other devices, which generate strong electrical or electromagnetic fields, near the device, as they may cause incorrect readings.
- 04.** Avoid speaking or moving body parts while testing.
- 05.** While testing, avoid strong electromagnetic interference such as microwave ovens and cell phones.
- 06.** Wait 3 minutes or longer before re-testing.
- 07.** Try to measure your blood pressure at the same time each day for consistency.
- 08.** Test comparisons should only be made when monitor is used on the same arm, in the same position, and at the same time of day.
- 09.** This blood pressure monitor is not recommended for people with severe arrhythmia.
- 10.** Do not use this blood pressure monitor if the device is damaged.

# QUICK START




Figure B

01. Install batteries. (See Figure A)
02. Insert cuff air plug into the left side of monitor unit. (See Figure B)
03. Remove thick clothing from the arm area.



Figure C

04. Rest for several minutes prior to testing. Sit down in a quiet place, preferably at a desk or table, with your arm resting on a firm surface, keep your legs uncrossed and your feet flat on the floor. (See Figure C)




Figure D

05. Apply cuff to your left arm and keep level with your heart. Bottom of cuff should be placed approximately 1-2cm (1/2") above elbow joint. (See Figures D&E)





Figure E

06. Press " ON/OFF " Button to start testing.

# UNIT OPERATION

## BATTERY INSTALLATION


- 01.** Slide battery cover off as indicated by arrow.
- 02.** Install 4 new AA alkaline batteries according to polarity.
- 03.** Close battery cover.



AC adapter jack is on the right side of the monitor. Medical AC adapter(DC 6.0 V, 600mA) can be used with the device (recommended, not provided).

The adapter connect pin should be positive inside and negative outside with a 2.1mm coaxial joint.

Do not use any other type of AC adapter as it may harm the unit.



AC adapter jack




Note: Power supply is specified as part of ME EQUIPMENT

## SYSTEM SETTINGS

With power off, press "SET" button to activate System Settings. The Memory Group icon flashes.

### SELECT MEMORY GROUP

- While in the System Setting mode, you may accumulate test results into 4 different groups. This allows multiple users to save individual test results (up to 30 memories per group.) Press "+" or "-" button to choose a group setting. Test results will automatically store in each selected group.



### TIME /DATE SETTING

- Press "SET" button again to set the Time/Date mode. Set the month first by adjusting the "+" or "-" button. Press "SET" button again to confirm current month. Continue setting the day, hour, and minute in the same way. Every time the "SET" button is pressed, it will lock in your selection and continue in succession (month, day, hour, minute.)



### SAVED SETTINGS


- While in any setting mode, press "ON/OFF" button to turn the unit off. All information will be saved.




If unit is left on and not in use for 3 minutes, it will automatically save all information and shut off.

## APPLYING THE ARM CUFF


- 01.** Firmly insert air plug into opening located on left side of monitor unit.



- 02.** With sticky nylon section facing outward, insert end of cuff underneath metal ring of cuff.



- 03.** Fasten cuff about 1-2cm (1/2") above the elbow joint. For best results apply cuff to bare arm and keep level with heart while testing.




Do not insert air plug into opening located on right side of monitor unit.


This opening is designed for an optional power supply only.

**01. TESTING**


Power On  
Press and hold ON/OFF button until a beep sounds. The LCD screen will appear for one second as unit performs a quick diagnosis. A long tone indicates device is ready for testing.



Unit will not function if residual air from previous testing is present in cuff.

The LCD will flash  until pressure is stabilized.

- 02. Pressurization**  
The unit will automatically inflate to the proper pressure value and stop inflating. During this time, please keep quiet.



Pressurization will gradually subside and ultimately stop when cuff is not properly applied to the arm. If this occurs, press "ON/OFF" button to turn the unit off.



**03. Testing**

After cuff inflation, air will slowly subside as indicated by the corresponding cuff pressure value. A flashing ♥ will appear simultaneously on screen signaling heart beat detection.



Keep relaxed during testing. Avoid speaking or moving body parts.

**04. Result Display**

Three short beeps sound when testing is complete. The screen will display measurements for systolic and diastolic blood pressure.



Refer to Page 22 for detail Blood Pressure Information.

**05. Deleting/Storing Test Results**

User may delete their current test result due to unfavorable testing conditions or for any other reason. To delete the last test result, press the "SET" button after result is displayed. If result is not deleted, it will automatically store by date within the previously configured Memory Group.



Be sure the appropriate Memory Group selection is made prior to testing.

If the number of tests surpasses the allotted 30 memories per group, the most recent tests will appear first, thus eliminating older readings


**POWER OFF**

Press "ON/OFF" button in any mode to turn off device. Unit automatically turns off after 3 minutes of inactivity.

**Safety Precaution:** If pressure in arm cuff becomes too extreme while testing, press the "ON/OFF" button to turn power off. The cuff pressure will rapidly dissipate once the unit is off.

## MEMORY CHECK


With power off, you may check past test results by using the “+” or “-” buttons. The most recent test result can be viewed by pressing and holding the “+” button. The oldest test result in memory can be viewed by pressing and holding the “-” button. Upon activating test results, you can press the “+” or “-” buttons to scroll through all test results stored in memory.



Previous test results will only be displayed from the most recently used memory group. To check previous test results in other memory groups, you must first select the desired group and then turn monitor off. (See “Select Memory Group” on Page 10.)

## MEMORY DELETION


Memory for a selected group may be deleted while in Memory Check mode. Press and hold the “SET” button for approximately 3 seconds to delete all memory records from the selected group. The monitor will beep indicating successful deletion and then transfer into testing mode. Press the “ON/OFF” button to turn the unit off.




Memory cannot be recovered once it has been deleted.

## LOW BATTERY INDICATOR

4 short warning beeps sound when battery life is depleting and unable to inflate cuff for testing.


The  appears simultaneously for approximately 5 seconds prior to shutting off. Replace batteries at this time. No memory loss will occur throughout this process.



## STATIC PRESSURE MEASUREMENT

In the power down state, press and hold the "ON/OFF" button, and then install the batteries.

Until the LCD screen is full, release the "ON/OFF" button. When the LCD screen displays the double zero, the blood pressure meter is in static state. Software version is displayed at the heart rate.



Only Service personnel permitted to access to this mode, the mode unavailable in normal use.

## TROUBLESHOOTING

Problem	Possible Cause	Solution
Blood pressure results are not within typical range	Cuff is too tight or not properly positioned on the arm	Firmly reposition cuff approximately 1-2cm above the elbow joint ( See Page 12)
	Inaccurate test results due to body movement or monitor movement	Sit in a relaxed position with arm placed near heart. Avoid speaking or moving body parts while testing. Make sure the monitor unit is placed in a stationary position throughout the testing period. ( See Page 7)
<b>Err</b> displayed	Cuff fails to inflate properly	Make sure hose is properly fastened to cuff and monitor unit
	Improper operation	Read user manual carefully and re-test properly.
	Pressurization is over cuff rated pressure 300mmHg	Read user manual carefully and re-test properly.


# BLOOD PRESSURE INFORMATION

## BLOOD PRESSURE

Blood pressure is the force of blood pushing against the walls of arteries. It is typically measured in millimeters of mercury (mmHg.) Systolic blood pressure is the maximum force exerted against blood vessel walls each time the heart beats. Diastolic blood pressure is the force exerted on blood vessels when the heart is resting between beats.

An individual's blood pressure frequently changes throughout the course of a day. Excitement and tension can cause blood pressure to rise, while drinking alcohol and bathing can lower blood pressure. Certain hormones like adrenaline (which your body releases under stress) can cause blood vessels to constrict, leading to a rise in blood pressure.

If these measuring numbers become too high, it means the heart is working harder than it should.




Example: fluctuation within a day (male, 35 years old)

## HEALTH REMINDER

Hypertension is a dangerous disease that can affect the quality of life. It can lead to a lot of problems including heart failure, kidney failure, and cerebral hemorrhaging.

By maintaining a healthy lifestyle and visiting your physician on a regular basis, hypertension and relative diseases are much easier to control when diagnosed in their early stages.



Do not be alarmed if an abnormal reading occurs. A better indication of an individual's blood pressure occurs after 2-3 readings are taken at the same time each day over an extended period of time. Consult your physician if test results remain abnormal.

# BLOOD PRESSURE

## Q&A



What is the difference between measuring blood pressure at home or at a professional healthcare clinic?



Blood pressure readings taken at home are now seen to give a more accurate account as they better reflect your daily life. Readings can be elevated when taken in a clinical or medical environment. This is known as White Coat Hypertension and may be caused by feeling anxious or nervous.



Abnormal test results may be caused by:

1. Improper cuff placement Make sure cuff is snug-not too tight or too loose. Make sure bottom of the cuff is approximately 1-2cm (1/2") above the elbow joint.
2. Improper body position Make sure to keep your body in an upright position.
3. Feeling anxious or nervous Take 2-3 deep breaths, wait a few minutes and resume testing.



What causes different readings?



Blood pressure varies throughout the course of a day. Many factors including diet, stress, cuff placement, etc. may affect an individual's blood pressure.



Should I apply the cuff to the left or right arm? What is the difference?



Either arm can be used when testing, however, when comparing results, the same arm should be used. Testing on your left arm may provide more accurate results as it is located closer to your heart.




What is the best time of day for testing?




Morning time or any time you feel relaxed and stress free.




# MAINTENANCE




- 01.** Avoid dropping, slamming, or throwing the unit.



- 02.** Avoid extreme temperatures. Do not expose unit directly under sunshine.




- 03.** When cleaning the unit, use a soft fabric and lightly wipe with mild detergent. Use a damp cloth to remove dirt and excess detergent.



- 04.** Cuff Cleaning: Do not soak cuff in water!

Apply a small amount of rubbing alcohol to a soft cloth to clean cuff's surface. Use a damp cloth (water-based) to wipe clean. Allow cuff to dry naturally at room temperature. The cuff must be cleaned and disinfected before use between different users.




- 05.** Do not use petrol, thinners or similar solvents.

- 06.** Remove batteries when not in operation for an extended period of time.

- 07.** Do not disassemble product.

- 08.** It is recommended the performance should be checked every 2 years.



- 09.** Expected service life: Approximately three years at 10 tests per day.

- 10.** No service and maintenance while it is in use and maintenance only be performed by service personnel. Service and maintenance require parts, repair, technical support will be provided.

# SPECIFICATIONS

Product Description	Arm-type Fully Automatic Digital Blood Pressure Monitor	
Model	DBP-1303	
Display	LCD Digital Display Size: 62.7mm×46.4mm (2.47" x1.83")	
Measurement Method	Oscillometric Method	
Measurement Range	Systolic Pressure	60mmHg-260 mmHg
	Diastolic Pressure	30mmHg-200 mmHg
	Pressure	0mmHg-300 mmHg
	Pressure	±3mmHg
	Pulse	30 ~ 180 Beats/Minute ±5%
	Pulse	±5%
Pressurization	Automatic Pressurization	
Memory	120 Memories in Four Groups with Date and Time	
Function	Low Battery Detection	
	Automatic Power-Off	
Power Source	4 AA batteries or Medical AC Adapter(DC6.0V, 600mA) (recommended, not provided)	
Battery Life	Approximately 2 months at 3 tests per day	
Unit Weight	Approx.395g (13.95 oz.) (excluding battery)	
Unit Dimensions	Approx.162 x 110 x 62.9mm (6.36" x 4.33" x 2.48" )(L x W x H)	
Cuff Circumference	Medium cuff: Fits arm circumference 22-36 cm	
Operating Environment	Temperature	10°C ~ 40°C (50°F-104°F)
	Humidity	15% ~93%RH
	Pressure	700hPa-1060hPa
Storage Environment	Temperature:	-25°C-70°C (-13°F-158°F)
	Humidity	≤93% RH
Classification:	Internal Powered Equipment,Type BF  ,Cuff is the Applied Part	
Ingress Protection Rating:	IP20, Indoor Use Only	

Specifications are subject to change without notice.

This Blood Pressure Monitor complies with the European regulations and bears the CE mark "CE 0197". This blood pressure monitor also complies with mainly following standards (included but not limited)

#### **SAFETY STANDARD:**

EN 60601-1 Medical electrical equipment part 1: General requirements for safety EMC standard:

EN 60601-1-2 Medical Electrical Equipment -- Part 1-2: General Requirements For Basic Safety And Essential Performance -- Collateral Standard: Electromagnetic Disturbances – Requirements And Tests.

#### **PERFORMANCE STANDARDS:**

IEC80601-2-30, Medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers.

EN 1060-3 Non-invasive sphygmomanometers – Supplementary requirements for electromechanical blood pressure measuring systems.

ISO 81060-2, non-invasive sphygmomanometers – part 2: clinical validation of automated measurement type.

# DISPOSAL



## USER INFORMATION

"Implementation of Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE)", pertaining to reduced use of hazardous substances in electrical and electronic equipment, as well as to waste disposal.

The symbol of the crossed-out wheeled bin on the equipment or on its packaging indicates that the product must be disposed of separately from other waste at the end of its service life.

The user must therefore take the dismissed equipment to suitable separate collection centres for electrical and electronic waste.

For more details, please contact the appropriate local authority.

Suitable segregated collection of the equipment for subsequent recycling, treatment or environmentally-friendly disposal helps prevent damage to the environment and to human health, and encourages the re-use and/or recycling of the materials that make up the equipment.

Abusive disposal of the product by the user shall result in the application of administrative fines in accordance with the laws in force.



## BATTERIES

Remove the batteries from the appliance before its disposal. Do not dispose the batteries with domestic waste. The batteries must be disposed of in the appropriate containers or in the specific collecting centre. Suitable segregated collection helps prevent damage to the environment and to human health.

To remove batteries, refer to Maintenance section.

## PACKAGING

In order to respect the environment, packaging material must be disposed of properly in accordance with separate collection.

Check local regulations.

# WARRANTY

## Conditions

The warranty is valid for 24 months from the date of purchase.

This warranty is valid only if it is correctly filled in and accompanied by the tax receipt proving the date of purchase.

The appliance must be delivered exclusively to our authorised Service Centre.

By warranty we mean the replacement or repair of appliance components that are faulty from the start due to manufacturing defects.

However, assistance (for a fee) is also guaranteed for products out of warranty.

The consumer is the owner of the rights applicable by the national legislation governing the sale of consumer goods; this warranty is without prejudice to those rights.

The manufacturer declines all responsibility for any damage to people, animals or things resulting from improper use of the appliance and failure to comply with the instructions indicated in the appropriate instruction booklet.

## Limitations

All warranty rights and all our responsibilities expire if the device has been:

- Tampered with by unauthorised personnel.
- Used, stored or transported in an inappropriate way.

In any case, loss of performance related to the aesthetics of the product or similar that does not compromise the essence of the functions is excluded from the warranty.

If, despite the careful selection of materials and our commitment to making the product you have just purchased, you find any defects, or if you need information, we recommend that you call your local dealer.

---

## CONTACT INFORMATION

Imported by:

Poly Pool S.p.A. Via Sottocorna, 21/B

24020 Parre (Bergamo) - Italy

Tel. +39 035 4104000 r.a. - Fax +39 035 702716

<http://www.ardes.it> - e-mail: [polypool@polypool.it](mailto:polypool@polypool.it)

MADE IN CHINA

# ELECTROMAGNETIC COMPATIBILITY INFORMATION

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

**Table 1**

Guidance and declaration of manufacturer-electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment guidance</b>
Radiated emission CISPR 11	Group 1, class B.	The device uses RF energy only for its internal function. Therefore, its emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted emission CISPR 11	Group 1, class B.	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	


# ELECTROMAGNETIC COMPATIBILITY INFORMATION

**Table 2**

<p>Guidance and declaration of manufacturer-electromagnetic immunity          The device is intended for use in the electromagnetic environment specified below.          The customer or the user of the device should assure that it is used in such an environment.</p>			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air	$\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient/burst IEC 61000-4-4	$\pm 2$ kV, 100kHz, for AC power port	$\pm 2$ kV, 100kHz, for AC power port	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 0.5$ kV, $\pm 1$ kV (differential mode)	$\pm 0.5$ kV, $\pm 1$ kV (differential mode)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°  0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°  0 % UT; 250/300 cycle	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°  0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0  % UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m; 50Hz or 60Hz	30 A/m; 50Hz or 60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.



**Table 3**

<p>Guidance and declaration of manufacturer-electromagnetic immunity  The device is intended for use in the electromagnetic environment specified below.  The customer or the user of the device should assure that it is used in such an environment.</p>			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Conducted RF IEC 61000-4-6	3V for 0.15-80MHz; 6V in ISM and amate -ur radio bands between 0.15-80MHz	3V for 0.15-80MHz; 6V in ISM and amate -ur radio bands between 0.15-80MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	385MHz, 27V/m	385MHz, 27V/m	Recommended separation distance
	450MHz, 28V/m	450MHz, 28V/m	$d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$
	710MHz, 745 MHz, 780MHz 9V/m	710MHz, 745 MHz, 780MHz 9V/m	80 MHz to 800 MHz
	810MHz, 870 MHz, 930MHz 28V/m	810MHz, 870 MHz, 930MHz 28V/m	$d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
	1720MHz, 1845 MHz, 1970MHz 28V/m	1720MHz, 1845 MHz, 1970MHz 28V/m	800 MHz to 2.7 GHz
	2450MHz, 28V/m	2450MHz, 28V/m	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).
	5240MHz, 5500 MHz, 5785MHz 9V/m	5240MHz, 5500 MHz, 5785MHz 9V/m	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:
			

**Table 4**

Recommended separation distances between portable and mobile RF communications equipment and the device.

The device is intended for use in an electromagnetic environment in which radiated therefore disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter	
	80 MHz to 800 MHz	800 MHz to 2.7 GHz
	$d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	$d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.