Healthcare-Manager.com

Questions or comments?

Please call toll-free:

1-855-822-6999 M-F 9 a.m. - 5 p.m. CST

E-mail: service@healthcare-manager.com

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Fabriqué en Chine

easy@Home®

Automatic Digital Blood Pressure Monitor

Model No.: EBP-073

V01:20230712

Contents

Safety Notice	01
Device Illustration	06
Important Testing Guidelines	08
Quick Start	09
Device Operation	10
Battery Installation	10
System Settings	11
Applying the Arm Cuff	13
Taking Reading	14
Power Off	16
Last 3 Test Average	16
Memory Check and Memory Deletion	17
Low Battery Indicator	18
Troubleshooting	20
Blood Pressure Information	21
Blood Pressure Q&A	24
Maintenance	25
Specifications	27
Warranty	29
FCC Statement.....	30
Electromagnetic Compatibility Information	30
Additional Notes	35

Safety Notice

Thank you for purchasing the Easy@Home EBP-073 Blood Pressure Monitor. The device is constructed using reliable circuitry and durable materials. Used properly, this device will provide years of satisfactory use.

Indications for use: The Automatic Digital Blood Pressure Monitor is intended to measure blood pressure (systolic and diastolic) and pulse rate of adults and adolescents over 12 years of age.

All functions can be used safely and values can be read out in one LCD DISPLAY. Measurement position is on adult upper arm only. The PATIENT is an intended OPERATOR.

Intended Use: The digital blood pressure monitor are reusable for clinical and home use and are non-invasive blood pressure measurement systems designed to measure the systolic and diastolic blood pressure and pulse rate of adolescents and adults individual by using a non-invasive technique, which is a well-known technique in the market called the "oscillometric method". It can measure the systolic blood pressure, diastolic blood pressure and pulse rate on up-arm, and the device is reusable for clinical or home use.

Readings taken by this device are equivalent to those obtained by a trained observer using the cuff/ stethoscope auscultation method, within the limits prescribed by the Recognized Consensus Standard (IEC 80601-2-30) for electronic sphygmomanometers.

Precautions to Ensure Safe, Reliable Operation

1. Do not drop the device. Protect it from sudden jars or shocks.
2. Do not insert foreign objects into any openings.
3. Do not attempt to disassemble the device.
4. Do not crush the pressure cuff.
5. If the device has been stored in an environment below 0 °C, leave it in a warm place for about 15 minutes before use it. Otherwise, the cuff may not inflate properly.
6. If the device has been stored in an environment above 40 °C, leave it in a cool place for about 15 minutes before use it. Otherwise, the cuff may not inflate properly.
7. Do not store the device in direct sunlight, high humidity or in a dusty area.
8. To avoid any possibility of accidental suffocation, keep this device away from children and do not wrap the air tube around anyone's neck.
9. Ensure that children do not use the instrument unsupervised; some parts are small enough to be swallowed.
10. Someone may experience skin irritation from the cuff after taking frequent readings over the course of the day, but this irritation typically goes away on its own after the monitor cuff is removed.

Important Instructions Before Use

1. Do not confuse self-monitoring with self-diagnosis. Blood pressure measurements should only be interpreted by a health professional who is familiar with your medical history.
2. Contact your physician if test results regularly indicate abnormal readings.

3. If you are taking medication, consult with your physician to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting with your physician.
4. Individuals with serious circulation problems may experience discomfort. Consult your physician before use.
5. For persons with irregular or unstable circulation resulting from diabetes, liver disease, arteriosclerosis or other medical conditions, there may be variations in blood pressure values measured at the wrist versus at the upper arm. Monitoring the trends in your blood pressure taken at either the arm or the wrist is nevertheless useful and important.
6. People suffering from vascular constriction, liver disorders or diabetes, people with cardiac pacemakers or a weak pulse, and women who are pregnant should consult their physician before measuring their blood pressure themselves. Different values may be obtained due to their condition.
7. People suffering from arrhythmias such as atrial or ventricular premature beats or atrial fibrillation should only use this blood pressure monitor in consultation with their physician. In certain cases, the oscillometric measurement method can produce incorrect readings.
8. Too frequent measurements can cause injury to the user due to blood flow interference.
9. The cuff should not be applied over a wound as this can cause further injury.
10. DO NOT attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially harming the patient.


11. The cuff should not be placed on the arm on the mastectomy side. In the case of bilateral mastectomy, use the least dominant arm.
12. Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring equipment on the same limb.
13. A compressed or kinked connection hose may cause constant cuff pressure resulting in blood flow interference and potentially harmful injury to the patient.
14. Check that operation of the device does not result in prolonged impairment of the circulation of the patient.
15. Product is designed for its intended use only. Do not misuse in any way.
16. Product is not intended for infants or individuals who cannot express their intentions.
17. Prolonged over-inflation of the bladder may cause ecchymoma of your arm.
18. Do not disassemble the device or arm cuff. Do not attempt to repair.
19. Use only the approved arm cuff for this device. Use of other arm cuffs may result in incorrect measurement results.
20. The system might produce incorrect readings if stored or used outside the manufacturer's specified temperature and humidity ranges. Make sure to store the blood pressure monitor properly, and make sure children, pets and pests are outside of accessible range.
21. Do not use the device near strong electrical or electromagnetic fields generated by cell phones or other devices, they may cause incorrect readings and interference or become interference source to the device.
22. Do not use new and old batteries simultaneously.
23. Replace batteries when Low Battery Indicator "  " appears on screen. Replace both batteries at the same time.

24. Do not mix battery types. Long-life alkaline batteries are recommended.
 25. Remove batteries from device when not in operation for more than 3 months to avoid corrosion and damage to the device.
 26. Dispose batteries properly; observe local laws and regulations.
 27. Only use a recommended class II AC Medical approved adaptor which comply with 2MOPP, and pass IEC 60601-1, IEC 60601-1-2 (or FCC) standard. An unauthorized adapter may cause fire and electric shock.
 28.  Always consult this instruction manual/ booklet for use information.
 29. Do not use the device in vehicles for influencing measurement accuracy.
 30. Contains small parts that may cause a choking hazard if swallowed by infants.
 31. Please align the polarities of each battery with the + and - signs imprinted on the battery compartment when you replace the batteries.
 32. Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
 33. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Protections include the following: IP21 indicate dustproof and waterproof grade. Number 2 indicate that Protected against solid foreign objects of 12.5mm and greater; Number 1 indicate that Protection against vertically falling water drops .


WARNING SIGNS AND SYMBOLS USED							
	Keep Dry		Refer to instruction manual		Keep away from sunlight		Type BF applied part
	Batch code		Manufacturing date		Recyclable		Serial number
	This device has not been tested for use in an MR environment, and should not be used exposed to MR environments while patients are wearing the device. Keep it outside the MPI scanner room.				IP21 indicate dust proof and waterproof grade. Number 2 indicate that Protected against solid foreign objects of 12.5 mm and greater; Number 1 indicate that Protection against vertically falling water drops .		
	WEEE (Waste Electrical and Electronic Equipment)						

Device Illustration

Monitor Device


Arm Cuff Large size cuff (fits arm circumference: 22.0 cm -42.0 cm).


If air is leaking from the arm cuff, replace the arm cuff with a new one. It is generally recommended to have the cuff replaced in a timely manner to ensure correct functioning and accuracy. Please consult your local authorized distributor.

Display

Contents

Important Testing Guidelines

1. Avoid eating, exercising, or bathing 30 minutes before the testing.
2. Sit in a calm environment for at least 5 minutes prior to testing.
3. Do not stand while testing. Sit in a relaxed position while keeping your arm level with your heart.
4. Avoid speaking or moving body parts while testing.
5. While testing, avoid strong electromagnetic interference such as microwave ovens and cell phones.
6. Wait 3 minutes or longer before re-testing.
7. Try to measure your blood pressure at the same time each day for consistency.
8. Test comparisons should only be made when monitor is used on the same arm, in the same position, and at the same time of day.

9. This blood pressure monitor is not recommended for people with severe arrhythmia.
10. Do not use this blood pressure monitor if the device is damaged.

Any blood pressure recording can be affected by the following factors:

1. The physical position of the subject and his or her physiologic condition;
 2. The performance and accuracy of the device;
 3. Cuff size: too small cuff will produce a higher blood pressure value than usual, too big cuff will produce a lower blood pressure value;
 4. Measuring position does not keep level with your heart;
 5. Speaking or moving body parts while testing;
- Not relaxing for about 5 minutes before taking the measurement.

Quick Start

1. Install batteries. (See Figure A)
2. Insert cuff air plug into the back of monitor device. (See Figure B)

Figure A

Figure B

3. Roll your sleeve on the arm area.
4. Rest for several minutes prior to testing. Sit down comfortably in a quiet place, back and arm supported on a chair and table with your legs uncrossed, your arm should be resting on a firm table or desk with your feet flat on the floor. (See Figure C)

Figure C

5. Apply cuff to your left arm and middle of the cuff at the level of your heart (unless any medical condition requires you to use your right arm). Bottom of cuff should be placed approximately 1-2cm (1/2") above elbow joint. (See Figures D&E)


Figure D

Figure E


6. Press " START/STOP " Button to start testing.

Device Operation

Battery Installation

Slide battery cover off as indicated by arrow.

Install 3 new AAA alkaline batteries according to polarity. Close battery cover.

This blood pressure provides an external power interface. The power supply port of the device is USB-C port, and the port applicable parameters is voltage 5v and current 1000mA. Do not use other types of AC adapters. Otherwise, the device may be damaged.

USB Adapter Jack

Note: Power supply is specified as part of ME EQUIPMENT.

System Settings

With power off, press "  " button to activate System Settings. The Memory Group icon flashes.

1. Select Memory Group


While in the System Setting mode, you may accumulate test results into 2 different groups. This allows multiple users to save individual test results (up to 60 memories per group.) Press " MEM " button to choose a group setting. Test results will automatically store in each selected group.

2. Time/Date setting

Press "  " button again to set the Time/Date mode. Set the year first by pressing the "MEM" button.

Press "  " button to confirm and lock in your selection for each. Continue to set the month, date, hour, minute, and 12/24 hour in the same way.

3. Time Format setting

Press "  " button again to set the time format setting mode. Set the time format by adjusting the "MEM" button. EU means European Time US means U.S Time.


4. Saved Settings

While in any setting mode, press " **START/STOP** " button to turn the device off. All information will be saved.

Note: If device is left on and not in use for 3 minutes, it will automatically save all information and shut off.


Applying the Arm Cuff

1. Firmly insert air plug into opening located on the back side of monitor device.

2. With the sticky nylon section facing outward, insert the end of cuff underneath the metal ring of cuff.

3. Fasten cuff about 1-2cm (0.4-0.8") above the elbow joint. For best results apply cuff to bare arm and keep level with heart while testing.

Note: Do not insert air plug into opening located on right side of monitor device.

This opening is designed for an optional power supply only.

Taking Reading

1. Power On

Press and hold "START/STOP" button to turn the device on. The LCD screen will appear for one second as device performs a quick diagnosis. A voice will indicate when device is ready for testing.

Note: Device will not function if residual air from previous testing is present in cuff.

The LCD will flash "↓" until pressure is stabilized.

2. Testing

After cuff inflation, air will slowly rise as indicated by the corresponding cuff pressure value. A flashing "♥" will appear simultaneously on-screen signaling heartbeat detection.

Note: Keep relaxed during testing. Avoid speaking or moving body parts.

3.Result Display

The screen will display measurements for systolic and diastolic blood pressure. An indicator representing the current measurement will appear next to the corresponding WHO Classification.

Note: Refer to Page 22~23 for detailed WHO Blood Pressure Classification Information.

Irregular Heartbeat Indicator

If the monitor detects an irregular heart rhythm two or more times during the measuring process, the Irregular Heartbeat Symbol "⊠" appears on screen along with measurement results. Irregular heartbeat rhythm is defined as rhythm that is either 25% slower or faster than the average rhythm detected while measuring systolic blood pressure and diastolic blood pressure. Consult your physician if the Irregular Heartbeat Symbol "⊠" frequently appears with your test results or if concerned.

Power Off

The "START/STOP" button can be pressed to turn off the device in any mode. The device automatically turns off the power itself about 3 minutes with no operation in any mode.

Safety Precaution: If pressure in arm cuff becomes too extreme while testing, press the "START/STOP" button to turn power off. The cuff pressure will rapidly dissipate once the device is off.

Last 3 Tests Average

With power off, press the "MEM" button to activate screen display. After the device performs a self-diagnosis, the screen will display the average test results from the last 3 readings of the last group used. The "AVG" symbol will appear along with the corresponding WHO Blood Pressure Indicator. The Memory check mode can be accessed by pressing "MEM" button. To check average results from other groups, select the desired group first prior to activating "SET" button in the off position. (See "Select Memory Group" on Page 11)

Memory Check

You may check previous test results by using the "MEM" button. The most recent test result and oldest test result in memory can be viewed by pressing and holding the "MEM" button. Upon activating test results, you can press the "MEM" button to scroll through all test results stored in memory.

Memory Deletion

Memory for a selected group may be deleted while in Memory Check mode. Press and hold the "⊠" button for approximately 3 seconds to delete all memory records from the selected group. With voice "Beep" and then the device switches to testing mode. Press the "START /STOP" button to turn the device off.

Note: Memory record cannot be recovered once it has been deleted.


Low Battery Indicator

The "⚡" appears for approximately 5 seconds prior to shutting off. Replace batteries at this time. No memory loss will occur throughout this process.

Static Pressure Measurement

Under the low battery state, press and hold the "START/STOP" button, and then install the batteries. Press the "START/STOP" button. Until all signs shows on the screen, release the button. When the LCD screen displays the double zero, the blood pressure meter is in the static state. Software version is displayed: 10 is a software version in the figure.

Note: Only professional service personnel permitted to access to this mode, the mode unavailable in normal use.

Arm Shake Indicator

If there is arm movement during the measurement, the "⚡" icon may flash. Indicates that the measurement results may be inaccurate, and the situation will be recorded at the end of the measurement as a reminder.

Cuff Loose Indicator

At the beginning of the measurement, when the cuff is properly worn, the "OK" icon will be displayed, and when the cuff is loose, The "⚡" icon will be displayed. When this happens, please turn off the machine, adjust the cuff and measure again. If the "⚡" icon is still displayed, please check if the cuff plug is connected to the device correctly and measure again. If it still cannot be resolved, please contact the seller or manufacturer.

Note that if another family member uses it:

- 1) Clean the cuff before use, and use it after it dries;
- 2) Press "SET" to set the memory group
- 3) Restart the measurement according to the operation steps.
- 4) If you want to view a past measured value, please refer to the memory view description on page 17.

Troubleshooting

Abnormal phenomenon	Cause analysis	Processing method
Abnormal sphygmomanometer	The arm cuff is tied too tight or too loose, Or the arm cuff is tied incorrectly.	Attach the arm cuff correctly
	Moved the arm during measurement or Electronic sphygmomanometer	Stay quiet, keep your arm steady, and do not move the monitor
	Speaking, nervous or emotional during measurement	Instead of talking, take deep breaths to calm your mood and relax your body
	Incorrect measurement posture	Adjust posture, see "Blood pressure gauge Wearing"
	There is interference in charging process or improper operation in measuring process.	See operation instructions.

The following table shows the error signs that may occur during measurement, possible causes and handling methods. Please measure again using the correct method

Error display	The cause of the problem	The solution
Er1	Can't detect high and low pressure	Please tighten the cuff before measuring
Er2	Cuff too loose or loose	Please tighten the cuff before measuring
Er3	Improper compression caused by arm or body movement	Hold the arm or body still and measure again
Er4	The pressure exceeds 300mmHg	Please tighten the cuff before measuring
Er5	The pressure exceeds 15mmHg for 3 minutes	Check whether the cuff is inflated or the vent valve is blocked. If the problem persists, contact the manufacturer
	Battery dead	Replace the battery or connect the power adapter (if any).

Note: If you cannot solve the abnormal situation by yourself, you can consult the manufacturer. It is forbidden to disassemble and repair without permission. If necessary, professional maintenance personnel can ask the manufacturer for the list of components and circuit schematic diagram.


Blood Pressure Information

Blood Pressure

Blood pressure is the force of blood pushing against the walls of arteries. It is typically measured in millimeters of mercury (mmHg.) Systolic blood pressure is the maximum force exerted against blood vessel walls each time the heart beats. Diastolic blood pressure is the force exerted on blood vessels when the heart is resting between beats.


An individual's blood pressure frequently changes throughout the course of a day. Excitement and tension can cause blood pressure to rise, while drinking alcohol and bathing can lower blood pressure. Certain hormones like adrenaline (which your body releases under stress) can cause blood vessels to constrict, leading to a rise in blood pressure.

If these measuring numbers become too high, it means the heart is working harder than it should.

WHO Blood Pressure Classification Indicator


The Easy@Home EBP-073 is equipped with a classification indicator based on established guidelines from the World Health Organization. The chart below (color coded on monitor device) indicates test results.

☐: Blood Pressure Classification Indicator

Health Reminder

Hypertension is a dangerous disease that can affect the quality of life. It can lead to a lot of problems including heart failure, kidney failure, and cerebral hemorrhaging. By maintaining a healthy lifestyle and visiting your physician on a regular basis, hypertension and relative diseases are much easier to control when diagnosed in their early stages.

Note: Do not be alarmed if an abnormal reading occurs. A better indication of an individual's blood pressure occurs after 2-3 readings are taken at the same time each day over an extended period of time. Consult your physician if test results remain abnormal.

Blood Pressure Q&A

Q: What is the difference between measuring blood pressure at home or at a professional healthcare clinic?

A: Blood pressure readings taken at home are now seen to give a more accurate account as they better reflect your daily life. Readings can be elevated when taken in a clinical or medical environment. This is known as White Coat Hypertension and may be caused by feeling anxious or nervous.

Note: Abnormal test results may be caused by:

1. Improper cuff placement

Make sure cuff is snug- not too tight nor too loose.

Make sure the bottom of the cuff is approximately 1-2cm (1/2") above the elbow joint.

2. Improper body position

Make sure to keep your body in an upright position.

3. Feeling anxious or nervous

Take 2-3 deep breaths, wait a few minutes and resume testing.

Q: What causes different readings?

A: Blood pressure varies throughout the course of a day. Many factors including diet, stress, cuff placement, etc. may affect an individual's blood pressure.

Q: Should I apply the cuff to the left or right arm? What is the difference?


A: Either arm can be used when testing; however, when comparing results, the same arm should be used. Testing on your left arm may provide more accurate results as it is located closer to your heart.

Q: What is the best time of day for testing?

A: Morning time or any time you feel relaxed and stress free.

Maintenance

1. Avoid dropping, slamming, or throwing the device.

2. Avoid extreme temperatures. Do not expose device directly under sunshine.

3. When cleaning the device, use a soft fabric and lightly wipe with mild detergent. Use a damp cloth to remove dirt and excess detergent. It is recommended to clean once in January to March under normal circumstances. If the previous cleaning instructions are not successful, repeat several times until the visible dirt is removed.

4. Cuff Cleaning and Disinfection:

A) Spread the cuff (skin-contact surface) upwards onto a clean table. Use a damp clean cloth (water-based) to wipe the skin-contact surface with a force.

B) Soak the cloth clean with drinking water and wring it dry. Repeat A) with the damp cloth (water-based) for 3 times.

C) Apply 70%-80% alcohol to a new cloth (or 75% alcohol cotton-ball), use it to wipe the skin-contact surface with a force. Then soak the cloth with the alcohol again (or change a new 75% alcohol cotton-ball), repeat the disinfection procedure for 3 times.

D) When the disinfection towards the skin-contact surface is finished, wipe the non-skin contact surface with a cloth (alcohol-based) or alcohol cotton-ball thoroughly for 3 times.

E) Leave the cuff to naturally dry, then it is ready for reuse.

Notice: Do not soak in water or splash water on it.

5. Do not use petrol, thinners or similar solvents.

6. Remove batteries when not in operation for an extended period to avoid corrosion and damage to the device.

7. Do not disassemble/ try to repair product.

8. It is recommended the performance should be checked every 2 years.

9. Expected service life: Approximately three years at 10 tests per day.

No service and maintenance while it is in use and maintenance only be performed by service personnel.

Specifications

Product Description	Automatic Digital Blood Pressure Monitor	
Model	EBP-073	
Display	LCD Digital Display Size:95.5mm*54.5mm (3.76" x 2.15")	
Measurement Method	Oscillometric Method	
Measurement Range	Systolic Pressure	60mmHg ~ 260mmHg
	Diastolic Pressure	30mmHg ~ 200mmHg
	Pressure	0mmHg ~ 290mmHg
	Pressure	±3mmHg
	Pulse	30 ~ 180 Beats/Minute
	Pulse	±5%
Pressurization	Automatic Pressurization	
Memory	2x60 Memories in Two Groups with Date and Time	

Function	Irregular Heartbeat Detection	
	WHO Classification Indicator	
	Last 3 Tests Average	
	Low Battery Detection	
	Automatic Power-Off	
Power Source	3 AAA batteries or Medical AC Adapter (DC5.0V, 1000mA) (recommended, not provided)	
Battery Life	Approximately 2 months at 3 tests per day	
Device Weight	Approx. 265g (9.35 oz.) (excluding battery)	
Device Dimensions	Approx. 142.5 x 107.2 x 44mm (5.61" x 4.22" x 1.73") (L x W x H)	
Cuff Circumference	Large cuff. Fits arm circumference 22-42 cm.	
Operating Environment	Temperature	10°C ~ 40°C (50°F ~ 104°F)
	Humidity	15% ~ 93% RH
	Pressure	800hPa ~ 1060hPa
Classification	Internal Powered Equipment, Type BF  . Cuff is the Applied Part	
Ingress Protection rating	IP 21, Indoor Used Only	
Battery Shelf life	60 months	
Battery Storage Temperature	-25°C ~ 55°C (-13°F ~ 131°F)	
Storage Environment	Temperature	-25°C ~ 55°C (-13°F ~ 131°F)
	Humidity	≤93% RH
Transport Environment	Temperature	-25°C ~ 55°C (-13°F ~ 131°F)
	Humidity	≤93% RH

Specifications are subject to change without notice.

Safety Standard (included but not limited):

1. IEC 80601-2-30, medical electrical equipment - part 2-30: particular requirements for the basic safety and essential performance of automated noninvasive sphygmomanometers. (Cardiovascular)
2. ISO 81060- 2, non- invasive sphygmomanometers - part 2: clinical validation of automated measurement type. (Cardiovascular)
3. AAMI / ANSI ES 60601- 1: 2005/(R) 2012 and C1: 2009/(R) 2012 and, a2: 2010/(r) 2012 (consolidated text) medical electrical equipment – part 1: general requirements for basic safety and essential performance
4. AAMI/ ANSI/ IEC 60601- 1- 2, Medical Electrical Equipment – Part 1- 2: General Requirements For Basic Safety And Essential Performance – Collateral Standard: Electromagnetic Disturbances – Requirements And Tests (General II (ES/ EMC)).
5. IEC 60601-1-11, medical electrical equipment - part 1-11: general requirements for basic safety and essential performance - collateral standard: requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment.

Warranty

The Blood Pressure Monitor is guaranteed for 1-year from the date of purchase. If the Blood Pressure Monitor does not function properly due to defective components or poor workmanship, we will repair or replace for no additional cost. The warranty does not cover damages to your Blood Pressure Monitor due to improper handling, loss, or theft. Please contact manufacturer for details.

FCC Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Electromagnetic Compatibility Information

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regards to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the device in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

Table 1

Guidance and declaration of manufacturer-electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment-guidance
Radiated emission CISPR 11	Group 1, class B.	The device uses RF energy only for its internal function. Therefore, its emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted emission CISPR 11	Group 1, class B.	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A.	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 2

Guidance and declaration of manufacturer-electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMDEVICEY test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient/burst IEC 61000-4-4	± 2 kV, 100kHz, for AC power port	± 2 kV, 100kHz, for AC power port	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5kV, ± 1kV (differential mode)	± 0.5kV, ± 1kV (differential mode)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT, 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	0 % UT, 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m; 50Hz or 60Hz	30 A/m; 50Hz or 60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 2(continued)

Guidance and declaration of manufacturer-electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMDEVICEY test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Radiated RF EM fields IEC 61000-4-3	10 Vm80MHz-2.7GHz 80%AM at 1kHz	10 Vm80MHz-2.7GHz 80%AM at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: 80 MHz to 800 MHz: 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V in 0.15 MHz-80 MHz 6 V in ISM and/or amateur radio bands between 0.15MHz and 80 MHz 80 % AM at 1kHz	3 V in 0.15 MHz-80 MHz 6 V in ISM and/or amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: 80 MHz to 800 MHz: 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:

Table 3

Guidance and declaration of manufacturer-electromagnetic immunity						
Nowadays, many RF wireless equipment have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. Arm-Type Fully Automatic Digital Blood Pressure Monitor has been tested with the immediacy test level in the below table and meet the related requirements of IEC 60601-1-2:2014. The customer and/or user should keep a minimum distance between RF wireless communications equipment and this medical equipment and/or systems as recommended below.						
Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immdevicey test level (V/m)
385	380-390	TETRA 400	Pulse modulation 16Hz	1.8	0.3	27
450	430-470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710 745 780	704-767	LTE Band 13, 17	Pulse modulation 217Hz	0.2	0.3	9
810 870 930	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation 16Hz	2	0.3	28
1720 1845 1970	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz	2	0.3	28
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217Hz	2	0.3	28
5240 5800	5100-5800	WLAN 802.11	Pulse modulation 217Hz	2	0.3	9

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device		
The device is intended for use in an electromagnetic environment in which radiated therefore disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.		
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m	
	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.7 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.		
NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

Additional Notes

1. Verify manometer pressure accuracy:

In the power down state, press and hold the " START/STOP" button, and then install the batteries. Until the LCD screen is full, release the " START/STOP" button. When the LCD screen displays the double zero, the blood pressure meter is in static state. At this point, 500ml gas capacity, calibrated standard pressure gauge and manual pressure device can be connected to the sphygmomanometer through the sleeve interface of the sphygmomanometer, and manual pressure can be applied to the effective display range of the sphygmomanometer, and then the difference between the reading of the sphygmomanometer and that of the standard pressure gauge can be compared. This mode can be used to verify manometer pressure accuracy.

2. The user of the device is the operator:

The user of the device is an intended OPERATOR.

The USER Do Not carry out other maintenance operations except to replace the battery.

3. WARNING:

Do not modify this equipment without authorization of the manufacturer.

4. ESSENTIAL PERFORMANCE Maintenance advice:

Pressure calibration will be carried out when this product leaves the factory. Patients can use the method described in the section "Verify Manometer Pressure Accuracy" to verify the accuracy. If the accuracy deviation is large, please contact the manufacturer to recalibration.

5. Mechanical strength and resistance to heat-the resistance to heat will be retained by device during the EXPECTED SERVICE LIFE of the ME EQUIPMENT.

6. Do not place the blood pressure monitor and cuff where children or pets can reach. It will cause asphyxiation if the child or pet swallows parts or can cause strangulation around his neck with the tube.

7. The cuff and the case of the blood pressure monitor have been tested for biocompatibility and do not contain allergenic or harmful materials. Please stop using it if allergy occurs during use.

8. Warning:

Non-professionals do not modify the equipment, otherwise it will make the measurement of the device inaccurate.

9. Warning:

Do not expose the equipment to extreme heat (such as sunlight) or cold for a long time, otherwise it will reduce the performance of the equipment.

10. Warning:

This device is not used for children and pets. Adults must keep all parts/device away from them.

11. Warning:

Do not use a damaged cuff for blood pressure measurement.

12. Warning:

When measuring with the cuff, if the tester feels seriously uncomfortable, press the button of the blood pressure monitor to deflate the cuff, or remove the cuff directly from the arm.

13. Warning:

If an unexpected reading occurs, the operator can take several more measurements and consult a doctor.

14. Warning:

If this device is not used in the prescribed environment, it may damage the equipment, or lead to inaccurate measurement.

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