# Dr Brown's

natural flow®

## Customflow™ **Double Electric Breast Pump** with SoftShape™ Silicone Shields

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Número de modelo BF100

## WARNING

- If you are a mother who is infected with Hepatitis B, Hepatitis C or Human Immunodeficiency Virus (HIV), pumping breast milk will not reduce or remove the risk of transmitting the virus to your baby through your breast milk.
- Powered breast pumps that are designed for single users should never be rented or shared.
- This breast pump should not be operated while traveling on an airplane.

## Contraindications

There are no known contraindications.

## Precautions

- If discomfort is felt, suction can be broken by inserting a finger between the breast and the breast shield.
- If irritation or discomfort continues to occur, discontinue use and contact your healthcare provider.
- Before each use, visually inspect the individual components for cracks, chips, tears, discoloration or deterioration. In the event that damage to the device is observed, please discontinue use until the parts have been replaced.
- Inspect the tubing. If you find condensation in the air tube, wash and dry it immediately or replace the tubing. If you see breast milk inside the tubing or tubing connector, do not attempt to wash, clean or sterilize the tubing - contact Dr. Brown's Customer Service. Outside of the US and Canada, contact your local distributor.

## Intended Use

The Dr. Brown's Natural Flow® Customflow™ Double Electric Breast Pump is a powered breast pump used by lactating women to express and collect milk from their breasts. This pump is intended for a single user.

## Product Description

The Dr. Brown's Natural Flow® Customflow™ Double Electric Breast Pump is a personal use electric breast pump capable of double pumping. The device is electrically powered from an external AC powered supply and is provided non-sterile.

The Double Electric Breast Pump consists of the pump unit, breast shields, tubing, bottle vessel, duckbill valve, and other related accessories. The materials used are silicone (SoftShape™ Silicone Shields, nipples, membranes, duckbill valves), polypropylene (pump body and bottle vessels) and acrylonitrile butadiene styrene (pump unit). The electric breast pump's drive unit employs a diaphragm type vacuum pump, powered by a DC motor, supervised by a microcontroller. The microcontroller provides control over motor speed (vacuum creation) and solenoid (vacuum release).

The pump is reusable and intended for direct user contact. The pump contains both backflow and overflow protection. The pump is not meant to be sterilized.

Additional parameters: Maximum suction pressure with no software control is -295mmHg. The shelf life of the pump is 2 years. This pump complies with safety standards IEC 60601-1 and IEC 60601-1-11.

## Model BF100 has 2 modes:

1. Letdown Mode – Used to stimulate the release of milk from the breast.
2. Expression Mode – Used after the breast is stimulated to release milk.

### **PLEASE DO NOT RETURN THIS PUMP TO THE RETAILER**

If you experience any issues with your Dr. Brown's Pump, please contact our dedicated Breast Pump Customer Service team by calling 1- 833-224-7675.

Outside of the US and Canada, contact your local distributor.

# Table of Contents

|                                    |    |
|------------------------------------|----|
| Important Safeguards.....          | 3  |
| Breast Pump Components .....       | 6  |
| Cleaning.....                      | 8  |
| Assembly .....                     | 10 |
| Using the Pump .....               | 14 |
| Pump Help .....                    | 17 |
| Storing & Feeding Breast Milk..... | 19 |
| Technical Data .....               | 21 |
| Electromagnetic Emissions.....     | 22 |
| Troubleshooting .....              | 24 |
| Important Symbols.....             | 25 |
| Disposal .....                     | 26 |
| Warranty .....                     | 27 |

# IMPORTANT SAFEGUARDS

**WARNING:** Read this ENTIRE Instruction Manual to become familiar with the features of the product before operating. Failure to operate the product correctly can result in damage to the product, personal property and cause serious injury. Failure to operate this product in a safe and responsible manner could result in injury or damage to the product or other property. This product is not intended for use by children. Do not use with incompatible components. This Manual contains instructions for safety, operation and maintenance. It is essential to read and follow all the instructions and warnings in the manual, prior to setup or use, in order to operate correctly and avoid damage or serious injury.

### **DANGER - To reduce the risk of electrocution:**

- Always unplug electrical product immediately after use.
- Do not use while bathing.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- DO NOT IMMERSE IN WATER.
- Do not reach for a product if it has fallen into water. Immediately unplug from electrical outlet.
- This appliance has a polarized plug (one blade is wider than the other). To reduce the risk of electric shock, this plug is intended to fit into a polarized outlet only one way. If the plug does not fit fully into the outlet, reverse the plug. If it still does not fit, contact a qualified electrician. Do not attempt to modify the plug in any way.

### **WARNING - To avoid fire, electrocution, or injury:**

- Danger of electrocution! Keep the device dry! Never immerse in water or other liquids!  
Do not run water over the pump.
- The breast pump and accessories are not heat-resistant: keep away from radiators, open flame and sunlight.
- Keep the power adapter away from heated surfaces.
- Do not reach for any electrical device if it has fallen into water. Unplug immediately.
- Do not use outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered.
- The breast pump should never be left unattended when plugged into a power source.

## Important Safeguards

---

### **WARNING (Cont):**

- Never operate an electrical device if it has a damaged cord, plug or battery pack, if it is not working properly, if it has been dropped or damaged, or dropped into water. If damage is found, immediately discontinue use of power adapter or battery pack and call Dr. Brown's Customer Service at 1-833-224-7675. Outside of the US and Canada, contact your local distributor.
- Close supervision is necessary when the breast pump is used in the vicinity of children to prevent strangulation due to power adapter cord or tubing.
- Follow the exact order of the steps in powering breast pump as listed on page 13 of this Instruction Manual.
- Only use the power adapter supplied with the breast pump.
- Make sure the voltage of the power adapter is compatible with the power source.
- Do not use an extension cord.

### **CAUTION - To avoid risk of minor injury:**

- Wash hands thoroughly with soap and water before touching breast pump, kit and breasts and avoid touching the inside of containers or lids.
- Separate and wash all parts that come in contact with breast milk immediately after use to avoid dry up of breast milk residue and prevent growth of bacteria. Always inspect SoftShape™ Silicone Shields, connectors, valves, membranes and tubing prior to use for cleanliness.
- Only use drinking-quality tap or bottled water for cleaning your breast pump and parts.
- Use the breast pump only for its intended use as described in this manual.
- Do not use the breast pump while sleeping or overly drowsy.
- This is a single user product. Use by more than one person may present a health risk.

# Important Safeguards

---

## **CAUTION (Cont):**

- Do not drive while pumping.
- Inspect all pump components before each use.
- Never use a damaged device. Replace damaged or worn parts.
- Do not store wet or damp parts as mold may develop.
- Do not run pump with wet tubing. Doing so will harm the pump motor.
- While some discomfort may be felt when first using a breast pump, using a breast pump should not cause pain.
- If you are experiencing discomfort at the base of the nipple due to rubbing of your breast tissue with the SoftShape™ Silicone Shields, discontinue use and contact your healthcare professional or lactation consultant.
- Contact your healthcare professional or lactation consultant if you can express only minimal or no milk or if expression is painful.
- Do not try to express with vacuum that is too high and uncomfortable (painful). The pain, along with potential breast and nipple trauma, may decrease milk output.
- Make sure tubing is not kinked or pinched while pumping.
- Do not hold the pump kit by the bottle. This can lead to blockage of the milk ducts and engorgement.
- Do not thaw frozen breast milk in a microwave or in a pan of boiling water.
- Do not microwave breast milk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving. (Microwaving can also change the composition of breast milk.)

**SAVE THESE INSTRUCTIONS**

Visit **[www.DrBrownsBaby.com](http://www.DrBrownsBaby.com)** for tips and help.

# Breast Pump Components

2x SoftShape™  
Silicone Shield, 25mm

2x Membrane Cap

2x Membrane

2x Tubing

2x Pump Connector

2x Duckbill Valve


2x Narrow Bottle Vessel

2 x Narrow Travel/Storage Cap


1x Wire Cleaning Brush

1x Slowest Flow/  
Premie Flow™ Nipple

2x Dr. Brown's® Options+™  
Bottle Components


2x Narrow to  
Wide Neck Adapter\*

\*Narrow to Wide Neck Adapter  
provided for use with Dr. Brown's  
Wide Neck bottles (bottles not included).


# Breast Pump Components

1x Customflow™ Double Electric Breast Pump

## BEFORE FIRST USE, AND AFTER EVERY USE

Clean and sterilize all pump parts that come in contact with milk:

## ALWAYS USE DRINKING QUALITY WATER FOR CLEANING!

### Cleaning Instructions for Milk-Contact Parts:


1. Disassemble pump following instructions on page 10.
2. Rinse milk residue from parts under cool running water.
3. Do not place parts directly in the sink. Soak in warm soapy water for several minutes.
4. Use a cloth or soft brush used only for cleaning infant feeding items to clean all parts.
5. Rinse parts under hot running water until all soap has been removed, being careful that temperature of water will not burn or irritate skin.
6. Alternately, wash all parts in the top rack of the dishwasher. Place small items in a closed-top basket or mesh bag. Run dishwasher using hot water and a heated drying/sanitizing cycle.
7. Rest parts on a clean paper towel or rack to air dry. All parts must be completely dry before storing to prevent bacteria and mold growth. Store clean, dry parts in a bag or container until next use.
8. Sanitize/sterilize before each use. (Microwave and electric steam sterilizers may be used to sterilize milk contact parts.)
9. For convenience, special wipes for cleaning breast pump parts may be used when away from home. However, all parts must be cleaned using soapy water or in the dishwasher before pumping again.

### NOTE

- Please clean all accessories according to cleaning guide every time before using.
- Let parts cool after washing in dishwasher or sterilizing. Handle carefully.
- Avoid using cloth towels to dry your pump parts because they can carry germs and bacteria that are harmful to your breast milk and your baby.
- Try not to touch the inside of any parts that will come in contact with your breast milk.
- See separate instructions for bottle cleaning & assembly.

## NON MILK-CONTACT PARTS:

The following parts do not come in contact milk and should only be cleaned as needed.

## Special Cleaning Instructions for Tubing:

1. Tubing does not come in contact with milk and should only be washed if dirty or condensation is present.
2. If condensation is from a previous washing process, there is no need to wash again.
3. Wash with warm soapy water and rinse clean with clear water.
4. Allow water to drain from tubing completely after washing.
5. Shake out any remaining water droplets and hang to dry.
6. Do not use tubes until completely dry.\*

\*A small amount of condensation is allowable; however, a tube with water droplets may damage the pump motor.

**DO NOT SUBMERGE THE  
PUMP MOTOR IN WATER!  
WIPE CLEAN WITH  
DAMP CLOTH ONLY  
AFTER EACH USE.**

Silicone stretch wall may be removed for cleaning, as necessary. Silicone stretch wall may simply be wiped clean with a damp cloth but is also dishwasher and sterilizer safe.

## DISASSEMBLY:

The SoftShape™ Silicone Shields/Milk Collection Kits come pre-assembled. All milk-contact components must be disassembled, washed and sterilized prior to first use.

Pull SoftShape™ Silicone Shield to remove from pump connector.

Twist membrane cap a  $\frac{1}{4}$  turn and lift to remove.

Lift membrane from edge to remove from pump connector.

Pull tab on duckbill valve to remove from pump connector.  
DO NOT PULL ON POINTED END OF VALVE.

Twist bottle to unthread from the pump connector.



## Wash hands before touching & assembling milk collection kit.

All milk-contact components must be disassembled, washed and sterilized prior to first use.


Push neck of SoftShape™ Silicone Shield onto pump connector. Ensure that both rings are attached.

Set membrane cap in place and twist a ¼ turn to attach.

Drop in membrane and press down edges to secure.


From the underside of the connector, push duckbill valve firmly onto the cylindrical piece inside the pump connector.

Twist bottle to thread on to pump connector.



## TUBING:

Male end of tubing connects to the membrane cap.


Female end of tubing connects to port on front of pump.


To use for single-sided pumping, remove the tube that is no longer needed and place the port cover in the empty port.

## STRETCH FIT WALL:


Wrap stretch wall around pump and slide each wall column into its holder from above.



## POWER ADAPTER:


To prevent damage to your power adapter cord, always grasp the body of the adapter as shown above. **Never pull by the cord to unplug.**

Always unplug the adapter from the pump before storage.

Fold the cord, do not wrap. **Wrapping the cord can break the wires.**


## Using the Pump

Before pumping, wash hands thoroughly with warm soapy water. Inspect all breast pump parts for cleanliness - IF UNWASHED, FOLLOW INSTRUCTIONS ON PAGE 8. Look closely at the silicone shield, membrane, and duckbill valve to inspect for any signs of wear or damage. Suction may decrease if parts are damaged or not properly assembled. Contact Dr. Brown's Customer Service to inquire about replacement parts. Outside of the US and Canada, contact your local distributor.

Do not try to express milk if the level of vacuum is painful, this may lead to nipple and/or breast trauma and cause difficulties with expression. Contact your medical professional or lactation consultant if you can only express minimal amounts of milk, or if expression is painful in any way.

### DISPLAY:

## LETDOWN and EXPRESSION MODE:

**Letdown Mode** is used to stimulate the release of milk from the breast. This mode mimics baby's way of rapidly sucking at the beginning of a breastfeeding session, prompting the reflex that releases breast milk.

**Expression Mode** is used once the breast milk has been released and is now expressing from the nipple. This mode mimics baby's more relaxed feeding that happens once milk is flowing and draws milk more deeply from the breast.

## Using Letdown Mode:

1. Press the POWER BUTTON to turn the unit on.
2. The pump will automatically display Letdown Mode at the last vacuum level used.
3. Position the assembled milk collection kit on breast as directed on page 16.
4. Press the START/PAUSE BUTTON to begin pumping. The timer will begin. Press the button again at any time you wish to pause pumping.
5. If needed, adjust vacuum level to your preferred starting level using the VACUUM UP and VACUUM DOWN BUTTONS. There are five user-selectable vacuum levels in the Letdown Mode: L1, L2, L3, L4, L5.

*Adjustments to vacuum level may be made at any time while pumping - adjustments cannot be made to the cycle level in Letdown Mode. It may take several pumping sessions to discover your most comfortable and productive vacuum level.*


## Using Expression Mode:

1. Once letdown has begun (milk is visibly expressing from nipples), switch to Expression Mode.
2. Press the MODE SELECTOR BUTTON. The LETDOWN MODE INDICATOR symbol will disappear when in Expression Mode.
3. The pump will automatically display Expression Mode at the last vacuum level used.
4. If needed, adjust cycle level (speed) to your preferred level using the CYCLE UP and CYCLE DOWN BUTTONS. There are five user-selectable cycle levels in the Expression Mode: 36, 40, 43, 46, 49.
5. If needed, adjust vacuum level to your preferred level using the VACUUM UP and VACUUM DOWN BUTTONS. There are twelve user-selectable vacuum levels in the Expression Mode: L1 through L12.
6. Once pumping is complete (breast milk is no longer expressing), simply press the POWER BUTTON to turn off the pump.

*Adjustments to vacuum & cycle levels may be made at any time while pumping in Expression Mode. It may take several pumping sessions to discover your most comfortable and productive cycle and vacuum levels.*


## Pumping with a Single Milk Collection Kit:

Remove the unneeded hose from the port on the front of the pump. Place the port cover in the empty port. Vacuum strength will intensify when switching from two to one port. Adjust vacuum levels to your preferred comfort level.

## Placement of the Milk Collection Kit:

Position the SoftShape™ Silicone Shield over your breast, centering the nipple in the neck of the shield. Your nipple should not touch or rub on the side walls of the shield. Check nipple positioning again once the pump is operating, to ensure that it remains away from the sides. If not, pause pump and try repositioning. Your breast should completely fill the SoftShape™ Silicone Shield so that a seal may form to achieve proper vacuum levels.

While you may feel some discomfort when first pumping, it should not cause pain. Adjust the cycle and vacuum levels to lower settings to achieve your preferred comfort level and milk flow.

## TIPS FOR EXPRESSING BREASTMILK

**Relax!** Create a comfortable and relaxing setting for your pumping sessions.

**Stay Hydrated.** Drink plenty of water throughout the day and keep a glass of water with you while you pump.

**Create a Routine.** Condition yourself to express milk at the same time each day by creating a routine around pumping. Try pumping each morning after you shower or pumping in the afternoon while watching your favorite television show.

**Think of Baby!** Simply imagining that you are breastfeeding may help induce letdown. If you are away from home, some moms find that a picture of baby or a recording of their cry helps.

**Massage.** A tried and true method - massaging your breasts before you begin pumping can help to quicken letdown. Massage again when you see your milk slowing to help release all milk from the breast.

## COMMONLY ASKED QUESTIONS

### **What is the difference between Letdown Mode and Expression Mode?**

Letdown Mode helps stimulate letdown (milk releasing from breast). Expression Mode is used after Letdown Mode - once the milk is visibly flowing from nipple - to draw milk from breast for collection.

### **What is the difference between cycle and vacuum levels?**

Vacuum levels represent the amount of suction. Cycle levels represent the speed of the on/off suction cycle.

### **How do I find my preferred cycle and vacuum levels?**

Every mom reacts differently to cycle and vacuum levels. Start low and work your way up to a comfortable and productive level. It may take a few pumping sessions to find your ideal settings, and your ideal settings may change over time. The Customflow™ settings are fully customizable so that you can have your perfect pumping experience every time!

### **What do I do if milk gets in the tubing?**

Wash the tubes as directed on page 9 of this instruction manual. Remember to dry the tubing completely before using - water in the tubes may damage the motor. If mold appears in tubing, discard and replace.

### **What do I do if there is low/no suction?**

Check that the tubing is securely connected. Disassemble and check parts for damage. Reassemble, ensuring that parts are connected properly, creating a good seal for the vacuum. While pumping, ensure that your breast is also creating a seal by fully filling the shield.

### **Why does my pump turn off after I have been pumping for a while?**

Your pump has a safety feature that shuts off the pump after 30 minutes.

### **If I have a problem, can I take my pump to a local repair center?**

No. Please contact Dr. Brown's dedicated Breast Pump Customer Service team at 1-833-224-7675 with any concerns about how your pump is operating. If outside of the US and Canada, please contact your local distributor.

**After every use**

- **Store milk safely.** Seal milk collection bottle or milk collection bag, label with date and time, and immediately place in a refrigerator, freezer, or cooler bag with ice packs. Refer to "Human Milk Storage Guidelines" on page 19 for more details.
- **Wipe down** breast pump with a clean soft towel, and disinfect pumping area.
- **Take apart** breast pump tubing and separate all parts that come in contact with breast/breast milk.
- **Rinse** breast pump parts that come into contact with breast/breast milk by holding under running water to remove remaining milk.
- **Clean** pump parts that come into contact with breast/breast milk as soon as possible after pumping. Follow the cleaning steps provided on page 8.

## Guidelines Are for Home Use and Not for Hospital Use.

Use the following guidelines for the storage of breast milk.

By following recommended storage and preparation techniques, breastfeeding mothers and caretakers of breastfed infants and children can maintain the safety and quality of expressed breast milk for the health of the baby. These are general guidelines for storing breast milk at different temperatures. Various factors (milk volume, room temperature when milk is expressed, temperature fluctuations in the refrigerator and freezer, and cleanliness of the environment) can affect how long breast milk can be stored safely. Always keep milk covered when storing.

| Human Milk Storage Guidelines                                                                                            |                                                           |                            |                                                                |
|--------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|----------------------------|----------------------------------------------------------------|
| Storage Locations and Temperatures                                                                                       |                                                           |                            |                                                                |
| Type of Breast Milk                                                                                                      | Countertop<br>77°F or colder (25°C)<br>(room temperature) | Refrigerator<br>40°F (4°C) | Freezer<br>0°F (-18°C)                                         |
| Freshly Expressed<br>or Pumped                                                                                           | Up to 4 hours                                             | Up to 4 days               | Within 6 months<br>is best<br>Up to 12 months<br>is acceptable |
| Thawed, Previously<br>or Pumped                                                                                          | 1-2 hours                                                 | Up to 1 day<br>(24 hours)  | Never<br>refreeze human milk<br>after it has<br>been thawed    |
| Leftover from<br>a Feeding<br>(baby did not finish bottle)                                                               | Use within 2 hours after<br>the baby is finished feeding  |                            |                                                                |
| Adapted from "ABM Clinical Protocol #8: Human Milk Storage Information for Home Use For Full-Term Infants," Revised 2017 |                                                           |                            |                                                                |

Breast Milk Storage Guidelines Center for Disease Control and Prevention ("CDC");  
[www.cdc.gov/breastfeeding/recommendations/handling\\_breastmilk.htm](http://www.cdc.gov/breastfeeding/recommendations/handling_breastmilk.htm)

# Storing & Feeding Breastmilk

---

## STORE

- Store breast milk in storage bags made for freezing human milk or clean glass or hard BPA-free plastic bottles with tight fitting lids. Do not use containers with the number 7 which may contain BPA. Do not use disposable bottle liners or other plastic bags to store breast milk.
- Label milk with date and baby's name if delivering to childcare.
- After each pumping session you can:
  - Keep milk at room temperature. Breast milk is safe for up to 4 hours at room temperature (up to 77°F) after pumping.
  - Refrigerate it. Breast milk is safe in the refrigerator for up to four days.
  - Place milk in the freezer. If you do not plan to use milk within four days of pumping, freeze right away.
  - Use cooler packs. You can store breast milk in an insulated cooler with frozen ice packs for up to 24 hours after pumping. Store in refrigerator or freezer after 24 hours.
  - Do not store milk on the door of the refrigerator or freezer.
- Freeze milk in small amounts of 2 to 4 ounces to avoid waste when warming.
- Breastmilk expands when frozen. Leave a little space at the top of your container. Wait to tighten caps until milk is completely frozen.

## THAW

- Always thaw and use the oldest milk first.
- Thaw milk under warm running water, in a container of warm water, in the refrigerator overnight, or in breast milk warming machine.
- Never thaw milk in the microwave. Microwaving destroys nutrients baby needs, and hot spots created in the microwave can burn baby's mouth.
- Use milk within 24 hours of thawing.
- Use thawed or warmed milk within two hours.
- Never refreeze thawed milk.


## FEED

- To heat milk, place the sealed container under warm running water, in a container of warm water, or in breast milk warming machine.
- Do not heat milk on the stove or in the microwave.
- Test the temperature before feeding it to your baby by putting a few drops on your wrist.
- Swirl the milk, mixing the fat back into the milk.
- Leftover milk should not be used after two hours.


## Double Electric Pump Motor

**Vacuum Range:** -35 to -230 mmHg    **Cycle Rate:** 36 – 63 cpm    **Mode of operation:** Continuous    **Weight:** 715 g

### Operation Conditions:

### Transport/Storage Environment:

### Power Supply

**AC Wall Plug Adapter:** Input: 100-240VAC 50/60 Hz 0.3A Max  
Output: 6 VDC, 2.0A

**AC Adapter Model:** YMK-12W060200

## EMF Compatibility

### Wireless Statement:

Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect the electric breast pump and should be kept at least a distance 1.0 m away from the equipment. Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007, Table 1)

### FCC Statement:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

### Stacking:

The electric breast pump should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the electric breast pump should be observed to verify normal operation in the configuration in which it will be used.

### Note:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced technician for help.

# Electromagnetic Emissions

The Dr. Brown's Natural Flow® Customflow™ Double Electric Breast Pump is intended for use in the electromagnetic environment specified below. The user of the Breast Pump should assure that it is used in such an environment.

## **Guidance and manufacturer's declaration – electromagnetic emissions**

| Emissions                                            | Compliance     | Electromagnetic environment – guidance                                                                                                                                                                                       |
|------------------------------------------------------|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emissions CISPR 11                                | Group 1        | The Breast Pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| RF emissions CISPR 11                                | Class B        |                                                                                                                                                                                                                              |
| Harmonic emissions IEC 61000-3-2                     | Class A        |                                                                                                                                                                                                                              |
| Voltage fluctuations/flicker emissions IEC 61000-3-3 | Not Applicable |                                                                                                                                                                                                                              |

## **Guidance and manufacturer's declaration – electromagnetic immunity**

| Immunity test                                                                                       | IEC 60601 test level                                                                                                                                           | Compliance level           | Electromagnetic environment – guidance                                                                                                                                                                                                                                                 |
|-----------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD) IEC 61000-4-2                                                         | ±6 kV contact<br>±8 kV air                                                                                                                                     | ±6 kV contact<br>±8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.                                                                                                                                         |
| Electrical fast transient/burst IEC 61000-4-4                                                       | ±2 kV for power supply lines<br>±1 kV for input/output lines                                                                                                   | Not Applicable             | Mains power quality should be that of a typical commercial or hospital environment                                                                                                                                                                                                     |
| Surge IEC 61000-4-5                                                                                 | ± 1 kV line(s) and neutral                                                                                                                                     | Not Applicable             | Mains power quality should be that of a typical commercial or hospital environment                                                                                                                                                                                                     |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <5 % UT (>95 % dip in UT) for 0,5 cycle<br>40 % UT (60 % dip in UT) for 5 cycles<br>70 % UT (30 % dip in UT) for 25 cycles<br><5 % UT (>95 % dip in UT) for 5s | Not Applicable             | Mains power quality should be that of a typical commercial or hospital environment. If the user of the Breast Pump requires continued operation during power mains interruptions, it is recommended that the Breast Pump be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8                                             | 3 A/m                                                                                                                                                          | 3 A/m                      | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                              |

NOTE: UT is the AC mains voltage prior to application of the test level

# Electromagnetic Emissions

## Guidance and manufacturer's declaration – electromagnetic immunity

| Immunity test                 | IEC 60601 test level       | Compliance level | Electromagnetic environment – guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|----------------------------|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted RF<br>IEC 61000-4-6 | 3V rms<br>150 kHz to 80MHz | Not Applicable   | <p>Portable and mobile RF communications equipment should be used no closer to any part of the Breast Pump, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 1,2\sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation Distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of Equipment marked with the following symbol:</p>  |
| Radiated RF<br>IEC 61000-4-3  | 3 V/m<br>80 MHz to 2.5GHz  | 3 V/m            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Breast Pump is used exceeds the applicable RF compliance level above, the Electric Breast Pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Breast Pump.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

| ISSUE       | POSSIBLE CAUSES            | CHECKS                                                                                                    | SOLUTION                                                                                                                                            |
|-------------|----------------------------|-----------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| No Power    | Power Adapter              | Check if the power cord is not connected or it is damaged.<br>Check that wall outlet is working properly. | Confirm connection of power cord; replace damaged power cord, move to a working wall outlet.                                                        |
| No Suction  | Tubing                     | Check connection to milk collection kit and to pump port(s) on the front of the unit.                     | Confirm connections.                                                                                                                                |
| Low Suction | Tubing                     | Inspect tubing for holes or cracks.                                                                       | Replace tubing if crack is in the middle of the length.<br>If cracked near tube ends, remove end and trim away damaged tube and reinstall tube end. |
|             | Assembly error             | Check for loose assembly of milk collection kit components.                                               | Disassemble and reassemble.<br>Take care to tighten the membrane cap.                                                                               |
|             | Duckbill valve             | Check for damage or wear from use.                                                                        | Replace valve if damaged or worn.<br>It is recommended to replace valves every 1-3 months depending on frequency of use.                            |
|             | Softshape™ Silicone Shield | Check for damage.<br>Softshape™ Silicone Shield is completely on the breast.                              | Replace if damaged.<br>Press more firmly against the breast.                                                                                        |
|             | Silicone membrane          | Check for damage or wear from use.                                                                        | Replace membrane if damaged or worn. It is recommended to replace membranes every 3-6 months depending on frequency of use.                         |
| Pain        | Too high suction           | Check the vacuum level.                                                                                   | Set to the lowest suction level or consult your lactation consultant or medical professional.                                                       |

# Important Symbols

|                                                                                     |                                                                                                                                |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------|
|    | This symbol indicates do not dispose the device together with unsorted municipal waste (in accordance with local regulations). |
|    | This symbol indicates compliance with international requirements for protection from electric shock. (Type BF applied parts)   |
|    | This symbol indicates that the power adapter is a class II device.                                                             |
|    | This symbol indicates that the power adapter is for indoor use only.                                                           |
|    | This symbol indicates the compliance with USA and Canada safety requirements.                                                  |
|    | This symbol indicates polarity of D.C. power connector.                                                                        |
|    | This symbol indicates alternating current. (The type from your wall outlet)                                                    |
|    | This symbol indicates direct current. (The type required to power the pump device)                                             |
|    | This symbol indicates that this device should be kept dry.                                                                     |
|    | This symbol indicates the need for the user to consult the instructions for use.                                               |
| <b>IP22</b>                                                                         | This symbol indicates the device meets IP22 classification.                                                                    |
|  | This symbol indicates when the device was manufactured.                                                                        |
|  | This symbol indicates where the device was manufactured.                                                                       |
| <b>SN</b>                                                                           | This symbol indicates the manufacturer's serial number that identifies the device.                                             |
| <b>LOT</b>                                                                          | This symbol indicates the manufacturer's batch code so that the batch or lot can be identified.                                |
| <b>REF</b>                                                                          | This symbol indicates the manufacturer's catalogue number.                                                                     |
|  | This symbol indicates the device should not be used if the package has been damaged or opened.                                 |

## Disposal of Breast Pump

The breast pump unit is made of a combination of plastic and metal. Be sure to dispose of it in accordance with local regulations as unsorted municipal waste. You may recycle it at your local PC or appliance recycling center.

## Disposal of Collection Kit Items

The hard, plastic pieces of the collection kit can be recycled per your local regulations. This includes the following parts, as outlined on page 6.

- Membrane Cap
- Pump Connector
- Bottle Vessel

The soft, silicone pieces of the collection kit cannot be recycled. This includes the following parts, as outlined on page 6.

- SoftShape™ Silicone Shields
- Membrane
- Tubing
- Duckbill Valve

## Disposal of Dr. Brown's® Options+™ Bottle and Accessories

The hard, plastic accessories and pieces of the bottle can be recycled per your local regulations. This includes the following parts, as outlined on page 6 and in the included separate bottle instructions.

- Caps
- Collars
- Reservoirs
- Storage/Travel Caps
- Narrow-to-Wide Neck Adapters

The soft, silicone pieces of the bottle and the bottle cleaning brush cannot be recycled. This includes the following parts, as outlined on page 6 and in the included separate bottle instructions.

- Nipples
- Inserts
- Bottle Cleaning Brush

# Warranty

---

## Warranty

This breast pump device is warranted against defective material or workmanship for a period of one year from the date of purchase. Any defective part will be repaired or replaced at no charge if it has not been tampered with and appliance has been used according to these printed instructions. This warranty gives you specific legal rights as well as other rights which vary from state to state. Warranty available to US and Canadian residents only.

This warranty does not cover a defect that has been caused by negligence, accidents or improper use of the appliance. To keep this appliance in good condition and retain your warranty, **USE ONLY AS DIRECTED.**

## Customer Service

If you experience unsatisfactory operation or have any questions/comments,  
call: **Handi-Craft's Customer Service** at 1-833-224-7675, or

write to: **Handi-Craft Company, Customer Service**, 4433 Fyler Avenue, St. Louis, Missouri 63116 USA.

*Outside of the US or Canada, please contact your local distributor.*

# SAVE THESE INSTRUCTIONS

Visit **[www.DrBrownsBaby.com](http://www.DrBrownsBaby.com)** for tips and help.

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©2021 Handi-Craft Company. All Rights Reserved. **F4**
DR. BROWN'S and all other trademarks are owned or used under
Guangzhou Yongyi Industrial Co., Ltd.