
S32

WEARABLE BREAST PUMP

User Manual


Visit www.nuliie.com to chat
with a live representative, or
email support@nuliie.com for
assistance and inquiries.


Contents

English	03 -10
Deutsch	11 -18
Français	19 -26
Italiano	27 -34
Español	35 -42
Nederlands	43 -50

Product Structure Diagram

Button Function and LED Display Description



- 1 Power/Mode Switch Button:** Press and hold for 2 seconds to turn the machine on or off. The product has an auto-shutdown feature after 30 minutes of continuous operation, ensuring energy efficiency and user convenience. Short press to cycle through stimulation mode, expression mode, bionic mode, and auto mode.
- 2 Pause/Start Button:** A single press initiates a pause, repress to start again.
- 3 Expression Mode:** This mode is suitable for efficient and rapid milk collection. It features 9 levels with a slow pumping frequency and powerful suction.
- 4 Stimulation Mode:** This mode effectively stimulates the mammary gland for lactation. It features 9 levels, with a rapid pumping frequency and gentle suction.

- 5 **Massage Mode:** This mode involves a gentle pulsating or rhythmic suction to simulate the initial stages of breastfeeding.. It also features 9 levels.
- 6 **Automatic Mode:** In this mode, expression and bionic modes operate simultaneously for a comprehensive pumping experience.
- 7 **Level Display:** Indicates the current suction level.
- 8 **Level-Up Button:** Increases suction level from level 1 to level 9.
- 9 **Level Down Button:** Decreases suction level.
- 10 **Battery Level:** Indicates three battery levels for easy monitoring.
- 11 **Running Time:** Indicates the current running time for user convenience.
- 12 **Unit for measuring running time.**

Battery Charging Procedure


This product comes with a built-in lithium battery. When the LED power display shows one red bar and flashes, it indicates the battery is low and needs recharging. Follow the steps below to recharge the battery:

- 1 Take the pump motor out from the milk collector to access the charging port (Refer to diagram 1).
- 2 Connect the DC Type-C port of the USB charging cable to the charging port. Then, plug the other end of the USB connector into a DC5V/1A adapter, power bank, or a socket with a USB port (Refer to diagram 2).

Note: Do not use the pump while charging. The display interface will only indicate the power level while charging.

Battery Status

LED Display	Battery Status
	Fully charged
	Medium battery level (41%-70%)
	Low on battery
& Flashing	Will power off soon
& Flashing	Charging

Assembly Procedure


Step1

Thoroughly wash your hands before assembling the cleaned components.


Step3

Insert the silicone suction bowl tightly inside the milk collector from the top. Press the edges with both hands to ensure a tight bond.


Step2

Install the valve to the bottom of the milk collector and install it in the correct direction.


Step4

Insert the silicone flange into the milk collector. Press the edges with both hands to ensure a tight bond.

Step 5

Install the pump motor onto the milk collector. The bottom part of the motor will go inside the suction bowl. Assemble it downward to fit with the milk collector and ensure it is appropriately affixed.


Usage Instructions

Before Pumping


- 1 Ensure the proper assembly of the pump by following the instructions above.
- 2 Warm the breasts for five minutes using a heated towel or a hot and cold compress to stimulate mammary gland growth.

Pumping Process


- 1 To apply suction with the breast pump, place the breast shield over your breast without any gaps and insert the nipple inside the flange. Keep your nipple upright and use your other hand to hold the breast, ensuring that no air enters. While doing this, avoid moving the device or your breasts.


- 2 Long press the power button "U" for two seconds to turn on the breast pump.

- 3 After turning on the breast pump, it will initiate the default stimulation mode. Short press the same button "U" to toggle between modes. You can choose from four modes: stimulation mode, expression mode, massage mode, and automatic mode.

- 4** Each mode features 9 levels of suction. Increase suction by pressing the "+" button and decrease it using the "-" button. Adjust the suction to your comfort level.


Warm Tips


1. This pump has a built-in memory function. When you turn on the device again, it automatically enters and shifts to the working mode and suction level you used last time.
2. The product has an auto-shutdown feature after 30 minutes of continuous operation, ensuring energy efficiency and user convenience.

After Pumping


- 1** Turn off the breast pump by long pressing the power button for two seconds.


- 2** Carefully lift the breast pump away from your breast, maintaining an upright position.

- 3** Afterwards, detach the motor from the milk collector.

- 4** To transfer the milk to the storage bottle, pour it through the outlet port. Afterwards, you can rotate the milk collector to extract any remaining milk.


Cleaning Washable Accessories

The washable accessories are the flange, silicone suction bowl, valve and milk collector.

Washing Process

- 1 Take apart all washable accessories from the breast pump, ensuring thorough disassembly.
- 2 Rinse each washable accessory under running water to remove any visible milk residue.
- 3 Submerge the accessories in warm, soapy water for 5 minutes to clean and disinfect them effectively.
- 4 Use a soft brush to scrub the components, paying attention to hard-to-reach areas gently.
- 5 Rinse all components thoroughly to eliminate any soap residue.

Disinfect

Allow the washed accessories to air dry entirely before reassembling or storing them. Alternatively, you can use a dishwasher or steam sterilization method to clean the accessories.

Tip

Sterilizing the breast pump accessories using a microwave oven or ultraviolet disinfection cabinet is not recommended, as it may reduce their durability.

Cautions


- 1 Nuliie breast pump is designed for individual use only. Sharing it with multiple users can compromise their health and safety.
- 2 Do not use the breast pump while charging.
- 3 Keep the motor away from high temperatures and ignition sources to avoid the risk of combustion and explosion.
- 4 If you experience little or no breast milk expression or pain while expressing breast milk, please consult a medical professional.
- 5 Do not submerge the motor in water for cleaning or disinfection, and avoid using steam sterilizers for disinfection.
- 6 During the usage and cleaning of the accessories, the parts undergo various temperature changes and exposure to chemical detergents. We recommend replacing the parts every 3 to 6 months, and if you notice any damage, please replace them promptly.

Troubleshooting

No.	Problem/Issue	Cause Analysis	Solution
1	Insufficient suction	Host failure	Test the suction power: Separate the main body from the cup, put the main body in the palm of your hand, and check whether there is suction in the palm of your hand. If there is suction, the machine is not the problem. Otherwise, the problem lies with the main unit.
		Improper installation	Check for installation issues: 1. To ensure proper assembly and a complete seal, follow the instructions provided to reassemble the fitting correctly. 2. Wash and dry the flange, look for damage, and see if disinfection has altered its appearance. 3. Verify whether the installation method is incorrect and whether the duckbill valve is damaged.
		Not tightly fitted to the breast	To ensure the flange and breast are properly fitted together during breastfeeding, follow the instructions for correct sitting posture provided.
2	Suction is too small	It's important to note that low suction can result in no milk being expressed during pumping. The baby's natural sucking pressure typically ranges from 10 to 45 kPa, and the international standard for breast pumps is a maximum suction pressure of 40 kPa. If this pressure is exceeded, it can damage the capillaries in a woman's breast. However, with Nuliie's breast pump, the suction pressure falls within the safe range of 10-45Kpa, which is similar to a baby's natural sucking pressure.	
		New machine	After ensuring that the new breast pump is functioning correctly, select the appropriate mode and gear. The atmosphere, your mood, and your physical condition may affect your milk flow. To increase milk production, consider using a warm compress with a towel. Generally, it doesn't matter if you utilize breast pumps from different brands.
		Suction power decreases with prolonged use	The primary causes of issues are deformation from cleaning and improperly installed accessories. Consumable components include silicone diaphragms and duckbill valves. It is advisable to change them every three to six months. You might try the following measures to see if the suction force decreases: 1. Examine the accessories for damage or significant deformation. 2. Assemble the accessories and ensure they are fully sealed following the provided instructions. 3. Ensure the silicone diaphragm is fully contracted and adjust to the breast pumping mode's highest gear.

3	Can't suck milk	Determine whether the breast is agglomerated	1. Gently massage the breast for five minutes with a heated towel or a hot and cold compress. 2. Start using the smallest "breast pumping" setting on the breast pump, gradually increasing the suction as you work to release any milk blockages.
		Determine whether the breast has any milk	If the mammary glands are inactive and not producing milk, nutrition therapy and training the body may help stimulate milk production. To promote milk production, consider consuming foods such as old hen soup, peanuts, trotters, soybeans, astragalus, angelica, soybean sprouts, tofu, enoki mushrooms, and other lactogenic foods.
4	Leaking milk	Improper sitting posture	Ensure that you are sitting in the correct position for breastfeeding. Lean forward at a 45-degree angle with your chin tucked into your chest. When using a bottle, hold it upright during and after the feeding.
		Damaged accessories, improper operation	1. Duckbill valves should be replaced immediately (the service life of the duckbill valve is about 3-6 months; contact customer service to purchase accessories). 2. Replace the silicone diaphragm immediately (the service life of the silicone diaphragm is about 3-6 months; contact customer service to purchase accessories). 3. Replace the flange immediately (yellow or deformed, needs to be replaced). 4. Use the designated milk pouring port to pour milk accurately. Avoid removing the flange while pouring milk, as this can cause milk to leak. 5. After pouring milk, a few drops of milk residue are typically visible on the linker.
5	The breast pump does not respond and cannot be turned on	The milk goes into the device	If the accessories (such as the silicone diaphragm) are damaged or if milk is accidentally poured into the main unit without removing the head, the breast pump may not function properly. Note: If any accessories are damaged, replace them promptly and be mindful of how you pour milk to prevent any issues.
		Low battery	Please use a certified 5V power adapter to charge the host. Please do not use an adapter with more than 5V, as it is easy to burn the machine.
6	Low noise	The breast pump creates suction by inhaling and deflating the pump body, resulting in a certain level of noise. However, this product has been designed to minimize noise, and the decibel level is maintained within the standard range of approximately 45 to 55 decibels.	
7	Damaged or missing components	Please contact after-sales customer service to purchase accessories.	

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FCC

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

CORRECT DISPOSAL OF THIS PRODUCT

(Waste Electrical & Electronic Equipment)

This marking shown on the product or its literature, indicate that it should not be disposed with other household wastes at the end of its working life.

To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources. Household user should contact either the retailer where they purchased this product or their local government office, for details of where and how they can take this item for environmentally safe recycling.

Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial wastes for disposal.

Manufacturer: Shenzhen TPH Technology Co., Ltd.

Address: 5th Floor, Building No.29 East Side, Lianchuang 2nd Technology Park, Bulan Road, Nanwan Community, Longgang District, Shenzhen, China

EC	REP
----	-----

UAB Qixiang Europe
Pranciškony g. 6-46, LT-01133 Vilnius, Lithuania
Technical Information Center
TGeschwister-Scholl-Str. 73, 60488 Frankfurt, Germany
Email: qixiang.europe@gmail.com
Phone: +49 6976 895074

UK	REP
----	-----

CET Product Service Ltd.
Beacon House Stokenchurch Business Park,
Ibstone Rd, Stokenchurch
High Wycombe HP14 3FE UK
Email: info.cetservice@gmail.com
Phone: +44741932 5266

Scan the QR Code
for FAQs