


Atellica® Solution

Troubleshooting Resources Virtual Training

CCE Team







How many of you have had the opportunity to Troubleshoot the Atellica Solution System?





Objectives:



At the end of this virtual Training Session you will be able to:

- 1. Identify System alerts on the Atellica Solution
- 2. Identify how to troubleshoot errors using the Operator Event Log.
- 3. Review the steps to create a Service Support Request.
- 4. Access PEP connect and Siemens Document Library.
- 5. Access Software Upgrade release notes.

Status Light





System Alerts

The Alert Module displays all active module and system alerts. The operator accesses the Alert Module using the colored tab above the Command bar. To filter alerts, select a Module or Category from the Alert Module.

The system alerts are designated as follows:

- A. Black tab with a red outline errors have been acknowledged but error conditions still exist
- Solid yellow tab there are one or more unacknowledged warnings
- C. Black tab with a yellow outline warnings have been acknowledged but warning conditions still exist
- Solid red tab there are one or more unacknowledged errors






Samples > Sample Handler



	<u> </u>	_ ^						A				
System	[UUU] Samples	Inventory	Calibration	QC A	Patient Orders	₩orklist	Maintena Maintena		ලාලි Setup			
Sample Handler	Sample List	intollisity	Cambradon	a.	r daoin ordere	, vortaist			Оскар			
SP00162: Standby Drawer: 2 Rack: Back												
003574						Racks		Sample Details				
A 1 Empty	В	C	D	E		Nacks		Sample Details				
	Empty	Empty	Empty	Empty								
2 Empty	Empty	Empty	Empty	Empty								
3 Empty	Empty	Empty	Empty	Empty								
						00000	△ ○○○○					
						00000	00000					
						Δ		∇	∇			
003664												
Α	В	C	D	E			00000	00000				
1 1111112 Intervention Needed	Empty	Empty	Empty	Empty			00000	00000				
2 Empty	Empty	Empty	Empty	Empty					00000			
3 Empty	Empty	Empty	Empty	Empty		Δ	\triangle	∇	\sqrt{00000}			
								V	V			
PEP P		abManager		品 🗘 😪				ay, November 18, 20				

Sample Details





Patient Sample Details





Worklist > Worklist Overview





Objectives:



- 1. Identify System alerts
- 2. Locate the Online Help in the Operator Event Log to investigate system errors.
- 3. Review the steps to create a Service Support Request.
- 4. Access PEP connect and Siemens Document Library.
- 5. Access Software Upgrade release notes.

System > Logs > Operator Event Log





System > Logs > Operator Event Log



Status	Logs	Operator Diagnostics	About	
Operator Event Log	System Message Log	LIS Logs	Audit Trail Log	
Filter	√ ×	Time: 24 Hours, 00180		Fi
Time 24 Hours	Sev	erity Module (SN)	Code Eve	nt
Module 00180 (IM0018	30)	△ 00180 (IM00180)		celed sample ID 83gdt, test CEA, replicate ID 1. Cannot comple order.
1		00180 (IM00180)		celed sample ID 83gdt, test CEA, replicate ID 1. Cannot comple order.
☐ ▲ Errors		00180 (IM00180)		me Check Error on reagent probe 1 at pack P04018110002608 test CEA
Unac	cknowledged	00180 (IM00180)		celed sample ID 83gdt, test CEA, replicate ID 1. Cannot comple order.
₹ 0 Corre	ective Action	00180 (IM00180)		me Check Error on reagent probe 1 at pack P04018110002608 test CEA
All Other Errors		00180 (IM00180)		celed sample ID 83gdt, test CEA, replicate ID 1. Cannot comple order.
Warnings		00180 (IM00180)		me Check Error on reagent probe 1 at pack P04018110002608 test CEA
		00180 (IM00180)	03 600 00 85 Ana	lyzer is inprocess.
\dagger 0 Unac	cknowledged	00180 (IM00180)	03 600 07 33 Rea	gent pack P04018110002608 is sufficiently mixed.
Pg 0 Corre	ective Action	△ 00180 (IM00180)		primary reagent for CEA, lot 181, is onboard in position 65, but vailable for tests .
All Other Warning	S	00180 (IM00180)	03 600 00 62 Ana	lyzer is ready.
Information		<u>↑</u> 00180 (IM00180)	03 400 09 21 Ren	noved primary reagent pack, barcode P04018110002608.
		00180 (IM00180)	03 600 00 47 Ana	lyzer is warming up.
Display Events		00180 (IM00180)		luzer is in standhu
Refresh Counters	Detail	s: Operating ModeInprocess, (CategoryDiagnostics, Code03 60	00 06 16, Volume Check Error on reagent probe 1 at pack P040
Advanced Filte	er			Count 61
				Count of

MATCH THE ICON ON THE LEFT TO THE DESCRIPTION ON THE RIGHT





ERRORS



CORRECTIVE ACTION




WARNING



UNACKNOWLEDGED

Filtering the Operator Event Log






Selecting the Event






Event Help






Event Help Continued






Following the blue hyperlink






Following the blue hyperlink





Inventory > IM Reagent Loader





Objectives:



- 1. Identify System alerts
- 2. Locate the Online Help in the Operator Event Log to investigate system errors.
- 3. Review the steps to create a Service Support Request.
- 4. Access PEP connect and Siemens Document Library.
- 5. Access Software Upgrade release notes.



Question:


How do you usually contact the Technical Solutions Center for System support?



How to create a Service Support Request

Service Support Request






























Question:

What other device will allow you to open a Service Support Request?

Tablet





Objectives:



- 1. Identify System alerts
- 2. Locate the Online Help in the Operator Event Log to investigate system errors.
- 3. Review the steps to create a Service Support Request.
- 4. Access PEP connect and Siemens Document Library.
- 5. Access Software Upgrade release notes.




Training Resources



Where can the Atellica Solution System Training Videos/ Learning Modules and Training support be found?

- a. Document Library
- b. Online Help (OLH)
- c. PEP Connect





Atellica® Solution Troubleshooting




Recover from a Slow User Interface

A slow user interface (UI) occurs when a window requires more than 5-10 seconds to load.

- For example, 200,000 patient results or 500,000 events take up to 10 seconds for the window to fully display
- This does not apply when data is loading (when features take up to 30 seconds to load) or progress window (indicates time in a progress bar)

To recover from a slow user interface:

- 1. Sign out and back into the system
 - A. On the Status bar, select
 - B. Select Sign Out > Yes
 - C. At the Sign In window, enter the Operator ID and Password
 - D. Select Sign In
 - *This procedure is continued on the next slide





Software Troubleshooting



On PEP Connect we have the PCC restart Job aid


Follow the steps to perform a successful PCC restart.



What types of documents can be found on the Siemens Document Library?

- a. Calibrator Cards
- b. Assay Instructions for Use
- c. Siemens Product Information
- d. All of the Above





Objectives:



- 1. Identify System alerts
- 2. Locate the Online Help in the Operator Event Log to investigate system errors.
- 3. Review the steps to create a Service Support Request.
- 4. Access PEP connect and Siemens Document Library.
- 5. Access Software Upgrade release notes.



Software Release Notes



Software Version 1.22 Highlights

- CH Open Channel Enhancements:
 - Smooth spline-interpolation calibration curve-fit
 - Updated high dose hook effect feature
- Streamlined QC control management with QC automated order creation by analyzer, assay, or onboard packs.

Reliability Enhancements:

- Improved Atellica Magline® Transport reliability by eliminating multiple causes of slow downs and sample overtravel
- Improved IM liquid waste management
- Resolved intermittent issues:
 - SH and IM analyzer maintenance and IMT advance clean complete more consistently
 - Full CH WBA, RPC4, and pretreatment packs are not ejected after maintenance
 - CH calibration orders complete more consistently
 - Streamlined reagent drawer alerts
 - Improved user interface response time when navigating between screens

Installing Software Updates



 Read the software release notes for installing software updates and Online Help using Smart Remote Services



Thank You!

We look forward to working with you again!