

# Chemgene MEDLAB surface disinfectant.

## PIPETTE DECONTAMINATION PROCEDURE

Recommended decontamination procedures for all models. Please note that these instructions are intended as a **GUIDE ONLY** and all pipettes requiring decontamination must be decontaminated according to your institution's own guidelines and the type of contaminant present.

### GENERAL DAILY CLEANING/DECONTAMINATION

- 1:10 dilution is generally recommended,
- Spray or wipe handle and shaft  
(insuring the tip ejector button and volumeter are also cleaned)
- Leave for a minimum contact time of 10 minutes

### ORDERING INFORMATION

					
1 Litre Self-dosing Container Clear, Unfragranced 100%	5 Litre Pump Container Clear, Unfragranced 100%	5 Litre Pump Container Blue, Eucalyptus 100%	750 ml Ready to Use Spray Clear, Unfragranced 10%	750 ml Ready to Use Spray Blue, Eucalyptus 10%	96-Wipe Resealable Flow Wrap 10%
<b>SKU074A</b> 6 x 1 Litre	<b>SKU075A</b> 2 x 5 Litres	<b>SKU077A</b> 2 x 5 Litres	<b>SKU093A</b> 6 x 750 ml	<b>SKU095A</b> 6 x 750 ml	<b>SKU101A</b> 6 x 96-Wipes

### COMPLETE DECONTAMINATION

Pipettes requiring full decontamination can be dismantled and cleaned:

- 1:10 dilution is generally recommended,
- Use a container that is deep enough to immerse the parts completely.  
(NOTE: Electronic pipettes must not be submerged in any liquids)
- Dismantle the pipette (as specified by the manufacturer) and **place only parts that can be immersed**, eg tip holder, piston, tip ejector etc completely into the solution for 30 minutes
- Remove parts and rinse with de-mineralised water
- Lie on a lint-free cloth and place in a dry oven at 60°C for at least one hour or until dry. Alternatively, place on the bench on a lint-free cloth and air dry for approximately 4 hours
- Reassemble the pipette following the manufacturer's instructions
- If complete decontamination is required prior to servicing, please ensure all parts are completely dry before reassembly and presenting for servicing

**NOTE: All pipettes for servicing must be accompanied by a completed Decontamination Declaration stating that the pipette is safe for us to work on, irrespective of whether the pipette required decontamination or not.**

### DECONTAMINATION DECLARATIONS

[starlab.click/decontamination](http://starlab.click/decontamination)