

User Manual &
Instruction Guide

Digivibe by sÜthe



Thanks for choosing Digivibe

by **süthe**[®]

All rights reserved. No part of this publication may be reproduced, distributed, or transmitted in any form or by any means, including photocopying, recording, or other electronic or mechanical methods, without the prior written permission of the publisher, except as permitted by U.S. copyright law.

Table of Contents

D-Link Symbols and Icons	4
The Diabetes Injection Kit Contents	5
Safety Information	6
Getting It Up Installation/Removal	8
Battery Installation/Removal	9
Loading Device Status	10
Removing the Used Insert	12
Cleaning your Loading Device	13
Taking a Glucose Test (Hand Held) – Patient, Parent	14
Taking a Glucose Test (Lance Free) – Patient, Parent	16
Abdominal Injections with D-Link	18
Care & Maintenance	20
Replacing & Disposal of Battery	21
Troubleshooting/FAQs	22
Limited 1-Year Warranty	23
User Profile	24
EMC Standard on Text Method	25

Digivibe Symbols and Icons



Consult User Manual



Manufacturer



Single person use only



Biohazard after use



Not waterproof



Warning




Pacemaker/Defibrillator Safe



European Conformity

THE DIGIVIBE INJECTION KIT CONTENTS



Safety Information



Before using the product carefully read the User Manual & Instruction Guide that comes with the Digivibe device



Lancets are for single use only.

Do not share the Digivibe lancing device or used lancet with anyone else, including family members.

Do not use the Digivibe lancing device on multiple patients.



After use and exposure to blood, all parts of this kit are considered bio-hazardous. A used kit may potentially transmit infectious diseases even after you have performed cleaning and disinfection.

Clean and disinfect the Digivibe device before initial use and after each use with warm soapy water and germicidal wipes.



Keep away from young children. Small parts may present a choking hazard.

Do not use in an oxygen-rich or explosive environment.

Do not expose the unit to high temperatures for an extended period of time (Example: Leaving it in a vehicle trunk on a hot day or in a sun-baked space). This may reduce the performance and/or cause a failure of the disposable lip.

For best performance and greatest comfort, if the device is stored at extreme temperatures such as -25°C to $+35^{\circ}\text{C}$, allow the device to come to room temperature or at least into the normal operational range of $+5^{\circ}\text{C}$ to -40°C before use. This should take about 2 hours when exposed to room temperature.



NOTE: Safe with pacemaker/defibrillator when used normally in accordance with proper intended use.



WARNING: Battery can explode or leak and cause burns if installed incorrectly, disassembled, or exposed to water, fire, or high temperature. No user-serviceable parts inside. Modification or disassembly not allowed.

Digivibe device setup



Before handling the Sonic Tip wash your hands thoroughly with warm soapy water. Rinse and dry hands completely.


Sonic Tip Installation / Removal



Install the Sonic tip onto the Digivibe device by aligning the dot located at the base of the Sonic tip with the triangle located on the Digivibe handle and push the tip down until you hear a click.

To remove, turn the Sonic tip $\frac{1}{4}$ turn counterclockwise and lift.

Battery Installation/Removal




Open the battery cover by pressing and lifting the latch located at the bottom rear of the device. After installing or removing the AA battery reinstall the battery cover properly.

Lancing device setup




Before handling your lancing device, wash your hands thoroughly with warm, soapy water. Rinse and dry hands completely.



1. Remove the lancing device cap by turning it counterclockwise and then pulling it straight up from the device.



2. Insert a sterile lancet into the grooves of the lancet holder. Push the lancet into the device until it is fully seated in the holder. Twist the protective cover until it separates and save it for the removal and disposal of the used lancet.



3. Place the lancing device cap back onto the device; gently turn clockwise to secure the cap.



Do not overtighten.



4. Adjust the depth setting of the lancing device by turning the wheel to the desired depth. Minimum starting from 1 and maximum ending at 5. A shallower puncture works for children and most adults.



5. Cock the lancing device trigger by sliding the device back until it clicks. If it does not click, it might have already been cocked. The device is now ready to use.

Removing the Used Lancet



1. Remove the cap by turning it counterclockwise and then pulling it straight up from the device.
2. Before removing the lancet, place the lancet protective cover on a hard surface, then push the lancet tip into the cupped side of the cover.
3. After securing the tip of the lancet with its cover, remove the lancet by gently pulling away from the lancet device.
4. Place the cap back onto the device; turn clockwise to secure the cap. Do not over-tighten; it is important to use a new lancet each time you obtain a blood sample. This will help prevent infections.
5. Discard the used lancet carefully after each use to avoid unintended lancet stick injuries. Used lancets may be considered bio-hazardous waste. Be sure to follow your healthcare professional's recommendations or local regulations for proper disposal. Wash hands thoroughly with warm soapy water after handling the Digivibe device, lancets, lancing device, and cap.

Cleaning your Lancing Device

The lancing device and cap should be disinfected after each use.

First, clean your lancing device and cap with warm soapy water prior to disinfecting.

Use a new germicidal wipe to rub or swab the outside of the lancing device until the surface is damp.

Wipe dry with clean, sterile gauze. Wash your hands thoroughly with warm soapy water after handling the Digivibe and lancing device.



Do not use alcohol or any other solvent.

Do not spray a cleaning solution on the Digivibe and lancing device.

Do not immerse the Digivibe device in any liquid, including water.

Taking a glucose test

(HAND HELD) PATIENT, PARENT

Purpose

The Digivibe device is designed to reduce the perception of pain during capillary blood sampling (fingerstick) through controlled vibration. This instruction outlines the proper use of the device prior to and during lancing for glucose monitoring.


Materials Needed

- Digivibe Device
- Lancing Device
- Glucose meter and test strip



1. Device Activation

Press the power button on the Digivibe to activate vibration.



2. Finger Positioning

Hold the device securely in one hand. Place the target finger either on top of or beneath the vibrating sonic tip, ensuring firm contact.



3. Pre-Lancing Vibration

Maintain finger contact on the vibrating tip for the recommended duration (see Vibration Chart below) prior to performing the fingerstick.

Vibration Timing Guidelines

Lancet Gauge	Vibration (Prior to Lancing)
28-35 gauge	15 seconds
23-27 gauge	20 seconds



4. Lancing procedure

While vibration is ongoing and the finger remains in place, apply the lancing device directly over the center of the vibrating tip. Ensure the lancet contacts the intended area of the finger and activate the device to obtain the blood sample.

Safety Notes

- The DigiVibe is a non-invasive Class I device.
- The device should not interfere with blood collection or glucose measurement accuracy.
- Do not use the device on broken skin or in the presence of infection.

Taking a glucose test

(HANDS FREE) PATIENT, PARENT

Purpose

The Digivibe device is designed to reduce the perception of pain during capillary blood sampling (fingerstick) through controlled vibration. This hands-free protocol uses the Digivibe stand for stability and ease of use.

Materials Needed

- Digivibe Device
- Digivibe Stand with suction base
- Flat, clean surface (e.g., table or counter)
- Lancing Device
- Glucose meter and test strip



1. Prepare the Surface

- Select a flat, clean surface to place the Digivibe stand.
- Secure the stand by pressing the suction activation button through the center of the stand using your thumb. Ensure the base is firmly suctioned to the surface.




2. Insert the Digivibe

Place the Digivibe device into the holder on the stand so that the vibrating tip is positioned horizontally and accessible for finger placement.

3. Activate the Device

Press the power button on the Digivibe to begin vibration.



4. Pre-Lancing Vibration

Place your target finger under the vibrating tip so that it makes firm, stable contact.

Maintain contact for the recommended vibration time based on needle gauge (see chart below).

Vibration Chart

Lancet Size	Vibration Time
28-30 gauge	15 seconds
23-27 gauge	20 seconds



5. Perform the Fingertick

While your finger remains in contact with the vibrating tip and vibration is active:

- Align the lancing device with the center of the vibrating tip.
- Apply the lancing device so it presses against the finger through the vibration zone.
- Activate the lancet to perform the fingertick.

Note: The Digivibe Stand is not included with the Digivibe by 50the Kit. It can be ordered separately by visiting digivibe.com

6. Post-Use

Remove your finger and wipe any residual blood.

Abdominal Injections with Digivibe


Purpose

The Digivibe device is designed to reduce the perception of pain during subcutaneous abdominal injections through the application of controlled vibration. This instruction outlines the proper use of the device to enhance comfort during self-administered or caregiver-assisted injections, such as insulin administration.



1. Prepare Injection Materials

Prepare your insulin or medication injection according to medical guidance. A 30-32 gauge needle with a length of 5-8 mm is recommended.



2. Clean the Injection Site

Use an alcohol wipe to clean the designated area of the abdomen where the injection will be administered.



3. Device Activation and Placement

Press the power button on the Digivibe to activate vibration. Place the vibrating sonic tip directly on the injection site.



4. Pre-Injection Vibration

Vibrate the injection site for up to 20 seconds prior to inserting the needle.

Vibration Chart

Needle Gauge	Vibration Time
30-32 gauge	20 seconds



5. Injection Procedure

While maintaining vibration on the injection site, insert the needle through the center area of the vibrating tip into the tissue.

6. Complete the Injection

Continue vibrating the site during the entire injection process until the needle is withdrawn.

Safety Notes

- The Digivibe is a non-invasive Class I device.
- The device should not interfere with medication absorption or injection accuracy.
- Do not use the device on broken skin or in the presence of infection.
- Follow standard hygiene protocols for all injection procedures.

Care & Maintenance

Cleaning your Digivibe device

Clean your Digivibe device by wiping it with a damp cloth and then with a germicidal wipe.

 DO NOT wet the Digivibe device.

Use germicidal disinfecting wipes. Clean the Digivibe sonic tip by first removing it from the device and then cleaning it with a mild liquid detergent and rinsing with warm water.

Technical Information

Permissible environmental conditions during normal use.

Temperature

Operating:	+5°C to +40°C	Storage:	-25°C to +55°C
------------	---------------	----------	----------------

Replacing & Disposal of Battery

Your Digivibe device uses one 1.5v AA lithium battery or one AA alkaline battery. A lithium battery will operate the Digivibe device for approximately 450 uses, assuming each use requires 15 seconds. If the device vibration is weak, replace the battery.



IMPORTANT: The recommended battery to use with this device is a 1.5v AA lithium battery, which will give the greatest performance and life of operation before needing to replace the battery.

Replacing the Battery

Start with the device turned off. Remove the battery cover by sliding it downward and pressing the release latch. Pull up on the battery to lift it out of the compartment. Remove the old battery.

Insert a new AA lithium battery or one AA alkaline battery into the battery compartment. Plus (+) and minus (-) signs will guide you in placing the battery. The battery needs to be the Energizer 1.5v AA LR1 or AA alkaline battery.

Disposal of Battery and Device

Dispose of battery and device according to your local environmental regulations.

Troubleshooting / FAQs

Digivibe device stand does not stick to the flat surface.

Detach the stand from the Digivibe device and wash the stand with a mild dishwashing detergent, then dry with a paper towel.

Digivibe device vibration is weak.

Replace the battery with a new lithium battery if the pain reduction effect is not felt. Do not use a rechargeable battery. We recommend a lithium battery. Alkaline batteries are not powerful enough to run the device at full strength.

What if my Digivibe device does not turn on?

Ensure that the battery is installed correctly.

What if my Digivibe device tip is cracked or damaged?

If your Digivibe device tip is cracked or damaged, contact Customer Service at 1-844-406-8425 or order a new Digivibe tip at digivibe.com.

Limited 1-Year Warranty

Ring Innovations, LLC warrants the Digivibe device and its contents to be free of defects in material and workmanship for a period of one year from the date of original purchase. If Digivibe exhibits such a defect, Ring Innovations, LLC will, at its option, repair or replace it without cost for parts or labor. The product must be shipped, prepaid, and insured (recommended) to Ring Innovations, LLC. Proof of date of original purchase is required.

The warranty does not cover damages resulting from accidents, misuse, water, tampering, unreasonable use, failure to provide reasonable and necessary maintenance, services performed or attempted by an unauthorized service agency, use of any tip other than those supplied by Ring Innovations, LLC.

The warranty extends only to the original purchaser and is not transferable.

User Profile

The DigiVibe device by eûthe is the latest product innovation from Bing Innovations, LLC. The DigiVibe device is designed to alleviate needle pain during injections and finger pricks or both. Its user profile varies from a wide range of users who regularly administer or receives injections. The DigiVibe device employs vibration technology making needle related procedures more comfortable. The DigiVibe device by eûthe is appropriate for use with children younger than 12 years of age but must be operated by an adult, for their safety.

IP rating: This device has an IP21 rating, which means that it is resistant to water being sprayed on it from the top only.



Do not service the DigiVibe device while equipment is in use.

Radio Frequency Emissions

Class B Group 1

Suitable for use in all environments, including those directly connected to a public-line-voltage power supply network

ELECTROSTATIC
DISCHARGE

IEC 61000-4-2

±8 kV contact
±2 kV, ±4 kV, ±6 kV, ±15 kV air

Immunity RF EM fields

IEC 61000-4-3
Table B Immunity Fields
from RF wireless
communications Equipment

10 V/m, 80 MHz ~ 2.7 GHz
80 % AM at 1 kHz

(SAR) 10 or 100g

Power Frequency
Magnetic Field

EN 61000-4-8

Parameter	IEC 61000-4 test method	Standard test limits	
		Professional healthcare facility (professional)	Other healthcare environment
Electromagnetic Immunity	IEC 61000-4-3	±1 kV contact ±2 kV, ±4 kV, ±6 kV, ±15 kV air	
Immunity to 50 Hz fields	IEC 61000-4-6	3 µT max 30 mG ~ ±1 kV/m ² 80 % AM at 1 kHz ¹	10 µT max 30 mG ~ ±1 kV/m ² 80 % AM at 1 kHz ¹
Immunity to 60 Hz AC power	IEC 61000-4-7	See 4.10	
Immunity to 50/60 Hz AC power	IEC 61000-4-8	See 4.10	
Immunity to 50/60 Hz AC power	IEC 61000-4-9	See 4.10	

¹ The difference between the maximum physiological magnetization of heart and the AC magnetization of any person shall be ignored within 0.1 m of the vertical plane of the conductive test area of each environment of the environment in use during.

² The maximum and minimum values that represent the maximum RF electromagnetic energy for the purpose of field exposure shall not be used for frequency of exposure. Testing shall be performed at other maximum frequencies identified by the test manufacturer. This test applies to the power supply and electrical equipment, and is not intended to be used as a general test for the equipment. It is recommended that the power supply be tested separately during the test.

³ Testing may be performed at other maximum frequencies identified by the test manufacturer.

⁴ Applies only to the maximum value of current with frequency between 50 Hz and 100 Hz.

or 100g

⁵ Before evaluation is applied.

Phenomenon	IEEE 61850-4 standard	Harmonics test levels	
		Professional hardware (factory environment)	Wide area (remote) environment
Electrical field (sensors) / cables	IEEE 61850-4-4	$\pm 1 \text{ kV}$ 100 kHz, repetition frequency	
Surge - 11 Line-to-line	IEEE 61850-4-5	$\pm 0.5 \text{ kV}$, $\pm 1 \text{ kV}$	
Surge - 16 Line-to-ground	IEEE 61850-4-5	$\pm 0.5 \text{ kV}$, $\pm 1 \text{ kV}$, $\pm 2 \text{ kV}$	
Conducted disturbances induced by RF fields 11-16	IEEE 61850-4-6	± 0 $0.15 \text{ MHz} - 50 \text{ MHz}$ $0.1 - 10 \text{ MHz}$ bands between 0.15 MHz and 50 MHz 50 Hz band at 1 kHz	± 0 $0.15 \text{ MHz} - 50 \text{ MHz}$ $0.1 - 10 \text{ MHz}$ and spectral radio bands between 0.15 MHz and 50 MHz 50 Hz band at 1 kHz
Voltage dips 11	IEEE 61850-4-11	0% to 0.1 cycle 40% , 45% , 50% , 100% , 150% , 200% , 250% and 300% 0% to 1 cycle and 10% to 2000 cycles Single phase at 0°	
Voltage interruptions 11	IEEE 61850-4-11	0% to 200000 cycles	

Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS	
		Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	
Radiated RF EM fields ⁽¹⁾	IEC 61000-4-3	3 V/m ⁽¹⁾ $50 \text{ MHz} - 2.7 \text{ GHz}$ ⁽¹⁾ $80 \% \text{ AM at } 1 \text{ kHz}$ ⁽¹⁾	10 V/m ⁽¹⁾ $50 \text{ MHz} - 2.7 \text{ GHz}$ ⁽¹⁾ $80 \% \text{ AM at } 1 \text{ kHz}$ ⁽¹⁾
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See 8.10.	
Rated power frequency magnetic fields ⁽¹⁾	IEC 61000-4-8	50 A/m $50 \text{ Hz or } 60 \text{ Hz}$	
Proximity magnetic fields	IEC 61000-4-38	See 8.11.	



Ring Innovations, LLC 1000 NW 6th Ave., Boca Raton FL 33437 1-844-303-8423 www.sütheinstruments.com

Digivibe by süthe



©2024 Ring Innovations, LLC. All rights reserved. Under the copyright laws, this manual may not be copied, in whole or in part, without the express written consent of Ring Innovations, LLC. Every effort has been made to ensure that the information in this manual is accurate. This product is covered by patents.

FDA Registered. Reverse engineering or disassembly is prohibited.

Digivibe Device:
Model Name: Digivibe by Süthe
Model Number: Dig001

Lancing Device:
Model Name: Lancing Device
Model Number: LDF

Digivibe Kit:
Model Name: Digivibe Device and Lancing Device Kit
Model Number: DigKit 001

Version: 20250101 Part: WDD-01 Rev 2

Manufactured in Mexico