

Stylish, Smart + **Simplified**



The F&P ICON+ **with added simplicity**

Featuring:

- + Simplified user experience including Simple Mode for patients
- + Stylish new color
- + Smart clinical technologies

We listened, we heard and we responded with an easy-to-access Clinician menu and a simple patient mode. The new patient mode offers easy access to the three most commonly used settings (on/off, ramp, humidity). These enhancements minimize your patient set-up and product instruction time, thus promoting a positive experience for you and your patients.



F&P ICON+




More humidity. More comfort.

ThermoSmart Technology is at the forefront of Fisher & Paykel Healthcare's humidification philosophy, and is compatible with each model of the F&P ICON+ family.

This Technology includes an advanced humidification algorithm and a heated breathing tube, which both work to provide optimal humidity and reduce condensation.

Heated Breathing Tube Humidification



Warm tube: Heat loss + heat gain from heated breathing tube = constant heat + reduced rain-out



Responsive pressure relief.

SensAwake promotes better overall sleep. We all experience subconscious waking through the night – at which time pressure intolerance is likely to occur. Whereas other pressure relief technologies provide partial relief during expiration, SensAwake provides a prompt and significant relief in pressure to the lowest most comfortable level upon waking. This eases the return to sleep and allows effective treatment to resume.

Simulates a technician's responsiveness to pressure intolerance:



Detects wakefulness

Promptly **relieves** pressure

Eases the **return** to sleep



More communication. More results.

InfoSmart Web is a flexible management system which delivers multiple data transfer options that helps to efficiently manage and report on patient adherence in the most cost-effective way.



F&P InfoSmart™ Web

F&P InfoGSM™

Data download
via GSM network




F&P InfoUSB™ *Via the Web*

Data download
via the Internet



F&P InfoUSB™ *Via the Mail*

Data download
via InfoUSB



F&P InfoView™

Data download via the display
screen of the F&P ICON+



Dashboard

InfoSmart Web gives the healthcare provider an instant overview of patient adherence on a real-time dashboard. Compliant and non-compliant patients are color-coded to indicate where early intervention is necessary. A single click of the mouse on a patient's report link can produce a compliance report.

Customizable Automated Patient Reminders

InfoSmart Web offers customized patient reminders which can be sent to patients at frequent intervals via telephone, text or email to remind them to upload their data.

PATIENTS CONTACTS DEVICES REPORTS							
Dashboard Events Reminders							
Current Filter: All regions		Update Filter		A B C D E F G H I J K L M N O P Q R S T U V W X Y Z All			
Patient Name	Device Serial Number	Initial Setup Date	Last Download	Compliant	Compliant Rate	Report	?
Kline, Mary	120801990038	05/31/2012	07/10/2012	Compliant	93% (30 Days)	Hide Report	What web browsers are compatible with F&P InfoSmart Web?
Theman, Carlos	110317087270	08/08/2012	08/17/2012	Non-Compliant	20% (10 Days)	Hide Report	How do I use the Dashboard?
Maxima, Edward	101213038887	02/09/2011	08/14/2012	Non-Compliant	33% (30 Days)	Hide Report	What should my data look like?
Hemsoth, Bob	120801990038	05/31/2012	08/20/2012	Tracking OK	90% (21 Days)	Hide Report	
Dem, Mark			No Data Available	N/A		Hide	

Product Specifications

Dimensions	6.3 x 6.7 x 8.7 in (160 x 170 x 220 mm)
Weight	4.8 lb (2.2 kg) Packed Weight: 8.7 lb (4.0 kg)
Performance	Pressure Range: 4 to 20 cmH ₂ O
Noise Level	< 30 dBA
Water Chamber Volume	420 ml up to the maximum fill line
Standards Compliance	IEC60601-1: 1998, AS/NZS 3200.1.0: 1998, EN60601-1: 1990 + A1 - A2
InfoUSB Data Recording	Auto and Premo The InfoUSB records up to 365 days of summary compliance and efficacy data, 7 days of detailed efficacy data (AHI, Pressure, Leak), as well as cumulative data from first use. Novo The InfoUSB records up to 365 days of compliance summary data and cumulative data from first use.

Electrical Ratings

Rated supply voltage	Rated current input	Rated supply frequency
100-115 V	1.27 A (1.43 A Max)	50-60 Hz/400 Hz
220-240 V	1.07 A (1.21 A Max)	50-60 Hz

F&P ICON+ Performance Features

Description	Auto	Premo	Novo
Fully integrated	•	•	•
ThermoSmart Technology*	•	•	•
Auto-Adjusting Pressure	•		
Efficacy Reporting	•	•	
Compliance Reporting	•	•	•
InfoUSB 2	•	•	•
SensAwake	•	•	
Proportional Ramp	•	•	•
Auto-Altitude Adjusting	Automatic	Automatic	Manual
Leak Compensation	•	•	
Clock and AlarmTones™	•	•	•
InfoSmart Software	•	•	•

F&P ICON+ Models

Plug type	Auto	Premo	Novo
North America (model numbers with HT include heated breathing tube)	ICONAAN/ICONAAN-HT	ICONPBN/ICONPBN-HT	ICONNAN/ICONNAN-HT

Included with the F&P ICON+

- 1 x CPAP
- 1 x Carry-bag
- 1 x Water Chamber
- 1 x Chamber Lid
- 1 x Breathing Tube*
- 1 x InfoUSB 2
- 1 x InfoUSB Mailer
- 1 x User Instructions
- 1 x Spare Air Filter

Replacement Parts

- 900ICON200 Water Chamber
- 900ICON202 InfoUSB 2
- 900ICON203 InfoUSB 2**
- 900ICON204 Elbow
- 900ICON205 Oxygen Elbow
- 900ICON206 Outlet Seal
- 900ICON208 ThermoSmart Breathing Tube
- 900HC221 Breathing Tube
- 900ICON503 Air Filter
- 900ICON218 Filter Grill
- 900ICON214 Silver Gloss Lid
- 900ICON216 Matt Gray Lid
- 900ICON217 Charcoal Lid
- 900ICON315 Soft Carry-case
- 185046221 User Instructions (English, French, Spanish)
- 185048028 User Instructions (English, Chinese, Traditional Chinese)
- 185047918 User Instructions (English, Japanese, Traditional Chinese)
- 185045354 User Instructions (Multi-Language)
- 900ICON110 Healthcare Provider CD
- 185047919 Healthcare Provider Manual (English, Japanese, Traditional Chinese)
- 900SW100 F&P InfoSmart CD
- 900ICON112 F&P Studio™ CD

Product Information

-  NZ Radio Interference C-tick mark
-  Latex free
-  Drip-proof
-  UL Classified
-  Class II Medical Electrical Equipment
-  Conforms with medical device directive 93/42/EEC
-  Device Disposal Instructions:
This device contains electronics. Please do not discard as regular rubbish. Dispose of electronics according to local guidelines.
-  Consumable Disposal Instructions:
Place the mask, Breathing Tube and Water Chamber in a rubbish bag at the end of use and discard with normal rubbish.

* In some countries the ThermoSmart Breathing Tube needs to be purchased as an accessory to activate ThermoSmart Technology.
** Part may vary depending on country.