

BARD

Site~Rite 5 Ultrasound System Manufactured by:

Bard Access Systems, Inc. Salt Lake City, UT 84116 U.S.A.

(801) 595-0700

Customer Service: (800) 545-0890 Technical/Clinical Support: (800) 443-3385

Fax: (801) 595-4948 www.bardaccess.com

An issued or revision date for these instructions is included for the users information. In the event two years have elapsed between this date and product use, the user should contact Bard Access Systems, Inc. to see if additional product information is available.

Revision date: May, 2006.

^{*} Bard and Site~Rite are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.







BAIRD



ULTRASOUND SYSTEM

TECHNICAL MANUAL





Do not operate in the presence of flammable anesthetics



Dot Markers Active





Warning: Refer to Manual Before Use



Power/Stand-by



BF Type Equipment



Medical Electrical Equipment Classified by ETL with respect to Electric Shock, Fire, and Mechanical Hazards only in accordance with UL60601-1 and CAN/CSA C22.2 No. 601.1



Do Not Dispose of Battery Pack In Fire



AC Adapter



Humidity Parameters



Class II Equipment






Table of Contents

Section

- 1 Warranty
- 2 Warnings, Precautions and Notes
- 3 Site~Rite* 5 Ultrasound System Description
 - 3.1 Power Source
 - 3.2 Scanner Assembly
 - 3.3 Transducer Assembly
 - 3.4 Front Panel Controls
 - 3.5 Display Screen Information
- 4 Installing Site~Rite* 5 Software
- 5 Calibrating the Internal Battery
- 6 Troubleshooting & Error Screens
- 7 Site~Rite* 5 System Specifications
 - 7.1 Operating and Storage Conditions
 - 7.2 Scanner Specifications
 - 7.3 Probe Specifications
 - 7.4 Power Supply Specifications
- 8 Standards Information
- 9 Disposal Information

1 Warranty

The manufacturer, Bard Access Systems, warrants this product against defects in material and work-manship for a period of one year from the date of original purchase, and agrees to repair, or at Bard Access Systems' discretion, replace any defective unit free of charge. The warranty on the repaired or replaced unit continues from the purchase date of the original unit. This warranty does not cover damages resulting from misuse, abuse, modification, or alteration of the Site~Rite* 5 System.

The following actions void the warranty of the Site~Rite* 5 System:

- · Opening or servicing the scanner or the probe housing.
- Removal of system labels by anyone other than by Bard Access Systems authorized service personnel.
- Opening or servicing the battery pack or the combination A/C adapter and battery charger by anyone other than Bard Access Systems authorized service personnel.
- Connecting the Site~Rite* 5 scanner to any power source other than the Site~Rite* 5 combination A/C adapter and battery pack.
- Connecting the Site~Rite* 5 scanner to any A/C adapter other than the one provided with the scanner.
- Connecting the Site~Rite* 5 System to any unauthorized accessory.



2 Warnings, Precautions and Notes

Warnings

Warning: This product should only be operated by qualified medical personnel.

Warning: Do not remove outer protective covers from the Site~Rite* 5 scanner. Hazardous

voltages exist at several points within the system.

Warning: Do not operate the Site~Rite* 5 Ultrasound System or the Site~Rite A/C adapter and

battery charger in the presence of flammable anesthetics or gases. Explosion

may result.

Warning: Do not use for ophthalmic indications. Ophthalmic use may cause patient injury.

Warning: Misuse of the Site~Rite* 5 System may result in damage to the equipment or

personal injury.

Warning: Use only the combination Site~Rite A/C adapter and battery charger to charge

Site~Rite* 5. Use of any other device to charge Site~Rite* 5 battery packs may

damage the battery packs and will void your warranty.

Warning: Only connect a Site~Rite* 5 combination A/C adapter and battery charger to the

Site~Rite* 5 System. Use of any other A/C adapter may cause intermittent or unpredictable operation, may damage the system and will void your warranty.

Warning: If a probe is damaged in any way, discontinue use immediately. Damage to the

scanner may occur.

Warning: Avoid subjecting the probe to excessive mechanical shock. Damage to the probe may

occur.

Warning: Use only Bard Access Systems probes with this system. Use of unapproved probes

may result in patient injury or equipment damage.

Warning: When using Site~Rite Needle Guides on the Site~Rite probe, use only sterile plastic

probe covers that are 1 mil (0.001 inch or 0.0254 mm) thick.

Warning: Do not allow liquid to enter the scanner, combination A/C adapter and battery

charger, probe connector or probe port. Damage to equipment may occur.

Warnings (continued)

Warning: Do not attempt to sterilize the Site~Rite* 5 scanner or probes with ethylene oxide or

heat sterilization methods. Damage to the equipment may occur.

Warning: Always properly dispose of dead battery packs in accordance with local regulations.

Improper disposal may present an environmental hazard.

Warning: Only qualified personnel should attempt to service this equipment. The Site~Rite* 5 contains static sensitive components and circuits. Failure to observe proper static

control procedures may result in damage to the system.

Warning: The following actions void the warranty of the Site~Rite* 5 System.

· Opening or servicing the scanner or the probe housing.

 Removal of system labels by anyone other than by Bard Access Systems authorized service personnel.

 Opening or servicing the battery pack or the combination A/C adapter and battery charger by anyone other than Bard Access Systems authorized service personnel.

Connecting the Site~Rite* 5 scanner to any power source other than the Site~Rite* 5 combination A/C adapter and battery pack.

 Connecting the Site~Rite* 5 scanner to any A/C adapter other than the one provided with the scanner.

Connecting the Site~Rite* 5 to any unauthorized accessory.

Warning: Inspect A/C adapter and battery cord for damage. If any of the prongs are damaged,

use battery power until replacement cord is obtained.

Warning: Verify that all accessories attached to the system comply to 60601 safety standards. Non-compliance may result in increased patient risk.

Warning: Use only IEC or ISO approved safety devices outside the patient environment. Failure to do so may damage the equipment.

Warning: Equipment that relies on basic insulation only shall not be used with this system. Failure to comply could result in increased patient risk.

Warning: Maximum shelf load on the VAD bedside roll stand is 22 lbs. Exceeding this weight may

damage the roll stand.



Warnings (continued)

Warning: Do not overtighten screws when attaching to the VESA roll stand mount. Doing so

may damage the scanner.

Warning: Use only screws provided in packaging. Ensure the unit is secure against the VESA

roll stand mount. Failure to do so may cause the scanner to disconnect from the VESA

roll stand mount.

Warning: Do not use the probe with high frequency surgical equipment. Doing so may damage the

equipment.

Warning: Do not pull on probe cable. Doing so may cause the system to tip.

Precautions

Caution: The adverse biological effects of ultrasound on tissue appear to be threshold effects.

When tissue is repeatedly exposed to ultrasound, with intervals in between, there will likely be no cumulative biological effect. If, however, a certain threshold has been passed biological effects may occur. While the Site~Rite* 5 acoustic output parameters fall well below all FDA thresholds for adverse biological effects, any given Ultrasound Procedure should be performed using the principle of ALARA (As Low As Reasonably Achievable). The licensed medical practitioner should limit the time of patient exposure to ultrasonic

radiation using the principle of ALARA.

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Caution: Do not pull the cable to disconnect the probe connector from the scanner. Pulling the

cable may damage the cable, cable connection or scanner.

Caution: Do not twist or bend the probe cable in excess of that required during normal use of

the probe. Excessive twisting or bending of the cable may cause failure, intermittent or

unpredictable operation.

Caution: When disinfecting the probe with a liquid disinfectant, do not soak the probe cable,

cable bend relief, probe connector or probe buttons. Doing so may damage the probe.

Caution: Only apply commercially available ultrasonic couplant, which has been specifically for-

mulated for use in medical applications, to the acoustic window (or face) of the probe.

Precautions (continued)

Caution: Use water or rubbing alcohol and a soft cloth to remove couplant from the acoustic

window (or face) of the probe. Failure to do so may scratch the acoustic window.

Caution: Do not to allow ultrasonic couplant to dry on the acoustic window (or face) of the

probe. If the couplant should dry, use water or rubbing alcohol and a soft cloth to remove it. Never use a tool of any kind to remove dry couplant from the acoustic

window (or face) of the probe.

Caution: Some commercially available probe covers contain latex. Natural rubber latex may

cause allergic reactions. Refer to the US FDA alert titled: "Medical Alert: Allergic

Reactions to Latex-Containing Medical Devices", issued March 29, 1991.

Bard Access Systems distributes sterile probe covers and needle quide kits that do

not contain latex

Caution: Do not force the probe connector. Damage to the connector and system could result.

Caution: Always snap the needle guides on to the probe hook. Do not slide the needle guide

on to the needle guide hook, as the sterile sheath may tear.

Caution: Do not subject the probe to excessive vibration. Vibration may dislodge sensitive

components and cause intermittent or unpredictable operation.

Caution: Prior to each use please inspect the integrity of all power cords and connectors as

well as the integrity of the unit itself. If any problems are found please discontinue use immediately and contact an authorized service representative. Use of a damaged

power cord could damage the machine.

Caution: Unapproved extension cords should not be used with this system. Doing so may

damage the system.

Caution: During use, the AC connector needs to be easily accessible. In case of emergency

remove the power cord as soon as possible.

Caution: To avoid unnecessary strain on the user, use the device in a comfortable manner.

Caution: Attach power source in such a way as to prevent damage. Improper installation may

damage power cords.

Caution: Inspect the probe prior to each use. If damage to the cable or transducer face is noted,

do not use the probe. Damage to the system may occur.



Precautions (continued)

Caution: Hot water (in excess of 113° F or 45° C) may damage the probe.

Caution: Use only Bard Access Systems cleaning and disinfection procedures. Failure to do so

may damage the device.

Notes

Note: When cleaning the system and components, it is important to remove all particles or

other matter from all surfaces and crevices.

Note: For 240 V applications use only center tapped 240 VAC single phase power.

3 Site~Rite* 5 Ultrasound System Description

The Site~Rite* 5 Ultrasound system contains three major components. (See Figure 1)

- Power Source (External AC Adapter or Base Unit Battery System/AC Adapter)
- Scanner Assembly
- · Transducer Assembly




Figure 1. Site~Rite* 5 System Block Diagram



3.1 Power Source

The power source supplies 15 volts direct current (VDC) nominal to the Scanner Assembly. Power source options include an External AC Adapter or a Base Unit Battery Assembly. The Base Unit Battery assembly contains a battery and built-in AC Adapter that permits portable system operation. Both power supply options provide power to the scanner and charge the internal battery simultaneously.





Figure 2. Site~Rite* 5 Power Architecture

3.2 Scanner Assembly

The Scanner Assembly transmits and receives electrical ultrasound data that is used to form a coherent image. A high-level block diagram of the system is shown in **Figure 2**.




NOTE: The Scanner Assembly is NOT field serviceable. If a problem is found with the Scanner Assembly please contact Bard Access Systems at (800) 443-3385.

3.3 Transducer Assembly


The Transducer Assembly converts electrical signals to acoustic pulses that are transmitted into the body. The pulses reflect off the target and return to the transducer, which converts the acoustic data back into electrical data. The electrical data is then transmitted to the scanner assembly where an ultrasonic image is generated from multiple pulse echo frames.

NOTE: The transducer assembly is NOT field-serviceable. If a problem is found with the transducer please contact Bard Access Systems at (800) 443-3385.

3.4 Front Panel Controls



3.5 Display Screen Information



Mechanical Index

The mechanical index is shown in increments of 0.1 with an accuracy of +/- 0.09.

··· Probe Orientation Marker

Time and Date

The time is shown in 24 hour "00:00" format and the date may be displayed in MM/DD/YY or DD.MM.YY format.

··· Image Depth Markers

When enabled, the depth markers are shown as green dots superimposed on the image at 1.0 cm pitch with smaller dash marks at 0.5 cm



Probe Depth Indicator Indicates current image depth.



Depth Marker Indicator

When depth markers are enabled, the indicator is highlighted in light green.



Gain Indicator Shows gain levels from 0-100%.



Run/Freeze Indicator

When the system is in freeze mode, the indicator is shown in red and blinks briefly every second.



Storage Device Indicator

When a storage device is connected, the indicator icon is shown in solid grayscale. When the storage device is full or has errors, the indicator is displayed with a red "!" overlay.





Printer Indicator

When a printer is connected, the indicator is shown in solid grayscale. If the printer has an error condition, the indicator is displayed with a red "!" overlay.



Display Screen Information (continued)



AC Power Indicator

The indicator is shown in solid grayscale when connected to AC wall power and is shown as a dark gray outline when disconnected.



Battery Indicator

The 5 charge levels for the batteries are displayed as shown. When the unit is recharging, the indicator is green for all 5 levels. When the unit is under battery power, the console battery indicator changes to yellow at 40% and red/blinking at 20% of remaining life. If a battery malfunction is detected, the system will show the battery with a red "!" overlay.

4 Installing Site~Rite* 5 Software

The Site~Rite* 5 Ultrasound System allows software installation through the USB connectors located on the rear of the scanner.

To install software:

- Press the Freeze button (► |) to pause the system.
- 2. Insert a USB drive into one of the USB connectors located on the rear of the Site~Rite 5 scanner.
- Wait until the storage device icon is illuminated before proceeding.
- Simultaneously press and hold the Gain buttons (▲▼) until the configuration screen appears.
- When the configuration screen appears, simultaneously press and hold the Depth button (cm) and Depth Marker button (/:\) until the software installation process begins.

Note: The Site~Rite 5 System will automatically upload the software from the attached USB drive.

Note: The screen will appear blank and/or inactive during the software installation process.

- 6. When a blinking cursor appears, power off the Site~Rite 5 scanner.
- 7. Disconnect the USB drive.
- 8. Power on the Site~Rite 5 scanner.
- 9. Verify that the correct software version appears on the upper left hand corner of the screen.
- Resume normal operation.

5 Calibrating the Internal Battery

The Site~Rite 5 internal battery may occasionally require calibration to ensure the battery power meter is accurate. The following icon indicates that the Site~Rite 5 battery requires calibration.

To calibrate the Site~Rite 5 Internal Battery:

- Disconnect all Site~Rite 5 power sources.
- Power on the Site~Rite 5 scanner and operate on internal battery power until the system powers off.
- Connect the Site~Rite 5 scanner to the Site~Rite 5 combination A/C adapter and battery charger to recharge the internal battery.

Note: At least six hours of charge time is recommended to fully charge the Site~Rite 5 internal battery.

- 4. Disconnect all Site~Rite 5 power sources.
- Power on the Site~Rite 5 Ultrasound System and operate on internal battery power until the system powers off.

Note: The internal battery is now calibrated.

Connect the Site~Rite 5 scanner to the Site~Rite 5 combination A/C adapter and battery charger to recharge the internal battery and continue normal use.





Troubleshooting & Error Screens

Solution:

Solution:

Missina/ **Invalid Probe**



Cause: Scanner does not recognize or identify a probe or probe not attached. Solution: Ensure that a Site-Rite 5 probe is properly connected to the system.

System Malfunction



Cause: Scanner is not operating within normal parameters.

Discontinue use immediately. Return to authorized repair facility.

Display Malfunction



Display malfunction.

Most display malfunctions can easily be corrected by resetting the system. To do so, power off the device, wait 60 seconds, then power the system back on. If the display malfunction is not resolved by resetting the system, discontinue use and

return to authorized repair facility.

Battery **Empty**



Cause: Battery empty.

Connect system to AC outlet for operation and battery recharge.

Battery Malfunction



Cause: Battery malfunction.

Solution: Send system to authorized repair facility for battery replacement.

Storage Device Indicator



Cause:

When a storage device is connected, the icon is shown in solid grayscale. If the storage device is full or has errors, the icon is displayed with a red "!" overlay.

Solution: Replace storage device.

Printer Indicator



Cause: When a printer is connected, the icon is shown in solid grayscale. If the printer has

an error condition, the icon is displayed with a red "!" overlay.

Solution: Check paper or service printer.

6 Troubleshooting & Error Screens Continued...

Troubleshooting Guide

Problem	Possible Cause	Solution	
The system will not power on. (Power LED is not illuminated, fan is not running, and no information is displayed.)	Did not press power button long enough	Press the power button for at least 1 second.	
	Battery is discharged	Plug the unit into a known working AC outlet to operate system and charge the battery.	
	AC Adapter is defective	Plug the unit into a known working AC outlet. Verify power LED on Scanner Assembly is blinking. If not, remove AC Adapter and measure voltage levels on the DIN connector. The DIN Shield is GND. Pin 1 & 3 should read ~15VDC.	
Gain up and gain down buttons on scanner or probe do not work	System is in freeze mode	Press the freeze button.	
	Button pad is defective	Use scanner buttons until unit can be returned for service.	
Depth buttons on scanner or	System is in freeze mode	Press the freeze button.	
probe do not work	Button pad is defective	Use scanner buttons until unit can be returned for service.	
Image is poor or penetration is inadequate	Not enough acoustic coupling gel is being used.	Use more coupling gel.	
	Fold or seam of sheath on the acoustic window of the probe.	Smooth out sheath over the acoustic window.	
Scanner Assembly powers on momentarily then shuts off	Battery is discharged	Plug the unit into a known working AC outlet to operate system and charge the battery.	



7 Site~Rite 5* System Specifications

7.1 Operating and Storage Conditions

Operating Temperature: 59°F to 100°F (15°C to 38°C) Storage Temperature: 50°F to 95°F (10°C to 35°C)

Operating Humidity: 5% to 85% Relative Humidity (non-condensing) Storage Humidity (packaged): 5% to 95% Relative Humidity (non-condensing) Storage Humidity (unpackaged): 5% to 85% Relative Humidity (non-condensing)

7.2 Scanner Specifications

Dimensions: 12" W x 12" H x 5" D

Weight: 10 lbs.

Power Sources: AC adapter, Internal and External DC Battery Pack

Power Consumption: 82 Watts Maximum Monitor Size: 82.1" diagonal

IEC 60601- 1: Class II, Type BF Applied Part, Continuous Operation, Internally

Powered Equipment, Not Category AP or APG Equipment,

Not protected against ingress of water.

7.3 Probe Specifications

L-VA: Linear Vascular Access Probe

 Frequency:
 5 -10 MHz

 Elevation Focus:
 1.8 cm

 Maximum Scan Depth:
 6.0 cm


 Scan Width:
 1.9 cm

Lateral Foci:

Image Depth	Focal Depth
1.5 cm	0.6 cm
3.0 cm	1.5 cm
4.5 cm	3.0 cm
6.0 cm	5.0 cm

7.4 Power Supply Specifications

Note: For 240 V applications use only center tapped 240 VAC single phase power as shown below.



A/C Adapter Specifications

Input Voltage:	100-240 VAC, 50/60 Hz.	Output Voltage:	15 VDC
Input Current (Max):	1 Amp	Output Current (Max):	6 Amps

Internal Battery Pack Specifications

Battery Chemistry:	Lithium Ion	Output Power (Full Charge):	52 Wh
Nominal Output Voltage:	10.8 VDC	System Run Time on Full Charge:	1.0 Hours
Output Current (Max):	6 Amps	Charge Time (Full):	1.75 Hours

Combination A/C Adapter Auxiliary Battery Specifications

Input Voltage:	100-240 VAC, 50/60 Hz.	Battery Output Current (Max):	6 Amps
Input Current (Max):	1 Amp	Battery Chemistry:	Lithium Ion
A/C Adapter Output Voltage:	15 VDC	Output Power (Full Charge):	95 Wh
A/C Adapter Output Current (Max):	7 Amps	System Run Time on Full Charge:	2.5 Hours
Nominal Battery Output Voltage:	10.8 VDC	Battery Charge Time (Full):	3 Hours



8 Standards Information

The Site~Rite* 5 Ultrasound System is designed to comply with applicable sections of the following International Standards:

- UL 60601-1: 2003, Medical Electrical Equipment, Part 1: General Requirements for Safety
- IEC 60601-1: 1988, Medical Electrical Equipment Part 1: General Requirements for Safety
- IEC 60601-1-1: 2000, Medical Electrical Equipment Part 1-1: General Requirements for Safety -Collateral Stanard: Safety Requirements for Medical Electrical Systems
- IEC 60601-1-2: 2004, Medical Electrical Equipment Part 1-2: General Requirements for Safety -Collateral Stanard: Electromagnetic Compatibility - Requirements and Tests
- IEC 60601-2-37:2005, Medical Electrical Equipment Part 2-37: Particular Requirements for the Safety of Ultrasonic Medical Diagnostic and Monitoring Equipment
- NEMA UD-2:2004, Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment
- NEMA UD-3:2004, Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment
- EN 55011:2000 Class A Industrial, Scientific, and Medical (ISM) Radio-Frequency Equipment -Radio Distrubance Characteristics-Limits and Methods of Measurement

9 Disposal Information

To return the Site~Rite* 5 System for end of life recycling, please contact your nearest Bard sales or distributor office in the country of purchase.

Warning: Always properly dispose of dead battery packs in accordance with local regulations. Improper disposal may present an environmental hazard.



The information contained in this document is subject to change without notice or obligation. This document contains proprietary information that is protected by copyright. No part of this document may be photocopied, reproduced, or translated without the expressed written consent of Bard Access Systems.

Site~Rite* 5 Ultrasound System Manufactured by:

Bard Access Systems, Inc. Salt Lake City, UT 84116

U.S.A.

(801) 595-0700

Customer Service: (800) 545-0890 Technical/Clinical Support: (800) 443-3385

Fax: (801) 595-4948

www.bardaccess.com

Site~Rite* Needle Guides Distributed by:

Bard Access Systems, Inc. Salt Lake City, UT 84116

U.S.A.

(801) 595-0700

Customer Service: (800) 545-0890 Clinical Information: (800) 443-3385

Fax: (801) 595-4948 www.bardaccess.com

