



CADScor<sup>®</sup> System

User Manual

Rx Only

## **1 Warranties and disclaimer**

The CADScor® System is covered by a general 1-year warranty from date of purchase, that covers faulty hardware, including the CADScor® Sensor, the CADScor® Docking station and the external power adaptor.

Furthermore, Acarix guarantees that the CADScor® System will be operational for at least 1000 recordings in the warranty period.

If, however the CADScor® System does not operate properly, please first consult the "Troubleshooting guide"-section 11 in this manual describing possible errors.

Any claims on warranty should be addressed to Acarix directly or to the Acarix distributor in your country, who will initiate replacement or repair of faulty equipment.

### *1.1 Disclaimer*

The information given in this document is adapted to an international market of healthcare professionals, based on illustrations, symbols and text. The document is giving no specific references to religious preference, ethnic origin, gender or political viewpoints. Any such relation is considered an interpretation of which Acarix cannot be held reliable. The information is intended to ensure the safe and accurate operation of the CADScor® System.

## **2 Copyright notice**

The contents of this document shall not be reproduced or communicated in any form to any third party without the prior written consent of Acarix.

### 3 Changes

This document is subject to change without notice and you are urged to contact Acarix to verify whether the document has been changed. The manual relevant to the specific equipment model can also be downloaded from the Acarix website. You can locate the specific model number (REF number) on the back of the aluminium frame of the CADScor® Docking station.

While every effort is made to ensure the correctness of the information provided in this printed manual Acarix disclaims any liability for errors and omissions herein.

The illustrations used in this manual may differ slightly from the appearance of the actual device, packaging materials or interface and reflects on-going efforts to improve safety, usability and clarify the overall instructions given.

### 4 Trademarks and Third-party software

CADScor® is a registered trademark owned by Acarix.

Third party software license agreements are listed in the back of this manual.

US-FDA revision 12.5, Mar. 20<sup>th</sup>, 2021.  
From software version 4.0 US-FDA.



Acarix  
Ryvangs Allé 81-83  
DK-2900 Hellerup  
Denmark  
[www.acarix.com](http://www.acarix.com)



## Table of content

|          |  |          |
|----------|--|----------|
| <b>1</b> | <b>Warranties and disclaimer .....</b>                         | <b>2</b> |
| 1.1      | <i>Disclaimer .....</i>  | 2        |
| <b>2</b> | <b>Copyright notice .....</b>                                  | <b>2</b> |
| <b>3</b> | <b>Changes .....</b>   | <b>3</b> |
| <b>4</b> | <b>Trademarks and Third-party software</b>                     | <b>3</b> |
| <b>5</b> | <b>Introduction .....</b>                                      | <b>6</b> |
| 5.1      | <i>About this manual .....</i>                                 | 6        |
| 5.2      | <i>The CADScor® System and the environment .....</i>           | 6        |
| 5.3      | <i>Product description .....</i>                               | 7        |
| 5.4      | <i>Intended use .....</i>                                      | 8        |
| 5.5      | <i>Use of the CAD-score .....</i>                              | 8        |
| 5.6      | <i>Risk-Benefit for use .....</i>                              | 10       |
| 5.7      | <i>Intended user profile .....</i>                             | 11       |
| 5.8      | <i>Intended patient population .....</i>                       | 11       |
| 5.9      | <i>Indications for use .....</i>                               | 12       |
| 5.10     | <i>Contra-indications for use .....</i>                        | 13       |
| 5.11     | <i>Intended conditions of use .....</i>                        | 14       |
| 5.12     | <i>Symbols used .....</i>                                      | 15       |
| 5.13     | <i>Symbols on the CADScor® System and in User manual .....</i> | 15       |
| 5.14     | <i>Warnings .....</i>  | 16       |
| 5.15     | <i>Precautions .....</i>                                       | 16       |
| 5.16     | <i>CADScor® System unboxing .....</i>                          | 18       |

|      |                                       |    |
|------|---------------------------------------|----|
| 5.17 | <i>The CADScor® Sensor .....</i>      | 20 |
| 5.18 | <i>CADScor® Docking station .....</i> | 22 |

|          |  |           |
|----------|--|-----------|
| <b>6</b> | <b>Installation of CADScor® System .....</b>                           | <b>23</b> |
| 6.1      | <i>Desktop model .....</i>   | 23        |
| 6.2      | <i>CADScor® System wall mounting ....</i>                              | 25        |
| 6.3      | <i>CADScor® Patch box opening .....</i>                                | 27        |
| 6.4      | <i>The CADScor® Patch .....</i>  | 28        |
| 6.5      | <i>The CADScor® Patch, assembly .....</i>                              | 29        |
| 6.6      | <i>Confirming correct sensor-patch assembly .....</i>                  | 30        |
| <b>7</b> | <b>Setting up the CADScor® System prior to use .....</b>               | <b>31</b> |
| 7.1      | <i>CADScor® Sensor .....</i>   | 31        |
| 7.2      | <i>First time configuration of the CADScor® Sensor .....</i>           | 31        |
| 7.3      | <i>Setting language .....</i>  | 34        |
| 7.4      | <i>Setting time .....</i>  | 35        |
| 7.5      | <i>Setting date .....</i>  | 36        |
| 7.6      | <i>Finishing Sensor SET-UP .....</i>                                   | 37        |
| 7.7      | <i>Universal access .....</i>  | 37        |
| 7.8      | <i>CADScor® Docking station .....</i>                                  | 39        |
| 7.9      | <i>Connecting/disconnecting power adapter to Docking station .....</i> | 40        |
| 7.10     | <i>Charging of the sensor battery .....</i>                            | 40        |

|          |   |           |
|----------|---|-----------|
| 7.11     | Docking station indication LED light .....                                | 41        |
| 7.12     | Qualification of sensor.....  | 42        |
| 7.13     | Forced qualification.....   | 43        |
| 7.14     | Re-setting sensor .....   | 44        |
| <b>8</b> | <b>Detection of heart sounds.....</b>                                     | <b>45</b> |
| 8.1      | Preparing the examination room ....                                       | 45        |
| 8.2      | Preparing the sensor.....   | 46        |
| 8.3      | Preparing the patient .....   | 47        |
| 8.4      | Identifying the fourth left Inter Costal space (IC4-L) .....              | 48        |
| 8.5      | Attachment of patch to sensor .....                                       | 51        |
| 8.6      | Explaining of the recording procedure to patient.....                     | 52        |
| 8.7      | Recording phase, patient risk factors . .....                             | 54        |
| 8.8      | Recording phase, Start .....  | 65        |
| 8.9      | Recording phase, Main recording.....                                      | 66        |
| 8.10     | Recording phase, cancelling .....   | 67        |
| 8.11     | Recording phase, CAD-scoring.....   | 68        |
| 8.12     | Short diastole CAD-score.....   | 69        |
| 8.13     | QR-code for printing, archiving and sending a CADScor patient report..... | 70        |
| 8.14     | CAD-score in the 30-39 year age group .....                               | 71        |

|      |   |    |
|------|---|----|
| 8.15 | Possible CADScor® System messages during or after recording ..... | 72 |
| 8.16 | Recalling a previous recording result .....                       | 73 |
| 8.17 | After the recording .....   | 75 |

|           |   |           |
|-----------|---|-----------|
| <b>9</b>  | <b>Sound guidance used in the interface . .....</b> | <b>76</b> |
| <b>10</b> | <b>Maintenance.....</b>                             | <b>77</b> |
| 10.1      | Disposal .....                                      | 78        |
| <b>11</b> | <b>Troubleshooting guide .....</b>                  | <b>79</b> |
| <b>12</b> | <b>CADScor® System requirements .....</b>           | <b>87</b> |
| 12.1      | System specification .....                          | 88        |
| 12.2      | Packaging materials .....                           | 89        |
| <b>13</b> | <b>Warranty .....</b>                               | <b>90</b> |
| <b>14</b> | <b>Approvals and EMC information.....</b>           | <b>91</b> |
| <b>15</b> | <b>Acronyms &amp; Definitions in text .....</b>     | <b>97</b> |
| <b>16</b> | <b>Licenses .....</b>                               | <b>98</b> |

## **5 Introduction**

### *5.1 About this manual*

This user manual is intended as a reference guide for the safe and correct use of the CADScor® System, after having been instructed in the proper operation of the CADScor® System.

This user manual contains both general and specific operating instructions, hereunder procedures for recording heart sounds, maintenance of the CADScor® System, troubleshooting, instructions and information of individual components.

To ensure optimal safety in operation and service of the CADScor® System, it is important to read this manual carefully and understand the use of the CADScor® System before starting to use the system professionally.

### *5.2 The CADScor® System and the environment*

The Acarix CADScor® System has been designed to minimize the environmental impact from fabrication, transport and use.

The System- and Patch boxing and packaging is made from lightweight recycled paper and cardboard, which can be disposed of as paper-waste and recycled.

The used patch and pouch can be disposed of as normal household waste.

The CADScor® Sensor, docking station and power adaptor contains electronic components, and should not be discarded in normal household waste, but returned for recycling at a regulated facility, local distributor or shipped back to Acarix for recycling.

### 5.3 Product description

The CADScor®System is a device for recording and quantifying acoustic noise arising from coronary artery stenosis micro-turbulence and myocardial movement. These noises are usually described as "coronary murmurs".

The CADScor®System calculates a patient specific CAD-score by computational processing of a recording obtained from the chest surface of the patient and the patient risk factors present.

The CADScor®System consists of two physical units; the CADScor®Sensor and the CADScor®Docking station for charging and qualification of the sensor. A specific power adaptor powers the Docking station. The patch for anchoring the sensor to the chest of the patient is a necessary accessory for proper functioning of the CADScor®System (all shown at right). The CADScor®System is operated by a graphical user touch-screen interface.

CADScor®System; Docking station holding Sensor.



CADScor®Patch and power adaptor below.



## 5.4 Intended use

The intended use of the CADScor®System is to record heart sounds, i.e. murmurs and vibration for calculation of a patient specific score, the CAD-score, indicating the risk of coronary stenosis, as an aid in cardiac analysis and diagnosis.

## 5.5 Use of the CAD-score

The CAD-score is a patient specific heart murmur score indicative of Coronary Artery Disease (CAD)/Chronic Coronary Syndrome (CCS) for immediate risk stratification, prior to potential secondary evaluation.

The CADScor®System risk stratification can be applied to patients with symptoms suggestive of CAD/CCS, in compliance with the indications for use.

Using the CAD-score to risk stratify patients prior to further testing will reduce un-necessary evaluation and risk.

The CAD-score can thus aid the decision to initiate additional evaluations or not, or to

observe the patient further prior to additional evaluations.

The presence of other patient risk factors or conditions may influence this decision.

### Definitions of CAD:

| Significant CAD  | Insignificant CAD  | Non-CAD   |
|--|--|---|
| ≥50% luminal diameter reduction by CAG.                      | < 50% luminal diameter reduction by CAG or CCT calcium > 0, or < 70% luminal area reduction. | Negative (0) calcium by CCT or no evidence of luminal stenosis. |
| CAG: Coronary Angiography; CCT: Coronary Computed Tomography |  |   |

The CAD-score is thus indicating risk of having significant CAD, defined as having ≥50% luminal diameter reduction.

Two risk categories are defined using the CADScor®System:

| CAD-score ≤20 | CAD-score > 20 |
|---------------|----------------|
| Low risk      | Elevated risk  |

The Negative (NPV) and Positive (PPV) predictive values vary with the prevalence



of significant CAD in the investigated population.

At lower significant CAD prevalence (e.g. 5-20%) in the investigated population, the NPV increases. The NPV of the CAD-score at or below 20, for ruling-out risk of significant CAD in the patient thus increases by falling prevalence, ranging from approximately 92.3-98.3% at 5-20% prevalence.

At higher significant CAD prevalence (e.g. 30-40%) in the investigated population, the NPV decreases. The NPV of the CAD-score for ruling out significant CAD thus decreases by increasing prevalence, ranging from approximately 87.5-81.8% at 30-40% prevalence.

A follow-up test, based on a CAD-score at or below 20, can be done after two to five years, if indicated.

At higher significant CAD prevalence in the investigated population, the PPV increases. The observed PPV of the CAD-score above

20 for predicting significant CAD is 14.4% (11.6-17.5%) in a 10.7% CAD prevalence population.

### **CADScor® System performance:**

Significant CAD vs Other (Significant CAD versus insignificant CAD and non-CAD combined).

#### Algorithm ver. 3.2-US:

Sensitivity: 87.5% (79.2-93.4%)

Specificity: 37.5% (34.2-41%)

Validation data CAD-prevalence: 10.7%

NPV: 96.2% (93.4-98%)

PPV: 14.4% (11.6-17.5%).

| CAD/CSS<br>Prevalence %  | NPV % *<br>CAD-score ≤20 |
|--|--------------------------|
| 5  | 98.3                     |
| 10   | 96.4                     |
| 15   | 94.4                     |
| 20   | 92.3                     |
| 30   | 87.5                     |
| 40   | 81.8                     |
| * modelled from Acarix clinical data for sensitivity/specificity above at threshold ≤20. |                          |

## 5.6 *Risk-Benefit for use*

By using the CADScor® System for evaluating patient risk of Coronary Artery Disease, a CADScor low risk of significant CAD can help to direct other cause evaluation decisions and also remove unnecessary anxiety.

Patients with a CADScor elevated risk of significant CAD may be earlier recognized and thus be diagnosed and enter treatment earlier.

By use of the CADScor® System no risks are associated directly from use, except in rare cases where a slight transient skin reaction towards the patch adhesive may be observed.

Inform patients who have CAD-score  $\leq 20$  to seek medical attention if symptoms persist or worsen after initial evaluation.


The risk of a false negative result for a true patient or a false positive result for a non-CAD patient should be taken into consideration when evaluating the patient.


### *5.7 Intended user profile*

The CADScor® System is intended to be operated only by registered nurses, clinical/medical laboratory technicians, medical doctors/physicians after having been instructed in the proper operation of the CADScor® System.

### *5.8 Intended patient population*

The CADScor® System is intended to be used in males and females (by born gender) above 40 years of age.

 Do not use the CADScor® System without the necessary qualifications, and only after instruction and having read and understood this User-Manual, due to risk of death or serious injury.

 The CAD-scores from patients in the 30-39-year group are currently outside the intended patient population. A warning triangle indicates the higher uncertainty of their CAD-scores (8.14).

### 5.9 Indications for use

The CADScor® System is indicated for use as a diagnostic aid in symptomatic patients suspected of stable Coronary Artery Disease/Chronic Coronary Syndrome.

Symptoms of CAD/CCS may manifest as\*

#### **Typical angina symptoms;**

(all three below)

- Sub-sternal chest discomfort of characteristic quality and duration;
- Provoked by exertion or emotional stress;
- Relieved by rest and/or nitrates within minutes.

#### **Atypical angina symptoms;**

(Two of the above)

#### **Non-anginal chest pain;**

Lacks or meets only one or none of the above characteristics.

#### **Dyspnea/ breathlessness**

---

\* as defined by ESC Guidelines 2019 & NICE UK Guidelines CG95, 2019.

### *5.10 Contra-indications for use*

Asymptomatic for angina or chest pain,  
Implanted donor heart,  
Previous Coronary Artery Bypass Graft  
(CABG),  
Previous Coronary stenting or known CAD  
Arrhythmia causing non-sinus rhythm,

Fragile or compromised skin,  
or abnormal anatomy or significant  
operation scars in the fourth left Inter  
Costal (IC4-L)-recording area.

Implanted mechanical heart or mechanical  
heart pump,  
Implanted Pacemaker or Cardioverter  
Defibrillator (ICD),  
Other implanted active electronics or active  
electronic support equipment closer than  
50 cm to the CADScor® System.

The CADScor® System has not been  
validated as a screening tool for Coronary  
Artery Disease/Chronic Coronary  
Syndrome in asymptomatic populations.



Do not use the CADScor® System on  
patients with implanted electronics like  
ICD, Pacemakers, heart-pumps or closer  
than 50 cm to similar active electronic  
support equipment, due to risk of  
equipment failure from  
CADScor® System electromagnetic RFID  
impulse.


### 5.11 Intended conditions of use

The CADScor®System is intended to be used at medical clinics and hospitals, at room temperature and below 2500 meter above sea level.

The patient shall be lying on the back on an examination bed (supine position) during identification of the fourth left Inter Costal (IC4-L) space and recording of heart sounds. Mark the IC4-L position after location, approximately two centimetres of the sternum bone for placement of patch.

The CADScor®Sensor must only be used together with a CADScor®Patch. Patches are single-use only and should not be moved on the body after first adherence. If a patch has been misplaced a new patch must be used. The applied time of the CADScor®Sensor (including patch) is normally less than 15 minutes.

Use the CADScor®System below 2500 meters.



- ⚠ Train to identify left IC4 location. Wrong position may influence result.
- ⚠ Do not attempt to re-use a CADScor®Patch. The patch is single use only for hygienic reasons.

### 5.12 Symbols used

#### Definitions:














The CADScor® System user manual contains symbols, operational warnings and cautions, which are important and should be read and understood carefully before performing the related procedures.

Symbols and warnings are also found on the CADScor® System, the CADScor® Patch and the system packaging boxes.

### 5.13 Symbols on the CADScor® System and in User manual

**Warnings:** Describes a condition or a situation where risk of death or serious injury can occur.

**Precautions:** Describes a condition or a situation where a non-serious injury can occur to patient or user or damage the equipment or property.

| Symbol  | Explanation   |
|---|---|
|   | Product code number.  |
|    | Serial number   |
|    | Warnings and precautions. See description below (5.13-5.15).  |
|    | Instruction to Read user manual.  |
|    | Instruction to Consult user manual.   |
| <b>Rx only</b>  | Caution: Federal law restricts this device to sale by or on the order of a physician.   |
|    | Magnetic Resonance Imaging room unsafe  |
| * /   | Indication of a short diastole (*) result or a notice for result (  )                        |
| IPXX  | Denotes liquid and dust ingress protection levels   |
|    | Defining the Sensor as the Applied part of CADScor® System  |
|    | Manufactured by   |
|  /  | Expiry date / single use of patches   |
|    | The CADScor® System, when reaching its end of life, must be collected and recycled separately from other waste according to national requirements (Directive 2012/19/EU, WEEE). |



## 5.14 Warnings



*Warnings are alerts described to alert and reduce the risk of death or serious injury from use of the CADScor® System.*

Warnings and precautions are labelled with a warning triangle and also refers to the relevant section in the User Manual.

### Patient and user risk warnings:

-  Do not use the CADScor® System without the necessary qualifications, and only after instruction and having read and understood this User-Manual, including, symbols, contra-indications, warnings and precautions, due to risk of death or serious injury.
-  Do not use the CADScor® System on patients with implanted electronics like ICD, Pacemakers, heart-pumps closer than 50 cm to similar active electronic support equipment, due to risk of equipment failure from

CADScor® System electromagnetic RFID impulse (5.10; 6.4).



Do not attempt to re-use a CADScor® Patch. The patch is single use only for hygienic reasons, to prevent possible cross contamination/infection between patients (5.11; 8.7).



The CADScor® System and patches are MRI-unsafe and should not be operated or placed in a Magnetic Resonance Imaging (MRI-) room, due to risk of skin burns and/or magnetic attraction and impact (7).



Never use other power adaptor than supplied with the CADScor® System, due to risk of electrical shock (5.17).

## 5.15 Precautions






*Precautions are alerts to reduce potentially hazards situations or risks of non-serious*




*injury or reduced effectiveness from use of the CADScor® System.*





*Risk of incorrect CAD-score:*

-  Incorrect placement of the CADScor® Sensor, outside patient IC4-L, may result in an incorrect CAD-score (8.4).
-  Inform patients who have CAD-score ≤20 to seek medical attention if symptoms persist or worsen after initial evaluation”, due to the risk of a false negative CADScor result (5.6).
-  The CAD-scores from patients in the 30-39-year group are currently outside the intended patient population. A warning triangle indicates the higher uncertainty of their CAD-score (8.14).

*Risk of damage to the CADScor® System:*

-  Do not touch the charge-point terminals on the CADScor® System, due to risk of permanent damage to the electronics

inside by electrostatic discharge (5.17-5.18)

-  Do not drop or exert excessive force to the CADScor® System since this may damage the CADScor® System permanently (6.1).
-  Do not expose the CADScor® System to liquids (water, oils, detergents or similar) or dust, since this may damage the CADScor® System permanently (10).
-  The CADScor® System cannot be heat- or radiation sterilized, or machine washed/cleaned, since this may damage the CADScor® System permanently (10).
-  Do not modify the CADScor® System or use or repair a defect CADScor® System, due to risk of malfunctioning. The CADScor® System must only be serviced by qualified Acarix personnel (10).

### 5.16 CADScor® System unboxing

Place the box on a table and open it as shown in figure below.

Establish the content of the CADScor® System box (one each):



- CADScor® Sensor.
- CADScor® Docking station.
- A box containing a power adaptor plus country specific plug.
- A box in the lid holding;
- a User Manual,
- a separator tool to separate the back plate from the top part of the Docking station and
- a drilling template for mounting the Docking station on the wall.




The CADScor® System is pre-assembled in desktop position.


Lift the CADScor® System from the system box by holding the system at the edges of the aluminium plate of the Docking station.

Place CADScor® System on steady surface. Inspect the CADScor® System for possible transport damages.

Take out the box containing the power adaptor and supplied plug. Assemble power adaptor unit by pressing the plug firmly onto the power adaptor until clicked securely into place. No voltage adjustments are needed for the power adaptor.

 Never use other power adaptor than supplied with the CADScor® System, due to risk of electrical shock.

Finally take out the lid box containing the user manual, separation key and drill template.



### 5.17 The CADScor® Sensor


The CADScor® Sensor is the part to record heart sounds at the fourth left Inter Costal region and to calculate the CAD-score based on the acoustic recording.

The sensor has one physical button, a touch display and a sensor head extending on a flexible arm from the sensor body (see figure at right).

The push button has several functions when pressed or pressed and held for extended time.

| Press       | Function                  | Beep        |
|-------------|---------------------------|-------------|
| 1x          | Turn sensor on            | -           |
| 2x          | Cancel on-going recording | Triple beep |
| 3x          | Enter Settings menu       | -           |
| Hold> 4 sec | Forced qualification      | Single beep |
| Hold> 8 sec | Sensor off and reset      | Double beep |

Main components of the CADScor® Sensor, top/side view.



The sensor LCD touch display shows a user interface to guide you through a recording session. Also, the sensor status and the CAD-score are shown on the display when it is calculated.

The flexible sensor arm connects the sensor body to the sensor head. Within the sensor head a heart microphone and an ambient microphone are placed.

Inside the sensor body is a built-in speaker that generates auditory cues for guiding the patient to pause breathing.

The sensor charge/contact point terminals are used to re-charge the sensor battery, when the sensor is placed in the Docking station (See figure at right).

Main components of the CADScor® Sensor, bottom-side view and labelling at right (SN refers to CADScor® Sensor serial number).




⚠ Do not touch the charge-point terminals on the CADScor® System, due to risk of permanent damage to the electronics inside by electrostatic discharge.

### 5.18 CADScor® Docking station

The Docking station serves as a "home" for the sensor to which the sensor is mated specifically from the factory.

In the Docking station the sensor is recharged and qualified. A specific sensor can only be qualified in its home Docking station.

 Do not touch the charge-point terminals on the CADScor® System, due to risk of permanent damage to the electronics inside by electrostatic discharge.

The Docking station LED shows the current status of the CADScor® System when the sensor is docked and powered (Colour codes in section 6.1).

The Docking station can be wall mounted for easy access and saving desk space (see section 6.2 below).

| Parts identification |                             |
|----------------------|-----------------------------|
| Top view             | Bottom view                 |
| Indicator LED        | Back plate                  |
| Sensor tester        | Product information marking |
| Contact terminals    |                             |
| DC-in                |                             |




## 6 Installation of CADScor® System

### 6.1 Desktop model

Place the Docking station on a stable surface (desktop/table). Desktop position is the pre-set operating mode of the Docking station.

Insert the CADScor® Sensor in the Docking station using both hands. Make sure that the sensor head is fully inserted into the sensor tester and that the red indicator line on the sensor head is completely hidden. A spring action mechanism will hold the sensor-body firmly in the Docking station.

Insert the small power plug to the Docking station on right hand side and the power adaptor into a 100-240 VAC power outlet. Turn on wall outlet power.











The LED on the Docking station lower right corner should light up.

Three different colours are used to indicate system status (table at right).

If no colour is lighting up the Docking station LED, wall outlet power has not been turned on or the power adaptor is faulty.

See troubleshooting guide section 11.

 Do not drop or exert excessive force to the CADScor® System since this may damage the CADScor® System permanently.

| LED indicator   | Sensor status  |
|---|--|
| Constant YELLOW<br> | Power on docking station   |
| Constant GREEN<br>  | Sensor fully charged; ready for use  |
| Flashing GREEN<br>  | Ongoing charging of sensor battery above minimum level                     |
| Constant RED<br>    | Error in Sensor or Docking station; Not ready for use (in sensor display). |
| Flashing RED<br>    | On-going qualification (in sensor display)                                 |
| Flashing RED<br>    | Low on Battery (in sensor display)   |
| The LED is off<br>  | Power not applied.   |




## 6.2 CADScor® System wall mounting

The CADScor® System can be wall mounted to save desktop space or for increased general visibility.

Determine suitability of the wall material and items in the wall before mounting.

To unlock the back plate from the top part of the Docking station, hold the Docking station by the aluminium plate, one end resting on a table. Place the separator tool into the two small holes in the edge and press gently using the index fingers of both hands and separate slowly (see figure at right). Do not apply strong force since slipping grip may happen and lead to injury.

Use drill template and a pencil to mark up the wall position. Ensure horizontal level using a spirit level or similar.



Drill two holes (Ø 4mm) through the centre of marks and mount back plate (thinner rim facing up) using appropriate screws and plugs.

Place the top part of the Docking station on the back plate, sensor tester placed at upper right corner, and click into locking position.

Ensure that top plate has engaged at all four corner positions.

Place CADScor® Sensor in the Docking station and apply power.



### 6.3 CADScor® Patch box opening

Place the CADScor® Patch box on a table and open along the perforated line at the box front and top sides (see figure at right).

Each CADScor® Patch box contains 20 patches individually packaged in peel-pouches. The pouch is opened at one-end, by holding and separating the two layers.

Additionally, an assembly tool is provided to ensure the correct alignment and assembly of the patch to the sensor head and body.

On the paper backing of the patch peel-pouch, instructions are given for the assembly. Follow these instructions carefully to obtain high quality heart recordings (see section 6.5-6.6).

Patch box with assembly tool and labeling:



#### 6.4 The CADScor® Patch


The disposable patch is used for anchoring the sensor onto the chest of the patient and provides optimal recording conditions (see figures at right).

The CADScor® System cannot be operated without a valid patch.

The CADScor® Patch has adhesive material on both sides and a coupling ring where the sensor head has to be inserted. On the patch is furthermore a numbering sequence to follow when applying the patch and sensor onto the patient's chest.

A Radio-frequency identification chip (RFID) can be seen under the top foil, used to validate the patch before use.

Upon removal of the sensor from the patch a pull-string is drawn to open the coupling ring for easier disassembly.



### 6.5 The CADScor® Patch, assembly

Place the cardboard assembly tool on an even hard surface (e.g. table top).

Align the un-packed patch on top of the assembly tool, so that the patch is within the outline on the tool. Make sure the coupling ring is situated correctly on top of the hole cut-out in the assembly tool (See figure at right).

Place the sensor on top of the patch, aligning to the axis of the patch. The sensor head shall be placed in the coupling ring of the patch.

Press the sensor head gently and evenly into the coupling ring, using the palm of your hand.

The sensor head will then correctly snap into the coupling ring.




## 6.6 *Confirming correct sensor-patch assembly*

The red colored assembly indicator ring on the sensor head must now be hidden in the coupling ring, when looked from the side.

Visually confirm the disappearance of the red colored indicator ring all around the sensor head (see figure at right).

From the patch underside a slight protrusion of the sensor head through the patch can be observed.



## 7 Setting up the CADScor®System prior to use

### 7.1 CADScor®Sensor

The CADScor®Sensor is the part of the CADScor®System to record heart sounds at the left IC4 region, and to display the calculated CAD-score.

### 7.2 First time configuration of the CADScor®Sensor

Connect the power adaptor to the Docking station and turn on supply voltage. Place the sensor in the Docking station, if not already docked, to charge the battery.

The sensor will automatically detect the Docking station and power on, showing the Acarix logo (as shown at right).



The CADScor®System and patches are MRI-unsafe and should not be operated or placed in a Magnetic Resonance Imaging (MRI-) room, due to risk of skin burns and/or magnetic attraction and impact.





As part of the power-on sequence an information display will show current software version and user relevant warnings and instructions (REF:101).

Since qualification of the sensor is needed every 23rd hours, a display message will show that qualification is in progress (REF:302).

Follow the instructions to keep quiet and not touching or removing the sensor from the Docking station.  
After qualification and first time the sensor is powered-on, configuration is needed.

While the sensor is "docked" in the Docking station, the "last-recording" display (REF:103) is shown.

Take out ("un-dock") the sensor from the Docking station to enable configuration.






The sensor will automatically go to the sensor "home" display (REF:102).


Press "Enter setup" on the "home" display to configure the sensor (REF:102) or triple-click the push button.

A sub-menu will appear, allowing adjustment of Date, Time and Language (REF:500).

Follow the instructions below to do the correct adjustments.

The SERVICE MODE option is intended for Acarix authorized personnel to e.g. install new software and is not accessible in normal use.

 Do not modify the CADScor® System or use or repair a defect CADScor® System, due to risk of malfunctioning. The CADScor® System must only be serviced by qualified Acarix personnel (10).



REF: 102

Last CAD-score: **43** Elevated risk


Record-ID: **3091**

@ time: **SEP 27, 2020, 10:35**

**New Examination**

**Enter setup**

Today: **SEP 28, 2020 09:18**



REF: 500

**Settings**

|            |              |
|------------|--------------|
| DATE       | LANGUAGE     |
| TIME       | SERVICE MODE |
| LOOK-UP ID | ABOUT        |



**DONE**

### 7.3 Setting language


Press LANGUAGE on the Settings display (REF:500) to change user interface language.

The default language is English.

In the "Set language" display (REF:501),

Press  or  to select a user interface language.

Press "DONE" to save language choice and exit "LANGUAGE" menu.





## 7.4 Setting time

Press "TIME" to set the correct time.

Time format is international 24 hours, and cannot be changed.

In the "Set time" display (REF:502),

Tap the hours (HH) or minutes (MM) fields individually to activate the wanted field.

Press  or  to change time (HH:MM)

Press "DONE" to save time setting and exit "TIME" menu.

When TIME has been changed, a re-qualification is automatically required.



## 7.5 Setting date



Press "DATE" to set the correct date.

Date format is Month: Day: Year, where month is given by three letter abbreviations.

The date format cannot be changed.


In the "Set date" display (REF:503),

Tap the month, day or year fields individually to activate field.

Press  or  to change date (M/D/Y).

Press "DONE" to save and exit "DATE" menu.

When DATE has been changed, a re-qualification is automatically required.



## 7.6 Finishing Sensor SET-UP



Press "DONE" to exit "Settings" display.  
The sensor will return to the "home" display (REF:102).

The sensor will turn off automatically after five minutes of inactivity.

## 7.7 Universal access

A functionality to enter the set-up menu directly is available by pushing the push button three times in succession. This cannot be done during a recording session. This functionality is also enabled while the sensor is in the Docking station.

In case an unknown language has been chosen, triple-push the button and select the upper right blue quadrant to enter language settings. Choose relevant language and press "DONE" to finish language selection.




While the sensor is docked, a pause display ("Last recording") will be shown (REF:103).

After a period of inactivity, the sensor will automatically power-down to save energy and the screen will go black.

The sensor will then only perform automatic qualification, every 23rd hour.

The sensor will power-on again after pressing the push button or taking the sensor out of the Docking station.

When powered on again, the sensor will return to its "last recording" display while in the Docking station (REF:103) or return to the "home" display outside the Docking station (REF:102).



### 7.8 CADScor® Docking station

The Docking station is used for:

- Holding the sensor and charge the sensor battery when docked.
- Qualify the sensor at regular intervals or on request.
- Indicate the CADScor® System status by LED light.
- Transfer of data from the sensor to a removable memory card (only available in clinical trials through arrangement with Acarix).
- Update the sensor software.



### *7.9 Connecting/disconnecting power adapter to Docking station*

The CADScor® Docking station has a dedicated external universal power adaptor (100-240 VAC/50-60Hz).

To supply the Docking station, connect the small power adaptor plug into the jack on the side of Docking station.

### *7.10 Charging of the sensor battery*

In the Docking station, charging of the sensor battery is initiated through the charging contact terminals of the Docking station.

Approximately 30 minutes are needed to charge the sensor battery allowing for approximately 5-10 standard recordings and CAD-score calculations.










Never use other power adaptor than supplied with the CADScor® System, due to risk of electrical shock (5.17).



### 7.11 Docking station indication LED light

When the sensor is placed in the Docking station, an LED on the aluminium front indicates the present status of the CADScor® System.

More information on the sensor status can be found in the display of the sensor if turned on in the Docking station.

| LED indicator   | Sensor status  |
|---|--|
| Constant YELLOW<br> | Power on docking station   |
| Constant GREEN<br>  | Sensor fully charged; ready for use  |
| Flashing GREEN<br>  | Ongoing charging of sensor battery above minimum level                     |
| Constant RED<br>    | Error in Sensor or Docking station; Not ready for use (in sensor display). |
| Flashing RED<br>    | On-going qualification (in sensor display)                                 |
| Flashing RED<br>    | Low on Battery (in sensor display)   |
| The LED is off<br>  | Power not applied.   |



### 7.12 Qualification of sensor

The CADScor® System is qualified at regular intervals to ensure safe operation. During qualification, a series of tones are played from the Docking station to the sensor.

The qualification procedure is automatically carried out every 23 hours, if the sensor is placed in the Docking station. If the sensor qualification has expired outside the Docking station, a display message will ask for sensor to be placed in the home Docking station (REF:300).

It is important that the sensor head is pressed firmly down into the built-in sensor tester.

The red indicator ring on the sensor head should be invisible when inserted correctly. When the qualification is on-going a display message is shown (REF:302).



### 7.13 Forced qualification



If you are bringing the CADScor® Sensor along and suspect qualification will expire before returning, a forced qualification can reset expiration clock by 23 hours.

Place the sensor in the home Docking station

The sensor will power-on automatically in the powered-on Docking station.

Press and hold the push button for 4 seconds (a single beep feedback will be given) to activate sensor qualification procedure (REF:302).

The sensor display will return to the "last recording" display when qualified (REF: 103).



### 7.14 Re-setting sensor

The CADScor® Sensor operating system can be reset in the event of unresponsive software failure. The re-setting will force the sensor to reload sensor-software and also re-qualify the sensor.

Procedure:

Press and hold the push button for 8 seconds (beyond the 4 seconds single beep, until a double beep signal feedback is given) to power off the sensor.

The sensor will automatically turn on again while in the powered Docking station or by button press.

The sensor will then reload operating software and qualify automatically.

The sensor will display the start-up display (REF:101) followed by the "last recording" display (REF:103) when the qualification has finished, and sensor is docked.



## **8 Detection of heart sounds**

### *8.1 Preparing the examination room*

The CADScor® System records very weak heart murmurs. Reducing external or environmental noise will increase the likelihood of a successful recording.

Close doors and windows to reduce noise from outside sources like traffic noise, talk, construction work and similar.

Turn off indoor noisy or unused electrical equipment (e.g. noisy PC, mobile telephones, ventilators, fans, air-condition).

Unplug power to patient bed if available.

Keep quiet (both user and patient) during the heart sound recording.

Reduce room noise!




The CADScor® System and patches are MRI-unsafe and should not be operated or placed in a Magnetic Resonance Imaging (MRI-) room, due to risk of skin burns and/or magnetic attraction and impact.







## 8.2 Preparing the sensor

Press the push button to turn on the sensor in the Docking station or take out the sensor from the Docking station. After power-on (and eventual qualification) the display will return to the "home"-display (REF:102) if sensor is ready, and outside Docking station.

The sensor battery status and the number of remaining recordings on current battery level is shown in the display, top right (20 REC). Battery level low: Charge the battery if the level is below minimum (red battery icon).

Qualification status is shown by circled blue check mark in upper display information line. Qualification expired: Qualify the sensor if the previous qualification has expired (red cross in circle).



| Icon  | Icon descriptions  |
|---|--|
| REF:  | Display reference number   |
| REC-ID:   | Recording identification number  |
|  | Sensor qualification status icon   |
|   |  Qualification valid.<br> Qualification expired. |
|  | Battery level icon   |
|   |  Full battery level<br> Low battery level.       |
| REC.  | Remaining capacity recording icon  |

### 8.3 *Preparing the patient*

Always use a paper/linen bedcover between patient and examination-bed to avoid rubbing noise. This will reduce situations potentially requiring a second recording.

#### **Hair removal**

Also, to obtain high-quality heart sound recordings and to reduce secondary recordings, body hair at the IC4-L region has to be removed as well as hair in the area of the applied patch.

The extended shaving area will ensure better anchoring of the CADScor® Sensor and easier removal of patch after recording.

Shaving can be done using surgical clipper or standard single use shaving tool.

Greasy or moist skin should be cleaned using standard alcohol skin wipes.

Always conduct proper skin cleaning and hair removal in patch area for optimal patch adhesion and high-quality heart sound recordings.

#### 8.4 Identifying the fourth left Inter Costal space (IC4-L)



Establishing the correct recording position on the patient is very important for obtaining a valid heart sound recording.

The fourth left Inter Costal space is most reliably identified using a simultaneous four-finger palpation on the patient's chest.

This method is also known as the "finger-method" for IC4-L identification, picture at right.

Un-dress the patient in the thorax area. It is important that clothing will not rub against the sensor during recording.

Place the patient flat on the examination bed in supine position (chest facing upwards).





Standing on the patients' right side, spread your left-hand fingers approximately 1.5 centimeters, and place your left-hand little finger up against the collar bone (clavicular). Just under the collarbone costa number one (C1) is found.

Press the other fingers against the chest left side and massage the fingers into the intercostal spaces at the sternum border.

A moderate to high palpation pressure is needed to recognize the costa and the intercostal spaces.

When each of the four fingers has found their intercostal spaces, mark up the fourth left Inter costal space under your index finger.

Mark the position of IC4-L, approximately two centimeters from the sternum border into the fourth Inter Costal space.



Use a clearly visible pen marker, to later reliably locate the established recording position to avoid repositioning the sensor-patch assembly and thereby reduce skin adherence.

The marked IC4-L spot indicates the centre of the recording site for the sensor head.

The sensor head should rest free of the sternum border.

Train to be able to identify the IC4-L position correctly.

The IC4-L position is equivalent to the electrode V2 position in a standard 12-lead ECG.

Train and counsel with a colleague to identify and establish the correct position for the IC4-L location, using the finger technique. Train on a number of patients until safe IC4-L locating experience has been obtained.

Train to obtain IC4-L locating experience.



Incorrect placement of the CADScor<sup>®</sup> Sensor, doing recording outside IC4-L on the patient, may result in an incorrect CAD-score.

### 8.5 Attachment of patch to sensor

An assembly tool is provided with the patches to mount the patch correctly to the sensor (see also section CADScor® Patch assembly).

Place the CADScor® Patch on top of the assembly tool as outlined on the tool (see figure at right).

Then place the CADScor® Sensor on top of the patch.

Press sensor head carefully into coupling ring by using the palm of your hand. Make sure that the sensor head is secured at the stop position and that the red indicator ring on the sensor head is now not visible. The sensor and patch assembly is now ready for attachment to the patient.



### 8.6 Explaining of the recording procedure to patient

Explain the recording procedure to patient. Being generally relaxed and pausing breathing at correct intervals is important for recording of high-quality heart sounds. Observe at least 5 minutes of patient resting to ensure hemodynamic balance before doing CAD-scoring and blood pressure readings.

The recording sequence is less than three minutes long and divided into a pre-recording and 4 recording loops.

During each recording loop, the patient should pause his breath for 8 seconds. A sound will be heard when it is time to pause breathing. The sound will be demonstrated to the patient before starting.

Instruct patient to relax respiration and do belly breathing instead of thorax breathing. Show belly breathing by patient hand on

belly, including inhalation, exhalation and natural breathing pause after exhalation.

| Procedure          | Step                 | Time (s) | Actions                        |
|--------------------|----------------------|----------|--------------------------------|
| Resting phase      |                      | > 300    | Patient lies in bed.           |
| Recording phase    | Pre-recording        | 30 + 30  | Analysing recording conditions |
|                    | Loop Nr. 1 Breathing | 18       | No recording                   |
|                    | Breathing pause      | 8        | Recording                      |
|                    | Loop Nr. 2 Breathing | 18       | No recording                   |
|                    | Breathing pause      | 8        | Recording                      |
|                    | Loop Nr.3 Breathing  | 18       | No recording                   |
|                    | Breathing pause      | 8        | Recording                      |
|                    | Loop Nr. 4 Breathing | 18       | No recording                   |
|                    | Breathing pause      | 8        | Recording                      |
| Result phase       | Post recording       | 60-120   | Calculation of CAD-score       |
| Approx. Total time |                      | 524-564  |                                |

### *Before Recording phase*

#### **Summary**

Check examination room for acoustic and electrical noise. Check sensor readiness.



Ask patient to remove clothing from the waist-up and lie on the back on the examination bed (supine position) and rest.

Identify the IC4-L region and do hair removal in IC4-L/patch region if needed.

Observe at least 5 minutes of patient resting to ensure hemodynamic balance before doing CAD-scoring and blood pressure readings.

Assemble new CADScor® Patch and sensor using assembly tool. Note absence of red indicator ring around sensor head in coupling ring of patch.

Instruct patient in recording sequence, breathing mode and instruction sounds.



### 8.7 Recording phase, patient risk factors

Initiating the recording sequence:

If not already powered-on press push button once to activate sensor. The Acarix start-up display (REF:101) will be shown during sensor self-test. Choose "New examination" on the "home"-display (REF:102), shown when sensor is undocked.

Whenever "New Examination" is chosen a new REC-ID will be assigned.

Four patient specific risk factors are required to initiate the recording and are included in the calculation of the CAD-score:


Patient Gender

Patient Age-group

Patient Symptom and

Patient hypertension status

The risk factors are entered sequentially and summarized for confirmation.



Patient birth Gender should be entered on display (REF:104).

If the entered data should be corrected press EXIT, to return to New Examination display (REF:102).

Patient Age-group, represented as 10-year intervals, should be registered as age on the day of the examination on display (REF:105).

If the entered data should be corrected press EXIT, to return to New Examination display (REF:102).

Patient symptom, as defined by ESC and NICE UK Guidelines (refer to section 5.9) should be entered in (REF:106).

If the entered data should be corrected press EXIT, to return to New Examination display (REF:102).

REF:104 REC-ID: 3091

### Gender of patient

Female patient

Male patient

EXIT

REF:105 REC-ID: 3091

### Age group of patient

30 - 39 40 - 49 50 - 59

60 - 69 70 - 79 80 -

EXIT

REF:106 REC-ID: 3091

### Symptom of patient

Typical angina Atypical angina

Non-anginal chest pain Dyspnea


EXIT

Presence (or not) of Patient hypertension should be entered in display (REF:107). Hypertension is defined as patient having systolic blood pressure equal or above 140 mmHg or being in medical treatment for systolic hypertension.

If systolic blood pressure is unknown and patient is not in medical treatment for hypertension, a standard blood pressure reading after the resting period and prior to the CADScor measurement, should be done to establish hypertension status. Observe at least 5 minutes patient rest to ensure hemodynamic balance before doing CAD-scoring and blood pressure readings.

If the entered data should be corrected press EXIT, to return to New Examination display (REF:102).

After entering hypertension status, a patient data summary display (REF:108) is shown. Confirm the patient data to continue or EXIT to re-enter data.



REF: 107 REC-ID: 3091 20 REC.

### Patient hypertension

Sys. BP  $\geq$  140 mmHg and/or in treatment with BP medication

No hypertension in patient

EXIT



REF: 108 REC-ID: 3091 20 REC.

### Patient summary

Female patient, 50-59  
Non-anginal chest pain  
Hypertensive

Confirm patient data

EXIT



If the patient age entered is within the patient age-group 30-39 years, a notice display (REF:219) is shown to indicate that the patient age is below the current intended patient population for the CADScor®System.

A warning-sign (⚠) behind the resulting CAD-score will indicate the higher uncertainty of the CAD-score from this group.

If the information of the patient was correct and approved, press Continue to acknowledge the information and proceed to Patch assembly and patient murmur recording or press EXIT to change the patient details or abort the current CADScor measurement.



## Attaching sensor/ patch assembly to patient

If sensor is ready to record, the sensor will ask for new patch (REF:110).

Mount a new patch as described above (section 8.5) or as described on the back of the patch pouch.


Inspect if sensor head is correctly seated in the coupling ring.

When patch is attached to sensor, press DONE.

The patch production data are checked by the sensor (REF:127).

Avoid placing the assembled sensor and patch on a metal-table during patch checking since this can interfere with the sensor reading the patch information.

*Always detect patch before patch is applied to patient.*




When the patch information has been verified and recognised a confirmation is given (REF:111).


Proceed to place sensor on IC4-L.

First remove the release liner marked “1” on the patch.

The part under the sensor microphone is taken off separately, divided from the remaining release liner by a pre-made cut.

Keep the remaining release liner part in place until sensor is placed at IC4-L.


 Do not attempt to re-use a CADScor® Patch. The patch is single use only for hygienic reasons, to prevent possible cross contamination/infection between patients.




Now, place the sensor head at the pre-marked IC4-L position, orienting the sensor body towards the patient head on the sternum.

Do not apply strong pressure to the sensor head when applying patch on the patient.

Ensure proper adherence of sensor head to chest by pressing on coupling ring around the sensor head.



Then, Peel-off release liner marked “**2**” on the patch and anchor the long part of the patch onto the chest.



Finally, remove release liner marked “**3**” on the patch and fix the sensor body to the top adhesive on the patch.

Do not pull the sensor to fix it on the patch, but let the sensor follow the curvature of the patient chest.



When sensor and patch assembly is securely positioned, press DONE to verify placement on IC4-L.



When the adhesive of the patch has been adhered to the patient, the Sensor and patch should not be moved.

Instruct patient to relax thorax by belly breathing and generally relaxing muscles: Relaxing while doing a CADScor® recording will increase the chances of a successful recording and subsequent CAD-scoring.

Find a suitable generally relaxing supine position, placing arms along the body and resting the head against the mattress. Offer a pillow for increased comfort or elevated headrest.

Ask and ensure patient is comfortable and resting without significant muscle tonus.

Observe at least 5 minutes patient rest to ensure hemodynamic balance before doing CAD-scoring and blood pressure readings.



Three options will then become available (in REF:112):

START SOUND DEMO to demonstrate recording sounds to patient



START RECORDING to initiate the actual recording, or

ROTATE DISPLAY to rotate the display by 180 degrees for better viewing comfort.

For the patient, it is important to become familiar with the sound instructions given during the recording phase.

Select START SOUND DEMO (REF:126), and demonstrate the PAUSE BREATHING (single beep), START BREATHING (double beep), RECORDING START (ascending tones) and RECORDING STOP sounds (descending tones).

Adjust speaker volume if needed.



Press DONE to exit SOUND DEMO after patient has become familiar with the sounds.


Again, instruct patient to relax respiration and do belly breathing instead of thorax breathing.

Belly breathing will ensure better recording results.

Explain belly breathing to patient.

Show e.g. belly breathing by patient hand on belly, including inhalation, exhalation and natural breathing pause after exhalation.

Explain to the patient that he or she will be guided by the operator as well, to pause or start belly breathing intervals.



Instruct belly breathing



## 8.8 Recording phase, Start




The CADScor® Sensor will evaluate current recording conditions during a Pre-recording, before continuing to the main heart sound recording.

Press START RECORDING (REF:112) to begin recording process.

During Pre-recording, no talking should take place (REF:113).

If external, internal and electrical noise conditions during Pre-recording are evaluated acceptable by the sensor, the sensor will proceed to heart sound recording automatically.

The Pre-recording and evaluation will take approximately one minute (REF:114).



### 8.9 Recording phase, Main recording


The Main recording has four recording loops of approximate 26 seconds each.

Each recording loop has two phases, 18 seconds of normal belly breathing and eight seconds of breathing pause.

During the Main recording, each recording phase will be indicated on the display in text as "Continue breathing" or "Pause breathing".

The time left in each phase of the loop is indicated by a count-down timer, from 18 to zero (breathing phase), or from eight to zero seconds (pause breathing).

The main heart sound recording phase will last approximately two minutes.



### 8.10 Recording phase, cancelling

The ongoing recording can be cancelled if conditions for recording heart sounds become unsuitable.

In case of increased external noise (like high traffic noise, helicopter landing, high sounding alarms) or internal noise (like patient coughing, stomach rumbling, increased hall-way talk), the actual recording can be retried without sacrificing a new patch.

To cancel the ongoing recording press the push button twice in quick succession. A triple beep warning will be heard.

A display message (REF:309) will appear asking if cancel is only concerning the current recording or will be to exit recording altogether. Exiting recording will invalidate the patch used.





### 8.11 Recording phase, CAD-scoring

After the heart sound recording has ended the sound-file will be evaluated for errors like excessive noise, and a CAD-score is calculated (REF:124).

The sensor/patch should not be removed from the patient at this moment, in case of a new recording has to be made due to recording errors.

The calculation time is approximately 1-2 minutes. The CAD-score will be shown on the display in a colored bar indicating the patient risk group, with the corresponding Record-ID. Note the CAD-score result and Record-ID in the patient's medical file.

The sensor and patch can now be gently removed by pulling patch ear (REF:125). The sensor is released by pulling the string of the coupling ring in the direction of the arrow.



### 8.12 Short diastole CAD-score

In a few cases a standard CAD-score cannot be calculated based on a minimum diastolic period length, however an estimated CAD-score will be shown. The short diastole CAD-score is indicated by an asterisk (\*) in front of the CAD-score.

When a short diastole-based CAD-score is presented, a possibility to redo the recording is offered. Try to relax (lower the heart rate) the patient further before attempting a new recording. A maximum of four recording attempts are allowed per patch. It is not recommended to pursue recording after four attempts on the same patient.

A short diastole-based CAD-score is not as precise as the standard CAD-score, which is why prompting for a new recording is encouraged (REF:128).



### *8.13 QR-code for printing, archiving and sending a CADScor patient report*

A functionality of the CADScor® System is the integration of a QR-code for printing, archiving and sending a patient report from a previous CAD-scoring, using a mobile device like a smart phone, or a QR code scanner set up for the purpose.

Currently, mobile devices running iOS system, at least version 13, with a camera can be used. The CADScor app can be downloaded from the Apple App Store and should be installed/configured prior to use.

On the result screen of the CADScor® Sensor, press the CAD-score result and the screen will convert to a QR-code. Point the activated CADScor App in “Scan Report” mode to transfer the data to the mobile device.

The QR-code functionality can be disabled/enabled by your local sales agent.



### 8.14 CAD-score in the 30-39 year age group

- ⚠ The CAD-scores from patients in the 30-39-year group are currently outside the intended patient population. A warning triangle indicates the higher uncertainty of their CAD-score.

REF: 125 REC-ID: 3091 20 REC.

CAD-score: 19 ⚠ Low risk

Record ID: 3091

Sensor: S15-07-1223

Sensor can be removed

DONE

### *8.15 Possible CADScor® System messages during or after recording*

The CADScor® Sensor will during a pre-recording, prior to the heart sound recording sequence, determine if ambient conditions are fulfilled.

If recording conditions are fulfilled the heart sound recording will continue.

In case of unsuitable recording conditions or other errors, a display message will be shown indicating the cause of error.

Depending on error or message type the CADScor® Sensor may prompt for an additional recording or aborting recording altogether.

In case of failure to obtain a CAD-score result, other patient evaluations should be pursued.

A comprehensive list of errors and messages is listed in the troubleshooting section 11.





### 8.16 Recalling a previous recording result

From the Settings display (REF:500) it is possible to select a LOOK-UP ID function to recall data from a previous recording session if the Record-ID (REC-ID) is known. Only REC-IDs with a resulting CAD-score can be found.

Within the LOOK-UP ID display (REF:504) enter the four-digit REC-ID number and press the OK button to continue.

Notice that the REC-ID is relating to the specific CADScor® Sensor that performed the recording and the subsequent CAD-score calculation.  
REC-IDs from cancelled recordings cannot be looked up.

Also notice that the REC-ID counter will reset to 1000 after REC-ID 9999.  
Previous recording REC-ID data will thus be overwritten with new data having the same REC-ID number.



Check in the LOOK-UP data that the LOOK-UP REC-ID recording date corresponds to the anticipated date for the recording that is being identified.


The other patients' parameters (e.g. gender and age-group), can also be used for establishing the relevance of the Look-up data to the patient (REF:506).

If the REC-ID is not found (REF:507), check if the entered REC-ID number or sensor number is correct.

If the LOOK-UP function is used on a different sensor than the original recording sensor, the LOOK-UP can display data relating to another patient, from a different recording session.

If the CADScor device is used in a blinded clinical investigation, the CAD-score will not be shown in the LOOK-UP display.

Press EXIT to return to settings menu.



### 8.17 After the recording

Clean the sensor using 70% ethanol wipes after use.


The sensor will automatically turn itself off.


Place the sensor in the Docking station for charging and automatic qualification until next recording session.

Do not use running tap- or splashing water to clean the sensor.

If the sensor has become accidentally wet, wipe it dry and leave to dry out for minimum 24 hours.

Do not attempt to dry in oven or using dry blower.

 Do not expose the CADScor® System to liquids (water, oils, detergents or similar) or dust, since this may damage the CADScor® System permanently.

 The CADScor® System cannot be heat- or radiation sterilized, or machine washed/cleaned, since this may damage the CADScor® System permanently.

## 9 Sound guidance used in the interface

During operation of the CADScor® System sounds are indicating the various events.

Some sounds are related to the timing of breathing/pause breathing intervals in the recording sequence and other sounds to confirm actions or as a warning of intended action.

| <b>Tone</b>                         | <b>Indication</b>  |
|-------------------------------------|--|
| Three tones in ascending frequency  | Recording sequence start.  |
| Single tone                         | Initiation of breathing pause.   |
| Double tone                         | Ending of breathing pause.   |
| Three tones in descending frequency | Recording sequence ended.  |
| 2x double tone                      | Completion of CAD-score calculation.   |
| Single beep                         | Forced qualification initiated after 4+ seconds push button activation.  |
| Double beep                         | Forced reset after 8+ seconds push button activation. Will turn sensor off, and require qualification after restart. |
| Triple beep                         | Warning after double-click in recording sequence to abort recording.   |




## 10 Maintenance

There are no user serviceable parts within the CADScor® System.

Maintenance is limited to cleaning the external surfaces of the sensor and the Docking station.

Only use 70% ethanol wipes for cleaning, do not use running tap water, splashing water or waterjets, or oils, detergents or similar substances.

The CADScor® System cannot be machine washed or autoclaved nor heat or radiation sterilized.

-  Do not expose the CADScor® System to liquids (water, oils, detergents or similar) or dust, since this may damage the CADScor® System permanently.
-  The CADScor® System cannot be heat- or radiation sterilized, or machine washed/cleaned, since this may damage the CADScor® System permanently.
-  Do not modify the CADScor® System or use or repair a defect CADScor® System, due to risk of malfunctioning. The CADScor® System must only be serviced by qualified Acarix personnel.

## *10.1 Disposal*

### **Patch**

The patch can be disposed of in ordinary waste.


### **Sensor and Docking station**

The CADScor® System, when reaching its end of life, must be collected and recycled separately from other waste according to national requirements. Please contact your local Acarix distributor for instructions.

Acarix and its distributors within the European Union and associated states have taken the necessary steps to comply with the directive, 2012/19/EU on waste electrical and electronic equipment (WEEE)

Waste electrical and Electronic Equipment environmental implications: WEEE contains materials that are potentially hazardous to the environment and to human health.

## 11 Troubleshooting guide

| Error message  | Display reference                                       | Possible cause  | Correction action  |
|--|---|---|--|
| # 1<br>The CAD-score has a *-symbol in front of it.  | REF:102,<br>REF:103,<br>REF:125,<br>REF:128,<br>REF:506 | The diastolic periods in the recorded acoustic patient file are shorter than optimal for precise CAD-score calculation, due to higher heart rate. See section 8.12. | Try to relax patient to lower heart rate below 80-90 bpm and retry recording.<br><br>Heart rates above 80-90 may result in error message "Irregular heart beat" see # 9 (REF:209; REF:210).  |
| # 2<br>The CAD-score has a warning symbol behind it.  | REF:102,<br>REF:103,<br>REF:125,<br>REF:128,<br>REF:506 | If a CAD-score is calculated from a patient below 40 years (outside intended patient group), the CAD-score is presented with a warning symbol.                      | Establish if the entered patient information was correct.<br><br>Be aware that the patient group below 40 years of age is outside current intended population group, and therefore conclude cautiously.                              |
| # 3<br>Sensor restarted due to battery supply error or Low battery power   | REF:306<br>REF:308                                      | Defect in battery.<br>The sensor battery power level is too low.<br>The sensor battery has been deep discharged.<br>See also # 20.                                  | Place the sensor in the Docking station for charging > 15 REC.<br>Set system for recharging, may take 1-2 days for recharging after deep discharge.<br>Contact local distributor/Acarix for available options if problem persisting. |

| <b>Error message</b>          | <b>Display reference</b> | <b>Possible cause</b>  | <b>Correction action</b>   |
|-------------------------------|--------------------------|--|--|
| # 4<br>Qualification expired  | REF:300                  | The previous qualification has expired.  | Qualify the sensor in the Docking station.   |
| # 5<br>Qualification failed   | REF:301                  | Incorrect Docking station.   | Place sensor in home Docking station.  |
|                               | REF:303                  | Leak between the sensor head and the Docking station sensor tester.<br><br>Error in the Docking station.<br>Error in the sensor. | Place the sensor head firmly into the sensor tester and observe that red indicator ring on the sensor head is not visible and retry qualification.<br>Contact local distributor/Acarix if sensor will not qualify.   |
| # 6<br>Ambient noise too high | REF:203<br>REF:204       | The acoustic noise level in the examination room was too high, above 65 dB SPL.  | Close doors and windows to shield from acoustic noise coming from outside recording room.<br>Keep silent during the recording.<br>Turn off running nearby fans, radios or similar.<br>Change examination room location.<br><br>Other examples of noise:<br>Construction/building machinery, road & rail traffic noise, emergency helicopter service and similar. |



| <b>Error message</b>             | <b>Display reference</b> | <b>Possible cause</b>   | <b>Correction action</b>  |
|----------------------------------|--------------------------|---|---|
| # 7<br>Electrical noise too high | REF:205<br>REF:206       | The electromagnetic field noise in the examination room was too high. | Unplug power (if any) to the examination bed.<br>Turn off unused electronic equipment in the examination room.<br>Turn off wireless electrical devices (e.g. cell phones).<br>Change examination room location. |
| # 8<br>Internal noise too high   | REF:207<br>REF:208       | Can be caused by intermittent stomach rumbling.                       | Typically stomach rumbling will occur in periods prior to eating as a consequence of air in the intestines. A minor meal prior to measurement may resolve frequent stomach rumble.                              |
| # 9<br>Irregular heartbeat       | REF:209<br>REF:210       | Patient heart rate was too irregular or too high.                     | Try to relax patient to lower heart rate below 80-90 bpm and retry recording.<br>Consider evaluating patient for arrhythmia.<br>Patient may already be diagnosed having known arrhythmia.                       |

| <b>Error message</b>                               | <b>Display reference</b>                            | <b>Possible cause</b>  | <b>Correction action</b>  |
|--|---|--|---|
| # 10<br>Heartbeat signal too low                   | REF:211<br>REF:212                                  | The sensor head was not anchored correctly at the left IC4 position.<br>The sensor head was incorrectly placed in coupling ring.<br>Heart action weak. | Re-identify left IC4 space and use a new patch for positioning correctly at left IC4.<br><br>Secure sensor appropriately in coupling ring using patch assembly tool.<br>Abort if recording is not possible. |
| # 11<br>The recording stopped after pre-recording. | REF:203<br>REF:205<br>REF:207<br>REF:209<br>REF:211 | Ambient noise or electrical noise was too high.<br>The heartbeat signal was too low.<br>The heartbeat signal was irregular.                            | See all individual causes above, # 6-10.<br><br>Reduce ambient- and electrical noise.<br><br>Guide the patient to relax to lower the heart rate, below 80-90 bpm.   |
| # 12<br>Inconsistent data analysis                 | REF:214<br>REF:215                                  | Data collected outside normal algorithm pattern.<br>Combination of heart related sounds not compatible with current algorithm.                         | Guide the patient to relax to lower the heart rate, below 80-90 bpm.  |

| <b>Error message</b>  | <b>Display reference</b>      | <b>Possible cause</b>  | <b>Correction action</b>  |
|---|-------------------------------|--|---|
| # 13<br>The patch was not accepted;<br>Patch is expired,<br>or has been used previously,<br>or<br>Patch RFID reading was disturbed. | REF:200<br>REF:201<br>REF:202 | The patch has past expiration date or was used previously.<br><br>Maximum storage period is 48 months from production date.  | Use a new non-expired or un-used patch.<br>Reading a patch RFID on a metal-table (or close to metal) may interfere with proper patch reading.<br><br>Do a retry patch RFID reading away from metal containing surfaces. |
| # 14<br>Temperature or pressure out of range  | REF:304<br>REF:313<br>REF:314 | The temperature in the room has exceeded the operating conditions for the CADScor®System.<br><br>The CADScor®System has been brought outside its height limitations. | Avoid exposing the CADScor®System to excessive heat, or direct sunlight.<br><br>Observe maximum operating height 2500 meters above sea level.   |
| # 15<br>Algorithm interaction failed  | REF:213<br>REF:218            | A technical problem inside the CADScor®Sensor has occurred.  | Retry recording.<br>If persistent failure contact local distributor/Acarix service.   |

| <b>Error message</b>                                 | <b>Display reference</b> | <b>Possible cause</b>   | <b>Correction action</b>  |
|--|--------------------------|---|---|
| # 16<br>Ambient microphone error                     | REF:216                  | The test signal received by the ambient microphone was lower than anticipated.<br>Error was triggered by liquid exposure to sensor. | Take care not to cover sensor or ambient microphone during recording.<br>Ensure CADScor®Sensor body is oriented along sternum towards chin.<br>Allow sensor to dry > 24 hours in Docking station.<br>If persistent failure contact local distributor/Acarix service.        |
| # 17<br>Self-test failed                             | REF:315                  | A technical problem inside the CADScor®System has occurred.   | Contact local distributor/Acarix service.   |
| # 18<br>Unexpected issue detected. Restart required. | REF:317                  | An internal error has occurred. Sensor needs to be restarted.<br>The USB connection between Docking station and Sensor failed.      | Place sensor in docking station and press "EXIT".<br>Wait for restart to finish.<br>Inspect charging pins in Docking station for free up and down movement. Exercise pin gently to free up stuck pin.<br>Clean charging pads on sensor using cotton swabs with 70% ethanol. |

| Error message            | Display reference | Possible cause                                  | Correction action  |
|--------------------------|-------------------|---|--|
| # 19<br>REC-ID not found | REF:507           | The patient REC-ID was not found on the Sensor. | <p>Check that the REC-ID number was typed in correctly.</p> <p>Check that the Look-up function was performed on the Sensor that made the original recording.</p> <p>Check if the last REC-ID number on the Sensor display is lower than the Look-up REC-ID. If the Sensor REC-ID number is lower, the Look-up use-data text-file may have been overwritten by new data.</p> <p>The sensor may have been memory reset as part of service; Recording ID has been erased.</p> <p>Notice that REC-ID's that have been cancelled do not appear in REC-ID log.</p> |

| <b>Error message</b>   | <b>Display reference</b> | <b>Possible cause</b>  | <b>Correction action</b>  |
|--|--------------------------|--|---|
| # 20<br>The LED indicator on Docking station is off, or constant yellow when Sensor is docked. |                          | The power adaptor has been turned off at mains.<br>The power jack has been disconnected.<br>Plug assembly not fully inserted into adaptor.<br>The power adaptor is faulty, or electrical supply is missing.<br>The Docking station contacts interface is not recognizing sensor. | Turn on power adaptor outlet.<br>Apply power jack supply to the docking station.<br>Establish full insertion of plug into power adaptor by pressing firmly.<br>Replace power supply unit if suspected malfunctioning.<br>Restore electrical supply.<br><br>Try restart procedure (7.14).<br>Contact local distributor/Acarix for available options if problem persists. |
| # 21<br>The display is black.  |                          | The sensor is in off mode.<br>The sensor battery is discharged.<br><br>See also # 3.   | Press push button to switch on sensor.<br>Charge the sensor battery in the Docking station.<br>For deep discharged batteries recharging may take 1-2 days.<br>It is recommended to always keep power on system, for ease of use and battery charging.   |
| # 22<br>Display cracked  |                          | The display is broken.   | Contact local distributor/Acarix for available options.   |
| # 23<br>Sensor unresponsive  |                          | Software error.  | Sensor needs restart. See section 7.14  |

## 12 CADScor® System requirements

### Normal operating environment:

Temperature +10 to +40°C.

Relative humidity 20-80% non-condensing.

Below 2500 m above sea level.

### Short term storage/ transportation environment (up to 96 hours):

*In packaging-box/shipper-box.*

Maximum -30 to +70°C.

Maximum 93% RH condensing, at +40°C.

Minimum ambient pressure 59,7 kPa (4.267 m above sea level).

### Normal storage environment:

Temperature +10 to +40°C.

Relative humidity 20-80% non-condensing.

### Ingress Protection (IP) ratings:

CADScor® Sensor IP44

CADScor® Docking station IP22

Power adaptors IP20

### Operating noise environment:

External noise level < 65dB SPL.

External electrical noise level < 65dB SPL.

### System altitude correction

Up to 2500 m above sea level

### System power source

Only the supplied power adaptor must be used, identified by a manufacturer name and type number:

XP-Power ACM06US05 (-XZ1505)

Rated input voltage: 100-240VAC, 50/60 Hz

Rated output voltage: 5VDC

Max. power: 6.0 W



Never use other power adaptor than supplied with the CADScor® System, due to risk of electrical shock.

## 12.1 System specification

Model: ACS-1401

Dimensions:

CADScor® Sensor: ASE-1401

75(W) x 160(L) x 26(H) mm

Approx. weight: 180 g

CADScor® Docking station: ADS-1401

193(W) x 109(D) x 53(H) mm

Approx. weight: 600 g

XP Power, Power adaptor: Excl. wire

42(W) x 70(D) x 69(H) mm

Approx. weight: 110 g (incl. wire, EU plug)

CADScor® Patch: ACP-1401

80(W) x 190(L) x 9(H) mm

Approx. weight: 10 g

### **RFID:**

Frequency band: 13.56 MHz ISM band

Communication standard: ISO14443

Tx power: Magnetic field short range device  
with < 0 dB  $\mu$ A/m @10m

Modulation: ASK

Bandwidth: 14 kHz



Do not use the CADScor® System on patients with implanted electronics like ICD, Pacemakers, heart-pumps or closer than 50 cm to similar active electronic support equipment, due to risk of equipment failure from CADScor® System electromagnetic RFID impulse.

### **Date of manufacture:**

The date of manufacture can be read from the serial [SN] numbers on the sensor or Docking station respectively, e.g: S14-02-0001 is manufactured Febr. 2014, and D14-04-0001 is manufactured April, 2014. The four last digits are lot numbers.



## *12.2 Packaging materials*

System Box: 1x

Outside box made of recyclable corrugated cardboard. Internal fixtures made from recyclable folded corrugated cardboard.

Internal boxes made from recyclable cardboard.

User manual made of 80g/m<sup>2</sup> paper.

Separator tool made from stainless steel.

Total approx. weight, incl. CADScor® System and power adaptor: 1500 g

System Shipper-box: 1x

Outside box made of recyclable corrugated cardboard.

Approx. weight: 120 g

System Shipper-box: 6x

Outside box made of recyclable corrugated cardboard.

Approx. weight: 580 g

Patch box: 1x

Outside box made of recyclable cardboard. Inside pouches made from Polyethylene and paper.

Assembly tool made from cardboard.

Total approx. weight, incl. patches: 240 g

Patch Shipper-box: 1x

Outside box made of recyclable corrugated cardboard.

Approx. weight: 100 g

Patch Shipper-box: 6x

Outside box made of recyclable corrugated cardboard.

Approx. weight: 580 g

## 13 Warranty

The CADScor® System is covered by a one (1) year warranty (subject to country regulations), replacing/correcting defective parts or system parts, or at least 1000 recordings in the warranty period.

Opening of the CADScor® Sensor or modifying the CADScor® System invalidates the warranty. The CADScor® Docking station can only be opened for wall mount installation. No serviceable parts inside.


The warranty does not cover any damages caused by improper operation.


This system comprises sensitive components and must be treated carefully to avoid strong vibrations or shock from dropping or handling.

Observe storage and operating conditions described in CADScor® System requirements, for temperature, altitude and humidity in section 12.

The manufacturer cannot be held liable for any damage caused by incorrect application or operation of CADScor® System.

The warranty does not extend to interpretation of results obtained using the CADScor® System or use of the system outside its intended use.

 Do not modify the CADScor® System or use or repair a defect CADScor® System, due to risk of malfunctioning. The CADScor® System must only be serviced by qualified Acarix personnel (10).

 Do not drop or exert excessive stress or force to the CADScor® System since this may damage the CADScor® System permanently.

## 14 Approvals and EMC information

The CADScor® System is a medical device class IIa, fulfilling the requirements in compliance with the Council Directive 93/42/EEC, Medical Device Directive (MDD).

The CADScor® System software is classified as safety class B, according to IEC/EN 62304, “Medical device software – Software life-cycle processes”

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

CADScor sensor is compliant with the requirement for RF exposure in US with less than

50 mm separation distance between the user and/or bystander of the device.

CADScor sensor is only compliant if no changes or modifications are made to the device.

|              |   |
|--------------|---|
| CE-mark      | Indicates compliance with the <i>Medical Device Directive</i> 93/42/EEC                         |
| IEC/EN 60601 | In compliance with IEC/EN 60601-1   |
| EMC-emission | The equipment complies with the emission requirements for Class B equipment in IEC/EN 60601-1-2 |
| EMC-immunity | The equipment complies with the immunity requirements in IEC/EN 60601-1-2                       |
| FCC ID       | 2AYXI-ASE1401   |

| <b>Table 1: Guidance and manufacturer's declaration – electromagnetic emissions</b>   |                   |  |
|---|-------------------|--|
| The CADScor system ACS-1401 is intended for use in the electromagnetic environment specified below. The customer or the user of the CADScor system ACS-1401 should assure that it is used in such an environment. |                   |  |
| <b>Emissions test</b>   | <b>Compliance</b> | <b>Electromagnetic environment – guidance</b>  |
| RF emissions<br>CISPR 11  | Group 1           | The CADScor system ACS-1401 must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.<br><br>The CADScor system ACS-1401 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| RF emissions<br>CISPR 11  | Class B           |  |
| Harmonic<br>emissions<br>IEC 61000-3-2  | Class A           |  |
| Voltage<br>fluctuations/<br>flicker emissions<br><br>IEC 61000-3-3  | Complies          |  |

| <b>Table 2: Guidance and manufacturer's declaration – electromagnetic immunity</b>  |  |   |  |
|---|--|---|--|
| The CADScor system ACS-1401 is intended for use in the electromagnetic environment specified below. The customer or the user of the CADScor system ACS-1401 should assure that it is used in such an environment. |  |   |  |
| <b>Immunity test</b>  | <b>IEC/ EN 60601-1-2 test level</b>                            | <b>Compliance level</b>                                       | <b>Electromagnetic environment – guidance</b>  |
| Electrostatic Discharge, (ESD)<br><br>IEC 61000-4-2   | 8 kV contact<br><br>15 kV air                                  | 8 kV contact<br><br>15 kV air                                 | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrical fast transient/burst<br><br>IEC 61000-4-4  | 2 kV for power supply lines<br><br>1 kV for input/output lines | 2 kV for power supply lines<br><br>N/A, no input/output lines | Mains power quality should be that of a typical commercial or hospital environment.  |
| Surge<br><br>IEC 61000-4-5  | 1 kV line(s) to line(s)<br><br>2 kV line(s) to earth           | 1 kV line(s) to line(s)<br><br>2 kV line(s) to earth          | Mains power quality should be that of a typical commercial or hospital environment.  |

| <b>Table 2, continued:</b><br><b>Guidance and manufacturer's declaration – electromagnetic immunity</b>    |  |  |   |
|--|--|--|---|
| Voltage dips, short interruptions and voltage variations on power supply input lines<br><br>IEC 61000-4-11 | 0% <i>UT</i><br>(100% dip in <i>UT</i> ) for 0.5 cycle<br><br>0% <i>UT</i><br>(100% dip in <i>UT</i> ) for 1 cycle<br><br>70 % <i>UT</i><br>(30 % dip in <i>UT</i> ) for 25 cycles<br><br>0% <i>UT</i><br>(100% dip in <i>UT</i> ) for 5 sec | 0% <i>UT</i><br>(100% dip in <i>UT</i> ) for 0.5 cycle<br><br>0% <i>UT</i><br>(100% dip in <i>UT</i> ) for 1 cycle<br><br>70 % <i>UT</i><br>(30 % dip in <i>UT</i> ) for 25 cycles<br><br>0% <i>UT</i><br>(100% dip in <i>UT</i> ) for 5 sec | Mains power quality should be that of a typical commercial or hospital environment.   |
| Power frequency (50/60 Hz) magnetic field<br><br>IEC 61000-4-8   | 30 A/m   | 30 A/m   | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| <b>NOTE</b> <i>UT</i> is the a.c. mains voltage prior to application of the test level.                    |  |  |   |

| Table 2, continued:<br>Guidance and manufacturer's declaration – electromagnetic immunity   |  |  |   |
|---|--|--|---|
| Conducted RF immunity<br><br>IEC 61000-4-6  | 3 Vrms<br>150 kHz to 80 MHz.<br>80 % AM at 1 kHz<br><br>6 Vrms in ISM bands<br>between 150 kHz and 80 MHz (see note) | 3 Vrms<br>150 kHz to 80 MHz.<br>80 % AM at 1 kHz<br><br>6 Vrms in ISM bands<br>between 150 kHz and 80 MHz (see note) | WARNING:<br><br>Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the CADScor system ACS-1401, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. |
| Radiated RF immunity<br><br>IEC 61000-4-3   | 3 V/m<br>80 MHz to 2.7 GHz<br>80 % AM at 1 kHz<br><br>Plus additional test levels according to the below table.      | 3 V/m<br>80 MHz to 2.7 GHz<br>80 % AM at 1 kHz<br><br>Plus, additional test levels according to the below table.     |   |
| NOTE:<br>The ISM (industrial, scientific and medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. |  |  |   |

**Table 3: Radiated Immunity Test Levels (Immunity to RF Wireless Communication Equipment).**

| Test frequency (MHz) | Band <sup>(1)</sup> (MHz) | Service <sup>(1)</sup>   | Modulation <sup>(1)</sup>                                | Maximum power (W) | Distance (m) | Immunity TEST LEVEL (V/m) |
|----------------------|---------------------------|--|--|-------------------|--------------|---------------------------|
| 385                  | 380 – 390                 | TETRA 400  | Pulse modulation <sup>(1)</sup><br>18 Hz                 | 1,8               | 0,3          | 27                        |
| 450                  | 430 – 470                 | GPRS 460,<br>FRS 460   | FM <sup>(1)</sup><br>$\pm 5$ kHz deviation<br>1 kHz sine | 2                 | 0,3          | 28                        |
| 710                  | 704 – 787                 | LTE Band 13,<br>17   | Pulse modulation <sup>(1)</sup><br>217 Hz                | 0,2               | 0,3          | 9                         |
| 745                  |                           |  |  |                   |              |                           |
| 780                  |                           |  |  |                   |              |                           |
| 810                  | 800 – 960                 | GSM 800/900,<br>TETRA 800,<br>iDEN 820,<br>CDMA 850,<br>LTE Band 5             | Pulse modulation <sup>(1)</sup><br>18 Hz                 | 2                 | 0,3          | 28                        |
| 870                  |                           |  |  |                   |              |                           |
| 930                  |                           |  |  |                   |              |                           |
| 1 720                | 1 780 –<br>1 990          | GSM 1800,<br>CDMA 1900,<br>GSM 1900,<br>DECT,<br>LTE Band 1, 3,<br>4, 25, UMTS | Pulse modulation <sup>(1)</sup><br>217 Hz                | 2                 | 0,3          | 28                        |
| 1 845                |                           |  |  |                   |              |                           |
| 1 970                |                           |  |  |                   |              |                           |
| 2 450                | 2 400 –<br>2 570          | Bluetooth,<br>WLAN,<br>802.11 b/g/n,<br>RFID 2450,<br>LTE Band 7               | Pulse modulation <sup>(1)</sup><br>217 Hz                | 2                 | 0,3          | 28                        |
| 5 240                | 5 100 –<br>5 800          | WLAN 802.11<br>a/n   | Pulse modulation <sup>(1)</sup><br>217 Hz                | 0,2               | 0,3          | 9                         |
| 5 500                |                           |  |  |                   |              |                           |
| 5 785                |                           |  |  |                   |              |                           |



## 15 Acronyms & Definitions in text

| Acronyms                             | Definitions   |
|--------------------------------------|---|
| CAD                                  | Coronary Artery Disease.  |
| CSS                                  | Chronic Coronary Syndrome   |
| CADScor <sup>®</sup> System          | Sensor, Docking station and Patch.  |
| CADScor <sup>®</sup> Sensor          | Acoustic recording sensor, part of the CADScor <sup>®</sup> System.   |
| CADScor <sup>®</sup> Docking station | Docking station for charging and qualification of the sensor, part of the CADScor <sup>®</sup> system.          |
| CADScor <sup>®</sup> Patch           | Patch for anchoring the sensor onto the chest of the patient, part of the CADScor <sup>®</sup> System           |
| Recording                            | Acoustic recording performed by the sensor.   |
| IC4-L                                | Left fourth intercostal space.  |
| Stenosis                             | Narrowing of the coronary arteries.   |
| Coronary murmurs                     | Heart sounds that are produced as a result of turbulent blood flow that is sufficient to produce audible noise. |
| CAD-score                            | Calculated heart murmur score (of coronary murmurs) on the basis of the acoustic recording.                     |
| RFID-chip                            | Radio-frequency identification chip   |
| Sensor qualification                 | A procedure to ensure sensor operational status.  |
| External noise                       | The ambient noise surrounding the heart sensor.   |
| Internal noise                       | The sounds arising from inside of human body (respiratory, muscle, bowel).                                      |
| Electromagnetic noise                | Non-acoustic noise generated by electromagnetic fields (EMC).   |

## 16 Licenses

The CADScor® System (r) makes use of several open source packages. Here is a list with package name, the name and version of the licenses used and, in some cases, other notes. Afterwards, the licenses will be written.

- alsa-lib, LGPL-2.1
- alsa-utils, GPL-2
- busybox, GPL-2
- dropbear, several licenses, see below
- freetype, GPL-2
- glib, LGPL-2.1
- i2c-tools, GPL-2
- jpeg, the jpeg library is work of the Independent JPEG Group.
- libpng, see <http://www.libpng.org>
- libsndfile, LGPL-2.1
- linux kernel, GPL-2
- memtester, GPL-2
- netbase, GPL-2
- pcre, BSD-3-Clause, see copyright holders below
- qt4-embedded, LGPL-2.1
- u-boot, GPL-2
- wireless-tools, GPL-2 and LGPL-2.1
- zlib, <http://zlib.net>

Some of the licenses grant the user full access to the source code for the specific package. To get this access, please contact Acarix at [info@acarix.com](mailto:info@acarix.com).

dropbear

The majority of code is written by Matt Johnston, under the MIT license below. Portions of the client-mode work are (c) 2004 Mihnea Stoenescu, under the same license

Copyright (c) 2002-2008 Matt Johnston

Portions copyright (c) 2004 Mihnea Stoenescu.

All rights reserved.

< The MIT License, see below>

=====

LibTomCrypt and LibTomMath are written by Tom St Denis, and are Public Domain.

=====

sshpty.c is taken from OpenSSH 3.5p1,

Copyright (c) 1995 Tatu Ylonen <ylo@cs.hut.fi>, Espoo, Finland

All rights reserved

"As far as I am concerned, the code I have written for this software can be used freely for any purpose. Any derived versions of this software must be clearly marked as such, and if the derived work is incompatible with the protocol description in the RFC file, it must be called by a name other than "ssh" or "Secure Shell". "

=====

loginrec.c, loginrec.h, atomicio.h, atomicio.c and strlcat() (included in util.c) are from OpenSSH 3.6.1p2, and are licensed under the 2 point BSD license (see below).

loginrec is written primarily by Andre Lucas, atomicio.c by Theo de Raadt.

strlcat() is (c) Todd C. Miller

=====

Import code in keyimport.c is modified from PuTTY's import.c, licensed as follows:

PuTTY is copyright 1997-2003 Simon Tatham.

The MIT License (MIT)

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

PCRE

PCRE is a library of functions to support regular expressions whose syntax and semantics are as close as possible to those of the Perl 5 language.

Release 8 of PCRE is distributed under the terms of the "BSD" licence, as specified below. The documentation for PCRE, supplied in the "doc" directory, is distributed under the same terms as the software itself.

The basic library functions are written in C and are freestanding. Also included in the distribution is a set of C++ wrapper functions, and a just-in-time compiler that can be used to optimize pattern matching. These are both optional features that can be omitted when the library is built.

#### THE BASIC LIBRARY FUNCTIONS

Written by: Philip Hazel

Email local part: ph10

Email domain: cam.ac.uk

University of Cambridge Computing Service, Cambridge, England.

Copyright (c) 1997-2011 University of Cambridge

All rights reserved.

#### PCRE JUST-IN-TIME COMPILATION SUPPORT

Written by: Zoltan Herczeg

Email local part: hzmester

Email domain: freemail.hu

Copyright(c) 2010-2011 Zoltan Herczeg

All rights reserved.

#### STACK-LESS JUST-IN-TIME COMPILER

Written by: Zoltan Herczeg

Email local part: hzmester

Email domain: freemail.hu

Copyright(c) 2009-2011 Zoltan Herczeg

All rights reserved.

#### THE C++ WRAPPER FUNCTIONS

Contributed by: Google Inc.

Copyright (c) 2007-2011, Google Inc.

All rights reserved.

Neither the name of the University of Cambridge nor the name of Google Inc. nor the names of their contributors may be used to endorse or promote products derived from this software without specific prior written permission.

<BSD 2-clause License, see below> (Note the above line + the BSD 2-Clause license make it a BSD 3-Clause license).

#### The BSD 2-Clause

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.

Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT OWNER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

#### LGPL-2.1

##### GNU LESSER GENERAL PUBLIC LICENSE

Version 2.1, February 1999

Copyright (C) 1991, 1999 Free Software Foundation, Inc. 51 Franklin Street, Fifth Floor, Boston, MA 02110-1301 USA

Everyone is permitted to copy and distribute verbatim copies of this license document, but changing it is not allowed.

[This is the first released version of the Lesser GPL. It also counts as the successor of the GNU Library Public License, version 2, hence the version number 2.1.]

#### Preamble

The licenses for most software are designed to take away your freedom to share and change it. By contrast, the GNU General Public Licenses are intended to guarantee your freedom to share and change free software--to make sure the software is free for all its users.

This license, the Lesser General Public License, applies to some specially designated software packages--typically libraries--of the Free Software Foundation and other authors who decide to use it. You can use it too, but we suggest you first think carefully about

whether this license or the ordinary General Public License is the better strategy to use in any particular case, based on the explanations below.

When we speak of free software, we are referring to freedom of use, not price. Our General Public Licenses are designed to make sure that you have the freedom to distribute copies of free software (and charge for this service if you wish); that you receive source code or can get it if you want it; that you can change the software and use pieces of it in new free programs; and that you are informed that you can do these things.

To protect your rights, we need to make restrictions that forbid distributors to deny you these rights or to ask you to surrender these rights. These restrictions translate to certain responsibilities for you if you distribute copies of the library or if you modify it. For example, if you distribute copies of the library, whether gratis or for a fee, you must give the recipients all the rights that we gave you. You must make sure that they, too, receive or can get the source code. If you link other code with the library, you must provide complete object files to the recipients, so that they can relink them with the library after making changes to the library and recompiling it. And you must show them these terms so they know their rights.

We protect your rights with a two-step method: (1) we copyright the library, and (2) we offer you this license, which gives you legal permission to copy, distribute and/or modify the library.

To protect each distributor, we want to make it very clear that there is no warranty for the free library. Also, if the library is modified by someone else and passed on, the recipients should know that what they have is not the original version, so that the original author's reputation will not be affected by problems that might be introduced by others.

Finally, software patents pose a constant threat to the existence of any free program. We wish to make sure that a company cannot effectively restrict the users of a free program by obtaining a restrictive license from a patent holder. Therefore, we insist that any patent license obtained for a version of the library must be consistent with the full freedom of use specified in this license. Most GNU software, including some libraries, is covered by the ordinary GNU General Public License. This license, the GNU Lesser General Public License, applies to certain designated libraries, and is quite different from the ordinary General Public License. We use this license for certain libraries in order to permit linking those libraries into non-free programs.

When a program is linked with a library, whether statically or using a shared library, the combination of the two is legally speaking a combined work, a derivative of the original library. The ordinary General Public License therefore permits such linking only if the entire combination fits its criteria of freedom. The Lesser General Public License permits more lax criteria for linking other code with the library.

We call this license the "Lesser" General Public License because it does Less to protect the user's freedom than the ordinary General Public License. It also provides other free software developers Less of an advantage over competing non-free programs. These disadvantages are the reason we use the ordinary General Public License for many libraries. However, the Lesser license provides advantages in certain special circumstances.

For example, on rare occasions, there may be a special need to encourage the widest possible use of a certain library, so that it becomes a de-facto standard. To achieve this, non-free programs must be allowed to use the library. A more frequent case is that

a free library does the same job as widely used non-free libraries. In this case, there is little to gain by limiting the free library to free software only, so we use the Lesser General Public License.

In other cases, permission to use a particular library in non-free programs enables a greater number of people to use a large body of free software. For example, permission to use the GNU C Library in non-free programs enables many more people to use the whole GNU operating system, as well as its variant, the GNU/Linux operating system.

Although the Lesser General Public License is less protective of the users' freedom, it does ensure that the user of a program that is linked with the Library has the freedom and the wherewithal to run that program using a modified version of the Library.

The precise terms and conditions for copying, distribution and modification follow. Pay close attention to the difference between a "work based on the library" and a "work that uses the library". The former contains code derived from the library, whereas the latter must be combined with the library in order to run.

#### TERMS AND CONDITIONS FOR COPYING, DISTRIBUTION AND MODIFICATION

**0.** This License Agreement applies to any software library or other program which contains a notice placed by the copyright holder or other authorized party saying it may be distributed under the terms of this Lesser General Public License (also called "this License"). Each licensee is addressed as "you".

A "library" means a collection of software functions and/or data prepared so as to be conveniently linked with application programs (which use some of those functions and data) to form executables.

The "Library", below, refers to any such software library or work which has been distributed under these terms. A "work based on the Library" means either the Library or any derivative work under copyright law: that is to say, a work containing the Library or a portion of it, either verbatim or with modifications and/or translated straightforwardly into another language. (Hereinafter, translation is included without limitation in the term "modification".)

"Source code" for a work means the preferred form of the work for making modifications to it. For a library, complete source code means all the source code for all modules it contains, plus any associated interface definition files, plus the scripts used to control compilation and installation of the library.

Activities other than copying, distribution and modification are not covered by this License; they are outside its scope. The act of running a program using the Library is not restricted, and output from such a program is covered only if its contents constitute a work based on the Library (independent of the use of the Library in a tool for writing it). Whether that is true depends on what the Library does and what the program that uses the Library does.

**1.** You may copy and distribute verbatim copies of the Library's complete source code as you receive it, in any medium, provided that you conspicuously and appropriately publish on each copy an appropriate copyright notice and disclaimer of warranty; keep intact all the notices that refer to this License and to the absence of any warranty; and distribute a copy of this License along with the Library.

You may charge a fee for the physical act of transferring a copy, and you may at your option offer warranty protection in exchange for a fee.

**2.** You may modify your copy or copies of the Library or any portion of it, thus forming a work based on the Library, and copy and distribute such modifications or work under the terms of Section 1 above, provided that you also meet all of these conditions:

**a)** The modified work must itself be a software library.

**b)** You must cause the files modified to carry prominent notices stating that you changed the files and the date of any change.

**c)** You must cause the whole of the work to be licensed at no charge to all third parties under the terms of this License.

**d)** If a facility in the modified Library refers to a function or a table of data to be supplied by an application program that uses the facility, other than as an argument passed when the facility is invoked, then you must make a good faith effort to ensure that, in the event an application does not supply such function or table, the facility still operates, and performs whatever part of its purpose remains meaningful. (For example, a function in a library to compute square roots has a purpose that is entirely well-defined independent of the application. Therefore, Subsection 2d requires that any application-supplied function or table used by this function must be optional: if the application does not supply it, the square root function must still compute square roots.)

These requirements apply to the modified work as a whole. If identifiable sections of that work are not derived from the Library, and can be reasonably considered independent and separate works in themselves, then this License, and its terms, do not apply to those sections when you distribute them as separate works. But when you distribute the same sections as part of a whole which is a work based on the Library, the distribution of the whole must be on the terms of this License, whose permissions for other licensees extend to the entire whole, and thus to each and every part regardless of who wrote it.

Thus, it is not the intent of this section to claim rights or contest your rights to work written entirely by you; rather, the intent is to exercise the right to control the distribution of derivative or collective works based on the Library.

In addition, mere aggregation of another work not based on the Library with the Library (or with a work based on the Library) on a volume of a storage or distribution medium does not bring the other work under the scope of this License.

**3.** You may opt to apply the terms of the ordinary GNU General Public License instead of this License to a given copy of the Library. To do this, you must alter all the notices that refer to this License, so that they refer to the ordinary GNU General Public License, version 2, instead of to this License. (If a newer version than version 2 of the ordinary GNU General Public License has appeared, then you can specify that version instead if you wish.) Do not make any other change in these notices.

Once this change is made in a given copy, it is irreversible for that copy, so the ordinary GNU General Public License applies to all subsequent copies and derivative works made from that copy.

This option is useful when you wish to copy part of the code of the Library into a program that is not a library.

**4.** You may copy and distribute the Library (or a portion or derivative of it, under Section 2) in object code or executable form under the terms of Sections 1 and 2 above provided that you accompany it with the complete corresponding machine-readable source code, which must be distributed under the terms of Sections 1 and 2 above on a medium customarily used for software interchange.



If distribution of object code is made by offering access to copy from a designated place, then offering equivalent access to copy the source code from the same place satisfies the requirement to distribute the source code, even though third parties are not compelled to copy the source along with the object code.

**5.** A program that contains no derivative of any portion of the Library, but is designed to work with the Library by being compiled or linked with it, is called a "work that uses the Library". Such a work, in isolation, is not a derivative work of the Library, and therefore falls outside the scope of this License.

However, linking a "work that uses the Library" with the Library creates an executable that is a derivative of the Library (because it contains portions of the Library), rather than a "work that uses the library". The executable is therefore covered by this License.

Section 6 states terms for distribution of such executables.

When a "work that uses the Library" uses material from a header file that is part of the Library, the object code for the work may be a derivative work of the Library even though the source code is not. Whether this is true is especially significant if the work can be linked without the Library, or if the work is itself a library. The threshold for this to be true is not precisely defined by law.

If such an object file uses only numerical parameters, data structure layouts and accessors, and small macros and small inline functions (ten lines or less in length), then the use of the object file is unrestricted, regardless of whether it is legally a derivative work. (Executables containing this object code plus portions of the Library will still fall under Section 6.)

Otherwise, if the work is a derivative of the Library, you may distribute the object code for the work under the terms of Section 6. Any executables containing that work also fall under Section 6, whether or not they are linked directly with the Library itself.

**6.** As an exception to the Sections above, you may also combine or link a "work that uses the Library" with the Library to produce a work containing portions of the Library, and distribute that work under terms of your choice, provided that the terms permit modification of the work for the customer's own use and reverse engineering for debugging such modifications.

You must give prominent notice with each copy of the work that the Library is used in it and that the Library and its use are covered by this License. You must supply a copy of this License. If the work during execution displays copyright notices, you must include the copyright notice for the Library among them, as well as a reference directing the user to the copy of this License. Also, you must do one of these things:

**a)** Accompany the work with the complete corresponding machine-readable source code for the Library including whatever changes were used in the work (which must be distributed under Sections 1 and 2 above); and, if the work is an executable linked with the Library, with the complete machine-readable "work that uses the Library", as object code and/or source code, so that the user can modify the Library and then relink to produce a modified executable containing the modified Library. (It is understood that the user who changes the contents of definitions files in the Library will not necessarily be able to recompile the application to use the modified definitions.)

**b)** Use a suitable shared library mechanism for linking with the Library. A suitable mechanism is one that (1) uses at run time a copy of the library already present on the user's computer system, rather than copying library functions into the executable, and

(2) will operate properly with a modified version of the library, if the user installs one, as long as the modified version is interface-compatible with the version that the work was made with.

**c)** Accompany the work with a written offer, valid for at least three years, to give the same user the materials specified in Subsection 6a, above, for a charge no more than the cost of performing this distribution.

**d)** If distribution of the work is made by offering access to copy from a designated place, offer equivalent access to copy the above specified materials from the same place.

**e)** Verify that the user has already received a copy of these materials or that you have already sent this user a copy.

For an executable, the required form of the "work that uses the Library" must include any data and utility programs needed for reproducing the executable from it. However, as a special exception, the materials to be distributed need not include anything that is normally distributed (in either source or binary form) with the major components (compiler, kernel, and so on) of the operating system on which the executable runs, unless that component itself accompanies the executable.

It may happen that this requirement contradicts the license restrictions of other proprietary libraries that do not normally accompany the operating system. Such a contradiction means you cannot use both them and the Library together in an executable that you distribute.

**7.** You may place library facilities that are a work based on the Library side-by-side in a single library together with other library facilities not covered by this License, and distribute such a combined library, provided that the separate distribution of the work based on the Library and of the other library facilities is otherwise permitted, and provided that you do these two things:

**a)** Accompany the combined library with a copy of the same work based on the Library, uncombined with any other library facilities. This must be distributed under the terms of the Sections above.

**b)** Give prominent notice with the combined library of the fact that part of it is a work based on the Library, and explaining where to find the accompanying uncombined form of the same work.

**8.** You may not copy, modify, sublicense, link with, or distribute the Library except as expressly provided under this License. Any attempt otherwise to copy, modify, sublicense, link with, or distribute the Library is void, and will automatically terminate your rights under this License. However, parties who have received copies, or rights, from you under this License will not have their licenses terminated so long as such parties remain in full compliance.

**9.** You are not required to accept this License, since you have not signed it. However, nothing else grants you permission to modify or distribute the Library or its derivative works. These actions are prohibited by law if you do not accept this License. Therefore, by modifying or distributing the Library (or any work based on the Library), you indicate your acceptance of this License to do so, and all its terms and conditions for copying, distributing or modifying the Library or works based on it.

**10.** Each time you redistribute the Library (or any work based on the Library), the recipient automatically receives a license from the original licensor to copy, distribute, link with or modify the Library subject to these terms and conditions. You may not impose any further restrictions on the recipients' exercise of the rights granted herein. You are not responsible for enforcing compliance by third parties with this License.

**11.** If, as a consequence of a court judgment or allegation of patent infringement or for any other reason (not limited to patent issues), conditions are imposed on you (whether by court order, agreement or otherwise) that contradict the conditions of this License, they do not excuse you from the conditions of this License. If you cannot distribute so as to satisfy simultaneously your obligations under this License and any other pertinent obligations, then as a consequence you may not distribute the Library at all. For example, if a patent license would not permit royalty-free redistribution of the Library by all those who receive copies directly or indirectly through you, then the only way you could satisfy both it and this License would be to refrain entirely from distribution of the Library.

If any portion of this section is held invalid or unenforceable under any particular circumstance, the balance of the section is intended to apply, and the section as a whole is intended to apply in other circumstances.

It is not the purpose of this section to induce you to infringe any patents or other property right claims or to contest validity of any such claims; this section has the sole purpose of protecting the integrity of the free software distribution system which is implemented by public license practices. Many people have made generous contributions to the wide range of software distributed through that system in reliance on consistent application of that system; it is up to the author/donor to decide if he or she is willing to distribute software through any other system and a licensee cannot impose that choice.

This section is intended to make thoroughly clear what is believed to be a consequence of the rest of this License.

**12.** If the distribution and/or use of the Library is restricted in certain countries either by patents or by copyrighted interfaces, the original copyright holder who places the Library under this License may add an explicit geographical distribution limitation excluding those countries, so that distribution is permitted only in or among countries not thus excluded. In such case, this License incorporates the limitation as if written in the body of this License.

**13.** The Free Software Foundation may publish revised and/or new versions of the Lesser General Public License from time to time. Such new versions will be similar in spirit to the present version, but may differ in detail to address new problems or concerns. Each version is given a distinguishing version number. If the Library specifies a version number of this License which applies to it and "any later version", you have the option of following the terms and conditions either of that version or of any later version published by the Free Software Foundation. If the Library does not specify a license version number, you may choose any version ever published by the Free Software Foundation.

**14.** If you wish to incorporate parts of the Library into other free programs whose distribution conditions are incompatible with these, write to the author to ask for permission. For software which is copyrighted by the Free Software Foundation, write to the Free Software Foundation; we sometimes make exceptions for this. Our decision will be guided by the two goals of preserving the free status of all derivatives of our free software and of promoting the sharing and reuse of software generally.

## **NO WARRANTY**

**15.** BECAUSE THE LIBRARY IS LICENSED FREE OF CHARGE, THERE IS NO WARRANTY FOR THE LIBRARY, TO THE EXTENT PERMITTED BY APPLICABLE LAW. EXCEPT WHEN OTHERWISE STATED IN WRITING THE COPYRIGHT HOLDERS AND/OR OTHER

PARTIES PROVIDE THE LIBRARY "AS IS" WITHOUT WARRANTY OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE ENTIRE RISK AS TO THE QUALITY AND PERFORMANCE OF THE LIBRARY IS WITH YOU. SHOULD THE LIBRARY PROVE DEFECTIVE, YOU ASSUME THE COST OF ALL NECESSARY SERVICING, REPAIR OR CORRECTION.

**16.** IN NO EVENT UNLESS REQUIRED BY APPLICABLE LAW OR AGREED TO IN WRITING WILL ANY COPYRIGHT HOLDER, OR ANY OTHER PARTY WHO MAY MODIFY AND/OR REDISTRIBUTE THE LIBRARY AS PERMITTED ABOVE, BE LIABLE TO YOU FOR DAMAGES, INCLUDING ANY GENERAL, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE USE OR INABILITY TO USE THE LIBRARY (INCLUDING BUT NOT LIMITED TO LOSS OF DATA OR DATA BEING RENDERED INACCURATE OR LOSSES SUSTAINED BY YOU OR THIRD PARTIES OR A FAILURE OF THE LIBRARY TO OPERATE WITH ANY OTHER SOFTWARE), EVEN IF SUCH HOLDER OR OTHER PARTY HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

## END OF TERMS AND CONDITIONS

### How to Apply These Terms to Your New Libraries

If you develop a new library, and you want it to be of the greatest possible use to the public, we recommend making it free software that everyone can redistribute and change. You can do so by permitting redistribution under these terms (or, alternatively, under the terms of the ordinary General Public License).

To apply these terms, attach the following notices to the library. It is safest to attach them to the start of each source file to most effectively convey the exclusion of warranty; and each file should have at least the "copyright" line and a pointer to where the full notice is found.

*one line to give the library's name and an idea of what it does.*

Copyright (C) *year* *name of author*

This library is free software; you can redistribute it and/or modify it under the terms of the GNU Lesser General Public License as published by the Free Software Foundation; either version 2.1 of the License, or (at your option) any later version.

This library is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU Lesser General Public License for more details.

You should have received a copy of the GNU Lesser General Public License along with this library; if not, write to the Free Software Foundation, Inc., 51 Franklin Street, Fifth Floor, Boston, MA 02110-1301 USA

Also add information on how to contact you by electronic and paper mail.

You should also get your employer (if you work as a programmer) or your school, if any, to sign a "copyright disclaimer" for the library, if necessary. Here is a sample; alter the names:

Yoyodyne, Inc., hereby disclaims all copyright interest in the library `Frob' (a library for tweaking knobs) written by James Random Hacker.

*signature of Ty Coon*, 1 April 1990

Ty Coon, President of Vice  
That's all there is to it!

## GPL-2

### GNU GENERAL PUBLIC LICENSE

Version 2, June 1991

Copyright (C) 1989, 1991 Free Software Foundation, Inc. 51 Franklin Street, Fifth Floor, Boston, MA 02110-1301, USA

Everyone is permitted to copy and distribute verbatim copies of this license document, but changing it is not allowed.

#### Preamble

The licenses for most software are designed to take away your freedom to share and change it. By contrast, the GNU General Public License is intended to guarantee your freedom to share and change free software--to make sure the software is free for all its users. This General Public License applies to most of the Free Software Foundation's software and to any other program whose authors commit to using it. (Some other Free Software Foundation software is covered by the GNU Lesser General Public License instead.) You can apply it to your programs, too.

When we speak of free software, we are referring to freedom, not price. Our General Public Licenses are designed to make sure that you have the freedom to distribute copies of free software (and charge for this service if you wish), that you receive source code or can get it if you want it, that you can change the software or use pieces of it in new free programs; and that you know you can do these things.

To protect your rights, we need to make restrictions that forbid anyone to deny you these rights or to ask you to surrender the rights. These restrictions translate to certain responsibilities for you if you distribute copies of the software, or if you modify it. For example, if you distribute copies of such a program, whether gratis or for a fee, you must give the recipients all the rights that you have. You must make sure that they, too, receive or can get the source code. And you must show them these terms so they know their rights.

We protect your rights with two steps: (1) copyright the software, and (2) offer you this license which gives you legal permission to copy, distribute and/or modify the software.

Also, for each author's protection and ours, we want to make certain that everyone understands that there is no warranty for this free software. If the software is modified by someone else and passed on, we want its recipients to know that what they have is not the original, so that any problems introduced by others will not reflect on the original authors' reputations.

Finally, any free program is threatened constantly by software patents. We wish to avoid the danger that redistributors of a free program will individually obtain patent licenses, in effect making the program proprietary. To prevent this, we have made it clear that any patent must be licensed for everyone's free use or not licensed at all.

The precise terms and conditions for copying, distribution and modification follow.

## TERMS AND CONDITIONS FOR COPYING, DISTRIBUTION AND MODIFICATION

**0.** This License applies to any program or other work which contains a notice placed by the copyright holder saying it may be distributed under the terms of this General Public License. The "Program", below, refers to any such program or work, and a "work based on the Program" means either the Program or any derivative work under copyright law: that is to say, a work containing the Program or a portion of it, either verbatim or with modifications and/or translated into another language. (Hereinafter, translation is included without limitation in the term "modification".) Each licensee is addressed as "you".

Activities other than copying, distribution and modification are not covered by this License; they are outside its scope. The act of running the Program is not restricted, and the output from the Program is covered only if its contents constitute a work based on the Program (independent of having been made by running the Program). Whether that is true depends on what the Program does.

**1.** You may copy and distribute verbatim copies of the Program's source code as you receive it, in any medium, provided that you conspicuously and appropriately publish on each copy an appropriate copyright notice and disclaimer of warranty; keep intact all the notices that refer to this License and to the absence of any warranty; and give any other recipients of the Program a copy of this License along with the Program.

You may charge a fee for the physical act of transferring a copy, and you may at your option offer warranty protection in exchange for a fee.

**2.** You may modify your copy or copies of the Program or any portion of it, thus forming a work based on the Program, and copy and distribute such modifications or work under the terms of Section 1 above, provided that you also meet all of these conditions:

**a)** You must cause the modified files to carry prominent notices stating that you changed the files and the date of any change.

**b)** You must cause any work that you distribute or publish, that in whole or in part contains or is derived from the Program or any part thereof, to be licensed as a whole at no charge to all third parties under the terms of this License.

**c)** If the modified program normally reads commands interactively when run, you must cause it, when started running for such interactive use in the most ordinary way, to print or display an announcement including an appropriate copyright notice and a notice that there is no warranty (or else, saying that you provide a warranty) and that users may redistribute the program under these conditions, and telling the user how to view a copy of this License. (Exception: if the Program itself is interactive but does not normally print such an announcement, your work based on the Program is not required to print an announcement.)

These requirements apply to the modified work as a whole. If identifiable sections of that work are not derived from the Program, and can be reasonably considered independent and separate works in themselves, then this License, and its terms, do not apply to those sections when you distribute them as separate works. But when you distribute the same sections as part of a whole which is a work based on the Program, the distribution of the whole must be on the terms of this License, whose permissions for other licensees extend to the entire whole, and thus to each and every part regardless of who wrote it.

Thus, it is not the intent of this section to claim rights or contest your rights to work written entirely by you; rather, the intent is to exercise the right to control the distribution of derivative or collective works based on the Program.

In addition, mere aggregation of another work not based on the Program with the Program (or with a work based on the Program) on a volume of a storage or distribution medium does not bring the other work under the scope of this License.

**3.** You may copy and distribute the Program (or a work based on it, under Section 2) in object code or executable form under the terms of Sections 1 and 2 above provided that you also do one of the following:

**a)** Accompany it with the complete corresponding machine-readable source code, which must be distributed under the terms of Sections 1 and 2 above on a medium customarily used for software interchange; or,

**b)** Accompany it with a written offer, valid for at least three years, to give any third party, for a charge no more than your cost of physically performing source distribution, a complete machine-readable copy of the corresponding source code, to be distributed under the terms of Sections 1 and 2 above on a medium customarily used for software interchange; or,

**c)** Accompany it with the information you received as to the offer to distribute corresponding source code. (This alternative is allowed only for noncommercial distribution and only if you received the program in object code or executable form with such an offer, in accord with Subsection b above.)

The source code for a work means the preferred form of the work for making modifications to it. For an executable work, complete source code means all the source code for all modules it contains, plus any associated interface definition files, plus the scripts used to control compilation and installation of the executable. However, as a special exception, the source code distributed need not include anything that is normally distributed (in either source or binary form) with the major components (compiler, kernel, and so on) of the operating system on which the executable runs, unless that component itself accompanies the executable.

If distribution of executable or object code is made by offering access to copy from a designated place, then offering equivalent access to copy the source code from the same place counts as distribution of the source code, even though third parties are not compelled to copy the source along with the object code.

**4.** You may not copy, modify, sublicense, or distribute the Program except as expressly provided under this License. Any attempt otherwise to copy, modify, sublicense or distribute the Program is void, and will automatically terminate your rights under this License. However, parties who have received copies, or rights, from you under this License will not have their licenses terminated so long as such parties remain in full compliance.

**5.** You are not required to accept this License, since you have not signed it. However, nothing else grants you permission to modify or distribute the Program or its derivative works. These actions are prohibited by law if you do not accept this License. Therefore, by modifying or distributing the Program (or any work based on the Program), you indicate your acceptance of this License to do so, and all its terms and conditions for copying, distributing or modifying the Program or works based on it.

**6.** Each time you redistribute the Program (or any work based on the Program), the recipient automatically receives a license from the original licensor to copy, distribute or modify the Program subject to these terms and conditions. You may not impose any

further restrictions on the recipients' exercise of the rights granted herein. You are not responsible for enforcing compliance by third parties to this License.

**7.** If, as a consequence of a court judgment or allegation of patent infringement or for any other reason (not limited to patent issues), conditions are imposed on you (whether by court order, agreement or otherwise) that contradict the conditions of this License, they do not excuse you from the conditions of this License. If you cannot distribute so as to satisfy simultaneously your obligations under this License and any other pertinent obligations, then as a consequence you may not distribute the Program at all. For example, if a patent license would not permit royalty-free redistribution of the Program by all those who receive copies directly or indirectly through you, then the only way you could satisfy both it and this License would be to refrain entirely from distribution of the Program.

If any portion of this section is held invalid or unenforceable under any particular circumstance, the balance of the section is intended to apply and the section as a whole is intended to apply in other circumstances.

It is not the purpose of this section to induce you to infringe any patents or other property right claims or to contest validity of any such claims; this section has the sole purpose of protecting the integrity of the free software distribution system, which is implemented by public license practices. Many people have made generous contributions to the wide range of software distributed through that system in reliance on consistent application of that system; it is up to the author/donor to decide if he or she is willing to distribute software through any other system and a licensee cannot impose that choice.

This section is intended to make thoroughly clear what is believed to be a consequence of the rest of this License.

**8.** If the distribution and/or use of the Program is restricted in certain countries either by patents or by copyrighted interfaces, the original copyright holder who places the Program under this License may add an explicit geographical distribution limitation excluding those countries, so that distribution is permitted only in or among countries not thus excluded. In such case, this License incorporates the limitation as if written in the body of this License.

**9.** The Free Software Foundation may publish revised and/or new versions of the General Public License from time to time. Such new versions will be similar in spirit to the present version, but may differ in detail to address new problems or concerns. Each version is given a distinguishing version number. If the Program specifies a version number of this License which applies to it and "any later version", you have the option of following the terms and conditions either of that version or of any later version published by the Free Software Foundation. If the Program does not specify a version number of this License, you may choose any version ever published by the Free Software Foundation.

**10.** If you wish to incorporate parts of the Program into other free programs whose distribution conditions are different, write to the author to ask for permission. For software which is copyrighted by the Free Software Foundation, write to the Free Software Foundation; we sometimes make exceptions for this. Our decision will be guided by the two goals of preserving the free status of all derivatives of our free software and of promoting the sharing and reuse of software generally.

## **NO WARRANTY**



11. BECAUSE THE PROGRAM IS LICENSED FREE OF CHARGE, THERE IS NO WARRANTY FOR THE PROGRAM, TO THE EXTENT PERMITTED BY APPLICABLE LAW. EXCEPT WHEN OTHERWISE STATED IN WRITING THE COPYRIGHT HOLDERS AND/OR OTHER PARTIES PROVIDE THE PROGRAM "AS IS" WITHOUT WARRANTY OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE ENTIRE RISK AS TO THE QUALITY AND PERFORMANCE OF THE PROGRAM IS WITH YOU. SHOULD THE PROGRAM PROVE DEFECTIVE, YOU ASSUME THE COST OF ALL NECESSARY SERVICING, REPAIR OR CORRECTION.

12. IN NO EVENT UNLESS REQUIRED BY APPLICABLE LAW OR AGREED TO IN WRITING WILL ANY COPYRIGHT HOLDER, OR ANY OTHER PARTY WHO MAY MODIFY AND/OR REDISTRIBUTE THE PROGRAM AS PERMITTED ABOVE, BE LIABLE TO YOU FOR DAMAGES, INCLUDING ANY GENERAL, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE USE OR INABILITY TO USE THE PROGRAM (INCLUDING BUT NOT LIMITED TO LOSS OF DATA OR DATA BEING RENDERED INACCURATE OR LOSSES SUSTAINED BY YOU OR THIRD PARTIES OR A FAILURE OF THE PROGRAM TO OPERATE WITH ANY OTHER PROGRAMS), EVEN IF SUCH HOLDER OR OTHER PARTY HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

END OF TERMS AND CONDITIONS

#### How to Apply These Terms to Your New Programs

If you develop a new program, and you want it to be of the greatest possible use to the public, the best way to achieve this is to make it free software which everyone can redistribute and change under these terms.

To do so, attach the following notices to the program. It is safest to attach them to the start of each source file to most effectively convey the exclusion of warranty; and each file should have at least the "copyright" line and a pointer to where the full notice is found.

*one line to give the program's name and an idea of what it does.*

Copyright (C) yyyy name of author

This program is free software; you can redistribute it and/or modify it under the terms of the GNU General Public License as published by the Free Software Foundation; either version 2 of the License, or (at your option) any later version.

This program is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU General Public License for more details.

You should have received a copy of the GNU General Public License along with this program; if not, write to the Free Software Foundation, Inc., 51 Franklin Street, Fifth Floor, Boston, MA 02110-1301, USA.

Also add information on how to contact you by electronic and paper mail.

If the program is interactive, make it output a short notice like this when it starts in an interactive mode:

Gnomovision version 69, Copyright (C) year name of author Gnomovision comes with ABSOLUTELY NO WARRANTY; for details type `show w'. This is free software, and you are welcome to redistribute it under certain conditions; type `show c' for details.

The hypothetical commands ``show w'` and ``show c'` should show the appropriate parts of the General Public License. Of course, the commands you use may be called something other than ``show w'` and ``show c'`; they could even be mouse-presss or menu items--whatever suits your program.

You should also get your employer (if you work as a programmer) or your school, if any, to sign a "copyright disclaimer" for the program, if necessary. Here is a sample; alter the names:

Yoyodyne, Inc., hereby disclaims all copyright interest in the program ``Gnomovision'` (which makes passes at compilers) written by James Hacker.

*signature of Ty Coon*, 1 April 1989

Ty Coon, President of Vice

This General Public License does not permit incorporating your program into proprietary programs. If your program is a subroutine library, you may consider it more useful to permit linking proprietary applications with the library. If this is what you want to do, use the GNU Lesser General Public License instead of this License.



US-FDA version.

CADScor®System User Manual, English

US-FDA revision 12.5, Mar. 20<sup>th</sup>, 2021.

From software version 4.0 US-FDA.