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Shenzhen Jamr Technology Co., Ltd.
A101-301, D101-201, Jamr Science & Technology Park, No. 2 Guiyuan Road,
Guixiang Community, Guanlan Street, Longhua District, Shenzhen 518100,
PEOPLE'S REPUBLIC OF CHINA

Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany.
shholding@hotmail.com

Company Name: SUNGO Certification Company Limited
Address: 3rd Floor, 70 Gracechurch street, London, EC3V 0HR, United Kingdom
Tel: 0044-20-3608 6516
Email: Jackie.zhu20@outlook.com/ ukrp@sungoglobal.com

Version Number:1.0

VIMIMK

Fully Automatic Upper Arm Blood Pressure Monitor

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CONTENTS

1.Introduction and Intended Use.....	2
2.Important Information on Blood Pressure and its Measurement.....	5
3.Components of Your Blood Pressure Monitor.....	6
4.Using Your Monitor for the First Time.....	8
5.Measurement Procedure.....	10
6.Care and Maintenance.....	17
7.Warranty/service.....	18
8.Certifications.....	18
9.Technical Specifications.....	18
10.EMC Declaration.....	19

1.Introduction and Intended Use

The Blood Pressure Monitor enables reliable measurement of systolic and diastolic blood pressure as well as pulse by the oscillometric method.

Before using, please read this instruction manual carefully and then keep it in a safe place for easy access.

1.1 Remember...

- Only a health-care professional is qualified to interpret blood pressure measurements.
- This device is NOT intended to replace regular medical checkups.
- Blood pressure readings obtained by this device should be verified before prescribing or making adjustments to any medications used to control hypertension.
- Under no circumstances should YOU alter the dosages of any drugs prescribed by your physician.
- This device is intended for adults only. Consult with a physician before using this device on a child.
- In cases of irregular heartbeat , measurements made with this device should only be evaluated after consultation with a physician.
- The main product, including its accessories, shall be processed in accordance with local regulations after reaching to the life cycle.

1.2 Warnings and Precautions

- Warning:** The use of other accessories other than those specified or provided by the equipment manufacturer may cause electromagnetic radiation to increase or decrease electromagnetic immunity resulting in operational failure
- Warning:** This system may fail to yield specified measurement accuracy if operated or stored in temperature or humidity conditions outside the limits stated in the specifications section of this manual.
- Warning:** The separate ac adapter which is intended to connect USB interface of Blood Pressure Monitor has not been evaluated according to IEC 60601-1. The safety of the product shall be reappraised when it power supply by a separate ac adapter.
- Warning:** The user must check that the device functions safely and see that it is in proper working condition before being used.
- Warning:** The device is not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.
- Warning:** The patient is an intended operator, the functions of monitoring blood pressure and pulse rate can be safely used by patient. The routine clean and changing batteries can be performed by the patient.
- Warning:** This device can not be used together with hf surgical equipment.
- Warning:** Use of power adapters
1. Adapter: input 100-240V, 50/60hz output DC 5V 1A
 2. Do not to position the device to make it difficult to operate the disconnection device while using adaptor.
 3. Do not be prone to water leakage, high temperature, moisture, direct sunlight and more or more corrosive gas environment. Do not use this device in the above environment.
- Warning:** Too frequent measurements can cause injury to the PATIENT due to blood flow interference.
- Warning:** Don't place the cuff over the wound part.
- Warning:** Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring ME EQUIPMENT on the same limb.
- Caution:** To avoid any possibility of accidental strangulation, keep this device away from children and do not drape tubing around your neck.

- Caution:** To avoid damaging the device, keep this device away from children and pets.
- Caution:** The standard material used for the bladder and tubing is latex-free.
- Attention:** Self-measurement means self-monitor/control instead of diagnosis or treatment. Unusual values must always be discussed with a physician. Under no circumstances should you alter the dosages of any drugs prescribed by your physician.
- Attention:** The pulse display is not suitable for checking the frequency of heart pacemakers!
- Attention:** In cases of irregular heartbeat, the measurements on this device should only be evaluated after consultation with a physician.
- Note:** To obtain the greatest accuracy from the device, it is recommended that the device should be used within the specified temperature and the relative humidity, please check the Technical Specification.
- Note:** The cuff is treated as the applied part. The user should contact the manufacturer for assistance, if needed, in setting up, using or maintaining the device.
- Note:** This device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens) during use. These can lead to erratic results.
- Note:** Do not attempt to service or repair this device by yourself. Should a malfunction occur, refer to local distributor or the manufacturer.

2.Important Information on Blood Pressure and its Measurement


2.1. How does high or low blood pressure arise?

Your level of blood pressure is determined in the circulatory center of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and speed of the heart (Pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is controlled by fine muscles in the blood vessel walls.

Your level of arterial blood pressure changes periodically during heart activity: During the "blood ejection" (Systole) the value is highest (systolic blood pressure value). At the end of the heart's "rest period" (Diastole) pressure is lowest (diastolic blood pressure value).

2.2. Which values are normal?

Please refer to the diagram below(Picture-01)


Picture-01


There are six grids in the display of device. Please refer to the picture-01-01. Different grids represent different interval scales of WHO.

Blood pressure value	WHO grids in device	WHO Classification
DIA<80 & SYS<120	1	Optimal blood pressure
DIA<85 & SYS<130	2	Normal blood pressure
DIA<90 & SYS<140	3	High normal value
DIA<100 & SYS<160	4	Mild hypertension
DIA<110 & SYS<180	5	Moderate hypertension
DIA>=110 or SYS>=180	6	Severe hypertension

picture-01-01


Picture-04


Picture-05

3. Components of your blood pressure monitor

3.1. Measuring unit

Picture-02

Picture-03

3.2 The symbols on the LED display

- 1.User A
- 2.User B
- 3.Irregular heartbeat symbol
- 4.Movement error symbol
- 5.Mute symbol
- 6.Cuff self-detecting function
- 7.Low battery symbol
- 8.WHO function symbol
- 9.Heartbeat symbol (Flashes during measurement)
- 10.Average value symbol
- 11.Memory symbol
- 12.Pulse symbol
- 13.Diastolic blood pressure
- 14.Systolic blood pressure

3.3 Features of Model C04

- 1.Cuff self-checking function
2. Double users: 2 x 120 sets memory
3. Average value function
4. Irregular heartbeat checking
5. WHO function
6. Low battery display
7. External power adapter support
8. Auto power-off

4. Using your Monitor for the First Time


4.1 Two Power Models for Choice

1. Battery Operated

Battery Installation/ Replacement

Use only 1.5V "AA" alkaline batteries with this device.

- Press the hook on the bottom of the battery cover and lift the cover off in the direction of the arrow (Picture-06).
- Install 3 "AA" size batteries and make sure the + (positive) and - (negative) polarities match the polarities of the battery compartment.
- Close the battery cover. Make sure that the battery cover is securely in position.

Picture-06

Low Battery Indicator

- When the Low Battery Indicator appears on the display, turn the device off and remove all the batteries. Replace with 3 new batteries at the same time. Long-life alkaline batteries are recommended.
- To prevent the damage of device from leaked battery fluid, please take out of battery if the device unused in a long time (generally more than 3 months.) If battery fluid get in your eyes, immediately rinse with plenty of clean water. Contact a physician immediately.
- Dispose of the device, components and optional accessories according to applicable local regulations. Unlawful disposal may cause environmental pollution.

2.AC Adapter Operated

You may also operate this device using the AC adapter (output 5V DC1A with Type-C connector).

Use only the approved AC adapter to avoid damaging the device.

- Ensure that the AC adapter and cable are not damaged.
- Plug the adapter cable into the AC adapter port on the right side of the device.
- Plug the adapter into your electrical outlet. When the AC adapter is connected, no battery current is consumed.

Note: No power is taken from the batteries while the AC adapter is connected to the device. If electrical power is interrupted, (e.g., by accidental removal of the AC adapter from the outlet) the device must be reset by removing the plug from the socket and reinserting the AC adapter connection.

4.2. System Settings

Setting the User:

After you load the battery or connect power for the device, press the SET button to select User B (IA). Then Press the START/ STOP button to confirm User B (IA).

Deleting the Record:

When you are checking the measurement memories of User A/B, long press the MEM button to delete all the existing measurement data of User A/B.

Note:

You can't delete a single measurement record. If Deleting Record is done, all the records of User A/B will be deleted. Please keep the records in another way in case that you may need it later.

5.Measurement Procedure

5.1. Before Measurement:


- Avoid eating and smoking as well as all forms of exertion directly before measurement. These factors will influence the measurement result.
- Relax by sitting in an armchair in a quiet atmosphere for about ten minutes before measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).

5.2. Fitting the Cuff


Please refer to picture-07

- a) Wrap the cuff around your upper left arm. To secure the cuff, wrap it around your arm and press the hook and loop closure together.
- b) The rubber tube should be on the inside of your arm extending downward to your hand. The on the edge of the cuff (Artery Mark) must lie over the artery which runs down the inner side of the arm.
- c) Make certain the cuff lies approximately 2 to 3 cm above the elbow. You should be able to fit 2 fingers between your arm and the cuff to keep some free space for inflating.
- d) Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked.

Note: Cuffs that don't fit properly will result in false measurement values. Measure your arm circumference if you are not sure of proper fit.

Picture-07

5.3. Measure Procedure


Refer to picture-08

1. Sit comfortably in a chair with your feet flat on the floor.
2. Select your User ID (A or B).
3. Stretch your arm forward on the desk and keep relaxing, make sure the palm is upturned.
4. Sit still and do not talk or move during the measurement to avoid body movement.
5. After the cuff has been appropriately positioned on the arm and connected to the device, the measurement will begin:
 - a) Press the START/STOP button. The pump begins to inflate the cuff. On the display, the increasing cuff pressure is continually displayed.
 - b) After automatically reaching to an individual pressure, the pump stops and the pressure slowly falls. The pressure is displayed during the measurement.
 - c) When the device detects your pulse, the heart symbol in the display begins to blink.
 - d) When the measurement has been concluded, the measured systolic and diastolic blood pressure values, as well as the pulse will be displayed.
 - e) The measurement results will be displayed until you turn the device off. If no operation for 60 seconds, the device will turn off automatically.
- f) **Cuff self-checking symbol** (☺☹)
The cuff correct symbol(☺) will be displayed if the cuff position is correct, otherwise the wrong symbol(☹) will be displayed.
Please check the cuff again if the cuff wrong symbol(☹) is displayed.
- g) **Movement error symbol** (🚫)
The Movement Error Symbol (🚫) will be displayed if you move your body during the measurement. Please remove the cuff, and wait 2-3 minutes. Reapply the cuff and take another measurement.

NOTE:

Correct Position:

- 1) Comfortably seated
- 2) Legs uncrossed
- 3) Feet flat on the floor
- 4) Back and arm supported
- 5) Middle of the CUFF at the level of the right atrium of the heart

5.4. Irregular Heartbeat Detector

The symbol  indicates that certain pulse irregularities were detected during the measurement. In this case, the result may deviate from your normal basal blood pressure, please repeat the measurement.

We recommend the patient to consult medical advice for the physician

- 1) If pulse irregularities occur during measurement, the irregular heartbeat symbol is displayed after the measurement.
- 2) If the symbol appears more frequently (e.g. several times per week on measurements performed daily) or if it suddenly appears more often than usual.

Note: The device is an oscillometric blood pressure monitor device that also analyzes pulse frequency during measurement. The device can not replace a cardiac examination, but it serves to detect pulse irregularities at an early stage.

5.5. Checking Memories

Each device stores 120 sets measurements for 2 users, totally 240 sets. (User A and B)

Viewing stored values

With the device off, press the MEM button. The first display shows "AVG" means an average of the latest 3 measurements stored in the device.

Then Press the MEM button to view the previous measurement records with Number of measurements. (M 1, M 2, M 3...).

The measurement records are sorted in chronological order and the latest measurement will be marked as M 1.

Note: Measurement records for User A/B are stored and averaged separately. Be certain that you are viewing the measurements for the correct user.

5.6. Discontinuing A Measurement

If it is necessary to interrupt the measuring for any reason (e.g the patient feels unwell), press the Start/Stop button at any time.

Then the device immediately lowers the cuff pressure automatically.

5.7. Error Indicators

The following symbol will appear on the display when measuring abnormal

SYMBOL	CAUSE	CORRECTION
No display appears	Weak battery or improper battery placement	Replace both batteries with new ones. Check the battery installation for proper placement of the battery polarities.
Er 1	Sensor abnormal	Check if the pump is working or not. If it is working, then the problem is sensor abnormal. Please send it to the local distributor.
Er 2	Monitor could not detect pulse wave or cannot calculate the measurement	Check if the air releasing is too slow or not. If it is too slow, please check if there is any dust in the tube plug of the cuff and the port in the device. If yes, please clean and start the measurement again. If no, please send the device back to the local distributor
Er 3	Measurement result is abnormal (SYS $\geq$ 35mmHg, DIA $\geq$ 23mmHg)	Occasionally- repeat measurement/ Always - send it to local distributor
Er 4	Too loose cuff or air leakage (Cannot inflate to 30mmHg within 15s)	Tie the cuff correctly and make sure the air plug is properly inserted in the device
Er 5	The air tube is kinked	position the tube correctly and repeat measurement
Er 6	The sensor is sensing great fluctuation in the pressure	Please keep quiet and don't move
Er 7	The pressure that the sensor sensing is over the limit	The sensor is sensing over-limited pressure
Er 8	The demarcation is incorrect or the device has not been demarcated	Please send back to the local distributor

Trouble removal

Problem	Check	Cause and solutions
No power	Check the battery power	Replace new one
	Check the polarity position	Installation for proper placement of the batteries polarities
No inflation	Whether the plug insert	Insert into the air socket tightly
	Whether the plug broken or leak	Change a new cuff
Err and stop working	Whether move the arm when inflate	Keep the body peaceful
	Check if chatting when measured	Keep quite when measure
Cuff leak	Whether the cuff wrap too loose	Wrap the cuff tightly
	Whether the cuff broken	Change a new cuff
 Please contact the distributor if you can't solve the problem, do not disassemble the unit by yourself!		

5.8. SYMBOL DESCRIPTIONS

The following symbols may appear in this manual, on the Digital Blood Pressure Monitor C04, or on its accessories. Some of the symbols represent standards and compliances associated with the Digital Blood Pressure Monitor C04 and its use.

	Authorized Representative in the European Community
	CE Mark: conforms to essential requirements of the Medical Device Directive 93/42/EEC.
	Date of manufacture.
	Manufacturer
SN	Specifies serial number
	Type BF applied part
	Direct current
	DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.
	Follow instructions for use

	Put up
	Fragile
	Afraid of the rain
	Fear of the sun
	Class II equipment
	Handle gently
	Temperature range
No Sterilize requirement	
Not category AP / APG equipment	
Mode of operation: continuous	

6. Care and Maintenance

Wash hands after each time measurement.

If one device is used by different patients, wash hands before and after each use.

- Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- The cuff contains a sensitive air-tight bubble. Handle this cuff carefully and avoid all types of stress through twisting or buckling.
- Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. The cuff with bladder must not be washed in a dishwasher, clothes washer, or submerged in water.
- Handle the tube carefully. Do not pull on it. Do not allow the tubing to kink and keep it away from sharp edges.
- Do not drop the monitor or treat it roughly in any way. Avoid strong vibrations.
- Never open the monitor! This invalidates the manufacturer's warranty.
- Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

6.1. Accuracy test

Sensitive measuring devices must be checked for accuracy from time to time. We recommend a periodical inspection of your unit by an authorized dealer every 1 year. Please turn to local distributor or the manufacturer.

7. Warranty/Service

Your blood pressure monitor is guaranteed for 1 year against manufacturers' defects for the original purchaser only, from date of purchase. The warranty does not apply to damage caused by improper handling, accidents, professional use, not following the operating instructions or alterations made to the instrument by third parties.

Warranty only applies to the main device and its cuff. All other accessories are not covered by warranty.

There are no user serviceable parts inside. Batteries or damage from old batteries is not covered by the warranty.

8. Certifications

Device standard:

This device is manufactured to meet the European blood pressure monitors: IEC 80601-2-30 / ISO81060-1 / IEC60601-1-11 / IEC60601-1

Electromagnetic compatibility:

Device fulfills the stipulations of the International standard IEC60601-1-2

9. Technical Specifications

Model: C04

Weight: 189g (Batteries and AC adapter are not included)

Display: 62.5*63.5mm **【2.46"x2.5"】** LED Digital Display

Size: 126 (W) x 85 (L) x 57 (H) mm **【4.96"(W)x3.35"(L)x2.24"(H)】**

Accessories: 1xMain Device, 1xCuff, 1xUsers manual, 1xWarranty card

Operating Conditions: Temperature: 5 °C to 40 °C; Humidity : 15% to 93% RH;
Storage And Shipping Conditions: Temperature: -25 °C to +70 °C ; Humidity: ≤ 93% RH;

Atmospheric pressure range: 70kPa–106kPa

Measuring method: Oscillometric

Pressure sensor: Resistive

Measuring range: DIA: 40-130mmHg; SYS: 60-230mmHg

Pulse: 40 to 170 per minute

Cuff pressure display range: <300mmHg

Memory: Automatically stores the last 120 measurements for 2 users (240 measurements in total)

Users: Adult

Measuring resolution: 1 mmHg

Accuracy: Pressure within ± 3 mmHg / pulse ± 5 % of the reading

Power source:

a) 3*AA batteries, 1.5 V

b) AC adapter INPUT: 100-240VAC 50/60HZ OUTPUT: 5V DC 1A

Accessories: Wide range rigid cuff 8.7" – 15.7" (22 - 40 cm)

Automatically power off : 60 seconds

Expected service life of the device and accessories: 5 years

Technical alterations reserved!

10. EMC Declaration

1) *This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.

2) * Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.

3) * Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!

4) * Caution: this machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used

Guidance and manufacture's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U _r ; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0 % U _r ; 1 cycle 70 % U _r ; 25/30 cycle 0% U _r ; 250/300 cycle	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30 A/m 50/60Hz	30 A/m 50/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U _r is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacture's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz 3 V RMS outside the ISM band, 6 V RMS in the ISM and amateur bands 80% AM at 1kHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance d=0.35√p d=1.2√p
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80% AM at 1kHz	10 V/m 80 MHz to 2.7 GHz 80% AM at 1kHz	80MHz to 800MHz: d=1.2√p 800MHzto 2.7GHz: d=2.3√p Where, P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance.  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. A Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. B Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Guidance and manufacture's declaration – electromagnetic emission		
The device is intended for use in the electromagnetic environment specified below. The customer of the user of the device should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The device use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Recommended separation distances between portable and mobile RF communications equipment and the device			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.7 GHz
	$d = 1,2\sqrt{P}$	$d = 1,2\sqrt{P}$	$d = 2,3\sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Guidance and manufacturer's declaration - electromagnetic immunity						
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device, should assure that it is used in such an environment.						
Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{a)}	Maximum power (w)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380-390	TETRA 400	Pulse Modulation ^{b)} 18 Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM ^{c)} ±5 kHz deviation 1 kHz sine	2	0.3	28
710 745 780	704-787	LTE Band 13, 17	Pulse Modulation ^{b)} 217 Hz	0.2	0.3	9
810 870 930	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse Modulation ^{b)} 18 Hz	2	0.3	28
1720 1845 1970	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3 4,25;UMTS	Pulse Modulation ^{b)} 217 Hz	2	0.3	28
2450	2400-2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse Modulation ^{b)} 217 Hz	2	0.3	28
5240 5500 5785	5100-5800	WLAN 802.11 a/n	Pulse Modulation ^{b)} 217 Hz	0.2	0.3	9
NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.						
a) For some services, only the uplink frequencies are included. b) The carrier shall be modulated using a 50% duty cycle square wave signal. c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.						
The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: $E = \frac{6}{d} \sqrt{P}$ Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.						

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