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


Introduction to Parts

How to charge?

Please take the band out in first charge, and charge by corresponding charger (charging method: as shown in below diagram, please place the band in charger, align two touch points of charger, to charge).

Wearing method

1. The watch is best worn behind the ulna styloid process;
2. Adjust its size according to adjustment holes, and fasten the wristband buckle;
3. Sensor should be close to skin, avoiding to be moved

Notices in monitoring heart rate:

- LED light of heart rate on the band back flashes;
- Sensor should be close to skin, avoiding the external light to affect the accuracy of test;
- Please keep clean in testing area, sweat or blot may affect test results.

On/off

1. Have a long press on band for 3 seconds in shutdown and charging state to power on;
2. Switch to shutdown interface in start-up state, and have a long press on shutdown interface for 3 seconds to power off.

Band Operation

1. Slide to light up screen or switch display in start-up state;
2. No operation default five-second screen, users can connect to the mobile APP to modify the screen time;
3. The factory default setting of the band (heart rate monitor for 24 hours). Users can turn on or off in the general settings of APP interface.

How to install APP

Scan the following QR code or download and install "GloryFit" in application markets.


Demands on Device: iOS 8.0 or above; Android 4.4 or above, Bluetooth 4.0 is supported

(APP)

Connect with device

It is required to correct by connecting with APP in first use. The band automatically synchronizes time after connecting successfully, otherwise, the step number and sleep data will be incorrect.

Bluetooth connection:After successful pairing, APP will automatically save band Bluetooth address. When APP is opened or operated in the background, it will automatically search and connect the band.

Data synchronization: manually drop APP main page down to synchronize data; the band can save offline date of 7 days. The more the data is, the longer the time to synchronize will be. It will take not more than 2 seconds to synchronize data. After synchronization, "synchronization finish" will be popped up.

APP functions and settings

Please set personal information after entering APP

My - Personal

- Settable profile photo - sex - age - height - weight, personal information can enhance the accuracy of data.

- Customize personal exercise goals. With reasonable sports planning, you can enhance your body quality

Reminder Function

It is required to open reminder switch in APP setting center, and keep the successful connection state between your mobile phone and the band. Notify by shaking when there is message.

Incoming reminder: the band will shake in case of incoming call, and the corresponding name and tel will be displayed. If the 'call rejection' function has been opened in APP, have a long press on it for 3 seconds to reject the call.

Message reminder: the band will remind you by shaking and displaying icon after receiving messages. Please click to check.

QQ reminder: the band will remind you by shaking and displaying icon after receiving QQ message. Please click to check.

WeChat reminder: the band will remind you by shaking and displaying icon after receiving WeChat message. Please click to check.

Smart alarm clock: you can set three alarm clocks. After successfully setting, it will synchronize to the band; smart alarm clock is silent, and reminds you by appropriate shaking. It can remind you in offline state.

Sedentary reminder: The default reminder interval is one hour; in the connected state, after the user is turned on, if the user does not exercise within one hour, the watch will vibrate to remind the user to exercise.

Other reminders: the band will remind you by shaking and displaying icon after receiving other messages such as emails. Please click to check.

Goals reminders: when the step number reaches the goals set in APP, the band will remind you by shaking and displaying icon.

Main function interface

Main interface

Smart Watch has 6 main interfaces. On the main interface, press and hold for two seconds. After shaking, you can switch the dial left and right. After the APP is paired for the first time and connected, the bracelet will synchronize the time, date and week of the phone.

Information interface

Slide up on the main interface to view the information. The bracelet can save the 8 pieces of information recently pushed by the APP (if there is no message, it will stop at the message interface), click to turn the page, slide to the right to exit to return to the upper interface, Automatically delete already viewed information, there are all delete options at the bottom, slide down Return to live interface.

More Functions

The main interface slides down to enter the function setting interface, and slide up to exit to return to the main interface.

Body temperature interface

Switch to the body temperature interface to start automatic body temperature measurement. Valid value measured after 5 minutes, or you can set the automatic interval measurement in the APP. By default, 37.3 °C is the alarm threshold. You can also set the alarm threshold in the APP.

Status interface

Switch to the status interface to view the current step count, distance, and calorie status. The distance and calories are calculated and displayed according to the current walking steps, the height and weight set by the APP personally.

Heart Rate

Switch to the heart rate interface and start to automatically detect the heart rate. After the test is completed, swipe to the right to exit the heart rate detection. You can also connect to the mobile APP to test on the heart rate interface.If the bracelet is not worn correctly or not worn, the heart rate detection is unsuccessful or does not work. In heart rate detection mode, the battery life of the bracelet will be reduced.

Sports Interface

Switch to the blood pressure interface to start automatic blood pressure measurement, switch the menu to re-measure.

Blood oxygen interface

Switch to the blood oxygen interface to start automatic measurement of blood oxygen, and switch the menu to re-measure.

Motion interface

Switch to the sports interface, click this interface to enter the specific sports interface, including outdoor running, cycling, skipping rope, badminton, table tennis, tennis, mountaineering, walking, volleyball, cricket, rugby, hockey, dancing, spinning, Yoga, sit-ups, indoor running, basketball, football, baseball, gymnastics, rowing, opening and closing jumps, and other 23 sports. Click on the sport item and start monitoring after three seconds. Swipe right and click "结束" to end, right Slide back to the previous menu.

Sleep detection function

Switch to the sleep interface and check the sleep status. When you fall asleep, the bracelet will automatically determine to enter the sleep monitoring mode, automatically detect the length of your deep sleep / light sleep / wake up all night, and calculate your sleep quality (Note: Only wearing the bracelet to fall asleep will have sleep data.)

Weather interface

Switch to the weather interface, you can view the weather and temperature conditions.


Stopwatch interface

Switch to the stopwatch interface, click the interface to enter the timing interface, and slide to the right to return to the upper-level interface when the timing is stopped.

Find phone interface

Switch to the search interface. When connected, press and hold the search phone interface, and the phone will send a ringtone reminder.

Set interface

Adjust brightness： Click to choose different brightness, and it has four-gear brightness.

Mobile phone search： Click the interface of mobile phone search in connection state, and the mobile phone will send out a bell warning

Restore: Have a long press the interface for 3 seconds to reset all data, and the data in that day in APP will be reset as well (other historical data will be kept).

Shut down:Have a long press on the interface for 3 seconds to power off.

Precautions

- Please prevent your bands from being shocked, and exposed under high temperature and blazing sun for a long time.
- Please don't disassemble, repair or change the device.
- The device should adopt 5V500MA to charge, it is forbidden to charge with over-voltage power supply.
- It should be used within the temperature of 0 degrees - 45 degrees. It is forbidden to throw into fire, avoiding explosion.
- Please charge after drying the band with soft cloth, otherwise, it will corrode charging touch points, leading to charge fault.
- Please don't touch chemical substance such as gasoline, cleaning solvent, propanol, alcohol or insect repellent.
- Please don't use the product in high voltage and high magnetic environment.
- In case of sensitive dry skin or too-tight wearing, you may be uncomfortable.
- Please don't take a hot bath by wearing it.
- Please dry sweat in your wrist timely for its possibility to lead to skin itch or allergy by long-time touch with soap, sweat, allergens or contaminated ingredients.
- In case of frequent use, it is suggested to clean weekly. Please wrap by wet clothes, and remove oil or dust by mild soap.
- Do not open more push without special requirements.

Common problem

How to do when it cannot connect with Bluetooth? (connection/failing re-connection or low speed while re-connection)
Please confirm whether your device is IOS8.0, Android 4.4 or above, and supports Bluetooth 4.0.

1. It takes time differently for each connection caused by the signal interference of Bluetooth wireless connection. If it fails to connect for a long time, please confirm that you use it in environment without magnetic field or interference of multiple Bluetooth device .
2. Turn on again after turning off the Bluetooth.
3. Turn off the background applications of mobile phones or restart your mobile phones.

4. Please don't make your mobile phones connect with other Bluetooth device or functions at the same time.

5. Please check whether the APP is operated normally in the background. Otherwise, it may fail to connect.

How to do when it fails to search the band?

1. Please make sure that the watch is powered and activated, and it is not bound to other mobile phones, bring the watch close to the phone.

If it still does not work, please turn off the phone Bluetooth and restart the phone Bluetooth after 20 seconds.

2. Make sure your phone turns on GPS.

Why is the wristband tight when measuring heart rate?

The band adopts light-reflection principle. It collects the signal reflected

FCC Caution:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.