USER INSTRUCTIONS

Dr Brown's

Natural Flow

Electric Breast Pump

Patented
NEW
Technology

**Pulsating honeycomb liner
uniquely stimulates breast
for increased milk supply**

**Patented One-Way airflow
to protect mom and baby
from bacteria**

***Clinically proven
to treat engorgement***

**IMPORTANT – Please carefully read the below.
Do not throw away these instructions.**

As a precaution against the risk of electrical shock please observe the following...

- Always unplug this product immediately after use.
- Do not use while bathing or showering.
- Do not place or store this product where it can be pulled or fall into a bath or a sink.
- Do not place or drop into water or other liquids.
- Do not reach for an electrical product that has fallen into water. Unplug the product from the wall socket immediately.

As a precaution against the risk of burns, electrocution, fire or injury to persons...

- This product should not be left unattended when plugged in.
- This product should not be used by children.
- The use of this product is only for its intended use as described in this manual. Only use the parts that come with this pump. Do not use attachments that were not specifically made for this pump.
- Do not operate this product if the cord or plug is damaged, if the unit is not working properly, if it has been dropped or damaged or immersed in water. Return the product to the manufacturer for examination and or repair.
- Ensure the electric cord is kept away from heated surfaces.
- Never use this product while sleeping or drowsy.
- Keep openings free of lint, hair or other foreign materials.
- Never insert or drop objects into openings or tubing.
- Do not use outdoors with AC adapter or operate where aerosol spray products are used or where Oxygen is being administered.
- Do not put the motor or power supply in a microwave.

Dr. Brown's Natural Flow Electric Breast Pump

Contents

Electric breast pump parts diagram	3
Disassembling the electric breast pump	4
Cleaning and sterilization	4
Reassembling the electric breast pump	5
How to use the electric breast pump	6
Expressed breast milk storage	6
Expressed breast milk defrosting and warming	7
Successful breast milk expression	7

Congratulations on the birth of your baby.

Thank you for choosing the Dr. Brown's Natural Flow Electric Breast Pump to help you provide breast milk to your baby. Whether returning to work or out for a few hours, your baby can still enjoy the benefits of your breast milk. Even if you can't be there to give it yourself.

Our pump stimulates the breast with the use of gentle pulsation and One-Way Airflow, which means faster let-down, more milk for your baby and unrivaled comfort for you.

Clinical trials show that stimulation is proven to help mothers initiate and maintain a long-term milk supply.

Many women do not express breast milk when they first try a new breast pump. Don't worry as this is common. Try not to give up, simply wait a few minutes and try again, ensuring that you are comfortable and relaxed when expressing.

Dr. Brown's Natural Flow Electric Breast Pump


Disassemble the Electric Breast Pump

Step 1. Unscrew the Bottle.

Step 2. Remove the honeycomb Liner (see photo).

Step 3. Hold Connector and turn the Motor Case to the unlocked position (see photo).

Twist motor case counter-clockwise to unlock.

Cleaning Instructions

Gently wash the Liner, Connector and Bottle in warm soapy water and rinse. Ensure all parts are rinsed thoroughly, then air or towel dry.


To clean the Motor Case, wipe down with a damp cloth. Be careful not to dislodge the purple rubber plugs. Do not remove purple rubber plugs.

Do not immerse the Motor Case in water or liquids.


Note: You may boil the Liner, Connector and Bottle...but normal hand washing is just fine. If you choose to boil these parts, no more than 3 minutes is needed. These parts are dishwasher safe, however dishwashers are harsh on plastic and will eventually reduce performance. Microwave sterilization is not recommend as it may reduce the performance of this product.

Reassemble the Electric Breast Pump

1. Wash your hands thoroughly before handling parts.
2. Insert the honeycomb Liner, insuring that the edges snap firmly into place.
3. Screw the Bottle in firmly.
4. Place the Motor Case onto the Connector top (see photo).
5. Align arrow with open padlock, press Motor Case firmly onto Connector, then turn clockwise to locked position (see photo).

Line up arrows in starting position

Twist the Motor Case clockwise to lock

Helpful Hints:

- Always examine the pump for firm and proper fit BEFORE USE. An airtight fit is essential for proper operation.
- At different vacuum settings, you may notice different sounds coming from the pump. For example, when the Vacuum Dial is at a lower setting, you may hear an airy, wheezing noise as the pump operates. This is perfectly normal. At higher settings on the Vacuum Adjustment Dial, you may notice that the air noises discontinue.
- The honeycomb Liner should be kept in its original shape. If folded or crimped for long periods, it will stay that way. If this occurs, submerge it in boiling water until the original shape returns.

- The Dr. Brown's pump feels like no other pump at the breast because it expresses milk differently than all other pumps. Some of our moms say it is a "fluttery" or "tickling" feeling that's lighter, more gentle on the breast. All other pumps tug and pull the nipple, which can be uncomfortable and is often painful.
- It is highly recommended that you drink plenty of water and remain well hydrated.

How to use the Electric Breast Pump

1. Plug the Power Cord into the wall socket before inserting other end into the Motor Case on the pump.
2. Make yourself comfortable, sitting in an upright position leaning slightly forward, well supported and relaxed.
3. Set the Vacuum Dial to minimum position before starting the pump and center the nipple into the honeycomb Liner. Turn the pump on.
4. After a few seconds, slowly increase the Vacuum Dial to a desired comfort level. Everyone is unique, but most first time users prefer a medium vacuum setting. Some moms, especially more experienced users, prefer a vacuum setting between medium and maximum.
5. If the vacuum is too low, simply adjust Vacuum Dial up. **If the vacuum ever becomes too high, press side of breast to break vacuum, adjust the vacuum setting lower, then center the nipple into the honeycomb Liner and resume pumping.**
6. When you have finished expressing, turn off the breast pump, press side of breast to break vacuum, remove the Power Cord from the Motor Case and unplug from wall socket.

NOTE: When the pump is in operation, do not lay it on its side.

Breast Milk Storage Guidelines

1. Storage containers should be clean before use.
2. Breast milk should be kept at room temperature for as short as time possible and refrigerated or frozen immediately after expression.
3. Breast milk can be stored at room temperature for up to six hours.
4. Breast milk stored in a refrigerator should be used within 5 days hours and stored at a temperature of 35-39°F (2-4°C).
5. Breast milk that has not been used within 5 days should be stored frozen at -4°F (-20°C) for a maximum of 6 months.

Breast Milk Defrosting and Warming

1. Wash your hands thoroughly before handling breast milk.
2. It is recommended to defrost frozen breast milk in the refrigerator.
If needed quickly -- defrost under cool, then warm running water.

DO NOT USE A MICROWAVE TO DEFROST OR WARM |BREAST MILK.

Successful Breast Milk Expression

The Dr. Brown's Natural Flow Electric Breast Pump is very simple and easy to use. Like all breast pumps, you may need to try a few times before you succeed initially. Don't give up – see below for some helpful tips:

1. Try not to express when you are rushed -- choose a time when you can relax, free of interruptions.
2. Let-down can be encouraged by having your baby (or your baby's photo) nearby.
3. Warmth and relaxation can encourage milk flow.
4. Try expressing after a bath or a shower.
5. Try expressing from one breast while your baby is feeding from the other, or immediately pumping after nursing your baby.

If you experience difficulties with breast milk expression, please consult a lactation consultant or your doctor.

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©2006 Handi-Craft Company