FORA® COMFORT pro GD40

Blood Glucose Test Strip

312-4395100-002 ver 2.0 2021/08

Warnings

- For ***in vitro*** diagnostic use (for use outside of the body only).
- For single use only.
- Healthcare professionals and other users testing multiple patients with this system should handle everything that has comes into contact with human blood carefully to prevent transmitting infectious diseases, including sanitized objects.
- Please read this sheet and your Blood Glucose Monitoring System Owner's Manual before you use this test strip. Use only FORA COMFORT pro GD40 Test Strips with FORA COMFORT pro GD40 Blood Glucose Monitoring System to obtain accurate results, and be covered by the manufacturer's warranty.
- Results may be inaccurate when testing on patients with abnormally low blood pressure or those who are in shock.
- Please do not use FORA COMFORT pro GD40 Blood Glucose Monitoring System on critically ill patients. While the blood glucose result is extremely hyperglycemia (over 600 mg/dL (33.3 mmol/L)), the collection of capillary blood from the approved sample sites is not advised when the peripheral circulation is impaired as the delivery of physiological blood glucose level might not be a true reflection. The following circumstances may apply: severe dehydration as a result of diabetic ketoacidosis or due to stress hyperglycemia, hyperosmolar non-ketotic coma, shock, decompensated heart failure NYHA Class IV or peripheral arterial occlusive disease.
- Keep test strips and lancets away from small children. If swallowed, consult a doctor immediately for advice.

Intended Use

FORA COMFORT pro GD40 test strips, when used together with FORA COMFORT pro GD40 Blood Glucose Monitoring System, allow your blood glucose levels to be measured by yourself at home or by healthcare professionals. It uses fresh whole blood samples from the finger, and the following areas: the palm, forearm and upper arm. This system is not intended for use in the diagnosis or screening of diabetes mellitus.

Professionals may test with capillary and venous blood sample; home use is limited to capillary whole blood testing. Use only heparin for anticoagulation of whole blood. Please do NOT use EDTA for anticoagulation.

Limitations

- Hematocrit: The hematocrit level is limited to between 0% and 70%. Please ask your healthcare professional if you do not know your hematocrit level.
- Neonatal Use: ***This test strip can be used for the testing of newborns.***
- This test strip is used for testing fresh capillary, venous, arterial and neonatal blood.
- Please see Appendix: Summary of substances and concentrations in excess of limitation with interference.
- Altitude Effects: Altitudes up to 3,275 m (10,742 ft) do not affect test results.
- Do not test blood glucose during or soon after a xylose absorption test. Xylose in the blood can produce elevated glucose results.

Storage and Handling

IMPORTANT: Do not use the test strips if they have expired.

- Test strips expire 6 months after first opening. Write the first opening date on the test strip vial when you first opened it.
- Store the test strips in a cool, dry place between 2°C and 32°C (35.6°F and 89.6°F) and below 85% relative humidity.
- Keep the test strips away from direct sunlight. Do not store the test strips in high humidity.
- Store the test strips in their original vial **ONLY**. Do not transfer them to a new vial or any other containers.
- Do not touch the test strips with wet hands.
- Use each test strip immediately after taking it out of the vial. Close the vial immediately after taking out a strip.
- Keep the vial closed at all times.
- Do not bend, cut, or alter the test strip.

Testing Your Blood Glucose

PLEASE WASH AND DRY YOUR HANDS BEFORE PERFORMING ANY TESTS.

Please refer to your Owner's Manual for more information.

The used lancet and test strip are potentially biohazardous. Please dispose of them carefully according to your local regulations.

Reading Your Result

Your blood glucose readings deliver plasma equivalent results and are displayed in mg/dL or mmol/L.

Reference values

Time of day	Normal plasma glucose range for people without diabetes
Fasting and before meal	< 100 mg/dL (5.6 mmol/L)
2 hours after meals	< 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association (2012). Clinical Practice Recommendations. Diabetes Care, 35 (Supplement 1): S1-100.

Please consult your doctor to determine a target range that works best for you.

Questionable or inconsistent results

If your test results are unusual or inconsistent with how you are feeling:

- Make sure the confirmation window of the test strip is completely filled with blood.
- Check the expiry date of the test strips.
- Check the performance of your meter and test strip with the control solutions.

Please Note: Unusually high or low blood glucose levels may be symptoms of a serious medical condition. If most of your results are unusually high or low, please contact your healthcare professional.

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* All concentrations in mg/dL unless otherwise stated.

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Substance	Limiting Concentration (mg/dL)	Therapeutic / Physiologic Concentration Range (or Upper Limit) (mg/dL)
Acetaminophen (Paracetamol)	> 6.25	0.45 – 3
Ascorbic acid	> 5.0	2.0
Bilirubin (Unconjugated)	> 20	0 – 2
Dopamine	> 1.25	0.03
Levo - Dopa	> 0.7	0.02 – 0.28
Methyl - Dopa	> 0.625	0.1 – 0.5
Reduced Glutathione	> 30	24.25 – 32.2
Tolazamide	> 6.25	2 – 2.5
Uric acid	> 10	2 – 8
Xylose	> 6.25	N/A
Lipemic Samples (Triglycerides)	> 3000	30 – 300

Table C Summary of anticoagulant and concentrations with interference. / **Tableau C** Résumé des anticoagulants et concentrations avec interférence. / **Tabela C** Riassunto degli anticoagulanti e delle concentrazioni con interferenza. / **Tabla C** Resumen de anticoagulantes y concentraciones que presentan interferencia. / **Tabla C** Resumo de anticoagulantes e concentrações com interferência.

Substance	Limiting Concentration (mg/dL)	Concentration of Blood Drawing Tube (mg/dL)
Na-Fluoride/K-Oxalate	< 250/200	250/200
K2EDTA	180	180
K3EDTA	175.5	175.5

Table D The following compounds within the therapeutic or physiologic concentration range may produce elevated glucose results. / **Tableau D** Les composés suivants dans la plage de concentration thérapeutique ou physiologique peuvent engendrer des résultats glycémiques élevés. / **Tabela D** I seguenti composti entro l'intervallo terapeutico o fisiologico possono causare risultati glicemici elevati. / **Tabla D** Los siguientes componentes dentro del intervalo de concentraciones terapéuticas o fisiológico pueden producir resultados de glucosa elevados. / **Tabla D** Os compostos seguintes dentro do intervalo de concentração terapêutica ou fisiológica podem produzir resultados elevados de glicose.

Substance	Concentration Tested (mg/dL)	Therapeutic / Physiologic Concentration Range (or Upper Limit) (mg/dL)
Pralidoxime iodide	> 5.0	~ 10 (i.v. Dose 500 mg)

Referencia / Référence / Referenza / Referência / Referencia

- [1] McEnroe, J Robert, et al. National Committee for Clinical Laboratory Standards. Interference testing in clinical chemistry; Approved Guideline – 2nd edition. NCCLS : 2005 - EP7-A2, volume 25, number 27.
- [2] EN ISO 15197: 2015 (E): In vitro diagnostic test systems – Requirements for blood-glucose monitoring systems for self-testing in managing diabetes mellitus, Second edition, section 6.4.4 Interference Testing, pp. 28-30.
- [3] Hardman JG, Limbird LE, Molinoff PB, et al. Goodman & Gilman's the Pharmacological Basis of Therapeutics. 9th ed. New York, NY: McGraw-Hill; 1996: 1570-1571.
- [4] Zaloga GP, Macgregor D. The Critical Care Drug Handbook. 2nd ed. St Louis, MO: Mosby-Year Book; 1997.
- [5] Kost GJ