

Victor Moonboot 3.0 Air Ankle Fitting Guide

Step 1.



1. Open all the straps.
2. Remove the liner from the boot.
3. Place your foot in the liner with your heel at the back.

Step 2.



Close the liner, starting at your toes, ensuring it fits snugly around your foot and leg.

Step 3.



Place your foot, with the liner, into the boot, ensuring your heel is fully at the back.

Step 4.



Fasten the straps from the toes upward.

- It's normal if your toes extend past the liner, but if they extend past the rigid shell, you may need a larger size. Ensure your heel is properly positioned at the back of the boot.
- If you experience numbness, tingling, or discomfort, loosen the strap until comfortable. If symptoms persist, contact your medical professional.
- Use the provided comfort pads to fill gaps and add extra padding as needed.

Step 5.



To inflate, press and release the blue air bulb until you get the desired support. The liner may expand initially, requiring additional pumps to maintain support.

Step 6.



To deflate, press the deflate valve pictured. If fully deflated, it may take up to 50 pumps to regain desired support.