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User Guide

Smart Heart Rate and Blood Pressure Wristband

User Guide

English

Deutsch

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Español

Français


Italiano

Polski

Hello, thank you for choosing our smart health bracelet, when using this device for the first time for your convenience, please read the instructions carefully and follow the steps in the instructions.

Hand ring quick use instructions

Component introduction

Instructions for Charging

For initial use, please take out the bracelet and charge it with the corresponding charger. When charging, the bracelet will turn on automatically. The charging method is shown in the following figure: take the bracelet body out of the wrist strap and plug the charging port into the standard USB port. The charging voltage is 5V. The battery grid displays in icons when charging and the grid will be full when it is fully charged.

Install bracelet APP

Scan the following QR codes or go to the major APP Stores to download and install the APP "GlobyFit".

Device requirements: iOS 9.0 and above, Android 4.4 and above, Bluetooth 4.0 is supported.


Open the mobile client to set up personal information (Switch to the "Device" page and click on the binding device)

Click on your device in the list of scanned devices

Binding completed

Main function interface

Main interface/dial Press and hold touch on the main interface to enter, click to switch, and press and hold ON. The dial has 5 styles.

Pedometer/Distance/Calories

View steps, distance, calories burn in real time. It can synchronize apps in time to view motion data.

Interface of body fat

Press and hold touch on the main interface to enter, click to switch, and press and hold ON. The dial has 5 styles.

Do not release your hand during the test. After the test is completed, the various health index data will be displayed.

Precautions for Body Fat Test

1. Before testing, please set the correct personal information with APP such as gender/weight, etc. If no personal information is set, there will be deviations from the system test.

2. During the test, gently press your index finger and middle finger on the electrodes at both ends of a bracelet.

3. During the test, do not move two fingers away from the position of the electrodes and do not have the left hand and right hand to contact each other. The two electrodes at the bottom should be completely attached to the wrist.

4. People with cardiac pacemakers and other implanted medical equipment are not recommended to use this product. For professional athletes and pregnant women or women in either of other special physical conditions, there may be deviation in composition test data.

5. Normally, the fat content of the upper body is higher than that of the lower body, so the body fat rate measured by the body fat bracelet is slightly higher than that measured by the body fat scale.

Heart rate monitoring

Switch to the heart rate monitor interface. Heart rate monitoring, in real-time. The measurement data can be synchronized with the APP in real time with a test report.

Blood pressure monitoring

Press and hold touch on the main interface to enter, click to switch, and press and hold ON. The dial has 5 styles.

Do not release your hand during the test. After the test is completed, the various health index data will be displayed.

Multiple sports modes

Long press to enter sport mode monitoring: Skipping, running, cycling, table tennis, badminton, tennis.

Press and hold touch on the main interface to enter, click to switch, and press and hold ON. The dial has 5 styles.

Notification Message push

On the message interface, long press three seconds to view the message. Click to turn the page, view the message to exit or grow Press three seconds to exit.

More features

On the more features page, long press three seconds to enter/leave menu page for more operations

About the interface

On the bracelet interface, after displaying the Bluetooth address Four digits and version number.

Restore interface

In the bracelet restore page, long press the touch key to restore wristband.

Find your phone

Keep the phone and the fitness tracker connected normally. Press and hold the touch button on your tracker. The phone will ring.

Brightness adjustment page

On the brightness page, long press the touch button to enter. The screen brightness can be chosen freely.

Turn off

On the shutdown page, long press three seconds. Vibration shutdown

Other feature reminders

Reminder function needs to set the opening reminder switch on the APP side, and keep the phone and bracelet in the state that the Bluetooth connection is successful.

Sleep monitoring function

When you fall asleep, the bracelet will automatically judge to enter the sleep detection mode and automatically detect your entire night. Deep sleep / light sleep / awake time, calculate your sleep quality, sleep data temporarily only View with the APP.

Note: Sleep data will only be available when you wear the bracelet to sleep.

Precautions

1. Do not use an adapter with a current over 2A. Charging time is about 2-3 hours.

2. Do not charge after the device was water damaged.

3. This is an electronic monitoring product, data received from it cannot be used as medical evidence.

4. Blood pressure testing tip: a) Keep your body relaxed and sit while testing. b) Keep device at the same height as your heart, and don't talk while you are testing.

How to make the body fat test more accurately?

1. Connect the APP and set the correct personal information, such as gender, age, weight, etc. During the test, press your index and middle fingers to the electrodes at both ends of the bracelet. During the test, keep the two fingers on the electrodes (try to make the contact area between the fingers and the electrodes as large as possible) and keep the left hand from touching the right hand, and the two electrodes at the bottom should be fully attached to the wrist. In addition, patients with cardiac pacemakers and other transplant medical devices are not recommended to use this product. For professional athletes and pregnant women or a small number of other people with special physical conditions, there may be deviation in the composition test data.


Note: If you turn on the functions such as incoming call and SMS reminder, you must keep Bluetooth connected, and the function of offline alarm clock reminder is supported.

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