

REF A0022583**Starrsed DL**

Supplier:
Mechatronics Instruments B.V.
De Corantijn 13
1689 AN Zwaag
The Netherlands
T +31 229 291 129

Intended purpose

Reagent for the automatic dilution of blood samples in Starrsed ESR analyzers.

Instructions For Use

The product is for use on the Starrsed analyzers only. See the instrument manual for detailed instructions.

1. Remove the cap from the new container, using the container cap tool A0033594.
2. Remove the cap with pickup tube from the empty reagent container.
3. Replace the empty container with the new one.
4. Place cap with pickup tube on the new container. Avoid touching the pickup tube which goes into the container.

Chemical/Physical data

Contains a slightly hypertonic amount of citrate (<3.5 % by weight) and preservative (4 ‰ by weight)

Packaging

Starrsed DL is supplied in 10 Liter jerrycans, LxWxH = approx.23x19x31 cm.

Safety

Any serious incident that has occurred in relation to this product must be reported to the distributor and/or manufacturer and the competent local authorities.

The latest version of the Safety Data Sheet is available on the website www.rrmechatronics.com

Storage and stability

Transport stability at 5 - 60 °C (maximum of 1 month)
Shelf life after production: 2 years, but no longer than the expiry date.
Store at 5 - 35 °C in dry conditions in well sealed containers.
In-use stability: Maximum of 3 months at advised temperature of 5 - 30°C.

Ecological information

Generally not hazardous for water.

Waste treatment and cleaning


Avoid discharge to lakes, streams, sewers, etc. Disposal must be made according to official regulations.
Wherever possible cleaning should be performed with normal cleaning agents. Avoid use of solvents.



RR Mechatronics Manufacturing B.V.
De Corantijn 13, 1689 AN Zwaag,
The Netherlands
www.rrmechatronics.com



Qserve Group UK Ltd
282 Farnborough Road,
Farnborough, GU14 7NA
United Kingdom



QUNIQUE GmbH
Bahnhofweg 17
5610 Wohlen
Switzerland

rev. A01
Initial release