HEALIX ADVANCE™ SP Biocomposite Anchor

en-cs-et-hu-it-lv-pl-sk

Issued 2024/04
© DePuy Synthes 2024.
All rights reserved.

Directory

en		Page 3
cs		Str. 12
et		Lehekülg 21
hu		30. oldal
lt		40 psl.
lv		49. lpp.
pl		Strona 59
sk		Strana 69

HEALIX ADVANCE™ SP Biocomposite Anchor

Some products may not be licensed in all jurisdictions.

DESCRIPTION

Contents are STERILE unless packaging is damaged or opened.

The HEALIX ADVANCE™ SP Biocomposite Anchor is a sterile implant designed to secure soft tissue to bone when used in conjunction with high-strength #2 suture.

PERFORMANCE CHARACTERISTICS

DePuy Synthes Mitek Sports Medicine has established the performance and safety of HEALIX ADVANCE SP Biocomposite Anchor and that it represents a state-of-the-art medical device and performs as intended for the function of soft tissue to bone reattachment when used according to the Instructions for Use and labeling.

MATERIALS

Anchor: BIOCRYL (Polylactic Acid (PLA) and Tricalcium Phosphate (TCP))

Dilator: PEEK (Polyetheretherketone)

Driver Shaft: Stainless steel

Driver Handle: Polycarbonate

Driver Handle Paddle: Acrylonitrile Butadiene Styrene (ABS) and stainless steel

Driver Collar: ABS

Threader Tab: ABS and stainless steel

Stay Suture: Ultra-High Molecular Weight Polyethylene (UHMWPE) and Green polyethylene terephthalate (PET), coated with poly {oxy-1,4 butanediyl oxy(1,6-dioxo-1,6 hexanediyl)} or Polybutilate.

The table below lists the materials and substances from which the implant components are made. Please contact the manufacturer for more information.

Implant Component	Implant Subcomponent and nominal weight %	Material Description	Material Composition (nominal, by weight)
Anchor (BIOCRYL)	Anchor Body (absorbable), 100%	PLLA/ (β -TCP) (Biocryl)	72% Poly (L-Lactide) (PLLA) 28% beta-tricalcium phosphate (β -TCP)
Dilator	Dilator (nonabsorbable), 100%	Poly-etheretherketone (PEEK)	100% PEEK

For indication that the device contains a restricted substance or material of animal origin see product label.

One or more nonimplantable components of this device such as Driver Shaft, Driver Awl Tip, Driver Hex Tip and Proximal Awl shaft contains the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight:

Cobalt; CAS No. 7440-48-4; EC No. 231-158-0

Current scientific evidence supports that medical device manufactured from cobalt alloys or stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

INTENDED USE

The HEALIX ADVANCE SP Biocomposite Anchor is intended for soft tissue to bone reattachment.

INDICATIONS

The HEALIX ADVANCE SP Biocomposite Anchor is indicated for use in the following procedures for reattachment of soft tissue to bone:

- **Shoulder:** Rotator cuff repair

CONTRAINDICATIONS

- Pathologic conditions of bone such as cystic changes or severe osteopenia that would impair its ability to securely fix the HEALIX ADVANCE SP Biocomposite Anchor.
- Pathologic changes in the soft tissues being fixated to bone that would prevent their secure fixation by the HEALIX ADVANCE SP Biocomposite Anchor.

- Comminuted bone surface that would prevent secure fixation of the HEALIX ADVANCE SP Biocomposite Anchor.
- Physical conditions that would eliminate or tend to eliminate adequate implant support or retard healing, e.g., blood supply limitation, previous infection, active infection, etc.
- Conditions which tend to limit the patient's ability to restrict activities or follow directions during the healing period.

INTENDED USER AND USER TRAINING

This device is intended to be used in a surgical setting by qualified health care professionals, e.g., surgeons, physicians, operating room staff, and individuals involved in preparation of the device. All personnel handling the device should be fully aware of the IFU and, if applicable, the surgical procedures.

PATIENT TARGET GROUP

For patients 18 years or older undergoing surgical procedures as per indications.

WARNINGS

- DePuy Mitek anchors are designed to lock into cortical or cancellous bone. Bone stock must be adequate to allow proper and secure anchor placement. Incomplete insertion or overtensioning or poor bone quality may result in anchor pullout.
- Immediate range of motion should be avoided to allow biological bony/soft tissue healing.
- This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the thoracic or lumbar spine.

- This product is for single use only. It has not been designed to be reused/resterilized. Reprocessing may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single-use instruments can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.
- When hard bone is encountered, it is recommended to use an awl, drill, or similar instrument to create a pilot hole in bone prior to anchor insertion. If insertion is attempted without a pilot hole and the anchor cannot be advanced into bone, the device should be discarded. The user should then create a pilot hole and insert a new device.

PRECAUTIONS

- The HEALIX ADVANCE SP Biocomposite Anchor is supplied STERILE and is intended for single use only. Do not resterilize. Do not use if sterile packaging appears to be damaged.
- The product should be stored under cool, dry conditions.
- Inspect all instruments for damage before use. Do not attempt to repair.
- A surgeon should not begin clinical use of the HEALIX ADVANCE SP Biocomposite Anchor without reviewing the Instructions for Use and practicing the procedure in a skills laboratory.
- Ensure that anchor is inserted axially to the bone hole (limit off-axis insertion).
- Ensure area of anchor insertion is free of soft tissue.
- Improper instrument use, axial misalignment or levering with the anchor upon insertion may result in anchor fracture or reduced performance.

ADVERSE EFFECTS

The following list of adverse effects may be associated with the use of this device:

- Adverse effects of implanted devices include mild inflammatory and foreign body reactions/adverse local tissue reaction
- Infection/sepsis
- Neurological complications
- Prolonged or additional surgical intervention
- Bone injury/damage
- Soft tissue injury
- Inadequate osteointegration
- Joint dislocation/joint disorder
- Pain, discomfort and/or abnormal sensations
- Treatment failure

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The Summary of Safety and Clinical Performance (SSCP) will be available upon launch of the European Medical Device Database (EUDAMED) at the following URL <https://ec.europa.eu/tools/eudamed>. The device tradename listed on the IFU cover can be used to search and find the intended SSCP in EUDAMED.

EXPECTED CLINICAL BENEFIT

Improvement in joint stability, pain, or movement.

MRI SAFETY INFORMATION

The HEALIX ADVANCE SP Biocomposite Anchor is MRI Safe.

INSTRUCTIONS FOR USE

1. Clean the bone surface of excess tissue around the desired point of insertion; do not de-corticate surface.
2. In hard bone, it may be necessary to create a pilot hole in bone using an instrument such as an awl or drill prior to anchor insertion.

NOTE: If insertion is attempted without a pilot hole and the anchor cannot be advanced into bone, the device should be discarded. The user should then create a pilot hole and insert a new device.

3. Load the desired sutures through the Threader Tab. The suture tails should extend approximately one inch (1")/2.5 cm through the Threader Tab kite opening.

NOTE: The device is intended to be used with high-strength #2 suture.

4. Pull the Threader Tab to load the sutures through the device. Discard the Threader Tab.
5. Hold slight tension on the sutures and slide the distal tip of the device toward the insertion site.
6. Use a mallet to insert the Dilator and Anchor Tip into bone, stopping as the first Anchor Thread enters bone.
7. Use downward force to hold the anchor nose in bone and apply desired tension to the sutures. Tensioned sutures may be placed in the slot on the driver handle for convenience.
8. To insert the anchor:

- a. Release tension on any sutures that are held by hand.

NOTE: Maintaining tension on the sutures exiting the driver by holding them while inserting the anchor will result in reduced fixation.

- b. Hold the Driver Paddle stationary and turn the Driver Handle clockwise, while applying downward force on the Driver and/or paddle, to insert the anchor. Stop insertion once the laser line is flush with bone.

9. Release the Stay Suture from the driver handle, release any sutures that have wrapped around the Driver during insertion, and remove the Driver from the anchor.
10. Remove the Stay Suture from the implant by pulling on one end of the suture and discard.

STORAGE

- Store below 25°C (77°F). Do not use after expiration date. Store in a cool, dry area.

DEVICE LIFETIME

The device lifetime of the implants includes the therapeutic lifetime period of twelve (12) weeks and the posttherapeutic lifetime.

The posttherapeutic lifetime of the Biocryl biocomposite Anchor, or implant, includes the absorption period that is generally defined to be about 5 years.

The posttherapeutic lifetime of the PEEK dilator, or implant, is the lifetime remaining of the patient as the implant is permanent and nonabsorbable.

DISPOSAL OF DEVICE

- This device should be handled and disposed of as healthcare medical waste in accordance with hospital procedures and applicable regulations.
- Any device that has been contaminated with potentially infectious substances of human origin (such as bodily fluids) should be handled according to hospital protocol for infectious medical waste.
- Any device that contains sharp edges should be discarded according to hospital protocol in the appropriate sharps container.

REPORTING OF SERIOUS INCIDENTS

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State or country in which the user and/or patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- a. the death of a patient, user or other person,
- b. the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health,
- c. a serious public health threat.

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