



ResMed

AcuCare™ HFNC

High flow nasal cannula

User Guide

English | Français | Português | Español | Deutsch | Italiano | Nederlands | Dansk | Norsk | Svenska |
Suomi | Polski | Ελληνικά | 日本語

Rx Only

Intended use

The AcuCare HFNC (high flow nasal cannula) system is intended to provide high flow therapy for patients >30 kg who would benefit from receiving high flow warmed and humidified respiratory gases. The AcuCare HFNC is intended for use in patients with acute respiratory failure or non-acute spontaneously breathing patients in the hospital/clinical environment, or for non-acute spontaneously breathing patients in the home. The AcuCare HFNC is designed to deliver air, with or without oxygen, up to a maximum flow of 60 l/min. The AcuCare HFNC is for single-patient use (maximum 7 days) in the hospital/clinical environment, and is for single-patient use (maximum 30 days) in the home environment.

Contraindications

High flow oxygen therapy may generate some positive airway pressure. Contraindications for positive airway pressure devices apply: pneumothorax, pathologically low blood pressure, cerebrospinal fluid leak, recent cranial surgery or trauma, severe bullous lung disease, dehydration.

⚠ GENERAL WARNINGS

- **Single patient reuse only.**
Maximum 7 days reuse in the hospital/clinical environment, 30 days reuse in the home environment if cleaning instructions are followed.
- Only use in approved configurations.
- Check for condensation regularly. Drain as required.

- Oxygen supports combustion. Do not use in the presence of an open flame.
- Do not use if components are damaged.
- Risk of minor nasal mucosa irritation and minor bleeding due to gas flow and mechanical contact.

⚠ CAUTION

- Before connecting the nasal cannula check for adequate gas flow, oxygen and temperature.
- Ensure the nasal cannula is sized correctly and does not seal nares.
- Do not crush or stretch tubing.
- Always follow cleaning instructions. Some cleaning products may damage the cannula, its parts and their function, or leave harmful residual vapours.

Cleaning your cannula at home

If there is any visible deterioration of the cannula (cracking, crazing, tears etc.) the product should be discarded and replaced.

Weekly:

1. Soak the cannula in warm water with a mild liquid detergent.
2. Hand wash the cannula with a soft bristle brush.
3. Thoroughly rinse the cannula under running water.
4. Leave the cannula to air dry out of direct sunlight.

If the cannula is not visibly clean, repeat the cleaning steps.

Technical specifications

Compatible devices: ResMed Lumis HFT

Setup: For use with high flow therapy humidification systems which deliver a warm,

humid air flow $\leq 37^{\circ}\text{C}$ ($\leq 98^{\circ}\text{F}$) and $\leq 100\%$ r.h. with a maximum oxygen content of $100\% \text{ O}_2$ and a flow of up to 60 L/min.

Connection: 22 mm Male Conical Connector (complies with ISO 5356-1).

Flow range: $\leq 60 \text{ L/min}$

Operating conditions: $+18^{\circ}\text{C}$ to $+35^{\circ}\text{C}$ ($+64^{\circ}\text{F}$ to $+95^{\circ}\text{F}$); $\leq 95\%$ r.h. non-condensing.

Storage conditions: -20°C to 60°C (-4°F to $+140^{\circ}\text{F}$); $\leq 95\%$ r.h. non-condensing.

Resistance to Flow at 60 l/min: Size Small (S) $29 \text{ cmH}_2\text{O}$; Medium (M) $17 \text{ cmH}_2\text{O}$; Large (L) $12 \text{ cmH}_2\text{O}$

Disposal

This nasal cannula does not contain any hazardous substances and may be disposed of in accordance with local environmental regulations.

Symbols

 Single patient reuse;  Not made with natural rubber latex;  Rx Only
Prescription only (In the US, Federal law restricts these devices to sale by or on the order of a physician);  Do not use if packaging is damaged

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See [ResMed.com](https://www.resmed.com) for other ResMed locations worldwide. AcuCare and Lumis are trademarks and/or registered trademarks of the ResMed family of companies. For patent and other intellectual property information, see [ResMed.com/ip](https://www.resmed.com/ip). © 2022 ResMed. 238011/1 2022-01
ResMed Pty Ltd