



# Q2ROLLER®

## Lateral Turning Device

---

# User Manual

Visit [www.HoverMatt.com](http://www.HoverMatt.com) for other languages

### TABLE OF CONTENTS

|                                                  |   |
|--------------------------------------------------|---|
| Intended Use and Precautions.....                | 2 |
| Part Identification - Q2Roller.....              | 2 |
| Part Identification - Air Supply .....           | 2 |
| Air Supply Keypad Functions.....                 | 3 |
| Instructions for Use .....                       | 3 |
| Product Specifications/Required Accessories..... | 3 |
| Cleaning & Maintenance .....                     | 4 |
| Returns and Repairs.....                         | 4 |

## Intended Use and Precautions

### INDICATIONS FOR USE

Patients who are unable to turn or laterally position themselves

### CONTRAINDICATIONS

Patients who are experiencing thoracic, cervical or lumbar fractures that are deemed unstable

### INTENDED CARE SETTINGS

Hospitals, long-term or extended care facilities


### PRECAUTIONS

- Ensure the Q2Roller meets the patient's needs before using.
  - Always use a minimum of two caregivers when operating the Q2Roller.
  - Ensure valves are not in contact with the patient while on the Q2Roller.
  - Product should only be used by trained personnel.
  - Only use attachments and/or accessories that are authorized by HoverTech International.
  - Use of this device with products or accessories not authorized by HoverTech International could result in injury or equipment malfunction and may void the Manufacturer's Warranty. HoverTech International will not be held responsible for any injuries or damages caused due to the improper use of this device.
- Use the Q2Roller according to instructions.
  - Only for use with the HT-Air air supply adjustable button.
  - Make sure patient is centered on the Q2Roller before inflating. Reference Center Line Indicator on the Q2Roller.
  - Caregivers must ensure patient is attended and stabilized during Q2Roller inflation, deflation or when inflated at angles greater than 30 degrees.
  - Do not launder.
  - Reference product-specific user manuals for additional operating instructions.


### PRECAUTIONS – HOVERTECH INTERNATIONAL AIR SUPPLY

- Route the power cord in a manner to ensure freedom from hazard.
- Avoid blocking the air intakes of the air supply.
- Use this product only for its intended purpose as described in this manual.
- CAUTION: Avoid electric shock. Do not open air supply.

## Part Identification - Q2Roller




## Part Identification - Air Supply



**WARNING:** The HT-Air is not compatible with DC power supplies.  
The HT-Air is not intended for use with the HoverJack Battery Cart.

**CAUTION:** NO USER SERVICEABLE PARTS.  
Only qualified service personnel shall perform repairs on the HoverTech International Air Supply.

Air Supply Keypad Functions



-  The ADJUSTABLE keypad function has four different settings. Each press of the button increases the air pressure and rate of inflation. The Green Flashing LED will indicate the inflation speed by the number of flashes (i.e. two flashes equals the second inflation speed).
-  **STANDBY:** Used to stop inflation/air flow (Amber LED indicates STANDBY mode).

Instructions for Use

1. Ensure the red deflation valves are capped tightly before beginning inflation.
2. Place the Q2Roller on the bed mattress, and center using the Center Line Indicator.
3. Tuck the securing apron of the Q2Roller under the bed mattress to keep it in place during use.
4. Remove the backing from the adhesive strips on the Q2R Pad, and center it on top of the Q2Roller (if the patient is already on the mattress, place the device and pad under the patient using a log-rolling technique).
5. With the patient centered on the device, ensure all bed/stretchers rails are up, or follow your facility's protocol.
6. Locate the inflation valve on the OPPOSITE side of the turn direction, and place the air supply hose over the valve. Press the ADJUSTABLE button up to 3 more times to increase air flow and rate of inflation.
7. When the turn is complete, remove the hose from the valve and press the STANDBY button on the air supply keypad to stop air flow. Chamber may not need to be fully inflated to achieve desired position.
8. To quickly deflate the Q2Roller, remove the red cap from the deflation valve. To slowly deflate or adjust the chamber, depress the center of the black inflation valve.

Product Specifications/Required Accessories

Q2ROLLER® LATERAL TURNING DEVICE

|               |                                                                            |
|---------------|----------------------------------------------------------------------------|
| Material:     | Polyurethane                                                               |
| Construction: | RF Welded                                                                  |
| Width:        | Chambers: 46.5" (118 cm)<br>(combined chamber width)<br>Apron: 72" (183cm) |
| Length:       | Chambers: 44.5" (113 cm)<br>Apron: 44" (112 cm)                            |

LATEX FREE

Model #: HTR-200  
Patient Weight Limit: 600 lbs. (272 kg)  
Single-Patient Use

Q2R® PAD

|               |                          |
|---------------|--------------------------|
| Material:     | Non-Woven, Breathable PE |
| Construction: | Thermal Bonding          |
| Width:        | 39" (99 cm)              |
| Length:       | 39" (99 cm)              |

Model #: HTR-CHX-5PACK (5 ea.)  
Model #: HTR-CHX-CASE (50 ea.)

REQUIRED ACCESSORY:

Model #: HTAIR1200 (North American Version) – 120V~, 60 Hz, 10A  
Model #: HTAIR2300 (European Version) – 230V~, 50 Hz, 6A  
Model #: HTAIR1000 (Japanese Version) – 100V~, 50/60 Hz, 12.5A  
Model #: HTAIR2356 (Korean Version) – 230V~, 50/60 Hz, 6A

## Cleaning & Maintenance

### CLEANING

The Q2Roller is a single-patient use device. If the device becomes soiled, completely wipe it down using a germicidal cleaner (phenolic disinfectant, quaternary solution or other intermediate level disinfectant according to facility procedure) to disinfect. For best results, follow the cleaner manufacturer's recommended dwell time and instructions for use. Apply germicidal cleaner directly into hard-to-reach areas. Let air dry. Do not launder.

### INFECTION CONTROL

The Q2R Pad is recommended to cover the Q2Roller to keep the device from getting soiled (also available for separate purchase). Other materials, such as pads or linens, may also be placed on top of the Q2Roller to keep it clean.

If the Q2Roller is used on an isolation patient, the hospital should employ the same protocols/procedures it utilizes for the bed mattress and/or for the linen in that patient room.

### AIR SUPPLY CLEANING AND MAINTENANCE

See air supply manual for reference.

NOTE: CHECK YOUR LOCAL/STATE/FEDERAL/INTERNATIONAL GUIDELINES BEFORE DISPOSAL.

### PREVENTIVE MAINTENANCE

Prior to use, a visual inspection should be performed on the Q2Roller to ensure that there is no visible damage that would render the Q2Roller unusable.

The Q2Roller should be periodically inspected to ensure the following:

- All inflation valves are self-sealing with no evidence of leakage
- There are no punctures or tears in the Q2Roller

If any damage is found that would cause the Q2Roller not to function as intended, the Q2Roller should be discarded.

---

## Returns and Repairs

All products being returned to HoverTech International (HTI) must have a Return Goods Authorization (RGA) number issued by the company. Please call 800-471-2776 and ask for a member of the RGA Team, who will issue you an RGA number. Any product returned without an RGA number will cause a delay in the repair time.

Returned products should be sent to:

HoverTech International  
4482 Innovation Way  
Allentown, PA 18109

Attn: RGA # \_\_\_\_\_

Phone: 800-471-2776

Fax: 610-694-9601

For European companies, send returned products to:

Attn: RGA # \_\_\_\_\_  
Kista Science Tower  
SE-164 51 Kista, Sweden  
[www.Etac.com](http://www.Etac.com)  
[OrderExport@Etac.com](mailto:OrderExport@Etac.com)



4482 Innovation Way  
Allentown, PA 18109

800.471.2776  
Fax 610.694.9601

[www.HoverMatt.com](http://www.HoverMatt.com)  
[Info@HoverMatt.com](mailto:Info@HoverMatt.com)

# HoverTech Symbols

---



CE MARKING OF CONFORMITY



SINGLE PATIENT – MULTIPLE USE



CAUTION



ELECTRICAL AND ELECTRONIC EQUIPMENT



OPERATING INSTRUCTIONS



UDI UNIQUE DEVICE IDENTIFIER



MEDICAL DEVICE




KEEP DRY



HUMIDITY LIMITATION



TEMPERATURE LIMITATION



AUTHORIZED REPRESENTATIVE



MANUFACTURER



SERIAL NUMBER




LOT NUMBER



HoverTech International  
4482 Innovation Way  
Allentown, PA 18109  
[www.HoverMatt.com](http://www.HoverMatt.com)  
[Info@hovermatt.com](mailto:Info@hovermatt.com)

*These products comply with the standards applicable for Class 1 products in the Medical Device Regulation (EU) 2017/745 on medical devices.*



CEpartner4U , ESDOORNLAAN 13,  
3951DB MAARN, THE NETHERLANDS.  
[www.cepartner4u.com](http://www.cepartner4u.com)

*In case an adverse event in relation to the device, incidents should be reported to our authorized representative, CEPartner4u. CEPartner4u will forward information to the manufacturer.*