



## maxtec R219P50 Muffled 2-in-1 Adapter Instruction Manual

[Home](#) » [maxtec](#) » maxtec R219P50 Muffled 2-in-1 Adapter Instruction Manual 

### maxtec R219P50 Muffled 2-in-1 Adapter



#### Contents

- [1 Muffled 2-in-1 Adapter Instructions For Use](#)
- [2 Documents / Resources](#)
- [3 Related Posts](#)

### Muffled 2-in-1 Adapter Instructions For Use

**INTENDED USE:** The Muffled 2-in-1 Adapter is intended to thread onto the outlet of a flowmeter for connection to a 22mm male or 15mm female tube. The device also provides muffling of the noise produced by flowmeters when in use.


**WARNING:** This device produces back-pressure of approximately 2.8psi at 70LPM.\* This device may cause inaccuracy in flow or oxygen readings if used with mixers or flow meters that are not pressure compensated.

**WARNING:** Do not attempt to clean the internal surfaces or muffler element of this device for any reason. If it is contaminated, discard immediately.

**INSTALLATION:** Thread the device onto the bottom of any flowmeter. Hand-tighten to a snug fit. Do not tighten with tools. Firmly push the mating hose onto the outlet of the device with a twisting motion to ensure a tight seal.

**CLEANING:** The outside surfaces of this device can be cleaned using standard hospital-grade cleaning solutions. Use caution to ensure no cleaning agents enter the device as this may result in contamination of delivered gas to a patient. Under normal circumstances, oxygen and air mixtures that pass through this single-direction device are medical-grade, therefore the need to clean internals of the Muffled 2-in-1 does not exist. If, for any reason, you expect that the internals of the device may be contaminated, discard the device immediately and replace it with a new one.

## Documents / Resources

	<a href="#">maxtec R219P50 Muffled 2-in-1 Adapter</a> [pdf] Instruction Manual R219P50, Muffled 2-in-1 Adapter
-----------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------