



# NeoDocs uACR Test App Instruction Manual


[Home](#) » [NeoDocs](#) » NeoDocs uACR Test App Instruction Manual 

## Contents

- [1 NeoDocs uACR Test App](#)
- [2 Product Information](#)
- [3 Download app & Signup](#)
- [4 How to take the test](#)
- [5 FAQ](#)
- [6 Documents / Resources](#)
  - [6.1 References](#)
- [7 Related Posts](#)



## NeoDocs uACR Test App



## Product Information

### Specifications:

- Product Name: uACR Test
- Application: Dr-Neodocs app
- Sample Volume: 30 ml
- Results Time: 30 seconds

### Download app & Signup

1. Download Dr-Neodocs app
- 2.




Enter phone number & OTP



3. Enter first & last name




4. Enter organization password



**Note:** If you don't have the password or face any difficulty, Please contact us at: +91 9987339111

### **Do's & Don'ts**

- Make sure you download the right app – Dr-Neodocs
- Make sure you select the correct organization




- Dip the bottom part of testcard completely
- 



Click photo immediately (Once the timer ends)


- 



Don't open the pouch before collecting the urine sample



- Don't re-use the testcard



- Always perform test on flat surface




- Make sure the room is brightly lit



If you face any issues during the process, please contact us at +91 9987339111 / +91 98336 94081

## How to take the test

1. Click on “Add New Test”



2. Enter patient details Name, Age, Gender Phone number & Case id (Optional)



3. Request the patient to collect a urine sample



Scan this QR to watch the video on “How to take uACR test”



4. Take out the testcard from pouch






5. Dip the bottom part of testcard in urine completely for 1-2 sec



6. Place testcard on control pad and start timer immediately



7. Take a clear picture of the testcard



8. Results in 30 seconds

<b>Name</b>	Dummy	<b>Report ID</b>	1a2b3c4d2e31e2e2e
<b>Age</b>	25 years	<b>Sample type</b>	Urine
<b>Gender</b>	Male	<b>Sample collected on</b>	07/06/2024, 12:23:51 M
<b>Reference</b>	Dr. Swastik Poojari	<b>Report generated on</b>	07/06/2024, 12:24:05 M
<b>Hospital</b>			

## Test Report

### Urine Examination CHEMICAL EXAMINATION

Investigation	Observed Value	Units	Biological reference interval
<b>ACR</b>	<b>A3 (2100)</b>	<b>mg/g</b>	<b>Normal: &lt;= 30 mg/g Slightly Abnormal: 30-300 mg/g Abnormal: &gt;300 mg/g</b>
<b>Microalbumin</b>	<b>210.0</b>	<b>mg/L</b>	<b>&lt;20 mg/L</b>
<b>Protein (Albumin)</b>	<b>300 (++)</b>	<b>mg/dL</b>	<b>&lt;=30 mg/dL</b>
Glucose	Negative	mg/dL	Negative
Creatinine	10	mg/dL	10-300 mg/dL
Specific gravity	1.01	-	1.005-1.025
Ketone	Negative	mg/dL	Negative
Nitrite	Negative	-	Negative
Leukocytes	Negative	Leu/ $\mu$ L	Negative
pH	6	-	5-8
Blood	0	Ery/ $\mu$ L	Negative

#### Note

- The urine routine is a screening test
- Pre-test conditions to be observed while testing- first void, mid-stream urine, collected in a clean, dry, sterile container is recommended for routine
- urine analysis, to avoid contamination with any discharge from the vaginal, and urethra
- During interpretation, points to be considered are Negative nitrite test does not exclude the presence of bacteria or urinary tract infections.
- Trace proteinuria can be seen with many physiological conditions like prolonged recumbency, exercise, high protein diet, etc
- False reactions for bile pigments, proteins, glucose, and nitrites can be caused by peroxidase-like activity by disinfectants, therapeutic dyes, ascorbic acid and certain drugs, etc.
- Physiological variations may affect the test results
- When trace results occur, it is recommended to retest using a fresh specimen from the same patient
- Ketones may occur in urine during fasting, pregnancy and frequent strenuous exercise
- Blood is often, but not invariably, found in the urine of menstruating females

Download the  
Neodocs app



2 of 2



#### Note

- The urine routine is a screening test
- Pre-test conditions to be observed while testing- first void, mid-stream urine, collected in a clean, dry, sterile

container is recommended for routine

- urine analysis, to avoid contamination with any discharge from the vaginal. and urethra
- During interpretation, points to be considered are Negative nitrite test does not exclude the presence of bacteria or urinary tract infections
- Trace proteinuria can be seen with many physiological conditions like prolonged recumbency, exercise, high protein diet, etc
- False reactions for bile pigments, proteins, glucose, and nitrites can be caused by peroxidase-like activity by disinfectants, therapeutic dyes, ascorbic acid and certain drugs, etc
- Physiological variations may affect the test results
- When trace results occur, it is recommended to retest using a fresh specimen from the same patient
- Ketones may occur in urine during fasting, pregnancy and frequent strenuous exercise
- Blood is often, but not invariably, found in the urine of menstruating females

**Note:** All these steps are available in app as well

[www.neodocs.in](http://www.neodocs.in)

## FAQ

- **Q: What should I do if I forget the organization password?**

A: If you don't have the password or face any difficulty, please contact us at: +91 9987339111.

- **Q: How long should I wait for the results?**

A: Results will be available in 30 seconds after performing the test following the instructions provided.

## Documents / Resources



[NeoDocs uACR Test App](#) [pdf] Instruction Manual  
uACR Test App, Test App, App

## References

- [NeoDocs](#)
- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.