



## 67-05 Redux Electrolyte Paste



# PARKER 67-05 Redux Electrolyte Paste Instruction Manual

[Home](#) » [Parker](#) » PARKER 67-05 Redux Electrolyte Paste Instruction Manual 

### Contents

- 1 [PARKER 67-05 Redux Electrolyte Paste](#)
- 2 [Product Information](#)
- 3 [Product Usage Instructions](#)
- 4 [Frequently Asked Questions](#)
- 5 [Redux® Creme](#)
- 6 [Redux® Gel](#)
- 7 [Redux® Paste](#)
- 8 [Documents / Resources](#)
  - 8.1 [References](#)



**PARKER 67-05 Redux Electrolyte Paste**



## Product Information

### Specifications:

- **Manufacturer:** Parker Laboratories, Inc.
- **Product Models:** REF 66-04, REF 65-04, REF 67-05
- **Recommended Use:** Redux Creme for routine ECG procedures, Redux Gel for routine ECG procedures and defibrillation, Redux Paste for short-term monitoring, ECG, and defibrillation
- **Storage:** Protect from direct sunlight
- **Disposal:** Follow federal, state, and local regulations for non-hazardous waste disposal

## Product Usage Instructions

### Redux Creme:

Redux Creme is recommended for routine ECG procedures. Follow these steps for usage:

1. Apply a thin layer of Redux Creme on the skin area where the ECG electrodes will be placed.
2. Ensure even coverage to promote good electrode contact.
3. Proceed with the ECG procedure as instructed by healthcare professionals.

### Redux Gel:

Redux Gel is recommended for routine ECG procedures and defibrillation. Follow these steps for usage:

1. Apply Redux Gel generously on the skin area where the electrodes will be placed.
2. Spread the gel evenly to ensure proper conduction during the procedure.
3. For defibrillation, follow standard protocols after applying the gel as directed.

**Redux Paste:**

Redux Paste is recommended for short-term monitoring, ECG, and defibrillation. Follow these steps for usage:

1. Apply a small amount of Redux Paste on the electrode area.
2. Spread the paste evenly to create a conductive surface.
3. Monitor the patient or proceed with the ECG/defibrillation procedure accordingly.

**Frequently Asked Questions****Q: How should I store the product?**

**A:** Protect the product from direct sunlight to maintain its quality and effectiveness.

**Q: What should I do with the product after use?**

**A:** Follow all federal, state, and local regulations for non-hazardous waste disposal when discarding the product.

**Redux® Creme****Description**

Redux Creme is a non-sterile, electrically conductive creme. It is bacteriostatic.

**Indications for Use**

Redux Creme is recommended for routine ECG procedures.

**Contraindications**

- Redux Creme is not indicated for defibrillation procedures.

**Precautions**

- For external use only on intact, unbroken skin.
- Discard product if package is damaged or evidence of contamination is present.
- To avoid contamination of product, the tube container tips should not come into contact with patient, staff or instruments.

**Redux® Gel****Description**

Redux Gel is a non-sterile, non-abrasive, water based, electrically conductive gel.

**Indications for Use**

Redux Gel is recommended for routine ECG procedures and defibrillation.

**Warnings**

- Do not use on patients with a known history of sensitivity to parabens.

**Precautions**

- For external use only on intact, unbroken skin.
- Discard product if package is damaged or evidence of contamination is present.
- To avoid contamination of product, the tube container tips should not come into contact with patient, staff or instruments.

## **Redux® Paste**

### **Description**

Redux Paste is a bacteriostatic, non-sterile, water based, abrasive, electrically conductive paste.

### **Indications for Use**

Redux Paste is recommended for short term monitoring, ECG, and defibrillation.

### **Precautions**

- For external use only on intact, unbroken skin.
- Discard product if package is damaged or evidence of contamination is present.
- To avoid contamination of product, the tube container tips should not come into contact with patient, staff or instruments.

## **Redux® Paste, Creme, Gel**

### **Directions for Use**

Refer to the Instructions For Use provided by the manufacturer of the electromedical device being used for information on how Redux should be used with their product. Redux can be removed using a moistened soft cloth. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

### **Handling**

No special handling required.

### **Storage**

Protect product from direct sunlight.

### **Disposal**

Follow all federal, state, and local regulations for non-hazardous waste disposal.




Catalogue number



Batch code



Manufacturer



Authorized representative in the European community, Representante



Non sterile



Caution



Use by date,



Protect from Sunlight,



Do not use if package damaged,



See Instructions for Use,



Medical Device



Authorized Representative for Switzerland

[parkerlabs.com](https://parkerlabs.com)



PARKER LABORATORIES, INC.

Parker Laboratories BV  
Bedrijvenpark Twente 305  
7602 KL Almelo  
The Netherlands

MDSS CH GmbH  
Laurenzenvorstadt 61  
5000 Aarau  
Switzerland

© 2022 Parker Laboratories, Inc.

**Documents / Resources**



[PARKER 67-05 Redux Electrolyte Paste](#) [pdf] Instruction Manual  
66-04, 65-04, 67-05, 67-05 Redux Electrolyte Paste, 67-05, Redux Electrolyte Paste,  
Electrolyte Paste, Paste

References

- [Home - ParkerLabs](#)
- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.